



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

CHEM.ARCH.7.4 TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

Test Name.	Sunquest Test Code	Units	Male		Female		Critical Values	AMR Analytical Measurement Range (Linearity)	CRR Clinical Reportable Range (Range with Dilutions) / Dilute with M (manual)
			Ranges		Ranges				
			Age	Normals	Age	Normals			
Acetaminophen	ACETO	ug/mL	0	10-30	0	10-30	> 150	3-377	3770 / Saline
Albumin	ALB	g/dL	15	3.5-5.0	15	3.5-5.0		0.4-8.0	8.0
			5D	3.8-5.4	5D	3.8-5.4			
			0	2.8-4.4	0	2.8-4.4			
Alkaline Phosphatase	ALKP	U/L	19	38-126	19	38-126		5-4555	9110 / Saline M 1:2
			16	65-260	16	50-130			
			14	130-525	14	70-230			
			12	200-495	12	105-420			
			10	135-530	10	130-560			
			7	175-420	7	175-420			
			4	150-380	4	150-380			
			5D	145-320	5D	145-320			
			0	110-300	0	110-300			



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Alanine Aminotransferase	ALT	U/L	19	11-58	19	9-52		6-4113	20565 / Saline
			16	10-40	16	6-35			
			14	10-45	14	6-30			
			12	10-55	12	10-30			
			10	10-35	10	10-30			
			7	10-35	7	10-35			
			4	10-25	4	10-25			
			1	6-45	1	6-45			
			0	6-50	0	6-50			
Alcohol	ETOH	mg/dL	21	0-10	21	0-10	> 300	10-500	1000 / ETOH NEG- CAL M 1:2
			0	0-10	0	0-10	> 10		
Ammonia	AMM	umol/L	0	11-72	0	11-72		5-600	1110 / Saline
Amylase	AMYL	U/L	5D	25-125	5D	25-125		3-6554	13108 / Saline
			0	5-65	0	5-65			
Aspartate Aminotransferase	AST	U/L	19	17-59	19	14-36		3-4202	8404 / Saline
			16	15-45	16	5-30			
			12	15-40	12	10-30			
			10	10-60	10	15-40			
			4	15-50	4	15-50			
			5D	20-60	5D	20-60			
			0	35-140	0	35-140			



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Beta Hydroxybutyrate	BETAH	mmol/L	0	<0.27	0	<0.27		0.02-4.50	18.00 / DI-H20 M 1:4
Bilirubin, Direct	DBIL	mg/dL	19	0.0-0.5	19	0.0-0.5		0.1-15.0	15.0
			15D	0.1-0.3	15D	0.1-0.3			
			0	0.3-0.7	0	0.3-0.7			
Bilirubin, Total	TBIL	mg/dL	2	0.1-1.2	2	0.1-1.2		0.1-25.0	125.0 / Saline
			16D	0.1-1.2	16D	0.1-1.2	>20.0		
			15D	0.1-12.0	15D	0.1-12.0	>20.0		
			3D	0.1-12.0	3D	0.1-12.0	>18.0		
			2D	0.1-10.0	2D	0.1-10.0	>15.0		
			1D	0.1-10.0	1D	0.1-10.0	>13.0		
			0	0.1-6.0	0	0.1-6.0	>12.0		
Bilirubin, Indirect Unconjugated	IBIL	mg/dL	16D	0.1-1.2	16D	0.1-1.2		0.1-25.0	
			3D	0.1-12.0	3D	0.1-12.0			
			1D	0.1-10.0	1D	0.1-10.0			
			0	0.1-6.0	0	0.1-6.0			
B Natriuretic Peptide	BNP	pg/ml	0	0-100	0	0-100		10-5000	5000



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Blood Urea Nitrogen	BUN	mg/dL	19	8-26	19	10-20		2-125	625 / Saline
			10	7-21	10	7-19			
			1	9-22	1	9-22			
			15D	3-17	15D	3-17			
			0	3-23	0	3-23			
Calcium	CA	mg/dL	19	8.4-10.5	19	8.4-10.5	<6.0 >13.0	2.0-20.0	20.0
			12	8.8-10.8	12	8.8-10.8			
			11D	9.0-11.0	11D	9.0-11.0			
			0	7.6-10.4	0	7.6-10.4			
Carbon Dioxide	CO2	mmol/L	19	20-29	19	20-29		5-50	100 / Saline
			7D	20-28	7D	20-28			
			0	13-27	0	13-27			
Cholesterol	CHOL	mg/dL	20	125-200	20	125-200		7-705	2820 / Saline
			0	125-170	0	125-170			
			Desirable		<200 mg/dL				
			Borderline		200-239 mg/dL				
Undesirable		>=240 mg/dL							
CHOL/HDL RATIO	CHR		0	<4.98	0	<4.45			
Chloride	CL	mmol/L	1M	98-110	0	98-110		50-150	150
			0	98-113	0	98-113			



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Chloride Urine Random	CLU	mmol/L						20-300	300
Cortisol Random	CORT	ug/dl	Therapeutic Range: AM DRAW: 3.7-19.4 ug/dl PM DRAW: 2.9-17.3 ug/dl					1.0-59.8	119.6 / Cortisol CAL-A M 1:2
Creatine Kinase	CPK	U/L	19	20-259	19	20-225		7-4267	42670 / Saline
			1	60-305	1	60-305			
CKMB MASS	CKMBM	ng/ml	0	0.1-6.6	0	0.1-6.6		0.1-300.0	300.0
Creatinine	CREAT	mg/dL	19	0.80-1.50	19	0.70-1.20		0.20-30.00	30.00
			17	0.69-1.10	17	0.60-0.88			
			15	0.65-1.04	15	0.59-0.86			
			12	0.57-0.80	12	0.57-0.80			
			7	0.52-0.69	7	0.52-0.69			
			4	0.44-0.65	4	0.44-0.65			
			1	0.39-0.55	1	0.39-0.55			
			15D	0.31-0.53	15D	0.31-0.53			
			0	0.42-1.05	0	0.42-1.05			
Creatinine Urine Random	UCRR	mg/dL						5.00-400.00	8000.00 / Saline



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Creatinine Urine Timed	UCR	g/24h						5.00-400.00	8000.00 / Saline
	CRUT								
	CRCK		0	0.60-2.50	0	0.60-2.50			
Creatinine Clearance	CRCL	mL/min/1.73 m2							
	CLR		0	66-163	0	66-165			
	CLR2		0	85-125	0	75-115			
Cardio CRP	CACRP (CCRP)	mg/L	<p>FOR AGES > 17 YEARS:</p> <p>RISK ACCORDING TO AHA/CDC CCRP mg/L GUIDELINES</p> <p>-----</p> <p><1.0 Low Cardiovascular Risk 1.0-3.0 Average Cardiovascular Risk 3.1-10.0 High Cardiovascular Risk >10.0 Persistent elevations may represent Non-Cardiovascular inflammation.</p> <p>-----</p>					0.1-160.0	1600.0 / Saline
C-Reactive Protein	CRP	mg/dL	15	0.1-0.8	15	0.1-0.8		0.1-32.0	160.0 / Saline
			11	0.1-0.8	11	0.1-0.8			
			4	0.1-0.8	4	0.1-1.0			
			13M	0.1-1.2	13M	0.1-0.8			
			91D	0.1-1.1	91D	0.1-0.8			
			0	0.1-1.6	0	0.1-1.6			
CSF Glucose	CSFGL	mg/dL	0	40-70	0	40-70	< 21	1-800	800



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
CSF Total Protein	CSFTP	mg/dL	1M	12-60	1M	12-60		7-200	400 / Saline
			0	20-80	0	20-80			
Cyclosporine PRE	CYCLB	ng/mL	<p>No definitive therapeutic or toxic ranges have been established.</p> <p>Optimal blood drug levels are influenced by type of transplant, patient response, time post-transplant, co-administration of other drugs, and drug formulation.</p> <p>The following trough ranges are suggested guidelines:</p> <p>Kidney Transplantation: 100-200ng/mL Other Organ Transplantation: 200-300ng/mL</p>					25-1500	3000 / CAL-A M 1:2
Cyclosporine POST	CYCLP	ng/mL	<p>No definitive therapeutic or toxic ranges have been established.</p> <p>Optimal blood drug levels are influenced by type of transplant, patient response, time post-transplant, co-administration of other drugs, and drug formulation.</p> <p>The following peak ranges are suggested guidelines:</p> <p>de novo: 1500-1700 ng/mL 1st 6 months: 800-1000 ng/mL 6-12 months: 600-800 ng/mL >1 year: 300-600 ng/mL</p>					25-1500	3000 / CAL-A M 1:2
Digoxin	DIG	ng/mL	0	0.8 - 2.0	0	0.8- 2.0	> 2.5	0.3-4.0	4.0



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR	
Estradiol	ESTD	pg/mL	0	0.0-39.8				10.0-1000.0	10000.0 / ARCH Estradiol Manual Diluent 1:10	
			ESTRADIOL REFERENCE RANGES IN pg/mL <hr/> MALES 0-39.8 MENSTRUATING FEMALES (BY DAYS FROM LH PEAK) -FOLLICULAR PHASE (-12 TO -4) 19.5-144.2 -MIDCYCLE (-3 TO +2) 63.9-356.7 -LUTEAL PHASE (+4 TO +12) 55.8-214.2 POSTMENOPAUSAL FEMALES 0-32.2							
Iron	FE	ug/dL	0	65-175	0	50-170		5-1000	1000	
Follicle Stimulating Hormone	FSH	mIU/mL	FEMALE: Normally Menstruating: Follicular Phase 3.0- 8.1 Mid-cycle Phase 2.5- 16.7 Luteal Phase 1.4- 5.5 Post-menopausal 26.7-133.4 MALE: 0.9-11.9						0.1-70.0	350.0 / ARCH Multi Assay Manual Diluent 1:5
Free T4	FT4	ng/dL	19	0.7-1.5	19	0.7-1.5		0.4-6.0	6.0	
			1	0.9-1.4	1	0.9-1.4				



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR		
GFR, GFRAA	GFR, GFRAA Interp and ram messages	mL/min/ 1.73 m2	GFRG GFRLT GFRUA GFRSX GFR ≥90 GFR 60-89 GFR 45-59 GFR 30-44 GFR 15-29 GFR,15	<p>Not calculated due to patient age <18.</p> <p>Not calculated due to creatinine <0.1.</p> <p>Not calculated due to unknown age.</p> <p>Not calculated due to unknown sex.</p> <p>An estimated GFR chronically in the range of >= 90 is categorized as normal or high, which corresponds to Stage G1 CKD. CKD-EPI equation used to estimate GFR.</p> <p>An estimated GFR chronically in the range of 60-89 is categorized as mildly decreased, which corresponds to Stage G2 CKD. CKD-EPI equation used to estimate GFR.</p> <p>An estimated GFR chronically in the range of 45-59 is categorized as mildly to moderately decreased, which corresponds to Stage G3a CKD. CKD-EPI equation used to estimate GFR.</p> <p>An estimated GFR chronically in the range of 30-44 is categorized as moderately to severely decreased, which corresponds to Stage G3b CKD. CKD-EPI equation used to estimate GFR.</p> <p>An estimated GFR chronically in the range of 15-29 is categorized as severely decreased, which corresponds to Stage G4 CKD. CKD-EPI equation used to estimate GFR.</p> <p>An estimated GFR chronically in the range of <15 is categorized as kidney failure, which corresponds to Stage 5 CKD. CKD-EPI equation used to estimate GFR.</p>							
Gamma GT	GGT	U/L	0	12-64	0	9-36		4-7000	7000		



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Glucose	GLU	mg/dL	29D	65-99	29D	65-99	<40>500	5-800	4000 / Saline
			0	65-99	0	65-99	<40>150		
			If result of random glucose > or = 200 or if result of fasting glucose is > 125 confirm Diabetes Mellitus diagnosis with second glucose on a different day.						
Glucose, Fasting	FBS	mg/dL	29D	65-99	29D	65-99	<40>500	5-800	4000 / Saline
			0	65-99	0	65-99	<40>150		
			If result of random glucose > or = 200 or if result of fasting glucose is > 125 confirm Diabetes Mellitus diagnosis with second glucose on a different day.						
Glucose 1H Gestational Screen	GLU1H	mg/dL			0	65-140	<40 >500	5-800	4000 / Saline
			If the glucose is >140, A three hour glucose tolerance gestational confirmation is recommended.						
Glucose Tolerance, Gest 2H	GTT2G	mg/dL					<40 >500	5-800	4000 / Saline
	GT0G			0	65-99				
	GT1G			0	65-180				
	GT2G			0	65-155				



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Glucose Tolerance Gest 3 H	GTT3G	mg/dL					<40 >500	5-800	4000 / Saline
	GT0G				0	65-99			
	GT1G				0	65-180			
	GT2G				0	65-155			
	GT3G				0	65-140			
	<p>Gestational Diabetes Mellitus if two or more of the following are true:</p> <p style="text-align: center;">FASTING: > 99 1 HOUR: > 180 2 HOUR: > 155 3 HOUR: > 140</p>								
Glucose 1HR post Pran	GLU1P	mg/dL	0		0		<40 >500	5-800	4000 / Saline
			Diabetes is diagnosed at blood glucose of greater than or equal to 200 mg/dl						
Glucose 2 Hour PP	GL2PP	mg/dL	0	65-139	0	65-139	<40 >500	5-800	4000 / Saline



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
HCG QUANT	HCGQ + XHCGQ	mIU/mL						2-15000	225000 / ARCH Multi Assay Manual Diluent
				Pregnancy ranges: 2-4 weeks 391-8,388 5-6 weeks 861-88,769 6-8 weeks 8,636-218,085 8-10 weeks 18,700-244,467 10-12 weeks 23,143-181,899 13-27 weeks 6,303-97,171 28-42 weeks 4,360-74,883 Nonpregnant or premenopausal range: <5 mIU/mL Postmenopausal range: <10 mIU/mL Reference Range for Males: <5 mIU/mL					
HGB A1C with EAG	A1C (HA1C)	%	0	<5.7	0	<5.7		4.0-14.0	14.0
				<5.7% Decreased risk of diabetes 5.7-6.4% Increased risk of diabetes > OR = 6.5% Consistent with diabetes These Reference Intervals are supported by the current "Standards of Medical Care in Diabetes" published in January 2013 in the Diabetes Care, the Journal of the American Diabetes Association.					



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
High D lipoprotein	HDL	mg/dL	21	>40	21	>46		5-139	139
			15	25-76	15	34-75			
			11	28-72	11	28-72			
			6	32-74	6	32-74			
			1	25-73	1	25-73			
			0	14-63	0	19-69			
Lactate	LAC	mmol/L	0	0.5-2.2	0	0.5-2.2	> 3.4	0.1-13.3	13.3
Lactate Dehydrogenase	LDH	U/L	19	125-243	19	125-243		30-1995	3325 / Saline
			15	130-250	15	130-250			
			10	170-283	10	157-272			
			1	192-321	1	192-321			
			15D	163-452	15D	163-452			
			0	309-1222	0	309-1222			
LDL (CALC)	LDL	mg/dL	Desirable <130 mg/dL Borderline 130-159 mg/dL Undesirable >=160 mg/dL Desirable Range <100 mg/dL for patients with diabetes or CHD <70 mg/dL for Diabetic patients with known heart disease						



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Luteinizing Hormone	LH	mIU/mL						0.1-250.0	500.0 / ARCH Multi Assay Manual Diluent 1:2
			AGE: >= 18 YEARS LH REFERENCE VALUES ===== FEMALE: Follicular Phase 1.9- 12.5 mIU/mL Mid-Cycle Peak 8.7- 76.3 mIU/mL Luteal Phase 0.5- 16.9 mIU/mL Post-Menopausal 5.0- 52.3 mIU/mL MALE: 1.5- 9.3 mIU/mL AGE: < 18 YEARS NO REFERENCE RANGE ESTABLISHED						
Lipase	LIPA	U/L	0	8-78	0	8-78		4-1200	12000 / Saline M 1:10
Lithium	LI	mmol/L	0	0.6 - 1.2	0	0.6 - 1.2	>1.5	0.1-3.5	14.0 / Saline
Magnesium	MG	mg/dL	19	1.6-2.6	19	1.6-2.6	<1.0 >7.0	0.7-7.5	15.0 / Saline M 1:2
			12	1.7-2.2	12	1.7-2.2			
			6	1.7-2.1	6	1.7-2.1			
			6M	1.7-2.3	6M	1.7-2.3			
			0	1.5-2.2	0	1.5-2.2			



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
NON HDL CHOL (CALC)	NHDLC		0	<160	0	<160			
Phenobarbital	PHENO	ug/mL	0	15.0-40.0	0	15.0-40.0	>60.0	1.1-80.0	80.0
Phosphorus	PHOS	mg/dL	19	2.5-4.7	19	2.5-4.7	<1.2 >8.9	0.7-20.0	20.0
			16	2.8-4.6	16	2.8-4.6			
			14	2.9-5.4	14	2.9-5.4			
			12	3.3-5.4	12	3.3-5.4			
			7	3.7-5.6	7	3.7-5.6			
			4	4.0-5.4	4	4.0-5.4			
			1	3.9-6.5	1	3.9-6.5			
			0	4.6-8.0	0	4.6-8.0			
Phenytoin/Dilantin	PTN	ug/mL	0	10.0-20.0	0	10.0-20.0	>40.0	0.5-40.0	40.0
Potassium -serum	K	mmol/L	29D	3.5-5.1	29D	3.5-5.1	<2.8>6.0	1.0-10.0	10.0
			5D	4.1-5.3	5D	4.1-5.3	<2.6>6.0		
			0	3.7-5.9	0	3.7-5.9	<2.6>6.0		
Potassium Urine Random	UKR	mmol/L						1.0-200.0	200.0
Prealbumin	PREAB	mg/dL	61	16-42	61	14-37		3-60	60
			13	18-45	13	16-38			
			2	11-34	2	20-30			
			0	7-25	0	8-25			



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Progesterone	PROG	ng/mL	PROGESTERONE REFERENCE RANGES IN ng/mL <hr/> MALE: <0.10 - 0.20 FEMALE: NORMAL MENSTRUATING FOLLICULAR <0.10 - 0.30 LUTEAL 1.20 - 15.90 POST MENOPAUSAL <0.10 - 0.20 PREGNANCY 1ST TRIMESTER 2.80 - 147.30 2ND TRIMESTER 22.50 - 95.30 3RD TRIMESTER 27.90 - 242.50					0.10-40.00	800.0 / PROG Manual Diluent 1:20
PROGRAF	PROGR	ng/mL	Normal Ranges: <6 mos post-transplant: 8-12 ng/mL >6 mos post-transplant: 5-10 ng/mL					3-30	60 / CAL-A M 1:2
Prolactin	PROL	ng/mL	0	3.5-19.4	0	5.2-26.5		0.6-150.0	1500.0 / ARCH Multi Assay manual Dil 1:10



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Rapamycin	RAPA	ng/mL						2-30	60 / CAL-A M 1:2
			No therapeutic or toxic ranges have been established and there are no suggested values for elapsed time post-transplant.						
Salicylate	SALI	mg/dL	0	2.0-20.0	0	2.0-20.0	>30	5.0-100.0	500.0 / Saline
Sodium	NA	mmol/L	5D	136-145	5D	136-145	<120>160	100-200	200
			0	133-146	0	133-146			
Sodium Urine Random	UNAR	mmol/L						20-400	400
Tegretol	TEG	ug/mL	0	4.0-12.0	0	4.0-12.0	>20.0	2.0-15.0	28.0 / CAL-A M 1:2
Theophylline	THEO	ug/mL	0	10.0-20.0	0	10.0-20.0	>20.0	0.1-40.0	40.0
Tobramycin PRE	TOBPR	ug/mL	0	0.0-0.9	0	0.0-0.9	>2.0	0.2-10.0	20.0 / CAL-A M 1:2
Tobramycin POST	TOBPO	ug/mL	0	5.0-8.0	0	5.0-8.0	>10.0	0.2-10.0	20.0 / CAL-A M 1:2
Tobramycin RANDOM	TOBR	ug/mL					>10.0	0.2-10.0	20.0 / CAL-A M 1:2



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Total Protein	TPROT	g/dL	19	6.0-8.3	19	6.0-8.3		0.8-12.0	12.0 / Saline
			10	6.3-8.6	10	6.3-8.6			
			7	6.2-8.1	7	6.2-8.1			
			4	5.9-7.8	4	5.9-7.8			
			1	5.9-7.0	1	5.9-7.0			
			0	3.8-7.0	0	3.8-7.0			
Protein Urine Random	UTP	mg/dL	0	1.00-14.00	0	1.00-14.00		6.80-200.00	2000.00/ Saline
PRO 24 HR W/O CRT	URTP	mg/dL						6.80-200.0	2000.0 / Saline
	PRUT		0	1.00-14.00	0	1.00-14.00			
	TPVOL		0	<300.00	0	<300.00			
Triglyceride	TRIG	mg/dL	0	15-150	0	15-150		7-1000	4000 / Saline
			Normal <150 mg/dL Borderline-High 150-199 mg/dL High 200-499 mg/dL Very High >500 mg/dL						
Triglyceride PEDS	NTRIG	mg/dL	0	15-150	0	15-150		7-1000	4000 / Saline



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Troponin I	TROPI	ng/mL	0	<0.03	0	<0.03		0.01-50.0	50.0 / (Auto 1:9) ARCH Cal A Manual Diluent 1:20
			<0.03 Normal 0.03-0.29 Indeterminate: >/=0.30 Consider Serial Testing Consistent with Myocardial Injury						
TSH, 3RD GEN Thyroid Stimulating Hormone	TSH	uIU/mL	19	0.35-4.94	19	0.35-4.94		0.05-85.00	425.0 / ARCH Multi Assay Manual Diluent 1:5
			14	0.47-3.41	14	0.47-3.41			
			6M	0.70-4.17	6M	0.70-4.17			
			4D	0.73-4.77	4D	0.73-4.77			
Urea Nitrogen Urine Random	UUNR	mg/dL						40-1200	1200
Uric Acid	URIC	mg/dL	19	3.5-8.5	35	3.0-7.5	> 15.0	1.0-25.0	25.0
			14	3.0-7.0	14	3.0-6.0			
			12	2.4-6.5	12	2.4-6.5			
			2	2.2-5.9	2	2.2-5.9			
			0	2.0-6.4	0	2.0-6.4			
Valproic Acid	VALP	ug/mL	0	50-100	0	50-100	>150	2-150	150
			M	NR	F	NR	CR	AMR	CRR
Vancomycin PRE	VANCB	ug/mL	0	10.0-20.0	0	10.0-20.0	>20.0	3.0-50.0	100.0 / CAL-A M 1:2



TEST CODES, AMR'S CRR'S CRITICAL AND NORMAL RANGES

			M	NR	F	NR	CR	AMR	CRR
Vancomycin POST	VANCP	ug/mL	19	20.0-40.0	19	20.0-40.0		3.0-50.0	100.0 / CAL-A M 1:2
			0	20.0-40.0	0	20.0-40.0	>45.0		
Vancomycin Random	VANCR	ug/mL						3.0-50.0	100.0 / CAL-A M 1:2