

## Anti-Xa

Test	Reagent	Preparation	Stability	Drugs	Action
<b>LMWH</b>	<b>Liquid Anti-Xa Chromogenic Substrate</b>	Invert to mix before use <b>*INITIAL/DATE PREP*</b>	30 DAYS 2°- 8° <b>*Keep in fridge when not in use</b> 7 DAYS O.B.	enoxaparin (Lovenox)	Run LMWH: Analyzer test = "Anti-Xa", choose <b>LMWH QC</b>
	<b>Liquid Anti-Xa Factor Xa Reagent</b>	Invert to mix before use <b>*INITIAL/DATE PREP*</b>		dalteparin (Fragmin)	
	<b>Clean B dil.</b>	2mL Cleaning Agent (Clean B) & 14mL dH <sub>2</sub> O in 20mL vial	*Daily Maintenance	Test will be frozen and run next M, W, F - Day shift	
<b>QC</b>	<b>LMWH Low LMWH High</b>	Pipette 1mL dH <sub>2</sub> O, invert & swirl. Let sit 30 mins, invert & swirl, <b>*INITIAL/DATE PREP*</b>	48 HOURS 2°- 8° <b>*Keep in fridge when not in use</b> 24 HOURS O.B.	If patient is on: fondaparinux (Arixtra)	
<b>UNFH</b>	<b>Same Reagent as LMWH, Liquid Anti-Xa Kit</b>			Run UNFH: Analyzer test = "Anti-Xa", choose <b>UNFH QC</b>	
<b>QC</b>	<b>UNFH Low UNFH High</b>	Pipette 1mL dH <sub>2</sub> O, invert & swirl. Let sit 30 mins, invert & swirl, <b>*INITIAL/DATE PREP*</b>	48 HOURS 2°- 8° <b>*Keep in fridge when not in use</b> 24 HOURS O.B.		
<b>Procedure:</b>					
<p>1) Use E-record to confirm which heparin patient is currently on</p> <p>2) As needed, change test ordered to correct test, log SPROB if changing test or .ND</p> <p>3) Liquid Anti-Xa reagents &amp; QC that are in use are found in the "specials tray" #7 SR800, in the respective boxes in SR800, or in the walk-in bins</p> <p>4) Determine which test to run and choose the corresponding Low &amp; High QC. Use reconstituted QC vials found in the tray if available. If not, reconstitute according to chart</p> <p>5) Load reagents in correct racks on STOP2 or STOP3</p> <p>6) Manually order QC: QC→Results List or Test Status List→Double click on any test name→Click the 🗑️ icon to clear any current selections. Choose appropriate _L and _H QC corresponding to LMWH or UNFH tests and click 🖱️</p> <p>7) Enter QC data into Soft QC:</p> <ul style="list-style-type: none"> <li>•LMWH QC will cross to the Instrument Menu for that analyzer. Click the Upload QC icon 📄 and properly handle any failures.</li> <li>•UNFH QC will <b>not</b> cross to SoftQC from the interface. It must be entered manually: SoftQC→"N"→'eNter by Tray Tube' screen→Instr: STOP2 or STOP3, Tray Tube: UNFH H_"XX" or UNFH L_"YY" or [F2], [Enter] to select correct QC_Lot, [F12] to accept→enter result and save. It does not need to be printed.</li> </ul>					
<b>Result Reporting:</b>				If a result holds on the STOP Interface:	
<b>Result</b>		<b>Action</b>			
<b>"Failed"</b>		Check sample integrity, repeat & follow lab protocol			
Remember, any <b>"HE"</b> or <b>"SE"</b> errors are <b>unacceptable</b>					

## Anti-Thrombin III

<b>AT3</b>	<b>Liquid Antithrombin Chromogenic Substrate</b>	Invert to mix before use <b>*INITIAL/DATE PREP*</b>	30 DAYS 2°- 8° <b>*Keep in fridge when not in use</b> 7 DAYS O.B.
	<b>Liquid Antithrombin Factor Xa Reagent</b>	Invert to mix before use <b>*INITIAL/DATE PREP*</b>	
<b>QC</b>	<b>SP TST2</b>	Pipette 1mL dH <sub>2</sub> O, invert & swirl. Let sit 30 mins, invert & swirl.	<b>8 HOURS O.B.</b>
	<b>NRM ASY</b>	<b>*INITIAL/DATE &amp; TIME PREP*</b>	24 HOURS O.B.