



TANGO infinity® Tips

USSD • February 2016



1

Table of Contents

Bromelin Preparation.....	4	Control Management.....	15
Steps to Getting Started.....	5	Flags Associated with Well/Sample Results.....	16
Powering the TANGO infinity® On/Off.....	6	Reagent Usage.....	18
Crossmatch Assay Procedure Notes.....	7	Required Reagents.....	21
Daily Care.....	8	Additional Supplies.....	25
Weekly Maintenance.....	9	Product Codes for Lot Reports.....	26
Monthly Maintenance.....	10		
Erytype Reaction Interpretation.....	11		
Solidscreen II Reaction Interpretation.....	12		
QC Requirements.....	13		
QC Sample Requirements.....	14		

2

TANGO infinity® Tips

3

Bromelin Preparation

Total Volume	Bromelin	Isoionic Saline	% Volume
65 mL	2.75 mL	62.25 mL	100
60 mL	2.50 mL	47.50 mL	90
45 mL	2.25 mL	42.75 mL	80
40 mL	2.00 mL	38.00 mL	70
35 mL	1.75 mL	33.25 mL	60
30 mL	1.50 mL	28.50 mL	50
25 mL	1.25 mL	23.75 mL	40
20 mL	1.00 mL	19.00 mL	30
15 mL	0.75 mL	14.25 mL	20
10 mL	0.50 mL	9.50 mL	10

4

Steps to Getting Started

1. Perform required Care or Maintenance.
Daily / Weekly / Monthly
2. Run and validate QC.
3. Load and run samples.

5

Powering the Tango infinity® On/Off

Power On Sequence

1. Power on the TANGO infinity®
2. Power on the PC
3. Enter password and the instrument will continue to initialization

Power Off Sequence

1. Exit TANGO infinity® Software
2. Shut down the computer if it does not automatically turn off
3. Turn off the TANGO infinity®

Note: Wait 30 seconds before powering the system on.

6

Crossmatch Assay Procedure Notes

Use of the Y rack*

- Not for use when TANGO **infinity**™ is interfaced to an LIS
- Recipient sample in position 1, position 2 is empty, donor red blood cells are placed into positions 3 to 12

Use of the N and K racks*

- Used with or without interface to an LIS
- N rack for patient samples, as with other assays
- K rack with spacers for donor red blood cells from segment
- Use TANGO **infinity**™ software to select which donors will be crossmatched to the selected recipient

*Crossmatch Assay not currently available in the USA.

7

Daily Care

• Clean, check and adjust SPOLV

Wipe pipettor needle with DI Water. Perform system rinse. Start automatic check and adjustment of SPOLV.

• Check Suspension Cup

Look for any residue in cup or red cells attached to the side of the cup, without removing it if possible. If large amount is present, replace the cup.

• Check the liquid supplies and liquid waste

Refill the System Liquid or Wash Solution bottle if volume is low. Empty the liquid waste if necessary.

• Check the strip waste

Empty or replace the small bag that catches the used strips. Alternately wash the strip waste container if a bag is not used.

8

Weekly Maintenance

• **Clean and disinfect surfaces**

Open front cover and remove all racks and acrylic glass cover from cooled reagent bay. Wipe splatter and spills with water or water based detergent. Disinfect by wiping with bleach or ethanol.

• **Clean Optical Filter**

Remove optical filter from front side of the Measurement Chamber and clean with canned air, lint free cloth or lens paper.

• **Check Pipettor Rinse Station**

Dry with a lint free wipe.

9

Monthly Maintenance

• **Disinfect Sample and Reagent Racks**

Remove all Sample and Reagent Racks. Clean racks by wiping with water or water based detergent. Disinfect racks with bleach or ethanol.

• **Replace Suspension Cup and Splash Bowl**

Remove the suspension cup and splash bowl. Properly seat new splash bowl and suspension cup. Be certain bowl and cup snap into place. Clean the gray plastic cover.

• **Disinfect Containers and**

Follow step by step directions located **BIO-RAD** BIO-RAD LABORATORIES
Microcide is used for disinfection. To c
750 mL DI water to make 1 L of disinfecting solution.

• **Disinfect Strip Waste Bin**

If not using a waste bag in the strip waste bin, or if container becomes dirty, clean the container with water. Then disinfect the container with bleach or ethanol. Rinse with water and wipe dry.

10

Erytype Reaction Interpretation



Negative

An even RBC suspension. No agglutinates detected.



+/-

Weak granularity in the RBC suspension. It is noted that there is a +/- on TRUQC with a 7 for the overall interpretation.



+

Many small agglutinates with a lot of free cells.



**

Large agglutinates in a sea of smaller clumps. No free RBCs.



Strong reaction. A number of large agglutinates.



A single agglutinate. No free RBCs detected.

11

Solidscreen II Reaction Interpretation



Negative

Compact cell button.



+/-

Cell button forming in the center of the well, some cells adhering to well bottom.



+

Some cells adhering to bottom of well, distinct settling of approximately 50% of the cells in the center of the well.



**

A layer of cells adhering to the bottom of the well, distinct settling of unattached cells toward center of well.



Homogenous layer of cells across the bottom, with some settling of unattached cells towards center of well.



Homogenous layer of cells across the bottom of the well.

12

QC Requirements

QC samples are tested:

- Every 24 hours or per facility SOP
- When a new lot number is used on the TANGO Infinity*
- When a new preparation of Bromelin is placed on the TANGO Infinity*
- After service or repair is performed on the TANGO Infinity*
- After Monthly Maintenance

13

QC Sample Requirements

ABO and Rh

- Hemo-QC commercial controls (Bio-Rad PN 816500100) or samples selected from previously tested samples up to 3 days old. (Samples up to 7 days may be used if necessary)
- Samples should be selected to provide a positive and negative reaction with every reagent.
For example:
AB Positive and O Negative
A Positive and B Positive and O Negative

Antibody Screen and Identification

- Hemo-QC Control Kit PN 816500100 may be used
 - WBCTL-01: A2B RhD+ Screen Negative
 - WBCTL-02: O RhD+ Screen Positive Anti-Fya
 - WBCTL-03: A1 RhD - Screen Positive Anti-D

Solidscreen® II Control is an anti-D for use as a control

Solidscreen II Control B is an anti-c for use as a control

Solidscreen II Negative control is available as a negative control

Previously tested patient samples may also be used as controls

Ensure target values for positive antibody controls match current lot of screening cells

14

Control Management

Loading & Running Controls

- Place properly barcoded control samples in "C" Rack. (Sample ID of a control begins with a lower case c). The first 3 digits of the sample barcode must match the designation in the Control Management Software of the TANGO.
- Ensure aliquot samples are placed into slots with spacers and full size tubes are placed into slots without spacers.
- Load the "C" rack into the sample area.
- Select the controls to be tested and press "Request Manually"
- Press the "Start" button on Main Menu.

15

Flags Associated with Well/Sample Results

- 6 Sample: Clot detected
- 7 Reagent Bay (cooled): Temperature exceeded
- 10 Suspension: No liquid detected
- 11 Suspension: Maximum immersion depth reached during aspiration
- 14 Washer: Insufficient dispense during washing process
- 15 Washer: Insufficient aspiration during washing process
- 16 TANGO **infinity**® paused for more than one minute
- 17 Sample rack removed too early
- 19 Incubator temperature out of range
- 21 Reagent: Not pipetted
- 22 Reagent: Maximum dive-in depth reached during aspiration
- 23 Reagent: No liquid
- 26 Sample: No Plasma / Serum
- 27 Sample: Not enough Plasma / Serum

Continued on the next page >

16

Flags Associated with Well/Sample Results

29 OLV: Wrong volume aspirated
 X Control Failed



Current result discrepant with previous result



Result edited manually

- () Empty well for Erytype assays
- [] Empty well for Solidscreen II assays
- ? Reaction strength can not be determined
- X Invalid results in compatibility testing

17

Reagent Usage

Reagent	Test	Vol. in Vial	Tests/Vial**	Vol./Test
Bromelin for Erytype*	ABO	50 mL	50 to 54	900 µL
Erytype cell A ₁ Cells	ABO	10 mL	100	50 µL
Erytype cell B Cells	ABO	10 mL	100	50 µL
MEB 2	ABS	50 mL	577	84 µL
AHG Anti-IgG Solidscreen I	ABS	55 mL	556	100 µL
Biobestcell Pool	ABS	10 mL	550	17 µL
Biobestcell 1 & 2, vial 1	ABS	10 mL	550	17 µL
Biobestcell 1 & 2, vial 2	ABS	10 mL	550	17 µL
Biobestcell 3, vial 1	ABS	10 mL	550	17 µL
Biobestcell 3, vial 2	ABS	10 mL	550	17 µL
Biobestcell 3, vial 3	ABS	10 mL	550	17 µL

Sample volume used for testing:

- 10 µL RBC for ABO/Rh
- 100 µL plasma/serum for ABO/Rh reverse (50 µL per well)
- 50 – 150 µL for antibody screen (50 µL per well)

* Diluted for use, as directed

** Accounts for Dead Volumes (1.5 mL in 50 – 55 mL bottle/0.5 mL in 10 mL bottle)

Continued on the next page >

18

Reagent Usage

Reagent	Test	Vol. in Vial	Tests/Vial**	Vol./Test
Soldscreen II Anti-D Blend	Weak D	4 mL	70	50 µL
MLB 2	Weak D	80 mL	80	500 µL
Altersers	Weak D	50 mL	258	198 µL
Anti-G, Anti-IgG Soldscreen II	Weak D	55 mL	535	100 µL per well
MLB 2	DAT	50 mL	50	500 µL
Altersers	DAT	50 mL	258	198 µL
Anti-G, Anti-IgG Soldscreen II	DAT	55 mL	535	100 µL per well
MLB 2	Crossmatch*Auto Control	80 mL	80	500 µL
Altersers	Crossmatch*Auto Control	50 mL	258	198 µL
Anti-G, Anti-IgG Soldscreen II	Crossmatch*Auto Control	55 mL	535	100 µL per well

Sample volume used for testing:

Weak D 50 µL of an approximate 1% cell suspension

DAT 50 µL of an approximate 1% cell suspension

Crossmatch* 50 µL of serum/plasma plus 50 µL of approximate 1% donor cell suspension

* Crossmatch is not FDA-approved

** Accounts for Dead Volume (1.5 mL in 50-55 mL bottle/0.5 mL for 5 mL bottle)

Continued on the next page >

19

Reagent Usage

Reagent	Test	Vol. in Vial	Tests/Vial**	Vol./Test
Biotestcell-B	0 Cell Antibody Identification	4 mL	70 panels	50 µL
MLB 2	8 Cell Antibody Identification	50 mL	670 wells, 72 panels	84 µL
Anti-G, Anti-IgG Soldscreen II	0 Cell Antibody Identification	55 mL	535 wells, 66 panels	100 µL
Biotestcell-111 and -111 Plus	11 Cell Antibody Identification	4 mL	70 panels	50 µL
MLB 2	11 Cell Antibody Identification	50 mL	672 wells, 52 panels	84 µL
Anti-G, Anti-IgG Soldscreen II	11 Cell Antibody Identification	55 mL	535 wells, 48 panels	100 µL

Sample volume used for testing:

Antibody Identification 50 µL per well

Rh+K Phenotype - 50 µL of an approximate 1% cell suspension

** Accounts for Dead Volumes (1.5 mL in 50-55 mL bottle/1 mL in 10 mL bottle)

20

Required Reagents

ABO/Rh		Weak D Testing	
Product Code	Product Name	Product Code	Product Name
806 127 100	Erytype S ABO+Rev A,B	806 521 100	Soldscreen II
818 956 100	Erytypecell A, & B Reagent Red Blood Cells (A ₁ & B cells)	805 500 100	MLB 2
806 210 100	Bromalin for Erytype	806 510 100	Wetvers Solution
		806 516 100	AHG, Anti-IgG Soldscreen II
		806 530 100	Soldscreen I Anti-D Bland

ABO Group		Rh Type	
Product Code	Product Name	Product Code	Product Name
806 130 100	Erytype S ABO Donor	806 190 100	Erytype S Rh Donor
806 210 100	Bromalin for Erytype	806 210 100	Bromalin for Erytype

Continued on the next page >

21

Required Reagents

Antibody Screen testing			
Product Code	Product Name	Product Code	Product Name
805 521 100	Soldscreen II	805 065 100	Biotechcell Pool Reagent Red Blood Cells
805 200 100	MLB 2	816 014 100	Biotechcell 1 & 2 Reagent Red Blood Cells
805 516 100	AHG, Anti-IgG Soldscreen II	816 086 100	Biotechcell 3 Reagent Red Blood Cells

Continued on the next page >

22

Required Reagents

Direct Antiglobulin Testing (DAT)		Rh+K Phenotyping	
Product Code	Product Name	Product Code	Product Name
805 521 100	Solidscreen II	806 105 100	Erytype SPH+K Type
805 200 100	MLB 2	806 210 100	Bromalin for Erytype
806 510 100	Alaways Solution		
806 516 100	AHG, Anti-IgG Solidscreen II		

Continued on the next page >

23

Required Reagents

Antibody Identification		Crossmatch Testing/Auto Control	
Product Code	Product Name	Product Code	Product Name
805 521 100	Solidscreen II	805 521 100	Solidscreen II
805 200 100	MLB 2	805 200 100	MLB 2
805 516 100	AHG, Anti-IgG Solidscreen II	805 510 100	Alaways Solution
816 020 100	Biotestcells -I6	805 166 100	AHG, Anti-IgG Solidscreen II
816 021 100	Biotestcells -I11		
816 022 100	Biotestcells -I11 Plus		

24

Additional Supplies

Product Code	Product Name	Package Size
848 000 091	Bio-Rad PBS Concentrate	6 x 1L liquid concentrate
848 000 010	Cellmixers	20 mixers
25099	Microcide-SQ™ Decontamination solution concentrate	250 mL
848 999 790	Cup and Bowl, Suspension	1 bowl, 1 cup
850 100 020	Cup and Bowl, Suspension	6 bowls, 6 cups
850 101 010	QC Barcode Set	14 sheets, 70 labels
850 800 090	Patent Sample Rack	1 Rack
850 800 100	Sample Rack Spacers	12 spacers per pack
848 000 190	50 mL bottles (empty for Bromelin)	2 Bottles
816 600 100	Hemo-QC	Set of 8 6 mL Vials for TANGO or manual reagent QC. Recommendation of 1 vial per week for TANGO.
806 514 100	Soldscreen II Control	4 mL positive control, Anti-D
806 519 100	Soldscreen Control B	4 mL positive control, Anti-c
806 509 100	Soldscreen Negative Control	4 mL negative control

25

Product Codes for Lot Reports

Pool cells	131	Bromelin	200
Cell 1 of 2 screen	212	USS-MLB2	201
Cell 2 of 2 screen	213	Cocoma	204
Cell 1 of 3 screen	214	ABO-Flav plates	013
Cell 2 of 3 screen	215	Soldscreen plates	060
Cell 3 of 3 screen	216	ABO donor plates	015
Erytype A1 cell	209	Flu donor plates	016
Erytype B cell	209	Soldscreen control	603
Alzovers	222	Soldscreen B ctrl	219
Anti-D blend	208	Soldscreen neg ctrl	220

26

Refer to the current version of the TANGO **infinity**[®] User's Guide for complete information.
