

 		
<b>University of Washington Medical Center</b> <b>1959 NE Pacific Street, Seattle, WA 98195</b> <b>Transfusion Services Laboratory</b> <b>Policies and Procedures Manual</b>	<b>Original Effective Date:</b> <b>03-11-16</b> <b>Revision Effective Date:</b>	<b>Number:</b> <b>PC-0017.02</b>
<b>TITLE: Returning Issued Blood Components to Inventory</b>		

**PURPOSE:**

To describe the process for evaluating blood components returned from an issued status prior to acceptance back into inventory or placing into quarantined status

**PRINCIPLE & CLINICAL SIGNIFICANCE:**

**Principle**

Blood components issued for transfusion and returned must be inspected to ensure all applicable Food and Drug Administration (FDA) and AABB requirements are met prior to acceptance back into inventory

**Clinical Significance**

Blood components issued for transfusions are susceptible to damage and increased risk for bacterial contamination if proper temperatures and handling practices are not maintained. Transfusion of components not maintained at the appropriate temperature may lead to reduced effectiveness of the product and can increase the risk of an adverse reaction in the recipient if transfused

**POLICIES:**

Blood components may be returned and accepted into inventory for reissue -only if the following conditions are met:

- Component is not spiked
- Appropriate temperature was maintained
  - Red Blood Cells & Plasma – 1-6°C if returned from issue in a monitored blood refrigerator
  - Red Blood Cells & Plasma - 1-10°C if returned from transport (not stored in a monitored blood refrigerator)
  - Platelets, cryoprecipitate and granulocytes – 20-24°C
- Red cell blood components have one integrally attached segment
- Component passes visual inspection and inspection is recorded
- Products stored outside of TSL monitored environment (returned > 30 minutes from issue time, in or out of temperature range), must be approved by the TSL MD prior to acceptance back into inventory.

**SPECIMEN REQUIREMENTS:**

NA

**REAGENTS/SUPPLIES/EQUIPMENT:**

Reagents:	Supplies:	Equipment:
NA	Calibrated thermometer "Quarantine" stickers	LIS

<b>TITLE: Returning Issued Blood Components to Inventory</b>	<b>Number: PC-0017.02</b>
--	---------------------------

**QUALITY CONTROL:**

NA

**INSTRUCTIONS:**

**TABLE of CONTENTS:**

[Verifying Blood Component Acceptability for Return to Inventory](#)  
[Returning Blood Component in LIS](#)

**Verifying Blood Component Acceptability for Return to Inventory**

STEP	ACTION							
1	<b>If</b>	<b>Then</b>						
	RBC or Thawed Plasma maintained in a monitored blood refrigerator	<ul style="list-style-type: none"> <li>Verify refrigerator maintained within acceptable temperature range of 1-6°C during storage</li> </ul> <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Acceptable</td> <td> <ul style="list-style-type: none"> <li>Go to step 5</li> </ul> </td> </tr> <tr> <td>Not acceptable</td> <td> <ul style="list-style-type: none"> <li>Quarantine Component</li> <li>Go to next section Returning Blood Components in LIS</li> </ul> </td> </tr> </tbody> </table>	If	Then	Acceptable	<ul style="list-style-type: none"> <li>Go to step 5</li> </ul>	Not acceptable	<ul style="list-style-type: none"> <li>Quarantine Component</li> <li>Go to next section Returning Blood Components in LIS</li> </ul>
	If	Then						
	Acceptable	<ul style="list-style-type: none"> <li>Go to step 5</li> </ul>						
Not acceptable	<ul style="list-style-type: none"> <li>Quarantine Component</li> <li>Go to next section Returning Blood Components in LIS</li> </ul>							
Product out of monitored storage	<ul style="list-style-type: none"> <li>Go to next section</li> </ul>							
2	<b>If out of monitored storage</b>	<b>Then</b>						
	< 30 minutes	Go to next step						
	≥ 30 minutes	<ul style="list-style-type: none"> <li>Quarantine Component</li> <li>Go to next section Returning Blood Components in LIS</li> </ul>						
3	Use a NIST calibrated thermometer to verify the temperature by placing the thermometer between two components (if possible) or sandwich the single product and read temperature after 3-5 minutes							
	<b>For</b>	<b>Acceptable Temperature Range</b>						
	RBC or Thawed Plasma issued for transfusion	1-10 °C						
Platelets, cryoprecipitate or granulocytes	20-24 °C							
4	<b>If components is</b>	<b>Then</b>						
	Within acceptable temperature range	Go to next step						
	Outside of acceptable temperature range	<ul style="list-style-type: none"> <li>Quarantine Component</li> <li>Go to next section Returning Blood Components in LIS</li> </ul>						
5	Perform visual inspection (refer to SOP <i>Visual Inspection</i> )							
	<b>If visual inspection</b>	<b>Then</b>						
	Passes	<ul style="list-style-type: none"> <li>Component is acceptable for return</li> <li>Go to next section</li> </ul>						

	Fails	<ul style="list-style-type: none"> <li>• Quarantine Component</li> <li>• Go to next section Returning Blood Components in LIS</li> </ul>
--	-------	--

**Returning Blood Component in LIS**

STEP	ACTION	
<b>1</b>	<b>If returning in SQ</b>	<b>Then</b>
	Immediately	Open "Blood Status Update" function and retain the default "Update Option" <Unit Update>
	Later	Record the return time on the Transfusion Record and place the component in the correct storage area until it can be returned in SQ
<p><b>NOTE:</b> If Sunquest is down the return is documented on the Downtime Issue Log <u>and must be entered into Sunquest as soon as feasible.</u></p>		
<b>2</b>	Scan the unit number and component type at the "Unit Selection" prompt <u>Note: The unit number may be manually entered if unable to scan</u>	
<b>3</b>	Tab through the current date and time or update if necessary when entering information from the downtime issue log	
<b>4</b>	Select the appropriate "New Status" <u>from the drop down box</u> and then <Tab>	
	<b>If unit is</b>	<b>Then select</b>
	Acceptable for return to inventory	<ul style="list-style-type: none"> <li>• Inventory (INV) and <u>clickpress &lt;tab&gt;, &lt;tab&gt; to bypass the Temperature field (Do not enter temperature date)</u></li> <li>• Perform <u>a visual inspection and answer the "Pass visual inspection&gt;"</u> in SQ by clicking the "Yes" or "No" box</li> </ul>
	Not acceptable for return to inventory	<ul style="list-style-type: none"> <li>• <del>Quarantine</del></li> <li>• Enter the reason for quarantine into Sunquest</li> </ul> <p><b>EXAMPLES:</b> QAPEAP = product expired, OMS30 =out of monitored storage greater than 30 minutes</p> <p><b>NOTE:</b> If an appropriate reason cannot be located use code BBR – Blood Bank Reason and enter a free text comment to explain the circumstances of the return</p>
<p><b>NOTE:</b> <u>Do not enter temperature data in this field. The system does not have logic to alert the user if the temperature is out of range. If there are concerns regarding product transport conditions, refer to SOP: Quarantine and Final Disposition of Blood Components.</u></p>		
<b>5</b>	Click <9.Unit Location>	
<b>6</b>	<u>Update the location if necessary and</u> Click <OK> in the "Location Update"	
<b>7</b>	Click <Save>	

**TITLE: Returning Issued Blood Components to Inventory**

**Number:  
PC-0017.02**

Select the appropriate response at the "Reallocation of unit" window	
If transfusion of the unit is	Then
8 Still anticipated	<ul style="list-style-type: none"><li>• Select 'Allocated' from the drop down menu</li><li>• Save</li><li>• Request reprint of Transfusion Record from the TSL MLS (refer to SOP <i>Sunquest Blood Order Processing Navigation</i>)</li><li>• Retag unit</li></ul>
Has been cancelled	<ul style="list-style-type: none"><li>• Select 'Released' from the drop down menu</li><li>• Save</li></ul>
9	<ul style="list-style-type: none"><li>• Return the unit to the appropriate storage location</li></ul>

**CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES**

NA

**CALIBRATION:**

NA

**NOTES AND LIMITATIONS:**

- If it is not possible to return a unit to the UWMC TSL due to breakage, notify a TSL Lead and complete a QI [Form](#) to be used for [tracking the event and](#) final disposition purposes

**REFERENCES:**

- Technical Manual. [Current edition](#). Bethesda, MD: AABB,
- Standards for Blood Banks and Transfusion Services. Bethesda, MD; AABB, current edition
- Sunquest users Guide version [8.1](#)

**RELATED DOCUMENTS:**

SOP: Visual Inspection of Blood Components  
SOP: Quarantine and Final Disposition of Blood Components Attaching the Transfusion  
SOP: Sunquest: Blood Order Process Navigation  
Form: Downtime Issue Log

**ADDENDUM:**

NA

**UWMC SOP Approval:**

<b>TITLE: Returning Issued Blood Components to Inventory</b>	<b>Number: PC-0017.02</b>
--	---------------------------

<b>Chief of Clinical Services (CLIA Medical Director)</b>	_____	Date _____
	Mark H. Wener, MD	
<b>Transfusion Service Manager</b>	_____	Date _____
	Deanne Stephens	
<b>Transfusion Service Compliance Analyst</b>	_____	Date _____
	Christine Clark	
<b>Transfusion Service Medical Director</b>	_____	Date _____
	Monica B. Pagano, MD	
<b>UWMC Biennial Review:</b>		
	_____	Date _____
	_____	Date _____

<b>SCCA SOP Approval:</b>		
<b>SCCA CLIA Medical Director</b>	_____	Date _____
	Brent L. Wood, MD	
<b>Director, Transfusion Services</b>	_____	Date _____
	Terry Gernsheimer, MD	
<b>Alliance Lab Manager</b>	_____	Date _____

**REVISION HISTORY:**  
04/22/2018: Updated to include changes due to Sunquest 8.1 upgrade.

