



<b>University of Washington Medical Center</b> <b>1959 NE Pacific Street. Seattle, WA 98195</b> <b>Transfusion Services Laboratory</b> <b>Policies and Procedures Manual</b>	<b>Original Effective Date:</b> <b>03-11-16</b>	<b>Number:</b> <b>PC-0036.01</b>
	<b>Revision Effective Date:</b> <b>05-09-18</b>	
<b>TITLE: Testing and Provision of Hemoglobin S Negative Blood</b>		

**PURPOSE:**

To provide instructions for performing sickle cell testing on red blood cell components using the SickLeDex kit and controls and provide hemoglobin S negative components for neonates and patients with Sickle Cell disease/trait and thalassemia

**PRINCIPLE & CLINICAL SIGNIFICANCE:**

**Principle**

SickleDex is a qualitative solubility test kit used to detect the presence of sickling hemoglobins (HbS) in blood or control material. *SickleDex* uses saponin to lyse the red blood cells. Sodium Hydrosulfite then reduces the released hemoglobin. Reduced HbS is insoluble in the concentrated phosphate buffer and forms a cloudy, turbid suspension.

**Clinical Significance**

Transfusion of units negative for hemoglobin S reduces tissue hypoxia in patients with Sickle Cells Disease (SCD) and reduces the risk of cellular breakdown from the hemoglobin S positive RBCs in neonates

**POLICIES:**

- HbS negative red blood cells components should be issued to the following:
  - Patients diagnosed with sickle cell disease, sickle cell trait and thalassemia
  - Children < 4 months old
- SCD patients must be flagged in Sunquest with the NHBS attribute to prevent release of untested units. Refer to SOP: *Blood Administrative Data Entry Updates*

**SPECIMEN REQUIREMENTS:**

Blood from donor segments stored at 1-6°C for up to 45 days

**REAGENTS/SUPPLIES/EQUIPMENT:**

Reagents:	Supplies:	Equipment:
<ul style="list-style-type: none"> <li>• SickLeDex Kit</li> <li>• SickLe-Chex Controls</li> </ul>	<ul style="list-style-type: none"> <li>• 12 x 75 mm test tubes</li> <li>• Adjustable pipette</li> <li>• Disposable pipette tips</li> <li>• Blood bank transfer pipettes</li> </ul>	Test tube rack with lines

**QUALITY CONTROL:**

Quality control is performed with each batch of testing

**INSTRUCTIONS:**

**TABLE OF CONTENTS:**

- [Working Solubility Buffer Preparation](#)
- [SickleDex Testing](#)

**Working Solubility Buffer Preparation**

STEP	ACTION
1	Bring buffer and reagent powder from the SickleDex Kit to room temperature before mixing
2	Add contents of one vial of SickleDex Reagent Powder to one bottle of SickleDex Solubility Buffer
3	Place a white dispenser cap on the bottle of working solubility buffer
4	Agitate the buffer vigorously until the reagent powder is dissolved
5	Label the buffer container with the following after reconstitution: <ul style="list-style-type: none"> <li>• Reconstitution date</li> <li>• 45-day expiration date</li> <li>• Tech ID</li> </ul>

**SickleDex Testing**

STEP	ACTION
1	Label 12 x 75 mm tubes for each donor unit to be tested, 1 positive control and 1 negative control
2	Allow controls and buffer to warm to room temperature (18°C to 30°C) for approximately 15 minutes before use
3	Place tubes in testing rack and fill each tube with working SickleDex Solubility Buffer to the red line on rack (approximately 2 mL) <b>NOTE:</b> Return buffer to refrigerator immediately after use
4	Mix controls by holding vertically between hands and rolling the vials back and forth for 20-30 seconds followed by inverting end over-end 20 times. <b>NOTE:</b> Visually inspect the bottom of the vial to ensure all cells are suspended
5	Add 1 drop of each control into the appropriate labeled control tube <b>NOTE:</b> The control MUST be inverted and held vertically directly over the test tube to ensure accurate delivery. Wipe threads on each control before returning cap, if necessary
6	Add 20uL of whole blood and swirl the contents to mix <b>NOTE:</b> If donor segment does not appear to contain whole blood (due to settling prior to sealing), use 10uL of packed red blood cells for testing
7	Allow tubes to stand in the testing rack at room temperature for at least 6 minutes and no longer than 60 minutes

<b>TITLE: Testing and Provision of Hemoglobin S Negative Blood</b>	<b>Number: PC-0036.02</b>
--	---------------------------

8	Read reaction macroscopically by looking through the test tubes at black lines on the back of the testing rack
---	--

**CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES**


**Interpretation:**

If black tube rack lines are	Interpret as
<b>CLEARLY VISIBLE</b> thru a transparent suspension	Negative
<b>NOT VISIBLE</b> thru a cloudy, turbid suspension	Positive

**Results Reporting in Sunquest**

STEP	ACTION												
1	Open 'Blood Order Processing' function												
2	Assign unit to patient using 'Blood Inventory Search' or by scanning unit barcode												
3	<p>QA Warning window will open with warning message about missing NHBS attribute on unit</p> <ul style="list-style-type: none"> <li>Check the box on the left to acknowledge QA Warning for each unit</li> </ul> <table border="1"> <thead> <tr> <th>Acknowledge</th> <th>Unit</th> <th>Comp Type</th> <th>Div #</th> <th>Warning message</th> <th>Have authority to override?</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/></td> <td>W141605200006</td> <td>E0382</td> <td>00</td> <td>Patient/unit attribute incompatibility Patient's Attribute(s): NHBS Attribute(s) Missing on Unit: NHBS</td> <td>Yes</td> </tr> </tbody> </table>	Acknowledge	Unit	Comp Type	Div #	Warning message	Have authority to override?	<input type="checkbox"/>	W141605200006	E0382	00	Patient/unit attribute incompatibility Patient's Attribute(s): NHBS Attribute(s) Missing on Unit: NHBS	Yes
Acknowledge	Unit	Comp Type	Div #	Warning message	Have authority to override?								
<input type="checkbox"/>	W141605200006	E0382	00	Patient/unit attribute incompatibility Patient's Attribute(s): NHBS Attribute(s) Missing on Unit: NHBS	Yes								
4	<ul style="list-style-type: none"> <li>Order a sickle cell test on the allocated unit</li> <li>Enter ;<b>SCKL</b> in the 'Add unit test' box at the bottom of the screen</li> </ul>												
5	<p>Enter the reactions according to the following tables:</p> <table border="1"> <thead> <tr> <th></th> <th>SKP</th> <th>SKN</th> <th>SKU</th> </tr> </thead> <tbody> <tr> <td>SCKL</td> <td>+</td> <td>-</td> <td>+ or -</td> </tr> </tbody> </table> <p>SKP = Sickle Positive Control; SKN = Sickle Negative Control; SKU = Unit Tested</p>		SKP	SKN	SKU	SCKL	+	-	+ or -				
	SKP	SKN	SKU										
SCKL	+	-	+ or -										
6	<p>Enter the Interpretation in the interpretation grid:</p> <table border="1"> <thead> <tr> <th>Interpretation</th> <th>SQ Result</th> <th>Sunquest Hot Key</th> </tr> </thead> <tbody> <tr> <td>Positive</td> <td>POSHBS</td> <td>s</td> </tr> <tr> <td>Negative</td> <td>NGHBS</td> <td>n</td> </tr> </tbody> </table>	Interpretation	SQ Result	Sunquest Hot Key	Positive	POSHBS	s	Negative	NGHBS	n			
Interpretation	SQ Result	Sunquest Hot Key											
Positive	POSHBS	s											
Negative	NGHBS	n											
7	<ul style="list-style-type: none"> <li>Click &lt;Save&gt;</li> <li>Override QA Failure by adding a 'TDONE' comment <ul style="list-style-type: none"> <li>In the free text box, add 'Sickle Cell Testing'</li> </ul> </li> <li>Click &lt;OK&gt; to acknowledge unit failure on the Electronic Crossmatch Eligibility Report</li> </ul>												
8	Open 'Blood Product Testing' function												
9	<ul style="list-style-type: none"> <li>Scan unit information</li> <li>Click &lt;Continue&gt;</li> </ul>												

<b>TITLE: Testing and Provision of Hemoglobin S Negative Blood</b>	<b>Number: PC-0036.02</b>
--	---------------------------

10	Order ‘;AO’ in the ‘Unit Testing’ field	
	If SCKL testing is	Then
	POSITIVE	Result with ; <b>POSHBS</b>
	NEGATIVE	Result with ; <b>NHBS</b>
11	Click <OK> to acknowledge the Blood Product Testing Comment	 <p>The screenshot shows a dialog box titled "Blood Product Testing" with a close button (X) in the top right corner. The text inside reads: "Results to AO will add to, but not delete from the unit's permanent antibody/antigen field. Use Blood Product Entry to delete any incorrect Ag/Ab code(s) previously entered." There is an "OK" button at the bottom right.</p>
12	Click <Save> and <Exit>	
13	Go to ‘Blood Order Processing’	
	If SCKL testing on unit is	Then
	NEGATIVE	Go to next step
	POSITIVE	In the Compatibility Testing area <ul style="list-style-type: none"> <li>• Uncheck ‘Use reaction grid results’ box</li> <li>• Enter ‘;ND’ in the ‘XM’ field</li> <li>• Enter ‘J’ in the ‘TS’ field for ‘NOK’</li> <li>• Go to step 15</li> </ul> <b>NOTE:</b> Unit will automatically be released from patient
14	Enter comment code on unit to display the status on the Transfusion Record <ul style="list-style-type: none"> <li>• Enter ;<b>CM</b> in the ‘Add unit test’ box at the bottom of the screen</li> <li>• Enter the following code in the CM box ;<b>NHBS</b></li> </ul>	
15	<b>If performing</b>	<b>Then</b>
	Electronic crossmatch	Click <Save>
	Serologic crossmatch (adult)	Enter results according to appropriate Crossmatch SOP then click <Save>
	TNRBC (Neonate crossmatch)	Answer ‘TS’ Box Click <Save>
16	Continue to Blood Product Issue if issuing unit	

**CALIBRATION:**

NA

**PROCEDURE NOTES AND LIMITATIONS:**

- Sickle –Chex controls expire 100 days after open or manufacturer’s expiration, whichever date is shorter
- Donors with Hgb S concentrations less than 30.9% may not be detected

**REFERENCES:**

SickleDex Manufacturer Insert. Omaha, Nebraska: Streck; current version  
 Sickle-Chek Manufacturer Insert. Omaha, Nebraska: Streck; current version

<b>TITLE: Testing and Provision of Hemoglobin S Negative Blood</b>	<b>Number: PC-0036.02</b>
--	-------------------------------

**RELATED DOCUMENTS:**

**APPENDIX:**  
NA

TRAINING

<b>TITLE: Testing and Provision of Hemoglobin S Negative Blood</b>	<b>Number: PC-0036.02</b>
--	---------------------------

<b>UWMC SOP Approval:</b>	
<b>UWMC CLIA Medical Director</b>	_____ Date _____
	Mark H. Wener, MD
<b>Transfusion Service Manager</b>	_____ Date _____
	Deanne Stephens
<b>Compliance Analyst</b>	_____ Date _____
	Christine Clark
<b>Transfusion Service Medical Director</b>	_____ Date _____
	Monica B. Pagano, MD
<b>UWMC Biennial Review:</b>	
	_____ Date _____
	_____ Date _____

**REVISION HISTORY:**  
04/16/18: Updated instructions for mixing controls to match revision to package insert.

