



University of Washington Medical Center
1959 NE Pacific Street. Seattle, WA 98195
Transfusion Services Laboratory
Policies and Procedures Manual

Original Effective Date:
10-28-2020

Revision Effective Date:

Number:
PC-0087.01

TITLE: Visual Inspection of Blood Components at Northwest Campus

PURPOSE:

To provide instructions for visual inspection of blood components

LOCATION:

Northwest Lab Transfusion Support Service (TSS)

PRINCIPLE & CLINICAL SIGNIFICANCE:

Principle

Visual inspection is a major control point of manufacturing to ensure the safety, purity and potency of blood components.

Clinical Significance

Blood components are susceptible to damage during shipment and components not shipped and stored at the proper temperatures are at increased risk for bacterial contamination, hemolysis and other deleterious effects or may otherwise not function as expected and should be discarded to protect the potential recipient.

POLICIES:

- Blood components and blood component labels are visually inspected for acceptability at critical points from receipt to final disposition and these points are defined in standard operating procedure specific for the procedure/process being performed.
- Visual inspection at some critical points is documented electronically or manually. When required, documentation includes the following:
 - Date and time of inspection
 - Donor Identification Number (DIN)
 - Outcome of inspection
 - Action taken (QI form)
 - Identity of person performing inspection
- Blood components found visually unacceptable are quarantined electronically and physically until final disposition is determined.
 - A Quality Improvement (QI) form is initiated with the reason for quarantine and faxed to Montlake TSL manager for determination of final disposition – refer to SOP ***Quarantine and Final Disposition of Blood Components at Northwest Campus***

Quarantine Terminology	Meaning
Electronic Quarantine	Component is quarantine in the laboratory information system (SQ)
Physical Quarantine	Place in an appropriate storage device in a location labeled as 'Quarantine' and segregated from non-quarantined blood components

TITLE: Visual Inspection of Blood Components at Northwest Campus	Number: PC-0087.01
---	-------------------------------

SPECIMEN REQUIREMENTS:

NA

REAGENTS/SUPPLIES/EQUIPMENT:

Reagents:	Supplies:	Equipment:
NA	NA	NA

QUALITY CONTROL:

NA

INSTRUCTIONS:

TABLE of CONTENTS

[Inspecting the Component Label and Bag](#)
[Verify Color](#) and Consistency of Contents

Inspecting the Component Label and Bag

STEP	ACTION	
1	Verify the component label (ISBT label) is legible and includes: <ul style="list-style-type: none"> • One Donor Identification Number (and aliquot number, if applicable) • In-date expiration date/time • Component name and E code • Attributes and/or processing if applicable • ABO/Rh label 	
2	Inspect the blood containers for leaks, breaks or tears NOTE: Light pressure may be applied to the unit.	
3	Inspect the following for presence of blood or plasma which may indicate inadequate closure or sealing <ul style="list-style-type: none"> • Bag seams • Port areas 	
4	If component is	Verify
	Thawed plasma or cryoprecipitate	Port areas are free of moisture
	RBC component	Presence of at least one integrally attached segment at the time of issue for transfusion
5	If component is	Verify
	Irradiated	Verify indicator window on the Rad-Sure indicator film is BLACK , obscuring the word "NOT" NOTE: RadSure labels may not be present on all irradiated components due to transfer of product from primary container following irradiation
6	If label and/or bag are	Go to
	Acceptable	Next section
	Unacceptable	Results Reporting in Sunquest

Verify Color and Consistency of Contents

STEP	ACTION	
1	Verify component is acceptable in appearance	
	Component	Unacceptable Appearance
	All components	<ul style="list-style-type: none"> • Blood or plasma observed in the ports at the sealing sites in tubing • Leaks or cracks in the bag or tubing • Murky, dark purple, brown in color (contamination should be expected if observed) • Grossly lipemic plasma (pronounced milky appearance) is usually considered unsuitable for transfusion • Foreign object
	Red cell components and segments	<ul style="list-style-type: none"> • Hemolysis: <ul style="list-style-type: none"> ○ Segments appear much lighter in color than red cells in bags ○ Zone of hemolysis observed just about the cell mass (less opaque, pink plasma) • Visible clots and/or aggregates: Dark purple to very dark burgundy masses that do not disperse easily by gently manipulation or change in temperature <p>EXCEPTIONS:</p> <ul style="list-style-type: none"> ○ A green hue from light-induced changes in bilirubin pigments is not cause for rejection ○ White particulate matter generally described as one of the following: white specks, fatty material, crystalline material waxy appearing globs - refer to Blood Component Visual Inspection Guide, AABB/AARC
Plasma, Cryoprecipitate Granulocytes	<ul style="list-style-type: none"> • Hemolysis: pink to red in color • Visible clots: a thick whitish, opaque mass that does not disperse easily by gentle manipulation or change in temperature <ul style="list-style-type: none"> ○ Opacity prevents inspection for clots or fibrin strands • Signs of unexpected thawing if components are frozen <p>EXCEPTIONS:</p> <ul style="list-style-type: none"> ○ Green hue due to light-induced changes in bilirubin or pale green color due to birth control pills are acceptable and not cause for rejection ○ Flocculent material that appears after thawing: describes as cloudy, fuzzy, or fluffy white precipitate that may have a tissue paper-like appearance. Can easily disperses by gentle manipulation or increase in temperature 	

TITLE: Visual Inspection of Blood Components at Northwest Campus	Number: PC-0087.01
---	---------------------------


STEP	ACTION	
	Platelets	<ul style="list-style-type: none"> • Hemolysis: pink to red in color • Visible clots and/or aggregates: a thick whitish, opaque mass that does not disperse easily by gentle manipulation or change in temperature • Gross lipemia or icterus • No swirl noted <p>EXCEPTIONS</p> <ul style="list-style-type: none"> ○ White particulate matter generally described as one of the following: white specks, fatty material, crystalline material waxy appearing globs - refer to Blood Component Visual Inspection Guide, AABB/AARC ○ Pathogen reduced platelets
2	Go to Results Reporting in Sunquest	





CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES

Interpretation

Any component not passing the visual inspection should not be issued for transfusion

Results Reporting in Sunquest

STEP	ACTION	
1	Record the result of the visual inspection in the LIS	
	If	Then record at the following step
	Issuing to a patient	
	If visual inspection is	Then answer the question "Select units-Visual inspection ok?"
ACCEPTABLE	<ul style="list-style-type: none"> • Yes • Continue following SOP <i>Issuing Blood Component at Northwest Campus</i> 	
UNACCEPTABLE	<ul style="list-style-type: none"> • No • Quarantine the component following SOP <i>Quarantine and Final Disposition of Blood Components at Northwest Campus</i> 	

STEP	ACTION							
								
	<p>Receiving , returning or shipping</p>	<div data-bbox="548 632 1068 772"> <p>Status entry</p> <p>Date: 09/28/2015 Time: 18:00 Destination: SCCA</p> <p>Pass visual inspection for all units <input type="checkbox"/> Yes <input type="checkbox"/> No</p> </div> <table border="1" data-bbox="548 783 1414 1713"> <thead> <tr> <th data-bbox="548 783 808 850">If visual inspection is</th> <th data-bbox="808 783 1414 850">Then answer the question "Pass Visual Inspection for all units"</th> </tr> </thead> <tbody> <tr> <td data-bbox="548 850 808 1224">ACCEPTABLE</td> <td data-bbox="808 850 1414 1224"> <ul style="list-style-type: none"> • Yes • Refer to the appropriate SOP <ul style="list-style-type: none"> ○ <i>Receiving Blood Components from Montlake at Northwest Campus</i> ○ <i>Returning Issued Blood Components to Inventory at Northwest Campus</i> ○ <i>Returning Blood Components to Montlake from Northwest Campus</i> </td> </tr> <tr> <td data-bbox="548 1224 808 1713">UNACCEPTABLE</td> <td data-bbox="808 1224 1414 1713"> <ul style="list-style-type: none"> • No <p>Quarantine the component following SOP <i>Quarantine and Final Disposition of Blood Components at Northwest Campus</i></p>  </td> </tr> </tbody> </table>	If visual inspection is	Then answer the question "Pass Visual Inspection for all units"	ACCEPTABLE	<ul style="list-style-type: none"> • Yes • Refer to the appropriate SOP <ul style="list-style-type: none"> ○ <i>Receiving Blood Components from Montlake at Northwest Campus</i> ○ <i>Returning Issued Blood Components to Inventory at Northwest Campus</i> ○ <i>Returning Blood Components to Montlake from Northwest Campus</i> 	UNACCEPTABLE	<ul style="list-style-type: none"> • No <p>Quarantine the component following SOP <i>Quarantine and Final Disposition of Blood Components at Northwest Campus</i></p> 
If visual inspection is	Then answer the question "Pass Visual Inspection for all units"							
ACCEPTABLE	<ul style="list-style-type: none"> • Yes • Refer to the appropriate SOP <ul style="list-style-type: none"> ○ <i>Receiving Blood Components from Montlake at Northwest Campus</i> ○ <i>Returning Issued Blood Components to Inventory at Northwest Campus</i> ○ <i>Returning Blood Components to Montlake from Northwest Campus</i> 							
UNACCEPTABLE	<ul style="list-style-type: none"> • No <p>Quarantine the component following SOP <i>Quarantine and Final Disposition of Blood Components at Northwest Campus</i></p> 							

PROCEDURE NOTES/LIMITATIONS

- It is not possible to determine all possible adverse conditions through a visual inspection of the component and where indicated, additional testing will be incorporated into the process.

- Blood labels are electronically checked for accuracy when blood component are modified and relabeled.

REFERENCES:

- Technical Manual. Bethesda, MD: AABB Press, current edition
- Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB Press, current edition

RELATED DOCUMENTS:

FORM *Quality Improvement Form*

SOP *Receiving Blood Components from Montlake at Northwest Campus L*

SOP *Returning Blood Components to Montlake from Northwest Campus*

SOP *Quarantine and Final Disposition of Blood Components at Northwest Campus*

SOP *Issuing Blood Components at Northwest Campus*

SOP *Returning Issued Blood Components to Inventory at Northwest Campus*

APPENDIX:

NA

TITLE: Visual Inspection of Blood Components at Northwest Campus

**Number:
PC-0087.01**

UWMC SOP Approval:

**UWMC CLIA
Medical Director** Mark Wener Date 10/20/20
Mark H. Wener, MD

**Transfusion
Service Manager** Nina Sen Date 10/16/20
Nina Sen

**Compliance
Analyst** Christine Clark Date 10-16-2020
Christine Clark

**Transfusion
Service
Medical Director** Monica Pagano Date 10-19-2020
Monica Pagano, MD

UWMC Biennial Review:

_____ Date _____

_____ Date _____