



University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA 98195 Transfusion Services Laboratory Policies and Procedures Manual	Original Effective Date: 02-11-16	Number: PC-0031.02
	Revision Effective Date:	
TITLE: LISS Indirect Antiglobulin Technique		

PURPOSE:

To provide instructions using tube Indirect Antiglobulin Test (IAT) technique with Low Ionic Strength Solution (LISS) as a secondary manual method to detect the presence of unexpected antibodies as a secondary method to a PeG Indirect Antiglobulin test

PRINCIPLE & CLINICAL SIGNIFICANCE:

An indirect antiglobulin test (IAT) demonstrates in-vitro reactions between red cells and antibodies, and is used in antibody detection, antibody identification, crossmatching, and blood group phenotyping.

Use of a LISS additive accelerates antibody binding to red cells

POLICIES:

- This method is used as an alternate method (see SOP *PeG Indirect Antiglobulin Technique* for primary method) for:
 - Detection and identification of unexpected antibodies (antibody screen or panel)
 - AHG crossmatching when the patient has a history of unexpected clinically significant antibodies

SPECIMEN REQUIREMENTS:

EDTA is preferred and if not tested soon after collection, should be stored at 1-6°C
 Red top tubes are acceptable
 Refer to SOP *Specimen Acceptability and Order Receipt*

REAGENTS/SUPPLIES/EQUIPMENT:

Reagents:	Supplies:	Equipment:
<ul style="list-style-type: none"> • LISS (Ortho AES) • Anti-IgG • Antibody Screen/Panel Cells • IgG coated control cells • Blood Bank Saline 	<ul style="list-style-type: none"> • 12 x 75 glass tubes • Blood bank transfer pipettes 	<ul style="list-style-type: none"> • Calibrated serologic centrifuge • Calibrated cell washer • 37°C Heat block • Agglutination viewer

QUALITY CONTROL:

Quality Control of LISS is performed on DAY OF USE. All other reagents are quality controlled daily

TITLE: LISS Indirect Antiglobulin Test	Number: PC-0031.02
---	-------------------------------

INSTRUCTIONS:

STEP	ACTION	
1	Verify quality control for LISS was performed	
	If	Then
	Performed	Go to next step
	Not performed	<ul style="list-style-type: none"> Perform QC as per SOP <i>Quality Control of Manual Testing Reagents</i> Go to next step
2	Label 12 X 75 mm tubes for each cell to be tested according to the <i>Labelling for Manual Testing SOP</i>	
	NOTE: Bio-Rad reagent red blood cells need to be room temperature before using	
3	Make an approximate 3-4% red cell suspension of patient or donor cells as per the following:	
	Test	3-4% Cell Suspension
	Antibody Panel	Patient cells
	Crossmatch	Donor cells
	NOTE: Reagent red blood cells can be used directly from the vial.	
4	Add 2 drops of patient plasma or serum to each tube	
5	Add 1 drop of red cells to the respectively labelled tubes and mix gently	
	Test	Red Cells
	Antibody screen	Screening cells
	Panel	Panel cells
	Auto control	Patient 3-4% red cell suspension
	Crossmatch	Donor 3-4% red cell suspension
6	Add 2 drops of LISS reagent to each tube	
7	Mix well and incubate all tubes at 37°C ±1 for 10-30 minutes	
8	Centrifuge according to calibration for LISS phase and observe for hemolysis	
9	Shake gently to resuspend the cell buttons and examine macroscopically for agglutination	
10	Read, grade and record reactions	
11	Wash the tubes 3 times with saline	
12	Add 2 drops of Anti-IgG and mix well and centrifuge according to calibration for AHG phase testing	
13	Shake gently to resuspend the cell buttons and examine macroscopically for agglutination	

TITLE: LISS Indirect Antiglobulin Test	Number: PC-0031.02
---	-------------------------------

STEP	ACTION
14	Read, grade and record reactions
15	Add 1 drop of IgG coated control cells to each tube with a negative antiglobulin test
16	Mix well and centrifuge according to calibration for AHG phase testing
17	Shake gently to resuspend the cell buttons, and examine for agglutination
18	Read, grade and record reactions – refer to Interpretation and Result Reporting below

CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES

Interpretation

If the following is observed	Interpret as
Agglutination of test cells after any phase of testing	Positive
No agglutination of test cells	Negative NOTE: Results are considered invalid and must be repeated if tube does not agglutinate with IgG coated control cells

If test performed is	Then						
Antibody Screen	Refer to SOP <i>Antibody Screen</i>						
Antibody Panel	Refer to SOP <i>Antibody Identification</i>						
Crossmatch	<table border="1"> <thead> <tr> <th>If</th> <th>Then RBC is</th> </tr> </thead> <tbody> <tr> <td>No hemolysis at 37°C AND No agglutination at any phase of testing.</td> <td> <ul style="list-style-type: none"> Compatible </td> </tr> <tr> <td>Hemolysis at 37°C AND/OR Agglutination at any phase of testing.</td> <td> <ul style="list-style-type: none"> Incompatible </td> </tr> </tbody> </table>	If	Then RBC is	No hemolysis at 37°C AND No agglutination at any phase of testing.	<ul style="list-style-type: none"> Compatible 	Hemolysis at 37°C AND/OR Agglutination at any phase of testing.	<ul style="list-style-type: none"> Incompatible
	If	Then RBC is					
No hemolysis at 37°C AND No agglutination at any phase of testing.	<ul style="list-style-type: none"> Compatible 						
Hemolysis at 37°C AND/OR Agglutination at any phase of testing.	<ul style="list-style-type: none"> Incompatible 						

Results Reporting:

STEP	ACTION							
1	Record reactions immediately upon reading in on appropriate form or Sunquest (SQ)							
2	If performing	Then						
	Antibody screen	Refer to SOP <i>Antibody Screen Testing</i>						
	Antibody ID	<ul style="list-style-type: none"> Record results on the appropriate antibody ID panel or Extended Testing Worksheet Refer to SOP <i>Antibody Identification</i> 						
	Crossmatch	Enter appropriate reaction in SQ Reaction Grid and interpretation fields per table below						
		XIS	XINC	XAHG	XCC	Interpretation	SQ Code	SQ Hot Key
		ND	0	0	+	LISS Compatible	CMP	C
		ND	0/+	0/+	ND	LISS Incompatible	ICMP	I

NOTES AND LIMITATIONS:

- Test method must be followed exactly to avoid changing the ionic strength of the mixture
- Weak reactions may be obtained if the tests are incubated less than 10 minutes for LISS
- Complement dependent antibodies may not be detected if plasma is used
- See current version of package insert
- 3-4% suspension of RBCs can be prepared by adding one drop of packed RBCs and approximately 1-1.5 mL of blood bank saline in an appropriately labeled tube. The suspension can be prepared using volume estimation with comparison to the reagent red cells for visual verification.

REFERENCES:

- Technical Manual. Bethesda, MD: AABB Press, current edition
- Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB Press, current edition
- ORTHO® Antibody Enhancement Solution Ortho-Clinical Diagnostics, Inc. Raritan, NJ 2011

RELATED DOCUMENTS:

- SOP *Antibody Identification*
- SOP *Antibody Screen Testing*
- SOP *Specimen Acceptability and Order Receipt*
- SOP *Grading Reactions*
- SOP *Quality Control for Manual Testing Reagents*
- SOP *Labelling Tubes and Gel Cards for Testing*

APPENDIX:

NA

TITLE: LISS Indirect Antiglobulin Test	Number: PC-0031.02
---	-------------------------------

UWMC SOP Approval:	
UWMC CLIA Medical Director	_____ Date _____
	Mark H. Wener, MD
Transfusion Service Manager	_____ Date _____
	Nina Sen
Compliance Analyst	_____ Date _____
	Christine Clark
Transfusion Service Medical Director	_____ Date _____
	Monica Pagano, MD
UWMC Biennial Review:	
	_____ Date _____
	_____ Date _____

REVISION:
03/30/2021: Instructions added for entering reactions and interpretation for crossmatches

