



<b>University of Washington Medical Center</b> <b>1959 NE Pacific Street. Seattle, WA 98195</b> <b>Transfusion Services Laboratory</b> <b>Policies and Procedures Manual</b>	<b>Original Effective Date:</b> <b>02-11-16</b>	<b>Number:</b> <b>PC-0032.02</b>
	<b>Revision Effective Date:</b>	
<b>TITLE: PeG Indirect Antiglobulin Technique</b>		

**PURPOSE:**

To provide instructions using tube Indirect Antiglobulin Test (IAT) technique with PEG to detect the presence of unexpected antibodies [and perform seriological crossmatches](#)

**PRINCIPLE & CLINICAL SIGNIFICANCE:**

An indirect antiglobulin test (IAT) demonstrates in-vitro reactions between red cells and antibodies, and is used in antibody detection, antibody identification, crossmatching, and blood group phenotyping.

PEG acts as an additive for tests to detect blood group antibodies, enhancing the sensitivity of the antibody detection by creating a low-ionic strength test environment that increases the rate of antibody uptake during incubation

**POLICIES:**

This is the primary manual method for:

- Detection and identification of unexpected antibodies (antibody screen or panel)
- AHG crossmatching when the patient has a history of unexpected clinically significant antibodies

**SPECIMEN REQUIREMENTS:**

EDTA is preferred and if not tested soon after collection, should be stored at 1-6°C.  
 Red top tubes are acceptable.  
 See SOP Specimen Acceptability [and Order Receipt](#)

**REAGENTS/SUPPLIES/EQUIPMENT:**

Reagents:	Supplies:	Equipment:
<ul style="list-style-type: none"> <li>• Gamma PeG</li> <li>• Anti-IgG</li> <li>• Antibody Screen/Panel Cells</li> <li>• IgG coated control cells</li> <li>• Blood Bank Saline</li> </ul>	<ul style="list-style-type: none"> <li>• 12 x 75 glass tubes</li> <li>• Blood bank transfer pipettes</li> </ul>	<ul style="list-style-type: none"> <li>• Calibrated serologic centrifuge</li> <li>• Calibrated cell washer</li> <li>• 37°C heat block</li> <li>• Agglutination viewer</li> </ul>

**QUALITY CONTROL:**

Quality Control is performed daily

**INSTRUCTIONS:**

STEP	ACTION										
1	Label 12 X 75 mm tubes for each cell to be tested according to the <i>Labelling for Manual Testing SOP</i>  <b>NOTE:</b> Bio-Rad reagent red blood cells need to be room temperature before using.										
2	Make an approximate 3-4% red cell suspension of patient or donor cells as per the following: <table border="1" data-bbox="313 541 1411 688"> <thead> <tr> <th data-bbox="313 541 613 590">Test</th> <th data-bbox="613 541 1411 590">3-4% Cell Suspension</th> </tr> </thead> <tbody> <tr> <td data-bbox="313 590 613 638">Antibody Panel</td> <td data-bbox="613 590 1411 638">Patient cells</td> </tr> <tr> <td data-bbox="313 638 613 688">Crossmatch</td> <td data-bbox="613 638 1411 688">Donor cells</td> </tr> </tbody> </table> <b>NOTE:</b> Reagent red blood cells can be used directly from the vial.	Test	3-4% Cell Suspension	Antibody Panel	Patient cells	Crossmatch	Donor cells				
Test	3-4% Cell Suspension										
Antibody Panel	Patient cells										
Crossmatch	Donor cells										
3	Add <b>2</b> drops of patient plasma or serum to each tube.										
4	Add <b>1</b> drop of red cells to the respectively labelled tubes and mix gently <table border="1" data-bbox="313 831 1411 1079"> <thead> <tr> <th data-bbox="313 831 565 879">Test</th> <th data-bbox="565 831 1411 879">Red Cells</th> </tr> </thead> <tbody> <tr> <td data-bbox="313 879 565 928">Antibody screen</td> <td data-bbox="565 879 1411 928">Screening cells</td> </tr> <tr> <td data-bbox="313 928 565 976">Panel</td> <td data-bbox="565 928 1411 976">Panel cells</td> </tr> <tr> <td data-bbox="313 976 565 1024">Auto control</td> <td data-bbox="565 976 1411 1024">Patient 3-4% red cell suspension</td> </tr> <tr> <td data-bbox="313 1024 565 1079">Crossmatch</td> <td data-bbox="565 1024 1411 1079">Donor 3-4% red cell suspension</td> </tr> </tbody> </table>	Test	Red Cells	Antibody screen	Screening cells	Panel	Panel cells	Auto control	Patient 3-4% red cell suspension	Crossmatch	Donor 3-4% red cell suspension
Test	Red Cells										
Antibody screen	Screening cells										
Panel	Panel cells										
Auto control	Patient 3-4% red cell suspension										
Crossmatch	Donor 3-4% red cell suspension										
5	<table border="1" data-bbox="313 1079 1411 1377"> <thead> <tr> <th data-bbox="313 1079 565 1127">If performing</th> <th data-bbox="565 1079 1411 1127">Then</th> </tr> </thead> <tbody> <tr> <td data-bbox="313 1127 565 1304">Crossmatch</td> <td data-bbox="565 1127 1411 1304"> <ul style="list-style-type: none"> <li>• Centrifuge and observe for hemolysis</li> <li>• Shake gently to resuspend the cell button and examine macroscopically for agglutination</li> <li>• Read, grade (refer to <i>SOP Grading Reactions</i>) and record reactions</li> </ul> </td> </tr> <tr> <td data-bbox="313 1304 565 1377">Antibody screen/panel</td> <td data-bbox="565 1304 1411 1377"> <ul style="list-style-type: none"> <li>• Go to next step (immediate spin phase is not required for antibody screen or panel)</li> </ul> </td> </tr> </tbody> </table>	If performing	Then	Crossmatch	<ul style="list-style-type: none"> <li>• Centrifuge and observe for hemolysis</li> <li>• Shake gently to resuspend the cell button and examine macroscopically for agglutination</li> <li>• Read, grade (refer to <i>SOP Grading Reactions</i>) and record reactions</li> </ul>	Antibody screen/panel	<ul style="list-style-type: none"> <li>• Go to next step (immediate spin phase is not required for antibody screen or panel)</li> </ul>				
If performing	Then										
Crossmatch	<ul style="list-style-type: none"> <li>• Centrifuge and observe for hemolysis</li> <li>• Shake gently to resuspend the cell button and examine macroscopically for agglutination</li> <li>• Read, grade (refer to <i>SOP Grading Reactions</i>) and record reactions</li> </ul>										
Antibody screen/panel	<ul style="list-style-type: none"> <li>• Go to next step (immediate spin phase is not required for antibody screen or panel)</li> </ul>										
6	Add 2 drops of PeG to each test tube										
7	Mix well and incubate at 37°C ±1 for 10-30 minutes										
8	Examine for hemolysis and document in the INC phase in the Sunquest testing grid										
9	Wash the tubes 3 times with saline										
10	Add 2 drops of Anti-IgG and mix well and centrifuge according to calibration for AHG phase testing										
11	Shake gently to resuspend the cell buttons, and examine macroscopically for agglutination										
12	Read, grade and record reactions										

STEP	ACTION
13	Add 1 drop of IgG coated control cells to each tube with a negative antiglobulin test
14	Mix gently and centrifuge according to calibration for AHG phase testing
15	Shake gently to resuspend the cell buttons, and examine for agglutination
16	Read, grade and record reactions
<u>17</u>	<u><a href="#">Go to section Interpretation and Results Reporting</a></u>

**CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES**

**Interpretation**

If Agglutination	Interpret as
Present after any phase testing	Positive
Not present	Negative
	<b>NOTE:</b> Results are considered invalid and must be repeated if tube does not agglutinate with IgG coated control cells

If test performed is	Then						
Antibody Screen	Refer to SOP <i>Antibody Screen Testing</i>						
Antibody Panel	Refer to SOP <i>Antibody Identification</i>						
Crossmatch	<table border="1"> <thead> <tr> <th>If</th> <th>Then RBC is</th> </tr> </thead> <tbody> <tr> <td>No hemolysis at 37°C <b>AND</b> No agglutination at any phase of testing</td> <td>Compatible <i>Refer to SOP Crossmatching</i></td> </tr> <tr> <td>Hemolysis at 37°C <b>AND/OR</b> Agglutination at any phase of testing</td> <td>Incompatible <i>Refer to SOP Crossmatching</i></td> </tr> </tbody> </table>	If	Then RBC is	No hemolysis at 37°C <b>AND</b> No agglutination at any phase of testing	Compatible <i>Refer to SOP Crossmatching</i>	Hemolysis at 37°C <b>AND/OR</b> Agglutination at any phase of testing	Incompatible <i>Refer to SOP Crossmatching</i>
	If	Then RBC is					
No hemolysis at 37°C <b>AND</b> No agglutination at any phase of testing	Compatible <i>Refer to SOP Crossmatching</i>						
Hemolysis at 37°C <b>AND/OR</b> Agglutination at any phase of testing	Incompatible <i>Refer to SOP Crossmatching</i>						

**Results Reporting:**

STEP	ACTION																																												
1	Record reactions immediately upon reading in on appropriate form or Sunquest																																												
2	<b>If performing</b>	<b>Then</b>																																											
	Antibody screen	Refer to SOP <i>Antibody Screen Testing</i>																																											
	Antibody ID	<ul style="list-style-type: none"> <li>• Record results on the appropriate antibody ID panel or Extended Testing Worksheet</li> <li>• Refer to SOP <i>Antibody Identification</i></li> </ul>																																											
	Crossmatch	Enter appropriate reaction and interpretation per table below																																											
	<table border="1" style="width: 100%; border-collapse: collapse; color: red;"> <thead> <tr style="background-color: #cccccc;"> <th colspan="4">Reaction Grid</th> <th rowspan="2">Interpretation</th> <th rowspan="2">SQ Code</th> <th rowspan="2">SQ Hot Key</th> </tr> <tr style="background-color: #cccccc;"> <th>XIS</th> <th>XINC</th> <th>XAHG</th> <th>XCC</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">0</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">0</td> <td style="text-align: center;">±</td> <td style="text-align: center;">PEG Compatible</td> <td style="text-align: center;">PCMP</td> <td style="text-align: center;">!</td> </tr> <tr> <td style="text-align: center;">0</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">±</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">PEG Incompatible</td> <td style="text-align: center;">PICMP</td> <td style="text-align: center;">%</td> </tr> <tr> <td style="text-align: center;">±</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">±</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">PEG Incompatible</td> <td style="text-align: center;">PICMP</td> <td style="text-align: center;">%</td> </tr> <tr> <td style="text-align: center;">±</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">0</td> <td style="text-align: center;">ND</td> <td style="text-align: center;">PEG Incompatible</td> <td style="text-align: center;">PICMP</td> <td style="text-align: center;">%</td> </tr> </tbody> </table>							Reaction Grid				Interpretation	SQ Code	SQ Hot Key	XIS	XINC	XAHG	XCC	0	ND	0	±	PEG Compatible	PCMP	!	0	ND	±	ND	PEG Incompatible	PICMP	%	±	ND	±	ND	PEG Incompatible	PICMP	%	±	ND	0	ND	PEG Incompatible	PICMP
Reaction Grid				Interpretation	SQ Code	SQ Hot Key																																							
XIS	XINC	XAHG	XCC																																										
0	ND	0	±	PEG Compatible	PCMP	!																																							
0	ND	±	ND	PEG Incompatible	PICMP	%																																							
±	ND	±	ND	PEG Incompatible	PICMP	%																																							
±	ND	0	ND	PEG Incompatible	PICMP	%																																							

**NOTES AND LIMITATIONS:**

- Polyethylene glycol ~~has a tendency to~~ tends to precipitate serum globulins. Accordingly, when using Gamma PeG to detect unexpected antibodies, it is especially important to assure that the red blood cells are thoroughly resuspended in each change of saline during the washing phases of the test. When testing samples containing elevated globulin levels, three washes may not be sufficient to remove unbound protein. If precipitated globulin remains enmeshed in the red blood cell button, it may neutralize the Anti-Human Globulin and cause a false negative test result. In the case of specimens having an exceptionally high level of globulin (as in multiple myeloma), the addition of polyethylene glycol may cause a gel to form. This will preclude the use of the polyethylene glycol test procedure when testing these specimens for unexpected antibodies.
- 
- Precipitation of fibrinogen may be observed when testing plasma samples. In such cases, as with elevated globulin levels, it may be necessary to wash the red blood cells more than three times to remove all unbound human protein.
- 
- Precipitation of serum proteins when PeG is added appears to be related to elevated serum globulin levels. The problem becomes apparent when the IgG-coated red cells are nonreactive or unexplained weak reactions are detected. At least four manual washes of the red cells at AHG phase, with agitation, will fully resuspend the red cells and usually prevent the problem from occurring.

- 
- Test method must be followed exactly to avoid changing the ionic strength of the mixture.
- 
- 3-4% suspension of RBCs can be prepared by adding one drop of packed RBCs and approximately 1-1.5 mL of blood bank saline in an appropriately labeled tube. The suspension can be prepared using volume estimation with comparison to the reagent red cells for visual verification.
- 
- It may be necessary to switch to an alternate test method such as LISS due to problems with protein precipitation that causes neutralization of the antiglobulin reagent causing false negative test results. In these cases, the IgG coated control cells will not react and the test will be invalid. (see SOP *LISS Indirect Antiglobulin Technique*)
- 
- Warm autoantibodies may react strongly in PeG test method and can mask the presence of underlying alloantibodies. In these cases, it may be desirable to switch to a less sensitive method such as LISS to attempt to detect and identify any underlying alloantibodies. (see SOP *LISS Indirect Antiglobulin Technique*)

**REFERENCES:**

- Technical Manual. Bethesda, MD: AABB Press, current edition
- Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB Press, current edition
- Gamma Peg - Polyethylene Glycol Additive for Antibody Detection Tests, Immucor, Norcross, GA 10/10.

**RELATED DOCUMENTS:**

SOP *Labelling Tubes and Gel Cards for Testing*  
SOP *Antibody Identification*  
SOP *Antibody Screen Testing*  
~~SOP *Crossmatching*~~  
SOP *Quality Control of Manual Testing Reagents*  
SOP *Grading Reactions*  
SOP *Specimen Acceptability and Order Receipt*  
SOP *LISS Indirect Antiglobulin Test*

**APPENDIX:**

NA

**UWMC SOP Approval:**

**UWMC CLIA  
Medical Director**

Mark H. Wener, MD

Date

**Transfusion  
Service Manager**

[Nina Sen](#)

Date

**Compliance  
Analyst**

Christine Clark

Date

**Transfusion  
Service  
Medical Director**

[Monica Pagano](#), MD

Date

**UWMC Biennial Review:**

Date

Date

[03/30/2021: Instructions added for entering reactions and interpretation for crossmatches](#)

TRAINING