



<b>University of Washington Medical Center</b> <b>1959 NE Pacific Street. Seattle, WA 98195</b> <b>Transfusion Services Laboratory</b> <b>Policies and Procedures Manual</b>	<b>Original Effective Date:</b> <b>03-11-16</b>	<b>Number:</b> <b>PC-0003.02</b>
	<b>Revision Effective Date:</b>	
<b>TITLE: Dividing Blood Components</b>		

**PURPOSE:**

To provide instructions for preparing blood splits using syringes or small bags with or without the sterile connecting device (SCD)

**PRINCIPLE & CLINICAL SIGNIFICANCE:**

NA

**POLICIES:**

- Product prepared by the TSL using a filter-syringe do not require further filtration at administration
- RBC's and FFP prepared in a syringe expire within 24 hours or the same time as the parent bag expiration, whichever is shorter
- Platelets and cryoprecipitate prepared in a syringe expire within 4 hours or the same time as the parent bag expiration, whichever is shorter
- Sterile connections:
  - Components remaining in the parent blood bag will maintain the original outdate when the secondary container is connected to the blood bag using a sterile method such as the sterile connecting device in a closed system
  - If a sterile connecting device (closed system) is not used or a weld fails, the parent container has the same outdate as the aliquots (24 hours for RBC and plasma, 4 hours for platelets and cryoprecipitate). The same as an open system.
- Splitting blood units using the SCD is the preferred method for preparing blood component aliquots, but is not required if the parent blood component will expire within 24 hours
- Expiration of divided products may not exceed the expiration for the parent blood component
- Blood Label Check must be performed following preparation of the aliquot on both the new product and the parent bag
- Lot number and expiration date of any additional containers added to the original product container must be tracked in Sunquest BOP or on the Downtime Component Prep Log along with a weld inspection if the SCD if used for attachment

**SPECIMEN REQUIREMENTS:**

NA

**REAGENTS/SUPPLIES/EQUIPMENT:**

Reagents:	Supplies:	Equipment:
None	<ul style="list-style-type: none"> <li>• Blood Component</li> <li>• Filter-Syringe Set (30mL or 60mL)</li> <li>• Syringe Tip Cap</li> <li>• Transfer Bag</li> <li>• Multi-bag set</li> <li>• Hemostat or clamp</li> </ul>	<ul style="list-style-type: none"> <li>• Tube Welder</li> <li>• Balance</li> <li>• Tube Sealer</li> </ul>

**QUALITY CONTROL:**

SCD weld inspection is performed and documented on the associated blood component with each use

**INSTRUCTIONS:**

STEP	ACTION							
1	Select the appropriate transfer bag or syringe filter set for the blood split to be performed							
2	Thoroughly mix the contents of the blood component							
3	<b>If</b>	<b>Then</b>						
	Dividing using the SCD	<table border="1"> <thead> <tr> <th>If using</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Syringe</td> <td> <ul style="list-style-type: none"> <li>• Tighten the tubing on the syringe to prevent leakage</li> <li>• Heat seal near the spike, remove and discard spike and blue clip</li> <li>• Attach the new syringe to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul> </td> </tr> <tr> <td>Bag or multi-bag set</td> <td> <ul style="list-style-type: none"> <li>• Place the empty bag on the scale and tare the scale (refer to SOP <i>Operation of Entris 2202-1S Balance</i>)</li> <li>• Attach the new bag to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul> </td> </tr> </tbody> </table>	If using	Then	Syringe	<ul style="list-style-type: none"> <li>• Tighten the tubing on the syringe to prevent leakage</li> <li>• Heat seal near the spike, remove and discard spike and blue clip</li> <li>• Attach the new syringe to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul>	Bag or multi-bag set	<ul style="list-style-type: none"> <li>• Place the empty bag on the scale and tare the scale (refer to SOP <i>Operation of Entris 2202-1S Balance</i>)</li> <li>• Attach the new bag to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul>
	If using	Then						
Syringe	<ul style="list-style-type: none"> <li>• Tighten the tubing on the syringe to prevent leakage</li> <li>• Heat seal near the spike, remove and discard spike and blue clip</li> <li>• Attach the new syringe to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul>							
Bag or multi-bag set	<ul style="list-style-type: none"> <li>• Place the empty bag on the scale and tare the scale (refer to SOP <i>Operation of Entris 2202-1S Balance</i>)</li> <li>• Attach the new bag to the original unit (refer to SOP <i>Sterile Welder Operation Cleaning Maintenance</i>)</li> <li>• Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul>							
	<ul style="list-style-type: none"> <li>• Express any residual air back into parent product container</li> <li>• Seal the tubing (refer to SOP <i>Tube Sealer: Operation and Cleaning</i>)</li> </ul>							

		<ul style="list-style-type: none"> <li>Do not separate the new component from the parent component</li> <li>Perform Blood Component Processing in Sunquest (refer to SOP <i>Blood Component Preparation in Sunquest</i>)</li> <li>Verify the expiration date/time of the products are correct and correct if needed</li> </ul>																
	<p>Dividing without using the SCD</p>	<table border="1" data-bbox="527 449 1386 1073"> <thead> <tr> <th data-bbox="527 449 669 485">If Using</th> <th data-bbox="669 449 1386 485">Then</th> </tr> </thead> <tbody> <tr> <td data-bbox="527 485 669 730">Syringe</td> <td data-bbox="669 485 1386 730"> <ul style="list-style-type: none"> <li>Tighten the tubing on the syringe to prevent leakage</li> <li>Sterilely insert the spike from the syringe into the blood bag port</li> <li>Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul> </td> </tr> <tr> <td data-bbox="527 730 669 1073">Bag or multi-bag set</td> <td data-bbox="669 730 1386 1073"> <ul style="list-style-type: none"> <li>Place the empty bag on the scale and tare the scale. See SOP <i>Operation of Entris 2202-1S Balance</i></li> <li>Clamp the bag tubing with a hemostat or clip</li> <li>Use the sterile method to insert the spike attached to the transfer bag into the blood bag port and open the clamp</li> <li>Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul> </td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>Express any residual air in the bag back into the parent product container Seal the tubing. See SOP <i>Tube Sealer: Operation and Cleaning</i></li> <li>Do not separate the new component from the parent component</li> <li>Perform Blood Component Processing in Sunquest (refer to SOP <i>Blood Component Preparation in Sunquest</i>)</li> <li>Verify the expiration date/time of the products are correct and correct if needed</li> </ul> <table border="1" data-bbox="574 1419 1386 1598"> <thead> <tr> <th data-bbox="574 1419 938 1455">Component</th> <th data-bbox="938 1419 1386 1455">Expiration Date/Time limit</th> </tr> </thead> <tbody> <tr> <td data-bbox="574 1455 938 1491">RBC</td> <td data-bbox="938 1455 1386 1491">24 Hours from "spiking" unit</td> </tr> <tr> <td data-bbox="574 1491 938 1526">FFP</td> <td data-bbox="938 1491 1386 1526">4 hours from "spiking" unit</td> </tr> <tr> <td data-bbox="574 1526 938 1562">Platelets</td> <td data-bbox="938 1526 1386 1562"></td> </tr> <tr> <td data-bbox="574 1562 938 1598">Cryoprecipitate</td> <td data-bbox="938 1562 1386 1598"></td> </tr> </tbody> </table> <p><b>NOTE:</b> Both the components will have the same expiration date/time</p>	If Using	Then	Syringe	<ul style="list-style-type: none"> <li>Tighten the tubing on the syringe to prevent leakage</li> <li>Sterilely insert the spike from the syringe into the blood bag port</li> <li>Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul>	Bag or multi-bag set	<ul style="list-style-type: none"> <li>Place the empty bag on the scale and tare the scale. See SOP <i>Operation of Entris 2202-1S Balance</i></li> <li>Clamp the bag tubing with a hemostat or clip</li> <li>Use the sterile method to insert the spike attached to the transfer bag into the blood bag port and open the clamp</li> <li>Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul>	Component	Expiration Date/Time limit	RBC	24 Hours from "spiking" unit	FFP	4 hours from "spiking" unit	Platelets		Cryoprecipitate	
If Using	Then																	
Syringe	<ul style="list-style-type: none"> <li>Tighten the tubing on the syringe to prevent leakage</li> <li>Sterilely insert the spike from the syringe into the blood bag port</li> <li>Draw the requested amount of blood product into the syringe allowing 3-5mL extra for the blood transfusion tubing</li> </ul>																	
Bag or multi-bag set	<ul style="list-style-type: none"> <li>Place the empty bag on the scale and tare the scale. See SOP <i>Operation of Entris 2202-1S Balance</i></li> <li>Clamp the bag tubing with a hemostat or clip</li> <li>Use the sterile method to insert the spike attached to the transfer bag into the blood bag port and open the clamp</li> <li>Allow the requested amount of blood to flow into the transfer bag allowing 10-15mL extra for the blood transfusion tubing</li> </ul>																	
Component	Expiration Date/Time limit																	
RBC	24 Hours from "spiking" unit																	
FFP	4 hours from "spiking" unit																	
Platelets																		
Cryoprecipitate																		
5	Separate the blood products by cutting the seal																	
6	<table border="1" data-bbox="310 1724 568 1759"> <thead> <tr> <th data-bbox="310 1724 568 1759">If container is a</th> <th data-bbox="568 1724 1427 1759">Then</th> </tr> </thead> <tbody> <tr> <td data-bbox="310 1759 568 1894">Syringe</td> <td data-bbox="568 1759 1427 1894"> <ul style="list-style-type: none"> <li>Remove the stub end from the syringe and replace with a sterile cap</li> <li>Division process is complete – continue with additional processing as necessary (such as irradiation)</li> </ul> </td> </tr> </tbody> </table>	If container is a	Then	Syringe	<ul style="list-style-type: none"> <li>Remove the stub end from the syringe and replace with a sterile cap</li> <li>Division process is complete – continue with additional processing as necessary (such as irradiation)</li> </ul>													
If container is a	Then																	
Syringe	<ul style="list-style-type: none"> <li>Remove the stub end from the syringe and replace with a sterile cap</li> <li>Division process is complete – continue with additional processing as necessary (such as irradiation)</li> </ul>																	

<b>TITLE: Dividing Blood Components</b>	<b>Number: PC-0003.02</b>
---	-------------------------------

	Bag	Process is complete – continue with additional processing as necessary (such as irradiation)
7	Return components to the appropriate storage location	

**CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES**

NA

**CALIBRATION:**

NA

**PROCEDURE NOTES AND LIMITATIONS:**

- Sunquest defaults the expiration date of ALIQ aliquots to 4 hours. If desired, it is acceptable to change the expiration of FFP to 24 hours from time of processing

**REFERENCES:**

- *Entris 2202-1S Balance Operation Manual*
- Standards for Blood Banks and Transfusion Services. Bethesda, MD; AABB, current edition

**RELATED DOCUMENTS:**

- FORM Component Prep Downtime Log
- SOP Blood Component Preparation
- SOP Operation of Entris 2202-1S Balance
- SOP Sterile Welder Operation Cleaning Maintenance

**APPENDIX**

NA

Trainline

**UWMC SOP Approval:**

**UWMC CLIA  
Medical Director**

Mark H. Wener, MD Date

**Transfusion  
Service Manager**

Nina Sen Date

**Compliance  
Analyst**

Christine Clark Date

**Transfusion  
Service  
Medical Director**

Monica Pagano, MD Date

**UWMC Biennial Review:**

Date

Date

**REVISIONS:**

07/12/2021: Add policy when sterile welder is not used or fails the component involved has the out-date of an open system.