



<b>University of Washington Medical Center</b> <b>1959 NE Pacific Street. Seattle, WA 98195</b> <b>Transfusion Services Laboratory</b> <b>Policies and Procedures Manual</b>	<b>Original Effective Date:</b> <b>10-28-20</b>	<b>Number:</b> <b>PC-0082.02</b>
	<b>Revision Effective Date:</b> <b>04-25-22</b>	
<b>TITLE: LIS Downtime Blood Component Labeling</b>		

**PURPOSE:**

To describe the procedure for printing ISBT labels and relabeling blood components following modifications during LIS Downtime using the HemaTrax Blood and Blood Component Stand-alone Labeling Application and Digi-Trax ZM400 printer

**LOCATION:**

Northwest Transfusion Support Service (TSS)  
Montlake Transfusion Service Laboratory (TSL)

**PRINCIPLE & CLINICAL SIGNIFICANCE:**

**Principle**

FDA standards for the labeling of blood products must continue to be met during LIS computer outages and other times when labels are printed using the non-interfaced stand-alone Digi-Trax ZM400 printer. The stand-alone HemaTrax system meets the United States Industry Consensus Standard for the Uniform Labeling of Blood and Blood Components using ISBT-128 format.

**Clinical Significance**

Consistency in labeling maintains regulatory compliance and helps to ensure the safety of the transfusion process by allowing the nursing bedside verification of the blood product label to remain unchanged.

**POLICIES:**

- A second person must verify accurate relabeling of blood components when Sunquest Blood Label Check is not available.
- Divided components requiring relabeling of the original parent product
- The original Donor Identification Number (DIN) must remain visible when relabeling the component. Never cover over the original DIN with a new label.

**SPECIMEN REQUIREMENTS:**

NA

**REAGENTS/SUPPLIES/EQUIPMENT:**

Reagents:	Supplies:	Equipment:
None	<ul style="list-style-type: none"> <li>• HemaTrax printer labels</li> <li>• Printer ribbons</li> </ul>	<ul style="list-style-type: none"> <li>• Dedicated PC with integrated stand-alone HemaTrax printer</li> <li>• Bar-code scanner</li> </ul>

**QUALITY CONTROL:**

NA

**INSTRUCTIONS:**

**TABLE OF CONTENTS:**

[Printing HemaTrax Labels](#)

[Blood Label Check](#)

[Appendix 1: Blood Component Prep E-codes](#)

**Printing HemaTrax Labels**

Step	Action						
<b>1</b>	Log into the HemaTrax computer using the following username and password <ul style="list-style-type: none"> <li>• Username: <b>hematrax</b></li> <li>• Password: <b>D1G1Trax</b></li> </ul>						
<b>2</b>	Click on <HemaTrax Unity Client> on the desktop using the following username and password <ul style="list-style-type: none"> <li>• Username: <b>hematrax</b></li> <li>• Password: <b>hematrax</b></li> </ul>						
<b>3</b>	Select the “Full Face (4x4)” format from the “Select label” drop down box						
<b>4</b>	Select the ‘Zebra’ printer from “Select printer” drop down box						
<b>5</b>	Check the “Set to Scan Mode” in the upper left						
<b>6</b>	<table border="1" style="width: 100%;"> <thead> <tr> <th style="background-color: #d3d3d3;">If component is</th> <th style="background-color: #d3d3d3;">Then</th> </tr> </thead> <tbody> <tr> <td>NOT divided</td> <td> <ul style="list-style-type: none"> <li>• Go to next step</li> </ul> </td> </tr> <tr> <td>Divided component</td> <td> <ul style="list-style-type: none"> <li>• Check the divided unit box</li> <li>• Enter aliquot letters for Division 1 and Division 2</li> </ul> </td> </tr> </tbody> </table>	If component is	Then	NOT divided	<ul style="list-style-type: none"> <li>• Go to next step</li> </ul>	Divided component	<ul style="list-style-type: none"> <li>• Check the divided unit box</li> <li>• Enter aliquot letters for Division 1 and Division 2</li> </ul>
If component is	Then						
NOT divided	<ul style="list-style-type: none"> <li>• Go to next step</li> </ul>						
Divided component	<ul style="list-style-type: none"> <li>• Check the divided unit box</li> <li>• Enter aliquot letters for Division 1 and Division 2</li> </ul>						
<b>7</b>	<table border="1" style="width: 100%;"> <thead> <tr> <th style="background-color: #d3d3d3;">If relabeling</th> <th style="background-color: #d3d3d3;">Then</th> </tr> </thead> <tbody> <tr> <td><b>Original container with a label</b></td> <td> <ul style="list-style-type: none"> <li>• <b>Cross out the supplier license number on the modified component</b></li> </ul> </td> </tr> <tr> <td>New container without a label</td> <td> <ul style="list-style-type: none"> <li>• Go to next step</li> </ul> </td> </tr> </tbody> </table>	If relabeling	Then	<b>Original container with a label</b>	<ul style="list-style-type: none"> <li>• <b>Cross out the supplier license number on the modified component</b></li> </ul>	New container without a label	<ul style="list-style-type: none"> <li>• Go to next step</li> </ul>
If relabeling	Then						
<b>Original container with a label</b>	<ul style="list-style-type: none"> <li>• <b>Cross out the supplier license number on the modified component</b></li> </ul>						
New container without a label	<ul style="list-style-type: none"> <li>• Go to next step</li> </ul>						
<b>8</b>	<p><b>Quadrant I: DIN</b></p> <ul style="list-style-type: none"> <li>• Click on the QI square</li> <li>• Click in the ‘Facility ID Number’ field and scan DIN barcode from the label</li> </ul> <table border="1" style="width: 100%;"> <thead> <tr> <th style="background-color: #d3d3d3;">If</th> <th style="background-color: #d3d3d3;">Then</th> </tr> </thead> <tbody> <tr> <td>Scanner fails</td> <td>                     Manually key in the following from the DIN                     <ul style="list-style-type: none"> <li>• Facility ID number (WXXXX)</li> <li>• Collection year (XX)</li> <li>• Serial number/unit number (XXXXXX)</li> </ul> </td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>• Check the “Include Check Digit” box</li> <li>• Click &lt;OK&gt;</li> </ul>	If	Then	Scanner fails	Manually key in the following from the DIN <ul style="list-style-type: none"> <li>• Facility ID number (WXXXX)</li> <li>• Collection year (XX)</li> <li>• Serial number/unit number (XXXXXX)</li> </ul>		
If	Then						
Scanner fails	Manually key in the following from the DIN <ul style="list-style-type: none"> <li>• Facility ID number (WXXXX)</li> <li>• Collection year (XX)</li> <li>• Serial number/unit number (XXXXXX)</li> </ul>						

Step	Action												
<b>9</b>	<p><b>Quadrant II: ABO Blood Group/Rh Factor</b></p> <ul style="list-style-type: none"> <li>Click on the QII square</li> <li>Select the blood type for new label</li> <li>Select the appropriate donation type from the 'Intended use' drop down box               <ul style="list-style-type: none"> <li>Volunteer Allogeneic Donation</li> <li>Autologous Use Only</li> </ul> </li> <li>Click &lt;OK&gt;</li> </ul>												
<b>10</b>	<p><b>Quadrant III: Product Code:</b></p> <ul style="list-style-type: none"> <li>Click on the QIII square</li> <li>Enter the component output E-code based on the initial product E-code and the modification performed - refer to <a href="#">Appendix 1: Blood Component Prep E-codes</a></li> <li>Click &lt;OK&gt;</li> </ul> <p><b>NOTE:</b> The output component Ecode can be found by entering information about the product category e.g. anticoagulant, draw volume, additive solution and clicking the search key.</p>												
<b>11</b>	<p><b>Quadrant VI: Expiration date and time:</b></p> <ul style="list-style-type: none"> <li>Click on the QIV square</li> <li>Select the correct Expiration option           <table border="1" style="margin-left: 20px; width: 80%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="padding: 2px;">If expiration time is</th> <th style="padding: 2px;">Then click on</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Midnight or 23:59</td> <td style="padding: 2px;">"Expiration date only"</td> </tr> <tr> <td style="padding: 2px;">NOT midnight</td> <td style="padding: 2px;">"Expiration date and Time"</td> </tr> </tbody> </table> </li> <li>Manually enter the correct new expiration date and time</li> </ul> <p><b>Special test field:</b></p> <ul style="list-style-type: none"> <li>Select any special test from the choices</li> <li>Click &lt;OK&gt;</li> </ul> <p><b>Processing Facility Box:</b></p> <table border="1" style="margin-left: 20px; width: 80%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="padding: 2px;">If</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">No additional processing</td> <td style="padding: 2px;">Go to next step</td> </tr> <tr> <td style="padding: 2px;">Addition processing performed</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> <li>Enter the facility code in the box under "Facility World Code", W2584 is for UWMC</li> <li>Select "Further Processing By" from "Processing Legend" drop down</li> <li>Check box to include FDA registration number</li> </ul> </td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>Click &lt;OK&gt;</li> </ul>	If expiration time is	Then click on	Midnight or 23:59	"Expiration date only"	NOT midnight	"Expiration date and Time"	If	Then	No additional processing	Go to next step	Addition processing performed	<ul style="list-style-type: none"> <li>Enter the facility code in the box under "Facility World Code", W2584 is for UWMC</li> <li>Select "Further Processing By" from "Processing Legend" drop down</li> <li>Check box to include FDA registration number</li> </ul>
If expiration time is	Then click on												
Midnight or 23:59	"Expiration date only"												
NOT midnight	"Expiration date and Time"												
If	Then												
No additional processing	Go to next step												
Addition processing performed	<ul style="list-style-type: none"> <li>Enter the facility code in the box under "Facility World Code", W2584 is for UWMC</li> <li>Select "Further Processing By" from "Processing Legend" drop down</li> <li>Check box to include FDA registration number</li> </ul>												
<b>12</b>	Click <PRINT> to print a single label												

Step	Action	
13	Verify the label accuracy by comparing printed label to original product label and verifying all information is correct (E-code, product type, expiration date/time, volume and division codes)	
	<b>If information is</b>	<b>Then</b>
	Correct	Go to next step
	Incorrect	Update any fields as necessary and reprint the label
14	<b>If</b>	<b>Then</b>
	No additional label is needed	<ul style="list-style-type: none"> <li>Exit the program and log out</li> </ul>
	Label for a different component type is needed	<ul style="list-style-type: none"> <li>Select "All Quadrants" in the "Clear Quadrants" box in the bottom left</li> <li>Click &lt;Clear&gt;</li> <li>Repeat steps 8 thru 14</li> </ul>
	Additional label of the same component type is needed	<ul style="list-style-type: none"> <li>Follow steps 7 thru 14 as required making any needed changes such as volume</li> </ul>
15	Go to section <a href="#">Blood Label Check</a> <b>IMPORTANT:</b> Do not relabel the component until a second tech verifies the label is accurate	

**[Blood Label Check](#)**

Step	Action
1	Obtain the <i>Component Prep Downtime Log</i> and ensure all information under the Input and Prep sections are complete  <b>Note:</b> UWMC NW campus only-The Prep section is not applicable. It can be left blank when using form for documentation of thawing components.
2	Ensure all required information in the Output section is complete except for "Label Verified By 2 <sup>nd</sup> Tech"
3	Ask a second trained staff member to perform and document a Blood Label Check following steps 4 thru 8 below
4	Verify the following information is identical on the original and new labels: <ul style="list-style-type: none"> <li>Unit Number</li> <li>ABO/Rh</li> </ul>
5	Verify that the following are accurate for the output product label based on the modification performed: <ul style="list-style-type: none"> <li>Product description and component code (Ecode) (refer to Appendix A: Blood Component Prep E-codes)</li> <li>Expiration date/time</li> <li>Product Volume</li> <li>Division code, when applicable</li> </ul>
6	Verify the collection facility and the modifying facility (Further Processed By:) are accurate and the license number of the collection facility is crossed out  <b>NOTE:</b> University of Washington Medical Center should always be selected as the

<b>TITLE: LIS Downtime Blood Component Labeling</b>	<b>Number: PC-0082.02</b>
---	-------------------------------

<b>Step</b>	<b>Action</b>
	modifying facility
<b>7</b>	Affix label to cover previous label without covering the original DIN
<b>8</b>	Have the 2 <sup>nd</sup> verifier initial the 'Verified By 2 <sup>nd</sup> Tech' box on the <i>Component Prep Downtime Log</i> to indicate the component is labeled correctly

**CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL**

NA

**VALUES/CRITICAL VALUES:**

NA

**CALIBRATION:**

NA

**NOTES AND LIMITATIONS:**

- HemaTrax tracks the activities of every user from the time they log-on to log-off
- The stand-alone HemaTrax database is populated with product labels defined in the International Council for Commonality in Blood Banking Automation (ICCBBA) product table and must be updated periodically to remain current
- The "clear quadrants" button can be used to facilitate removal of prior unit information prior to creating additional labels

**REFERENCES:**

- Hematrax Blood and Blood Component Stand Alone Labeling Application Version 6.2.1
- Standards for Blood Banks and Transfusion Services, American Association of Blood Banks, Bethesda, MD. Current Edition.

**RELATED DOCUMENTS:**

FORM *Component Prep Downtime Log*  
FORM *Irradiation Downtime Log*

**APPENDIX:**

Appendix 1: Blood Component Prep E-codes (see attached separate document)

<b>TITLE: LIS Downtime Blood Component Labeling</b>	<b>Number: PC-0082.02</b>
---	-------------------------------

<b>UWMC SOP Approval:</b>	
<b>UWMC CLIA Medical Director</b>	_____ Date _____
	Andrew Bryan, MD
<b>Transfusion Service Manager</b>	_____ Date _____
	Nina Sen
<b>QA Manager</b>	_____ Date _____
	Taylor Reeves
<b>Transfusion Service Medical Director</b>	_____ Date _____
	Monica Pagano, MD
<b>UWMC Biennial Review:</b>	
	_____ Date _____
	_____ Date _____

**REVISIONS:**

4/19/22- Step 1 of Blood Label Check section revised for leaving prep section blank when thawing components