



University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA, 98195 Transfusion Services Laboratory Policies and Procedures Manual	Original Effective Date: 02-11-2016	Number: PC-0038.02
	Revision Effective Date: 05-23-2022	
TITLE: Weak D Manual Tube Testing		

PURPOSE:

To provide instructions for performing Weak D testing using a manual tube technique

PRINCIPLE & CLINICAL SIGNIFICANCE:

Principle

Red blood cells from patients are tested with commercial anti-D reagent to determine if the D antigen is present on the patient's red blood cells. In the presence of D antigen, Anti-D reagents bind to the antigen on sample red blood cells and cause an antigen-antibody reaction visible as red blood cell agglutination. RBCs with weakened D antigen may require incubation with the anti-D reagent or testing through the antiglobulin phase for detection.

Clinical Significance

In some recipients determined to be weak D, the possibility exists that pregnancy or transfusion with D+ RBC containing components could result in the production of an allo anti-D.

POLICY

- Weak D testing is performed on the following specimens when the initial Rh result is negative:
 - Cord blood specimens for determination of maternal Rh immune globulin eligibility of a Rh-negative mother
 - Specimens from donors or potential donors of blood or stem cell components (does not include red blood cell component type confirmation testing)
- Weak D testing may also be utilized for ABO/Rh type discrepancy resolution
- For newly identified weak D positive perinatal females:
 - Send the sample to reference lab for weak D type 1,2,3, genotype evaluation
 - Notify the TSL MD on-call 8 am to 5 pm to determine the need for Rh immune globulin.

SPECIMEN REQUIREMENTS:

EDTA is preferred and if not tested soon after collection, should be stored at 1-6°C

Red top clotted blood samples are also acceptable.

See SOP *Specimen Acceptability and Order Receipt*

REAGENTS/SUPPLIES/EQUIPMENT:

Reagents:	Supplies:	Equipment:
<ul style="list-style-type: none"> • Anti-D • ABO + Rh Control • Blood Bank Saline • Anti-IgG • IgG coated control cells 	<ul style="list-style-type: none"> • 12 x 75 glass tubes • Blood bank transfer pipettes 	<ul style="list-style-type: none"> • Calibrated Serologic Centrifuge • Calibrated Cell washer • 37 °C Heat block • Agglutination viewer

QUALITY CONTROL:

Quality Control is performed each day of use.

INSTRUCTIONS:

STEP	ACTION						
1	Label tube(s) per SOP <i>Labeling for Manual Testing</i>						
	<table border="1"> <thead> <tr> <th data-bbox="313 457 800 506">If initial ABO/Rh test tubes are</th> <th data-bbox="800 457 1412 506">Then</th> </tr> </thead> <tbody> <tr> <td data-bbox="313 506 800 648">Available</td> <td data-bbox="800 506 1412 648"> Use the D tube and patient cell suspension from ABO/Rh typing to continue testing <ul style="list-style-type: none"> • Label one tube for control • Go to next step </td> </tr> <tr> <td data-bbox="313 648 800 957">Not available</td> <td data-bbox="800 648 1412 957"> <ul style="list-style-type: none"> • Label 3 tubes: one for patient red cell suspension, one for control and one for anti-D • Prepare an approximate 3-4% patient cell suspension • Add the following to the anti-D tube <ul style="list-style-type: none"> ○ 1 drop of Anti-D ○ Add 1 drop of the patient's cell suspension </td> </tr> </tbody> </table>	If initial ABO/Rh test tubes are	Then	Available	Use the D tube and patient cell suspension from ABO/Rh typing to continue testing <ul style="list-style-type: none"> • Label one tube for control • Go to next step 	Not available	<ul style="list-style-type: none"> • Label 3 tubes: one for patient red cell suspension, one for control and one for anti-D • Prepare an approximate 3-4% patient cell suspension • Add the following to the anti-D tube <ul style="list-style-type: none"> ○ 1 drop of Anti-D ○ Add 1 drop of the patient's cell suspension
	If initial ABO/Rh test tubes are	Then					
Available	Use the D tube and patient cell suspension from ABO/Rh typing to continue testing <ul style="list-style-type: none"> • Label one tube for control • Go to next step 						
Not available	<ul style="list-style-type: none"> • Label 3 tubes: one for patient red cell suspension, one for control and one for anti-D • Prepare an approximate 3-4% patient cell suspension • Add the following to the anti-D tube <ul style="list-style-type: none"> ○ 1 drop of Anti-D ○ Add 1 drop of the patient's cell suspension 						
2	Add 1 drop of the ABO/Rh Control reagent to the control tube						
3	Add 1 drop of the 3-4% patient cell suspension to the control tube						
4	Mix and incubate the anti-D and control tubes at 37°C ± 1 for 15-30 minutes						
5	Wash tubes 3 times with blood bank saline						
6	Add 2 drops of anti-IgG to each tube						
7	Mix and centrifuge						
8	Resuspend the cell buttons gently, and examine for agglutination						
9	Read, grade, and record the reactions						
10	Add 1 drop of IgG Coated Control Cells to all negative tests						
11	Mix and centrifuge						
12	Read, grade, and record the results						
13	Go to section "Interpreting & Reporting Results"						

CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES

Results Reporting in Sunquest

Step	Action																													
1	Select the appropriate patient specimen in "Blood Order Processing" (BOP)																													
2	Add; DU (Weak D test) in the <u>A</u> dd Spec. Test window, if not previously ordered																													
3	<p>Enter the results and interpretation</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="3" style="background-color: #cccccc;">If</th> <th style="background-color: #cccccc;">Interpretation</th> <th style="background-color: #cccccc;">Sunquest Key</th> </tr> <tr> <th style="background-color: #cccccc;">Field</th> <th style="background-color: #cccccc;">DAHG</th> <th style="background-color: #cccccc;">DUCC</th> <td rowspan="3" style="text-align: center;">Negative</td> <td rowspan="3" style="text-align: center;">N</td> </tr> </thead> <tbody> <tr> <td>PT</td> <td style="text-align: center;">0</td> <td style="text-align: center;">+</td> </tr> <tr> <td>C</td> <td style="text-align: center;">0</td> <td style="text-align: center;">+</td> </tr> <tr> <th colspan="3" style="background-color: #cccccc;">Field</th> <th style="background-color: #cccccc;">Interpretation</th> <th style="background-color: #cccccc;">Sunquest Key</th> </tr> <tr> <td>PT</td> <td style="text-align: center;">+</td> <td style="text-align: center;">ND</td> <td rowspan="2" style="text-align: center;">Positive</td> <td rowspan="2" style="text-align: center;">P</td> </tr> <tr> <td>C</td> <td style="text-align: center;">0</td> <td style="text-align: center;">+</td> </tr> </tbody> </table>	If			Interpretation	Sunquest Key	Field	DAHG	DUCC	Negative	N	PT	0	+	C	0	+	Field			Interpretation	Sunquest Key	PT	+	ND	Positive	P	C	0	+
If			Interpretation	Sunquest Key																										
Field	DAHG	DUCC	Negative	N																										
PT	0	+																												
C	0	+																												
Field			Interpretation	Sunquest Key																										
PT	+	ND	Positive	P																										
C	0	+																												
4	<p>Add the Antigen/Antibody code corresponding to Du test interpretation to flag the patient's historical record:</p> <ul style="list-style-type: none"> Enter ;PB (Patient Problem Info) in the <u>A</u>dd Spec. Test window Enter the code below: <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="background-color: #cccccc;">If</th> <th style="background-color: #cccccc;">Enter</th> </tr> </thead> <tbody> <tr> <td>Weak D Positive</td> <td>;WKDP</td> </tr> <tr> <td>Weak D Negative</td> <td>;WKDN</td> </tr> </tbody> </table>	If	Enter	Weak D Positive	;WKDP	Weak D Negative	;WKDN																							
If	Enter																													
Weak D Positive	;WKDP																													
Weak D Negative	;WKDN																													
5	<p>Select the initial ABO/Rh test and interpret the Rh as follows:</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="background-color: #cccccc;">Weak D Interpretation is</th> <th style="background-color: #cccccc;">And patient is</th> <th style="background-color: #cccccc;">Rh Interpretation</th> </tr> </thead> <tbody> <tr> <td>Negative</td> <td></td> <td> <ul style="list-style-type: none"> Rh Negative </td> </tr> <tr> <td>Positive</td> <td> <ul style="list-style-type: none"> Female >50 years old Male non-SCCA </td> <td> <ul style="list-style-type: none"> Rh Positive </td> </tr> <tr> <td>Positive</td> <td> <ul style="list-style-type: none"> SCCA patients Female ≤ 50 years old who are not pregnant </td> <td> <ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive </td> </tr> <tr> <td>Positive</td> <td> <ul style="list-style-type: none"> Prenatal </td> <td> <ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive Add a BBC comment: RHREC = Further D antigen characterization should be done through Rh D genotyping to determine the need for Rhogam. If genotyping is not done or not available, patient is a candidate for Rh Immune globulin. </td> </tr> </tbody> </table>	Weak D Interpretation is	And patient is	Rh Interpretation	Negative		<ul style="list-style-type: none"> Rh Negative 	Positive	<ul style="list-style-type: none"> Female >50 years old Male non-SCCA 	<ul style="list-style-type: none"> Rh Positive 	Positive	<ul style="list-style-type: none"> SCCA patients Female ≤ 50 years old who are not pregnant 	<ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive 	Positive	<ul style="list-style-type: none"> Prenatal 	<ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive Add a BBC comment: RHREC = Further D antigen characterization should be done through Rh D genotyping to determine the need for Rhogam. If genotyping is not done or not available, patient is a candidate for Rh Immune globulin. 														
Weak D Interpretation is	And patient is	Rh Interpretation																												
Negative		<ul style="list-style-type: none"> Rh Negative 																												
Positive	<ul style="list-style-type: none"> Female >50 years old Male non-SCCA 	<ul style="list-style-type: none"> Rh Positive 																												
Positive	<ul style="list-style-type: none"> SCCA patients Female ≤ 50 years old who are not pregnant 	<ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive 																												
Positive	<ul style="list-style-type: none"> Prenatal 	<ul style="list-style-type: none"> Rh Negative Add a PB comment: WKDP = Patient is weak D positive Add a BBC comment: RHREC = Further D antigen characterization should be done through Rh D genotyping to determine the need for Rhogam. If genotyping is not done or not available, patient is a candidate for Rh Immune globulin. 																												

Step	Action		
			<ul style="list-style-type: none"> • Send the sample to reference lab for weak D type 1,2,3, genotype evaluation • Notify the TSL MD on-call 8 am to 5 pm (to determine the need for Rh immune globulin. <ul style="list-style-type: none"> ○ Document the following in SQ BOP as a BBCS comment: <ul style="list-style-type: none"> ▪ Date and time of notification ▪ Name of TSL MD notified
	Positive	Cord Blood - male	<ul style="list-style-type: none"> • Rh Positive
	Positive	Cord Blood - female	<ul style="list-style-type: none"> • Rh Positive • Add a PB comment indicating the following “Consider Rh negative for transfusion purposes • Notify the TSL MD on-call 8 am to 5 pm (to determine the need for Rh immune globulin. <ul style="list-style-type: none"> ○ Document the following in SQ BOP as a BBCS comment: <ul style="list-style-type: none"> ▪ Date and time of notification ▪ Name of TSL MD notified

CALIBRATION:
NA

PROCEDURE NOTES AND LIMITATIONS:

- Calls from patient providers regarding Rh genotyping should be referred to the TSL MD on-call.
- Resulting the Rh type as positive in a weak D positive patient may result in a QA failure that can be overridden using BBR along with a free text comment indicating the “Patient is weak D positive”.
- Mixed-field agglutination should be investigated for possible cause – refer to *SOP ABO/Rh Discrepancy Resolution*
- 3-4% suspension of RBCs may be prepared by one of the following:
 - Using volume estimation with comparison to the reagent red cells for visual verification
 - By adding one drop of packed RBCs to approximately 1-1.5 mL of blood bank saline
- Results are considered invalid and must be repeated if negative reactions do not produce agglutination following the addition of IgG Coated Control Cells

REFERENCES:

- Technical Manual. Bethesda, MD: AABB Press, current edition
- Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB Press, current edition

RELATED DOCUMENTS:

SOP Specimen Acceptability and Order Receipt
SOP Quality Control for Manual Testing Reagents
SOP Labeling Tubes and Gel Cards for Testing
SOP Grading Reactions
SOP ABO/Rh Manual Tube Testing
SOP ABO/Rh Discrepancy Resolution

APPENDICES:

UWMC SOP Approval:					
UWMC CLIA Medical Director	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%; padding: 2px 5px;">Andrew Bryan, MD</td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> </table>	Andrew Bryan, MD	Date		
Andrew Bryan, MD	Date				
Transfusion Service Manager	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%; padding: 2px 5px;">Nina Sen</td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> </table>	Nina Sen	Date		
Nina Sen	Date				
QA Manager	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%; padding: 2px 5px;">Tayler Reeves</td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> </table>	Tayler Reeves	Date		
Tayler Reeves	Date				
Transfusion Service Medical Director	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%; padding: 2px 5px;">Monica Pagano, MD</td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> </table>	Monica Pagano, MD	Date		
Monica Pagano, MD	Date				
UWMC Biennial Review:					
	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; padding: 2px 5px;">Date</td> </tr> </table>		Date		Date
	Date				
	Date				

REVISIONS:
07/21/2021: Clarification when to notify the TSL MD for consideration of genotyping, Rhogam eligibility, and adding “consider genotyping” comment to test results. Updated related SOPs.