



University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA 98195 Transfusion Services Laboratory Policies and Procedures Manual	Original Effective Date: 01-30-2017	Number: PC-0060.02
	Revision Effective Date: 05-23-2022	
TITLE: ABO/Rh Discrepancy Resolution		

PURPOSE:

To provide guidance for resolving ABO and Rh testing discrepancies. A discrepancy exists when:

- Results of the forward group do not agree with the reverse group
- Unexplained weak or mixed field reactivity
- Current testing is not in agreement with previous records

PRINCIPLE & CLINICAL SIGNIFICANCE:

The test methods employed depend on the type of discrepancy detected. Failure to resolve a blood typing discrepancy accurately can lead to transfusion of incompatible blood components which can result in hemolytic transfusion reactions.

- Discrepancies may arise from patient identification or labeling errors, intrinsic problems with red cell (forward type), serum (back type), patient disease, treatment, or technical errors in test performance.
- Mixed field (MF) reaction and other discrepancies may arise due to transfusion of ABO compatible but non-identical blood components and or allogeneic BMT
- Weak D positive donors considered Rh positive can lead to Rh discrepancies between the patient and transfused circulating red blood cells.

POLICIES:

- Difficulty in typing and resolution must be documented in the LIS as a BBCS comment.
- Mixed field reactions seen with cord blood samples may be due to contamination from maternal blood, Warton’s jelly and weak expression of RBC antigens. Follow the instructions in section [Cord Blood Mixed Field Reaction](#) for resolution and action to take if not resolved.

SPECIMEN REQUIREMENTS:

EDTA is preferred and if not tested soon after collection, should be stored at 1-6°C. Red top tubes are also acceptable. See SOP *Specimen Acceptability and Order Receipt*

REAGENTS/SUPPLIES/EQUIPMENT:

Refer to the SOP specific for the test being performed.

QUALITY CONTROL:

Refer to the SOP specific for the test being performed.

INSTRUCTIONS:

TABLE OF CONTENTS:

- [All Discrepancies](#)
- [Unexpected positive reactions in the forward type \(patient cells\)](#)
- [Weak or unexpected negative reaction in the forward type \(patient cells\)](#)
- [Unexpected positive reactions in reverse type \(patient serum\)](#)
- [Weak or unexpected negative reaction in the reverse type \(patient serum\)](#)
- [Rh Discrepancy](#)
- [Cord Blood Mixed Field Reaction](#)

All Discrepancies

STEP	ACTION						
1	Review the patient’s history in the Laboratory Information System (LIS) and/or Electronic Medical Record (eMR) for indications of the cause of the discrepancy						
	NOTE: It may be necessary to contact other facilities where the patient may have been treated to obtain transfusion and/or transplant history.						
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 30%;">If</th> <th>Then</th> </tr> </thead> </table>	If	Then				
	If	Then					
	Patient has received a BMT or other type of HSCT	<ul style="list-style-type: none"> If the discrepancy is explained by a mixed population of donor and recipient cells and/or antibodies -report as No Type Determined (NTD). Go to section Reporting Results 					
	Patient received out of group blood	<ul style="list-style-type: none"> If the discrepancy is explained by a mixed population including the donor RBCs, report as the patient’s original blood type Go to section Reporting Results 					
Current does not match historical	Request a new sample and repeat testing. <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 30%;">If new sample</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td style="background-color: #cccccc;">Agrees with historical</td> <td> <ul style="list-style-type: none"> Proceed to interpretation and reporting results. Reject first sample as MISLABEL Look for a reciprocal mislabeled specimen Notify a manager or lead Write a QI </td> </tr> <tr> <td style="background-color: #cccccc;">Disagrees with historical</td> <td> <ul style="list-style-type: none"> Examine historical for misinterpretation Ask RN to verify photo-ID Notify department management and obtain Manager/Lead approval before issuing any products to patient if current and new samples both disagree with historical type </td> </tr> </tbody> </table>	If new sample	Then	Agrees with historical	<ul style="list-style-type: none"> Proceed to interpretation and reporting results. Reject first sample as MISLABEL Look for a reciprocal mislabeled specimen Notify a manager or lead Write a QI 	Disagrees with historical	<ul style="list-style-type: none"> Examine historical for misinterpretation Ask RN to verify photo-ID Notify department management and obtain Manager/Lead approval before issuing any products to patient if current and new samples both disagree with historical type
If new sample	Then						
Agrees with historical	<ul style="list-style-type: none"> Proceed to interpretation and reporting results. Reject first sample as MISLABEL Look for a reciprocal mislabeled specimen Notify a manager or lead Write a QI 						
Disagrees with historical	<ul style="list-style-type: none"> Examine historical for misinterpretation Ask RN to verify photo-ID Notify department management and obtain Manager/Lead approval before issuing any products to patient if current and new samples both disagree with historical type 						
Discrepancy not resolved	Go to next step						

STEP	ACTION	
2	Prepare a new 3-4% patient cell suspension and repeat the ABO/Rh testing using the manual test method	
	If discrepancy is	Then
	Resolved	Go to section Results Reporting
	NOT resolved <ul style="list-style-type: none"> Save reaction results, but do not enter an interpretation in Sunquest Go to next step 	
3	Review patient's transfusion/transplant/age history for cause of the discrepancy. The following sources may be used to obtain relevant patient information: <ul style="list-style-type: none"> Sunquest BBI Epic Contact the patient, patient's physician or caregiver Contact other medical facilities where the patient was provided care 	
4	Select the type of discrepancy below and follow the suggested technique as appropriate to resolve the discrepancy <ul style="list-style-type: none"> Unexpected positive reactions in the forward type (patient cells) Weak or unexpected negative reaction in the forward type (patient cells) Unexpected positive reactions in reverse type (patient serum) Weak or unexpected negative reaction in the reverse type (patient serum) Rh Discrepancy Cord Blood Mixed Field Reaction 	

Unexpected positive reactions in the forward type (patient cells)

STEP	ACTION	
1	Perform an IgG DAT	
	If	Then
	Positive	Treat with EGA and repeat testing with EGA treated cells
	Negative	Wash patient cells with warm saline 4X and repeat testing
2	If resolved	Then
	Yes	Go to section Results Reporting
	No	Consult with TSL management or submit to IRL

Weak or unexpected negative reaction in the forward type (patient cells)

STEP	ACTION	
1	If	Then
	Weak subgroup of A or B or depression of antigen expression	Try the following in order and discontinue additional testing once the discrepancy is resolved: <ul style="list-style-type: none"> Incubate patient's cells with anti-A, anti-B, 7% albumin and anti-A,B (if applicable) for 15-20 minutes at RT and repeat testing Incubate the tubes from repeat ABO/Rh testing (along with 1 drop group O screen cells and two drops of patient's plasma as a control) at 4°C for 10 min Try using a different manufacturer of antisera if available
	Suspect neutralization of blood group antisera	Wash patient cells 3-4 times in blood bank saline and repeat testing
2	If resolved	Then
	Yes	Go to section Results Reporting
	No	Consult with TSL management or submit to IRL

Unexpected positive reactions in reverse type (patient serum)

STEP	ACTION	
1	Perform a microscopic evaluation of the reaction to look for rouleaux	
	If rouleaux is	Then
	Observed (retractile stacking of red cells resembling a stack of coins)	Perform Saline Replacement and repeat testing
	Not Observed	Go to next step
2	Perform an antibody screen at IS phase	
	If Antibody screen is	Then
	Positive	Suspect possible cold antibody: <ul style="list-style-type: none"> Perform a cold panel to detect cold-reacting antibodies – If identified repeat testing with antigen negative reverse cells Prewarm reverse typing for cold agglutinin
Negative	Suspect isoantibodies such as anti-A ₁ <ul style="list-style-type: none"> Type patient cells with A₁ lectin and test patient plasma with A₂ cells 	
3	If resolved	Then
	Yes	Go to section Results Reporting
	No	Consult with TSL management or submit to IRL

Weak or unexpected negative reaction in the reverse type (patient serum)

STEP	ACTION	
1	Suspect depressed antibody production and try the following in the order listed and stop once discrepancy is resolved: <ul style="list-style-type: none"> Incubate patient's plasma with A₁ and B cells, screen cells and an auto control for 15-30 minutes at room temperature (RT) Incubate at 4°C for at least 10 minutes Increase the serum to cell ratio by using 3-4 drops of patient's plasma and one drop of reagent A₁ and B cells and repeat testing 	
2	If resolved	Then
	Yes	Go to section Results Reporting
	No	Consult with TSL management or submit to IRL

Rh Discrepancy

STEP	ACTION	
1	Repeat ABO/Rh using alternate method - refer to SOPs <i>ABO/Rh Manual Tube Method</i> or <i>Ortho Vision Patient and Donor Testing</i> and the appropriate package insert	
	If	Then
	Resolved	Go to section Results Reporting
	Not Resolved	Go to next step
2	Perform Weak D testing per SOP <i>Weak D Tube Testing Procedure</i>	
	If	Then
	Resolved	Report results according to SOP <i>Weak D Tube Testing Procedure</i>
	Not Resolved	Consult with TSL management

Cord Blood Mixed Field Reactions

STEP	ACTION	
1	Wash cord blood specimens 3 times with saline	
2	Repeat ABO/Rh front type	
3	If mixed field is resolved	Then
	Yes	Go to section Results Reporting below
	No	Go to next step
4	If cord blood ABO/Rh is	Then
	For maternal Rh immunoglobulin evaluation (mother is Rh negative)	<ul style="list-style-type: none"> Reject specimen Notify the RN a new specimen is required

STEP	ACTION							
		<ul style="list-style-type: none"> Cancel the BBCORD test in SQ and create a QI <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 30%; padding: 5px;">If</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">First collected specimen</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Reject specimen Notify the RN a new specimen is required Cancel the BBCORD test in SQ and create a QI </td> </tr> <tr> <td style="padding: 5px;">Repeat specimen</td> <td style="padding: 5px;">Go to section Results Reporting below</td> </tr> </tbody> </table>	If	Then	First collected specimen	<ul style="list-style-type: none"> Reject specimen Notify the RN a new specimen is required Cancel the BBCORD test in SQ and create a QI 	Repeat specimen	Go to section Results Reporting below
If	Then							
First collected specimen	<ul style="list-style-type: none"> Reject specimen Notify the RN a new specimen is required Cancel the BBCORD test in SQ and create a QI 							
Repeat specimen	Go to section Results Reporting below							
	NOT for Maternal Rh immunoglobulin evaluation	Go to section Results Reporting below						

CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES

Results Reporting

STEP	ACTION	
	Record all reactions immediately in Sunquest or the appropriate manual result form.	
	If discrepancy is	Then
	Resolved	Refer to interpretation and reporting results according to SOP <i>ABO/Rh Testing by Tube Method</i>
1	Due to a bone marrow transplant (BMT) with discrepant forward/reverse reactions or in between typing	Report the ABO/Rh interpretation as NTD NOTE: In some cases, the patient may not make the new isoantibody following conversion. Testing may still be reported with the new blood group with TSL MD approval on file.
	Due to recent transfusion of blood components that were not identical to the patient's own ABO/Rh	Result the patient's actual blood type and add a BBBS comment explaining the cause of the mixed cell population
	Resolved by an Immunohematology Reference Lab (IRL)	TSL manager or designee will enter the results

STEP	ACTION						
	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 35%; padding: 5px;">Unresolved cord blood specimen mixed field where the source is not contaminated</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> • Report the ABO/Rh interpretation as NTD • Document the following on the cord blood log <ul style="list-style-type: none"> ○ ABO: as “NTD” ○ Rh: record as POS or NEG based on anti-D reaction results – refer to SOP <i>ABO/Rh Manual Tube Method</i> </td> </tr> </table>	Unresolved cord blood specimen mixed field where the source is not contaminated	<ul style="list-style-type: none"> • Report the ABO/Rh interpretation as NTD • Document the following on the cord blood log <ul style="list-style-type: none"> ○ ABO: as “NTD” ○ Rh: record as POS or NEG based on anti-D reaction results – refer to SOP <i>ABO/Rh Manual Tube Method</i> 				
Unresolved cord blood specimen mixed field where the source is not contaminated	<ul style="list-style-type: none"> • Report the ABO/Rh interpretation as NTD • Document the following on the cord blood log <ul style="list-style-type: none"> ○ ABO: as “NTD” ○ Rh: record as POS or NEG based on anti-D reaction results – refer to SOP <i>ABO/Rh Manual Tube Method</i> 						
2	<p>Add a PB comment indicating the cause of the discrepancy in SQ for future reference</p> <p>NOTE: If the cause of the discrepancy is already noted in the BAD file, the PB comment is not required</p>						
3	<p>Determine if the new blood type is discrepant with the type of record in the BAD file</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 50%; padding: 5px;">If</th> <th style="width: 50%; padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Does match the historical type or the historical type is NTD</td> <td style="padding: 5px;">No action needed</td> </tr> <tr> <td style="padding: 5px;">Does NOT match the historical type</td> <td style="padding: 5px;"> <p>Refer the discrepancy to TSL management for evaluation for record correction.</p> <p>NOTE: BMT patient’s historical record should not be updated following type conversion until the conversion is verified to be complete and approved by the TSL MD</p> </td> </tr> </tbody> </table>	If	Then	Does match the historical type or the historical type is NTD	No action needed	Does NOT match the historical type	<p>Refer the discrepancy to TSL management for evaluation for record correction.</p> <p>NOTE: BMT patient’s historical record should not be updated following type conversion until the conversion is verified to be complete and approved by the TSL MD</p>
If	Then						
Does match the historical type or the historical type is NTD	No action needed						
Does NOT match the historical type	<p>Refer the discrepancy to TSL management for evaluation for record correction.</p> <p>NOTE: BMT patient’s historical record should not be updated following type conversion until the conversion is verified to be complete and approved by the TSL MD</p>						

CALIBRATION:

NA

PROCEDURE NOTES AND LIMITATIONS:

- Refer to reagent manufacturer’s package insert for instructions and limitations for use of reagent

REFERENCES:

- Standards for Blood Banks and Transfusion Services, Bethesda, MD; AABB, current edition
- Technical Manual, Bethesda, MD; AABB, current edition

RELATED DOCUMENTS:

- FORM *Extended Workup*
- FORM *Cord Blood Log*
- SOP *Specimen Acceptability and Order Receipt*
- SOP *Quality Control for Manual Testing Reagents*
- SOP *Labelling for Manual Testing*
- SOP *ABO/Rh Manual Tube Method*
- SOP *Antibody Identification*
- SOP *Grading Reactions*
- SOP *Saline Replacement*
- SOP *Weak D Manual Tube Testing*
- SOP *EGA Treatment of Red Blood Cells*
- SOP *DAT (Direct Antiglobulin Test) by Tube Method*
- SOP *Cold Panel*

APPENDIX: Possible Causes of ABO Typing Discrepancies

TYPE OF DISCREPANCY	POSSIBLE CAUSES
Mixed-field forward type	<ul style="list-style-type: none"> • Recent transfusion • Transplantation • Fetomaternal hemorrhage • Twin or dispermic (tetragametic) chimerism
Extra reaction in the forward type (includes positive control)	<ul style="list-style-type: none"> • Autoagglutinins/excess protein coating red cells • Unwashed red cells: <ul style="list-style-type: none"> ○ plasma proteins ○ antibody in patient's serum to reagent constituent • Transplantation • Acquired B antigen • B(A) phenomenon • Out-of-group transfusion
Negative and/or weak forward type	<ul style="list-style-type: none"> • ABO subgroup • Leukemia/malignancy • Recent transfusion • Intrauterine fetal transfusion • Transplantation • Excessive soluble blood group substance
Extra reaction in the reverse type	<ul style="list-style-type: none"> • Cold autoantibody • Cold alloantibody • Serum antibody to reagent constituent • Excess serum protein • Recent transfusion of plasma components • Transplantation • Infusion of intravenous immune globulin
Negative, weak reverse type	<ul style="list-style-type: none"> • Age related (<4-6 months old, elderly) • ABO subgroup • Hypogammaglobulinemia • Transplantation
Rh does not match historical Type	<ul style="list-style-type: none"> • Weak D antigen • Differences in antisera specificity • Recent transfusion • Transplantation

UWMC SOP Approval:					
UWMC CLIA Medical Director	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> <tr> <td style="padding-left: 20px;">Andrew Bryan, MD</td> <td></td> </tr> </table>		Date	Andrew Bryan, MD	
	Date				
Andrew Bryan, MD					
Transfusion Service Manager	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> <tr> <td style="padding-left: 20px;">Nina Sen</td> <td></td> </tr> </table>		Date	Nina Sen	
	Date				
Nina Sen					
QA Manager	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> <tr> <td style="padding-left: 20px;">Tayler Reeves</td> <td></td> </tr> </table>		Date	Tayler Reeves	
	Date				
Tayler Reeves					
Transfusion Service Medical Director	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> <tr> <td style="padding-left: 20px;">Monica Pagano, MD</td> <td></td> </tr> </table>		Date	Monica Pagano, MD	
	Date				
Monica Pagano, MD					
UWMC Biennial Review:					
	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> <tr> <td style="border-bottom: 1px solid black; width: 60%;"></td> <td style="border-bottom: 1px solid black; width: 40%; text-align: right;">Date</td> </tr> </table>		Date		Date
	Date				
	Date				

REVISIONS:

07/21/2020: Add instructions for resolution of mixed field reactions from cord blood and neonatal specimens