



| | | |
|--|--|-------------------------------------|
| University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA, 98195 Transfusion Services Laboratory Policies and Procedures Manual | Original Effective Date: 06-13-2022 | Number: EQ-0018.01 |
| | Revision Effective Date: | |
| TITLE: Helmer Freezer Operation & Maintenance | | |

PURPOSE:

To provide instructions on the operations and maintenance of the Helmer Freezer

PRINCIPLE & CLINICAL SIGNIFICANCE:

Plasma and cryoprecipitate are stored frozen at -18°C or colder to help maintain factor activity and provide an extended shelf life. Blood storage devices should be continuously monitored to ensure required temperatures are maintained during product storage. In the absence of continuous monitoring, temperatures should be recorded at least every four hours

Clinical Significance

Blood Components must be stored within strict temperature guidelines to ensure the safety, purity and potency of the product

POLICIES:

- Monitored 24/7 by TempTrak temperature monitoring system
- Quarterly and annual maintenance performed by UWMC Montlake Refrigeration shop
- For any maintenance and service issues call Operations and Maintenance
- Daily temperatures are recorded on the **Blood Component Storage Temperatures Form**
- Freezers are set to alarm at the following temperatures

| Alarm | Temperature |
|------------|-------------|
| High Alarm | -20.0°C |
| Low Alarm | -35.0°C |

SPECIMEN REQUIREMENTS:

NA

REAGENTS/SUPPLIES/EQUIPMENT:

| Reagents | Supplies | Equipment |
|--|---|---|
| <ul style="list-style-type: none"> • Propylene Glycol | <ul style="list-style-type: none"> • Calibrated Thermometer • Recording Charts • Distilled Water | <ul style="list-style-type: none"> • Blood Bank Freezers • Remote Alarms and Pagers |

QUALITY CONTROL:

Temperatures are verified daily to ensure the acceptable range $\leq -18^{\circ}\text{C}$ for product storage is maintained.

INSTRUCTIONS:

TABLE of CONTENTS

| | |
|---|-----------------------|
| TITLE: Helmer Freezer Operation and Maintenance | Number: EQ-0018.01 |
|---|-----------------------|

[Periodic Maintenance](#)
[Installing Temperature Chart Paper](#)

Periodic Maintenance

| STEP | ACTION | | |
|-----------|--|--|--|
| 1 | Perform the following maintenance tasks at the required frequency: | | |
| | If performing maintenance | Then | |
| | Quarterly | <ul style="list-style-type: none"> • Check temperature calibration on the monitor. Change if necessary • Clean the condenser grill, and | <ul style="list-style-type: none"> • High and low temperature alarms • Power failure alarms (as required by UWMC's protocols) |
| | | Test | |
| | | Inspect | <ul style="list-style-type: none"> • Electrical components and wiring terminals in the electrical box for discoloration. If discoloration is found, contact Helmer Technical Service • All wiring for terminals for secure connection and tighten as necessary |
| | Annually | Examine the probe bottles. Clean or replace if necessary | |
| As needed | Clean | <ul style="list-style-type: none"> • Freezer interior and exterior • Door gaskets | |
| | Test | <ul style="list-style-type: none"> • Door alarm (as required by UWMC protocols) • Ground fault circuit interrupter on the internal outlet (if applicable) | |
| | Check | <ul style="list-style-type: none"> • Backup battery for the chart recorder after a prolonged power failure and change it if necessary or change the battery if it has been in service for one year. Refer to the Temperature Chart Recorder Operation and Service Manual • Level of the solution in the probe bottles containing a volume of solution no larger than the smallest unit stored in the device. Freezers: Propylene glycol. | |

| | |
|---|-----------------------|
| TITLE: Helmer Freezer Operation and Maintenance | Number: EQ-0018.01 |
|---|-----------------------|

Installing Temperature Chart Paper

| STEP | ACTION |
|------|--|
| 1 | Press and hold "C" button. |
| 2 | Release button when stylus begins to move to the left. NOTE: LED flashes to indicate current temperature range. |
| 3 | Remove chart knob when stylus stops moving. Move the knob up and away |
| 4 | Place new chart paper on the chart recorder |
| 5 | Lift the stylus gently and rotate chart paper so that the current time line corresponds to the time line groove |
| 6 | Hold the chart paper and reinstall the chart knob, verify time line after the chart knob is tightened |
| 7 | Confirm temperature range is set to the correct value |
| 8 | Press and hold "C", releasing the button when the stylus begins to move right |
| 9 | Confirm stylus is marking temperature correctly |

CALCULATIONS/INTERPRETATIONS/RESULTS REPORTING/NORMAL VALUES/CRITICAL VALUES

**CALIBRATION:
Calibration of the Circular graph**

| STEP | ACTION |
|------|--|
| 1 | Determine the direction the pen needs to be moved (left or right) to match the internal temperature. |
| 2 | Press and hold either the left or right arrow buttons until the stylus begins to move. |
| 3 | Release the button when the pen is recording the correct temperature. |

NOTES AND LIMITATIONS:

- Replace rechargeable battery if it has been in service for two years

REFERENCES:

- Helmer Scientific Freezer Operation Manual
- Temperature Chart Recorder Operation and Service Manual
- Technical Manual. Bethesda, MD: AABB Press, current edition

RELATED DOCUMENTS:

FORM Blood Component Storage and Temperatures

APPENDICES:

NA

| | |
|---|-----------------------|
| TITLE: Helmer Freezer Operation and Maintenance | Number: EQ-0018.01 |
|---|-----------------------|

| | | | | | |
|---|---|-------|------------|-------------------|--|
| UWMC SOP Approval: | | | | | |
| UWMC CLIA Medical Director | <table border="0"> <tr> <td data-bbox="454 321 1039 430">_____</td> <td data-bbox="1039 321 1451 430">Date _____</td> </tr> <tr> <td data-bbox="454 430 1039 493">Andrew Bryan, MD</td> <td></td> </tr> </table> | _____ | Date _____ | Andrew Bryan, MD | |
| _____ | Date _____ | | | | |
| Andrew Bryan, MD | | | | | |
| Transfusion Service Manager | <table border="0"> <tr> <td data-bbox="454 493 1039 577">_____</td> <td data-bbox="1039 493 1451 577">Date _____</td> </tr> <tr> <td data-bbox="454 577 1039 640">Nina Sen</td> <td></td> </tr> </table> | _____ | Date _____ | Nina Sen | |
| _____ | Date _____ | | | | |
| Nina Sen | | | | | |
| QA Manager | <table border="0"> <tr> <td data-bbox="454 640 1039 724">_____</td> <td data-bbox="1039 640 1451 724">Date _____</td> </tr> <tr> <td data-bbox="454 724 1039 787">Tayler Reeves</td> <td></td> </tr> </table> | _____ | Date _____ | Tayler Reeves | |
| _____ | Date _____ | | | | |
| Tayler Reeves | | | | | |
| Transfusion Service Medical Director | <table border="0"> <tr> <td data-bbox="454 724 1039 913">_____</td> <td data-bbox="1039 724 1451 913">Date _____</td> </tr> <tr> <td data-bbox="454 913 1039 976">Monica Pagano, MD</td> <td></td> </tr> </table> | _____ | Date _____ | Monica Pagano, MD | |
| _____ | Date _____ | | | | |
| Monica Pagano, MD | | | | | |
| UWMC Biennial Review: | | | | | |
| | <table border="0"> <tr> <td data-bbox="454 955 1039 1050">_____</td> <td data-bbox="1039 955 1451 1050">Date _____</td> </tr> </table> | _____ | Date _____ | | |
| _____ | Date _____ | | | | |
| | <table border="0"> <tr> <td data-bbox="454 1050 1039 1150">_____</td> <td data-bbox="1039 1050 1451 1150">Date _____</td> </tr> </table> | _____ | Date _____ | | |
| _____ | Date _____ | | | | |