



University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA, 98195 Transfusion Services Laboratory Policies and Procedures Manual	Original Effective Date: 8/15/2022	Number: PE-0001.00
	Revision Effective Date:	
TITLE: Training and Competency Assessment		

PURPOSE:

To describe the process for initial training, ongoing training, and assessment of competency of both testing personnel and personnel performing critical tasks. Prior to hire, employee's qualifications and credentials will be assessed by University of Washington (UW) Human Resources.

POLICIES:

Training:

- Initial training is conducted prior to a new employee reporting patient test results and, where applicable, encompasses pre-analytical (specimen collection, receipt, and processing), analytical (test performance) and post-analytical (result evaluation and reporting) phases of testing.
- Training may include but is not limited to the following:
 1. Reading relevant policies and procedures
 2. Observation of proper collection, labeling and processing of specimens
 3. Observation of test and/or task performance
 4. Hands-on training on the use of the computer system (proper accessioning, resulting and inquiry)
 5. Hands-on training on the use of equipment
 6. Hands-on training with patient samples, quality control samples, proficiency samples or blind samples
- As new instruments and/or analytical procedures are introduced, personnel will receive documented training prior to performing critical tasks, and/or tests. All ongoing training will be documented in each employee's training folder.
- Additionally, training will be provided and documented as part of the corrective action process in the event of failed proficiency testing or if deemed necessary after a non-conforming event.

Competency Assessment:

- After testing personnel have been trained and are performing tasks independently, competency is assessed semi-annually i.e., at 6-months and 1-year after training is completed. After the first year, a reassessment of competency is annual and includes all testing personnel.
 - Critical Task competency will be assessed before the employee performs tasks independently and will be reassessed annually.
 - Semi-annual assessment is not required for current employees trained on new methods or instruments.
- Retraining and reassessment of competency must occur if the competency is not satisfactory. Additionally, if testing personnel fail to demonstrate satisfactory performance on the competency assessment, a corrective action plan will be used to retrain and reassess competency.

- Competency assessments for moderate and high complexity tests will include six assessment elements for each test system and for each employee.
- For waived test systems and critical tasks, only the competency assessment elements applicable to the process will be used.
- Competency assessment elements include, but are not limited to:
 1. Direct observation (DO) of routine patient test performance
 2. Monitoring the recording and reporting of test results.
 3. Review of intermediate test results/worksheets, quality control records, proficiency testing results and preventive maintenance records.
 4. Direct observations of performance of instrument maintenance and function checks.
 5. Assessment of test performance through testing previously analyzed specimens, internal blind testing samples, or external proficiency testing samples.
 6. Assessment of problem-solving skills.
- **Direct Observation:**
 - A Direct Observation task list for each test system is outlined in Appendix A, if a step is not performed or is performed incorrectly during a direct observation, the step number designated in this policy will be documented on the appropriate area of the competency assessment form.
 - Personnel performing direct observations for high complexity testing must meet the CAP (College of American Pathology) requirements for general supervisor and are delegated in writing by the section director (per CAP GEN.55500).

Performance Assessment of Managers:

- Responsibilities of Managers (Technical Supervisors, Technical consultants, General supervisors) are delegated in writing per Policy ***TSL Medical Director Designation of Record Review***
- Performance assessment will occur annually in line with UW performance evaluation season.
- Unsatisfactory performance will be addressed in a corrective action plan or performance improvement plan.
- Individuals performing non-waived patient testing will complete competency assessment requirements for testing personnel.

INSTRUCTIONS:

Training:

- New employees are to complete the Initial Training Checklist appropriate for their role.
- The date documented for “Date Trained” on the Initial Training Checklist should indicate the date that training is completed, and the employee is ready for Competency Assessment.
- If an employee will not be trained on an item listed on the training checklist, **NPBE** (Not Performed by Employee) will be documented in the “Date Trained” field for that task.
- It is the trainee’s responsibility to communicate which tasks remain on their checklist with the trainer to ensure that all checklist items are completed.

Competency Assessment:

- See **Appendix A** for Direct Observation Critical Step Lists

- If a task/test will not be performed by an employee, it will be documented on the competency assessment form by answering “No” to the question “Is this task performed by this employee?” AND indicating “C” in the “Competency Level” box.
- **Initial Competency Assessment** will require the minimum number of blind testing samples outlined below:

Tests	Minimum # Blind Samples	Comments
Ortho Vision	1	
ABO/Rh	5	
Antibody Screen	5	
Crossmatch	5	(At least 1 of each of the following: <ul style="list-style-type: none"> • Electronic Crossmatch • Immediate Spin Crossmatch • LISS Crossmatch • PEG Crossmatch
Direct Antiglobulin Test (DAT)	5	At least 2 positives
Antibody Identification	5	
Antigen Typing	1	Patient or Donor sample
Fetal Maternal Hemorrhage RapidScreen	1	
Sicklelex Hemoglobin S	1	
Testing with DTT treated cells	1	
EGA	1	
Eluate	1	

- **Subsequent (Semi-annual and Annual) Competency Assessment** will require a single blind sample for each test listed above. CAP survey specimens or selected samples may be used.

REFERENCES:

AABB Standards for Blood Banks and Transfusion Services, AABB Press
 CAP Checklists- General, All Common, Transfusion Medicine
 AABB Technical Manual, AABB Press

RELATED DOCUMENTS:

- SOP TSL Medical Director Designation of Record Review (OR-0002.02)*
- Form Initial Training Checklist - Clinical Laboratory Technician (TC-0006.00)*
- Form Initial Training Checklist - Medical Laboratory Scientist (TC-0007.00)*
- Form Transfusion Services Critical Task Competency Assessment (TC-0004.00)*
- Form Transfusion Services Testing Personnel Competency Assessment (TC-0005.00)*
- Form ABORh and Antibody Screen Initial Competency Assessment Blind Testing (TC-0001.00)*
- Form Crossmatch Initial Competency Assessment Blind Testing (TC-0002.00)*
- Form DAT Initial Competency Assessment Blind Testing (TC-0003.00)*

UWMC SOP Approval:	
UWMC CLIA	
Medical Director	Date
Andrew Bryan, MD	
Transfusion Service Manager	Date
Nina Sen	
QA Manager	Date
Tayler Reeves	
Transfusion Service Medical Director	Date
Monica Pagano, MD	
UWMC Biennial Review:	
_____	Date _____
_____	Date _____

Appendix A – Direct Observation Critical Steps**Table of Contents:****Critical Task Direct Observations**[Specimen Receipt and Processing – Interfaced Orders](#)[Issuing Blood Components](#)[Irradiating Blood Components](#)[Thawing Frozen Blood Components](#)[Combining Double Bagged Apheresis Platelets](#)[Dividing Blood Components](#)[Volume Reduced Platelets](#)[Washing Components \(Platelet or RBC\)](#)[Preparing RBCs for Intrauterine Transfusion \(IUT\)](#)**Test System Direct Observations**[Ortho Vision Patient Testing](#)[Ortho Vision Donor Confirmation Testing](#)[Ortho Vision Daily Maintenance](#)[ABO/RH Tube Testing](#)[Antibody Screen \(PEG\)](#)[Antibody Screen \(LISS\)](#)[Heat Block Daily Maintenance](#)[Cell Washer Daily Maintenance](#)[DAT in Tube](#)[Electronic Crossmatch](#)[Immediate Spin Crossmatch](#)[AHG Crossmatch \(PEG or LISS\)](#)[Antigen Typing](#)[Fetal Maternal hemorrhage RapidScreen](#)[SickleDex Hemoglobin S](#)[Testing with DTT Treated Cells](#)[EGA Testing](#)[Eluate Testing](#)

Specimen Receipt and Processing – Interfaced Orders							
1	Confirm specimen is: <ul style="list-style-type: none"> Collected in correct container Labeled with Name, MRN, draw date, draw time, and 2 signatures 						
2	Open Blood Bank Inquiry						
3	Select Lookup by 'Patient ID' and enter the patient medical record number (MRN)						
4	Determine if the patient has multiple HID <table border="1" style="width: 100%; margin-top: 5px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 30%; padding: 2px;">If the Patient has</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Only a "U" HID</td> <td style="padding: 2px;">Go to the next step</td> </tr> <tr> <td style="padding: 2px;">"H" and "U" HID</td> <td style="padding: 2px;">Select the U HID</td> </tr> </tbody> </table>	If the Patient has	Then	Only a "U" HID	Go to the next step	"H" and "U" HID	Select the U HID
If the Patient has	Then						
Only a "U" HID	Go to the next step						
"H" and "U" HID	Select the U HID						
5	Open General Laboratory in Sunquest						
6	Click on 'Orders' and select <Orders Receipt/Modify>						
7	Enter/Scan the CID or MRN from the specimen						
8	Click <Get Patient>						
9	Click <Display Orders> NOTE: It might be necessary to change the Day(s) of Activity						
10	Update/enter the collection information <table border="1" style="width: 100%; margin-top: 5px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 40%; padding: 2px;">If</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Collection time on the tube is different from SQ</td> <td style="padding: 2px;">Update SQ to match the tube</td> </tr> <tr> <td style="padding: 2px;">There is a Tech ID for the Phlebotomist</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Enter: VP in the Order Workload Code field Enter the Tech ID in the Phlebotomist Code field </td> </tr> </tbody> </table>	If	Then	Collection time on the tube is different from SQ	Update SQ to match the tube	There is a Tech ID for the Phlebotomist	<ul style="list-style-type: none"> Enter: VP in the Order Workload Code field Enter the Tech ID in the Phlebotomist Code field
If	Then						
Collection time on the tube is different from SQ	Update SQ to match the tube						
There is a Tech ID for the Phlebotomist	<ul style="list-style-type: none"> Enter: VP in the Order Workload Code field Enter the Tech ID in the Phlebotomist Code field 						
11	<table border="1" style="width: 100%; margin-top: 5px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 25%; padding: 2px;">If the sample has a(n)</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">SQ label</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Click <Receive> Verify checkmark X is next to the CID number that is the same as the labeled specimen Click <SAVE> </td> </tr> <tr> <td style="padding: 2px;">EPIC ADT label</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Update the container type, if necessary Click <Route> Click <SAVE> when Result Entry Box appears. </td> </tr> </tbody> </table>	If the sample has a(n)	Then	SQ label	<ul style="list-style-type: none"> Click <Receive> Verify checkmark X is next to the CID number that is the same as the labeled specimen Click <SAVE> 	EPIC ADT label	<ul style="list-style-type: none"> Update the container type, if necessary Click <Route> Click <SAVE> when Result Entry Box appears.
If the sample has a(n)	Then						
SQ label	<ul style="list-style-type: none"> Click <Receive> Verify checkmark X is next to the CID number that is the same as the labeled specimen Click <SAVE> 						
EPIC ADT label	<ul style="list-style-type: none"> Update the container type, if necessary Click <Route> Click <SAVE> when Result Entry Box appears. 						
12	<table border="1" style="width: 100%; margin-top: 5px;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 45%; padding: 2px;">If the Result Entry Box</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Has a units ordered field</td> <td style="padding: 2px;">Enter the units ordered</td> </tr> <tr> <td style="padding: 2px;">Does not have a units ordered field</td> <td style="padding: 2px;">Continue to the Next Step</td> </tr> </tbody> </table>	If the Result Entry Box	Then	Has a units ordered field	Enter the units ordered	Does not have a units ordered field	Continue to the Next Step
If the Result Entry Box	Then						
Has a units ordered field	Enter the units ordered						
Does not have a units ordered field	Continue to the Next Step						
13	Click <Save>						

[Return to Table of Contents](#)

Issuing Blood Components							
1	Timestamp the Blood Product Release Form						
2	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 35%; padding: 5px;">If issuing a unit that was</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Previously allocated</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Open Blood Product Issue Scan or enter the MRN from the Blood Product Release Form Input component group based on Blood Product Release Form Choose the appropriate unit to issue from storage location </td> </tr> <tr> <td style="padding: 5px;">Not Previously allocated</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Open Blood Order Processing Scan or enter the MRN from the Blood Product Release Form Allocate a unit Click <SAVE> Click <ISSUE> </td> </tr> </tbody> </table>	If issuing a unit that was	Then	Previously allocated	<ul style="list-style-type: none"> Open Blood Product Issue Scan or enter the MRN from the Blood Product Release Form Input component group based on Blood Product Release Form Choose the appropriate unit to issue from storage location 	Not Previously allocated	<ul style="list-style-type: none"> Open Blood Order Processing Scan or enter the MRN from the Blood Product Release Form Allocate a unit Click <SAVE> Click <ISSUE>
If issuing a unit that was	Then						
Previously allocated	<ul style="list-style-type: none"> Open Blood Product Issue Scan or enter the MRN from the Blood Product Release Form Input component group based on Blood Product Release Form Choose the appropriate unit to issue from storage location 						
Not Previously allocated	<ul style="list-style-type: none"> Open Blood Order Processing Scan or enter the MRN from the Blood Product Release Form Allocate a unit Click <SAVE> Click <ISSUE> 						
3	Time stamp or manually record date and time the Transfusion Record						
4	Verify unit meets all the patient transfusion requirements Note: All red cells, granulocytes and non pathogen reduced platelets must be irradiated						
5	Scan unit and product code in Sunquest						
6	Perform visual inspection for the following: <ul style="list-style-type: none"> Expiration date has not passed Correct Labeling Intact Container No Clots, turbidity, hemolysis or other abnormal appearance of the component						
7	Document in Sunquest						
8	Verify patient and component information matches on the transfusion record, Blood Product Release Form, ISBT label, and Sunquest						
9	Click <Continue>						
10	Complete issue in Sunquest						
11	Click <Save> and click <Cancel> in the Add Billing window						
12	Give unit to second staff member to perform clerical check of Unit label and Transfusion Record						

[Return to Table of Contents](#)

Irradiating Blood Components											
1	Date and initial a Rad-Sure indicator label and adhere to component without obliterating other required labeling elements										
2	Verify key is in “cycle” mode on the irradiator										
3	Open the RSTScan program										
4	Scan the following into the appropriate field on the tablet <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 40%; padding: 5px;">Field</th> <th style="padding: 5px;">Scan</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">User ID</td> <td style="padding: 5px;">User ID badge (this is only scanned once for each batch of components)</td> </tr> <tr> <td style="padding: 5px;">Indicator Batch ID</td> <td style="padding: 5px;">Rad-Sure XR 25 Gy Indicator label Batch ID (Lot number)</td> </tr> <tr> <td style="padding: 5px;">Product Code</td> <td style="padding: 5px;">Product Code</td> </tr> <tr> <td style="padding: 5px;">Donor ID</td> <td style="padding: 5px;">Donor Unit Number</td> </tr> </tbody> </table>	Field	Scan	User ID	User ID badge (this is only scanned once for each batch of components)	Indicator Batch ID	Rad-Sure XR 25 Gy Indicator label Batch ID (Lot number)	Product Code	Product Code	Donor ID	Donor Unit Number
Field	Scan										
User ID	User ID badge (this is only scanned once for each batch of components)										
Indicator Batch ID	Rad-Sure XR 25 Gy Indicator label Batch ID (Lot number)										
Product Code	Product Code										
Donor ID	Donor Unit Number										
5	Select ‘Add’ to populate the chart below the data fields for each unit added to the batch										
6	Scan all units that will be irradiated in one batch										
7	Load the unit into the canister										
8	Load canister into canister holder with lid facing out										
9	Close the chamber door										
10	Start the cycle by pressing the “Start” button										
11	Press “Door Release” when buzzer sounds										
12	Remove product and close chamber door										
13	Enter Rad-Sure indicator result into RSTScan program										
14	Save data										
15	Use Blood Component Preparation in Sunquest and choose correct function based on the product’s E code										
16	Enter the correct new expiration date/time in Sunquest (28 days from date of Irradiation or original expiration if shorter)										
17	Place new label on product										
18	Verify that the US license number is obscured or mark out the license number										
19	Perform Blood Label Check in Sunquest										
20	Place product in correct storage location in a timely manor										

[Return to Table of Contents](#)

Thawing Frozen Blood Components											
1	Select an appropriate unit to thaw										
2	Examine the product for defects and ensure unit is not expired										
3	Place the product in an overwrap bag, ports up										
4	Check water level and temperature is between 30-37°C										
5	Place the bag in the thaw bath basket attaching the overwrap bag to the basket										
6	<div style="border: 1px solid black; padding: 5px;"> Select an appropriate thawing time for the selected product following these guidelines: <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 40%; padding: 5px;">If Component is</th> <th style="padding: 5px;">Median Thaw Time (minutes) is</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">≈ 10-15 mL cryoprecipitate</td> <td style="text-align: center; padding: 5px;">5</td> </tr> <tr> <td style="padding: 5px;">Pooled cryoprecipitate ≈ 250 mL</td> <td style="text-align: center; padding: 5px;">8</td> </tr> <tr> <td style="padding: 5px;">≈ 250 mL plasma</td> <td style="text-align: center; padding: 5px;">10</td> </tr> <tr> <td style="padding: 5px;">≈ 300 mL plasma</td> <td style="text-align: center; padding: 5px;">14</td> </tr> </tbody> </table> </div>	If Component is	Median Thaw Time (minutes) is	≈ 10-15 mL cryoprecipitate	5	Pooled cryoprecipitate ≈ 250 mL	8	≈ 250 mL plasma	10	≈ 300 mL plasma	14
If Component is	Median Thaw Time (minutes) is										
≈ 10-15 mL cryoprecipitate	5										
Pooled cryoprecipitate ≈ 250 mL	8										
≈ 250 mL plasma	10										
≈ 300 mL plasma	14										
7	Start cycle										
8	Examine product at the end of the cycle, verify no leaks and unit thawed correctly										
9	Use Blood Component Preparation in Sunquest and choose correct function based on the product's E code										
10	Verify the expiration date and time is accurate										
11	Verify the label information is correct and place new label on the unit										
12	Perform Blood Label Check in Sunquest										

Combining Double Bagged Apheresis Platelets	
1	Select the bag with the complete labeling (ISBT label including attached tags and additional stickers if applicable) as the final product container NOTE: Labeling of double bag products may vary by blood supplier
2	Hold or hang the bag to be discarded higher than the bag selected as the final product container
3	Open any slider valve or clamp between the bags
4	Allow the upper bag to drain into the final container by gravity
5	Heat seal the tubing between the two bags
6	Modify the expiration date/time to 24 hours or the original expiration date/time whichever is shorter in the BB LIS

[Return to Table of Contents](#)

Dividing Blood Components							
1	Mix the blood product prior to dividing						
2	Tighten the tubing on the syringe, remove spike and clip						
3	Check the weld count on the Sterile Welder						
4	Follow directions on Sterile Welder						
5	Remove tubing and discard stubs						
6	Inspect Weld for leaks and check integrity of weld by gently pulling and squeezing tubing						
7	Press Reset on the sterile welder						
8	Draw required volume plus 3-5 ml of blood into the syringe						
9	Seal the tubing						
10	Perform modification in Sunquest						
11	Verify new label matches the original unit expiration and product code and adhere label to the unit						
12	Verify large label unit number matches the original unit, the unit expiration date and product code are correct and place on syringe						
13	Perform Blood Label Check in Sunquest						
14	Cut the syringe from the bag carefully						
15	Aseptically remove the stub end from the syringe and replace with a sterile cap						
16	Irradiate the blood product						
17	Perform Irradiation in Sunquest- Label and Label Check unit						
18	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 25%; padding: 5px;">If employee is</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">MLS</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Allocate unit to correct patient Add LOTNO and WELD to unit testing Result testing </td> </tr> <tr> <td style="padding: 5px;">CLT</td> <td style="padding: 5px;">Continue to the next step</td> </tr> </tbody> </table>	If employee is	Then	MLS	<ul style="list-style-type: none"> Allocate unit to correct patient Add LOTNO and WELD to unit testing Result testing 	CLT	Continue to the next step
If employee is	Then						
MLS	<ul style="list-style-type: none"> Allocate unit to correct patient Add LOTNO and WELD to unit testing Result testing 						
CLT	Continue to the next step						
19	Replace units in correct storage locations						

[Return to Table of Contents](#)

Volume Reduced Platelets							
1	Use the SCD Rapidweld to add a transfer pack of the appropriate size to the unit prior to centrifugation (if necessary)						
2	Inspect Weld for leaks and check integrity of weld by gently pulling and squeezing tubing						
3	Clamp the tubing between the product bag and transfer pack prior to centrifugation						
4	Balance the centrifuge cups prior to loading the centrifuge						
5	Load the product with ports up and product label facing outward						
6	Select the correct program for the modification being performed and verify temperature requirements are met						
7	Document centrifuge QC on the Refrigerated Centrifuge QC log						
8	Remove product without disturbing the platelet pellet						
9	Hang or insert Platelet unit into plasma extractor						
10	Tare scale with empty bag and remove all but 100 ml from the component or alternate volume if specifically requested						
11	Heat seal off waste bag and perform modification in Sunquest						
12	Use Blood Component Preparation in Sunquest and choose correct function based on the product's E code						
13	Use the date and time that plasma expression began as the process date and time						
14	Enter the new volume and verify the expiration date and time is accurate						
15	Verify the label information is correct and place new label on the unit						
16	Perform Blood Label Check in Sunquest						
17	Discard waste bag						
18	Allow component to rest undisturbed for 20 minutes and resuspend the component						
19	<div style="border: 1px solid black; padding: 5px;"> Place component on platelet agitator for 20 minutes and inspect <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="text-align: left; padding: 2px;">If component is</th> <th style="text-align: left; padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Acceptable</td> <td style="padding: 2px;">Continue to the next step</td> </tr> <tr> <td style="padding: 2px;">Not acceptable</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Break up large aggregates with gentle finger pressure Agitate for an additional 20 minutes </td> </tr> </tbody> </table> </div>	If component is	Then	Acceptable	Continue to the next step	Not acceptable	<ul style="list-style-type: none"> Break up large aggregates with gentle finger pressure Agitate for an additional 20 minutes
If component is	Then						
Acceptable	Continue to the next step						
Not acceptable	<ul style="list-style-type: none"> Break up large aggregates with gentle finger pressure Agitate for an additional 20 minutes 						
20	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="text-align: left; padding: 2px;">If the employee is</th> <th style="text-align: left; padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">MLS</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add 'WELD' and 'LOTNO' to the unit testing Enter the test results </td> </tr> <tr> <td style="padding: 2px;">CLT</td> <td style="padding: 2px;">Take unit to MLS for further processing in Sunquest</td> </tr> </tbody> </table>	If the employee is	Then	MLS	<ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add 'WELD' and 'LOTNO' to the unit testing Enter the test results 	CLT	Take unit to MLS for further processing in Sunquest
If the employee is	Then						
MLS	<ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add 'WELD' and 'LOTNO' to the unit testing Enter the test results 						
CLT	Take unit to MLS for further processing in Sunquest						

[Return to Table of Contents](#)

Washing Components (Platelet or RBC)							
Check Prime							
1	Install priming cushion						
2	Replace centrifuge bowl cover						
3	Ready light illuminates when powered on						
4	Verify Excessive Pressure light illuminates and COBE alarms within 15 seconds						
Installing Cell Processing Set							
5	Inspect cell processing set for damage, kinks, or missing caps						
6	Inspect washing solution for leaks or open ports						
	If washing a	Then inspects					
	Platelet	Plasma-lyte					
	RBC	Saline					
7	Place unit number on the hexagonal seal						
8	Place junction manifold						
	If washing a	Then place the Junction Manifold about 1 inch above the Red Cell Detector (RCD) and the tubing is:					
	Platelet	<ul style="list-style-type: none"> • Left outside of the RCD 					
	RBC	<ul style="list-style-type: none"> • Inserted fully into the back of the slot in the RCD 					
9	Place all tubing from cell processing set into the correct valve						
	Tubing	Valve					
	Red striped	V1					
	Purple striped	SOV					
	Green striped	2					
	Yellow striped	3					
10	Blue and yellow striped tubing are either clamped or heat sealed						
11	Install cell processing set bag flat with spike at an angle						
12	Position white alignment blocks flat side up and replace centrifuge bowl cover						
13	Align rotating seal so that a point of rotating seal is pointed to the front of the machine and is at the rear of the seal weight						
Attaching Component and Washing Solution							
14	Clamp hemostat on clear tubing between RCD and rotating seal						
15	Use aseptic technique to remove caps from spikes and insert: <ul style="list-style-type: none"> • green tubing into the washing solution red tubing into the component						
16	Record new expiration time of washed component						
Predilution and Loading of Component							
17	If washing a	Then set the LED program board to:					
	Platelet	2					
	RBC	1					
18	Verify controls on main control panel are in correct position and adjust if necessary						
	Component	Centrifuge Speed	Super Out Rate	Min. Agitate Time	SUPER OUT Volume	Auto/Manual	Collect/SOV
	Platelet	3000 rpm	100 mL/min	30 sec	Platelet Vol. +50mL	Auto	SOV
	RBC	3000 rpm	450 mL/min	60 sec	600 mL	Auto	SOV

CONTINUED - Predilution and Loading of Component									
19	Allow approximately 100mls of washing solution into the component								
20	Unclamp hemostat between RCD and hexagonal seal.								
21	Press <BLOOD IN> and then <AIR OUT> to remove air bubbles from tubing								
22	Press <BLOOD IN> and then <STOP RESET> before all the platelets drain into the hexagonal seal								
Processing and Removal of Processing Set									
23	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="2" style="padding: 5px;">Press <START SPIN> to start wash process</td> </tr> <tr style="background-color: #e0e0e0;"> <th style="width: 25%; padding: 5px;">If washing a</th> <th style="padding: 5px;">Then</th> </tr> <tr> <td style="padding: 5px;">Platelet</td> <td style="padding: 5px;">Change SUPER OUT volume to 400 after first cycle</td> </tr> <tr> <td style="padding: 5px;">RBC</td> <td style="padding: 5px;">Continue to next step</td> </tr> </table>	Press <START SPIN> to start wash process		If washing a	Then	Platelet	Change SUPER OUT volume to 400 after first cycle	RBC	Continue to next step
Press <START SPIN> to start wash process									
If washing a	Then								
Platelet	Change SUPER OUT volume to 400 after first cycle								
RBC	Continue to next step								
24	Press <STOP RESET> when the audible alarm sounds								
25	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr style="background-color: #e0e0e0;"> <th style="width: 25%; padding: 5px;">If washing a</th> <th style="padding: 5px;">Then</th> </tr> <tr> <td style="padding: 5px;">Platelet</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Press PREDILUTE and then STOP to allow 50-100 mL of plasma-lyte to drain into donut Allow platelet to rest undisturbed for 30 minutes </td> </tr> <tr> <td style="padding: 5px;">RBC</td> <td style="padding: 5px;">Perform visual inspection of waste line for hemolysis and take necessary steps if hemolyzed</td> </tr> </table>	If washing a	Then	Platelet	<ul style="list-style-type: none"> Press PREDILUTE and then STOP to allow 50-100 mL of plasma-lyte to drain into donut Allow platelet to rest undisturbed for 30 minutes 	RBC	Perform visual inspection of waste line for hemolysis and take necessary steps if hemolyzed		
If washing a	Then								
Platelet	<ul style="list-style-type: none"> Press PREDILUTE and then STOP to allow 50-100 mL of plasma-lyte to drain into donut Allow platelet to rest undisturbed for 30 minutes 								
RBC	Perform visual inspection of waste line for hemolysis and take necessary steps if hemolyzed								
26	Replace hemostat in clear tubing near RCD								
27	Make 3 segments with tube sealer								
28	Remove cell processing set from COBE disposing of waste								
29	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr style="background-color: #e0e0e0;"> <th style="width: 25%; padding: 5px;">If washing a</th> <th style="padding: 5px;">Then</th> </tr> <tr> <td style="padding: 5px;">Platelet</td> <td style="padding: 5px;">Gently massage platelet back into suspension and then allow it to rock for at least 20 minutes</td> </tr> <tr> <td style="padding: 5px;">RBC</td> <td style="padding: 5px;">Continue to next step</td> </tr> </table>	If washing a	Then	Platelet	Gently massage platelet back into suspension and then allow it to rock for at least 20 minutes	RBC	Continue to next step		
If washing a	Then								
Platelet	Gently massage platelet back into suspension and then allow it to rock for at least 20 minutes								
RBC	Continue to next step								
30	Perform electronic processing in Sunquest								

[Return to Table of Contents](#)

Preparing RBCs for Intrauterine Transfusion (IUT)							
1	Use the SCD Rapidweld to add a transfer pack to the unit prior to centrifugation						
2	Inspect Weld for leaks and check integrity of weld by gently pulling or squeezing tubing						
3	Clamp the tubing between the product bag and transfer pack prior to centrifugation						
4	Balance the centrifuge cups prior to loading the centrifuge						
5	Load the product with ports up and product label facing outward						
6	Select the correct program for the modification being performed and verify temperature requirements are met						
7	Document centrifuge QC on the Refrigerated Centrifuge QC log						
8	Remove product without disturbing the RBCs						
9	Hang the RBC unit on the plasma extractor						
10	Release the plasma extractor plate handle and unclamp the tubing allowing the supernatant to express into the transfer pack						
11	Stop expression when the RBCs reach the top of the bag						
12	Heat seal off waste bag and discard leaving sufficient tubing attached for HCT sampling						
13	Strip and mix the tubing 4 times						
14	Seal the stripped tubing and remove a segment for HCT verification by Hematology						
15	Label test tube and submit sample for testing						
16	Perform processing in Sunquest using the date and time that supernatant expression began as the process date and time						
17	Perform Blood Label Check in Sunquest						
18	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 25%; padding: 5px;">If the employee is a</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">MLS</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add the WELD and LOTNO test in the unit testing Enter the test results </td> </tr> <tr> <td style="padding: 5px;">CLT</td> <td style="padding: 5px;">Take unit to MLS for further processing in Sunquest</td> </tr> </tbody> </table>	If the employee is a	Then	MLS	<ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add the WELD and LOTNO test in the unit testing Enter the test results 	CLT	Take unit to MLS for further processing in Sunquest
If the employee is a	Then						
MLS	<ul style="list-style-type: none"> Allocate unit to the patient in Blood Order Processing Add the WELD and LOTNO test in the unit testing Enter the test results 						
CLT	Take unit to MLS for further processing in Sunquest						

[Return to Table of Contents](#)

Ortho Vision Patient Testing															
1	Select a ring position into which to load samples														
2	Touch <Load/Unload> and open the door														
3	<table border="1"> <tr> <td colspan="2">Load the sample rack</td> </tr> <tr> <td style="text-align: center;">If loading</td> <td style="text-align: center;">Then</td> </tr> <tr> <td>A single rack</td> <td>Proceed to next step</td> </tr> <tr> <td>Multiple racks</td> <td> <ul style="list-style-type: none"> • Select additional ring positions to load • Load sample rack </td> </tr> </table>	Load the sample rack		If loading	Then	A single rack	Proceed to next step	Multiple racks	<ul style="list-style-type: none"> • Select additional ring positions to load • Load sample rack 						
Load the sample rack															
If loading	Then														
A single rack	Proceed to next step														
Multiple racks	<ul style="list-style-type: none"> • Select additional ring positions to load • Load sample rack 														
4	<table border="1"> <tr> <td colspan="2">Close the Load Station Door</td> </tr> <tr> <td style="text-align: center;">If system</td> <td style="text-align: center;">Then</td> </tr> <tr> <td>Downloads orders from the LIS</td> <td>The Vision automatically scans sample ID and starts running the assay</td> </tr> <tr> <td rowspan="2" style="vertical-align: top;">Does not download from LIS</td> <td> <table border="1"> <tr> <td style="text-align: center;">If creating:</td> <td style="text-align: center;">Then:</td> </tr> <tr> <td>Order for a single sample</td> <td> <ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> </td> </tr> <tr> <td>Order with the same profile for multiple samples</td> <td> <ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start> </td> </tr> </table> </td> </tr> </table>	Close the Load Station Door		If system	Then	Downloads orders from the LIS	The Vision automatically scans sample ID and starts running the assay	Does not download from LIS	<table border="1"> <tr> <td style="text-align: center;">If creating:</td> <td style="text-align: center;">Then:</td> </tr> <tr> <td>Order for a single sample</td> <td> <ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> </td> </tr> <tr> <td>Order with the same profile for multiple samples</td> <td> <ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start> </td> </tr> </table>	If creating:	Then:	Order for a single sample	<ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> 	Order with the same profile for multiple samples	<ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start>
Close the Load Station Door															
If system	Then														
Downloads orders from the LIS	The Vision automatically scans sample ID and starts running the assay														
Does not download from LIS	<table border="1"> <tr> <td style="text-align: center;">If creating:</td> <td style="text-align: center;">Then:</td> </tr> <tr> <td>Order for a single sample</td> <td> <ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> </td> </tr> <tr> <td>Order with the same profile for multiple samples</td> <td> <ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start> </td> </tr> </table>	If creating:	Then:	Order for a single sample	<ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> 	Order with the same profile for multiple samples	<ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start> 								
	If creating:	Then:													
Order for a single sample	<ul style="list-style-type: none"> • Touch the sample in yellow • Touch <Create Order> • Fill in the required details • Touch <Save and Start> 														
Order with the same profile for multiple samples	<ul style="list-style-type: none"> • Touch the samples • Touch <Batch Order> • Touch <Sample ID> and select sample IDs • Fill in the required details for the assay to be run • Touch <Save and Start> 														
5	<table border="1"> <tr> <td colspan="2">Touch Results menu button</td> </tr> <tr> <td style="text-align: center;">If</td> <td style="text-align: center;">Then</td> </tr> <tr> <td>Results have been completed and sent to LIS automatically</td> <td>Go to next step</td> </tr> <tr> <td>Results have been flagged for review</td> <td> <ul style="list-style-type: none"> • Retrieve card if available • Select the test result from the Results menu • Touch Show details action button • Edit gel card column grade and/or test interpretation • Accept or Reject test result as applicable • Send Accepted results to LIS </td> </tr> </table>	Touch Results menu button		If	Then	Results have been completed and sent to LIS automatically	Go to next step	Results have been flagged for review	<ul style="list-style-type: none"> • Retrieve card if available • Select the test result from the Results menu • Touch Show details action button • Edit gel card column grade and/or test interpretation • Accept or Reject test result as applicable • Send Accepted results to LIS 						
Touch Results menu button															
If	Then														
Results have been completed and sent to LIS automatically	Go to next step														
Results have been flagged for review	<ul style="list-style-type: none"> • Retrieve card if available • Select the test result from the Results menu • Touch Show details action button • Edit gel card column grade and/or test interpretation • Accept or Reject test result as applicable • Send Accepted results to LIS 														
6	Open Blood Order Processing in Sunquest														
7	Scan the CID of the sample to be result														
8	Load on the 'On-line Results Available' window														
9	Save all appropriate results - Click <Save>														

[Return to Table of Contents](#)

Ortho Vision Donor Confirmation Testing							
1		Drain segments into labeled 10X75mm tube					
2		Load samples into appropriate rack with barcode labels facing out					
3		Load sample rack onto the vision					
4		Touch <Samples>					
5		Touch <Batch Order>					
6		Select samples IDs					
7		Select Donor Rh Pos or Donor Rh Neg test profile as appropriate					
8		Fill in the required details					
9		Touch <Save and Start>					
10		Open 'Blood Bank Instruments' in Sunquest					
11		Select the appropriate analyzer and click ok					
12	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 25%; padding: 5px;">If result is</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Not discrepant</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Check the box in the Release column for each component with acceptable results Select "Release batch" button Select OK </td> </tr> <tr> <td style="padding: 5px;">Discrepant</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Do not check the box Do not release the unit into available inventory Clear the interface results </td> </tr> </tbody> </table>	If result is	Then	Not discrepant	<ul style="list-style-type: none"> Check the box in the Release column for each component with acceptable results Select "Release batch" button Select OK 	Discrepant	<ul style="list-style-type: none"> Do not check the box Do not release the unit into available inventory Clear the interface results
If result is	Then						
Not discrepant	<ul style="list-style-type: none"> Check the box in the Release column for each component with acceptable results Select "Release batch" button Select OK 						
Discrepant	<ul style="list-style-type: none"> Do not check the box Do not release the unit into available inventory Clear the interface results 						

Ortho Vision Daily Maintenance							
1		Execute Daily Probe Maintenance and follow the on-screen prompts					
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 50%; padding: 5px;">If on screen prompt to maintain waste/saline/DI</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Appears</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Empty waste Refill saline Refill DI Water </td> </tr> <tr> <td style="padding: 5px;">Does not appear</td> <td style="padding: 5px;">Continue to next step</td> </tr> </tbody> </table>	If on screen prompt to maintain waste/saline/DI	Then	Appears	<ul style="list-style-type: none"> Empty waste Refill saline Refill DI Water 	Does not appear	Continue to next step
If on screen prompt to maintain waste/saline/DI	Then						
Appears	<ul style="list-style-type: none"> Empty waste Refill saline Refill DI Water 						
Does not appear	Continue to next step						
2		Add 5 mL of 0.1 NaOH to a 10 mL vial with a supported barcode and place into position 3 of a Diluent Rack					
3		Place a new 5 mL vial of 7% BSA into position 2					
4		Load the Diluent Rack and follow the on-screen prompts					
5		Open the maintenance door and follow the on-screen prompts					
6		Clean the probe with Kimwipe moistened with 70% Isopropyl alcohol and follow the on-screen prompts					
7		Close the maintenance door and follow the on-screen prompts					
8		Open Load Station Door and follow the on-screen prompts					
9		Remove Diluent Rack and follow the on-screen prompts					
10		Close Load Station Door and follow the on-screen prompts					
11		Print Maintenance History Report and Complete Ortho Vision Maintenance Log					

[Return to Table of Contents](#)

ABO/RH Tube Testing	
1	Label tubes per SOP Labeling Tubes
2	Add 1 drop of reagent anti-A, -B, -D and ABO/Rh control in respective labeled tubes
3	Add 2 drops of patient plasma in respective labeled tubes
4	Prepare an approximate 3-4% patient cell suspension in respective labeled tube
5	Add 1 drop of the 3-4% patient's red cell suspension in respective tubes
6	Add 1 drop of reagent A1 and B cells in respective labeled tubes
7	Mix gently and centrifuge for time posted
8	Resuspend the cell button and examine for agglutination and/or evidence of hemolysis
9	Read, grade, and record the reactions per SOP Grading Reactions in Sunquest
10	Interpret ABO/Rh correctly

Antibody Screen (PEG)							
1	Label tubes per SOP Labeling Tubes						
2	Add 2 drops of patient plasma in respective labeled tubes						
3	Add 1 drop of the appropriate reagent red cell in respective tubes						
4	Add 2 drops of PeG in respective labeled tubes						
5	Mix gently and incubates at 37°C ±1 for 10-30 minutes						
6	Examine for evidence of hemolysis						
7	Wash the tubes 3 times with saline						
8	Add 2 drops of Anti-IgG, mixes well and centrifuges according to calibration for AHG phase testing						
9	Shake gently to resuspend the cell buttons, and examine macroscopically for agglutination						
10	Read, grade, and record the reaction per SOP <i>Grading Reactions</i> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 30%; padding: 2px;">If the reaction is</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Negative</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination </td> </tr> <tr> <td style="padding: 2px;">Positive</td> <td style="padding: 2px;">Continue to the next step</td> </tr> </tbody> </table>	If the reaction is	Then	Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination 	Positive	Continue to the next step
If the reaction is	Then						
Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination 						
Positive	Continue to the next step						
11	Interpret the Antibody Screen correctly						

[Return to Table of Contents](#)

Antibody Screen (LISS)							
1	Labels tube per SOP Labeling Tubes						
2	Add 2 drops of patient plasma in respective labeled tubes						
3	Add 1 drop of the appropriate reagent red cell in respective tubes						
4	Add 2 drops of LISS in respective labeled tubes						
5	Mix gently and incubate at 37°C ±1 for 10-30 minutes						
6	Centrifuge according to calibration for LISS phase and observe for hemolysis						
7	Shake gently to resuspend the cell buttons and examine macroscopically for agglutination						
8	Read, grade, and record reactions per SOP Grading Reactions						
9	Wash the tubes 3 times with saline						
10	Add 2 drops of Anti-IgG and mix well and centrifuge according to calibration for AHG phase testing						
11	Shake gently to resuspend the cell buttons and examine macroscopically for agglutination						
12	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 30%; padding: 5px;">If the reaction is</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Negative</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination </td> </tr> <tr> <td style="padding: 5px;">Positive</td> <td style="padding: 5px;">Continue to the next step</td> </tr> </tbody> </table>	If the reaction is	Then	Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination 	Positive	Continue to the next step
If the reaction is	Then						
Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative reaction Mix gently and centrifuge on AHG program Shake gently to resuspend the cell button and examine for agglutination 						
Positive	Continue to the next step						
13	Interpret the Antibody Screen correctly						

Heat Block Daily Maintenance							
1	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #cccccc;"> <th style="width: 40%; padding: 5px;">If the thermometer bulb is</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Covered in saline</td> <td style="padding: 5px;">Read the temperature</td> </tr> <tr> <td style="padding: 5px;">Not covered in saline</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Fill with saline until the bulb is covered Allow the temperature to equilibrate Read the temperature </td> </tr> </tbody> </table>	If the thermometer bulb is	Then	Covered in saline	Read the temperature	Not covered in saline	<ul style="list-style-type: none"> Fill with saline until the bulb is covered Allow the temperature to equilibrate Read the temperature
If the thermometer bulb is	Then						
Covered in saline	Read the temperature						
Not covered in saline	<ul style="list-style-type: none"> Fill with saline until the bulb is covered Allow the temperature to equilibrate Read the temperature 						
2	Read temperature, verify acceptability, and document on the Thaw Bath & Heat Block QC form.						

Cell Washer Daily Maintenance	
1	Inspect tubing and drain to ensure they are: <ul style="list-style-type: none"> Clear of obstructions Tubing connections are secure
2	Press and hold SALINE until the "Calibrate 56.4mL" message displays
3	Hold a graduated cylinder under the saline nozzle; press CHECK button and collect the dispensed saline
4	Verify saline volume is between 54 and 59 mL. If not, adjust volume and prime until dispensed saline volume is between 54 and 59mL
5	Clean and dry the interior
6	Document saline volume and tubing inspection on the Cell Washer & Scale QC form

[Return to Table of Contents](#)

DAT in Tube		
Polyspecific		
1	Label tube(s) per SOP Labeling Tubes	
2	Prepare an approximate 3-4% cell suspension of patient red cells	
3	Add 1 drop of patient red cell suspension labeled tube	
4	Wash tube(s) at least 3 times in saline	
5	Add 2 drops of Anti-IgG/C3d Polyspecific AHG to the dry cell button in the tube(s)	
6	Mix gently and centrifuge for time posted	
7	Resuspend the cell button and examine for agglutination	
8	Read, grade, and record reactions per SOP Grading Reactions	
	If DAT is	Then
	Positive	Indicate IgG and C3 testing is required
	Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to the tube Mix tube(s) gently and centrifuge for time posted Resuspend the cell button and examine for agglutination
9	Interpret results correctly	
IgG & C₃		
14	Label 3 tubes for reagents IgG, C3 and saline control per SOP Labeling Tubes	
15	Prepare an approximate 3-4% cell suspension of patient red cells	
16	Add 1 drop of patient red cell suspension to appropriately labeled tube(s)	
17	Wash tube(s) at least 3 times in saline	
18	Add 2 drops of appropriate reagents or saline to the dry cell button in each tube	
19	Mix Anti-IgG and saline control tubes gently and centrifuges for time posted	
20	Read, grade, and record the reactions per SOP Grading Reactions	
21	Incubate Anti-C3 tube for 5 minutes	
22	Centrifuge Anti-C3 tube for time posted	
23	Resuspend the cell button and examine for agglutination	
24	Read, grade, and record the reactions per SOP Grading Reactions	
	If Interpretation is	Then
	Positive	No further action
	Negative	<ul style="list-style-type: none"> Add 1 drop of the appropriate control cells to the tube Mix tube(s) gently and centrifuge for time posted Resuspend the cell button and examine for agglutination
25	Interpret IgG and C3 AHG correctly	

[Return to Table of Contents](#)

Electronic Crossmatch

1	Open the patient order in Blood Order Processing module
2	Select the appropriate crossmatch eligible battery
3	Review patient history and document in Sunquest if not previously completed
4	Select and allocate the appropriate red cell component that meets all special requirements
5	Click save
6	Review the "Electronic Crossmatch Eligibility Report" Click OK or Cancel, as appropriate
7	Click the button in the "Call BPI" pop up box to issue or not

Immediate Spin Crossmatch

1	Label tubes for each crossmatch to be tested per SOP Labeling Tubes
2	Prepare an approximate 3-4% cell suspension for each donor cell to be tested
3	Add 2 drops of patient plasma/serum to tubes labeled for each donor cell
4	Add 1 drop of 3-4% donor cells to labeled tubes
5	Mix gently and centrifuge for time specified
6	Read, grade, and record results, per SOP Grading Reactions

[Return to Table of Contents](#)

AHG Crossmatch (PEG or LISS)							
1	Label tubes for each crossmatch to be tested per SOP Labeling for Manual Testing						
2	Prepare an approximate 3-4% cell suspension for each donor cell to be tested						
3	Add 2 drops of patient plasma/serum to tubes labeled for each donor cell						
4	Add 1 drop of 3-4% donor cells to labeled tubes						
5	Mix gently and centrifuge for time specified						
6	Read, grade, and record results per SOP Grading Reactions and enters results in Sunquest						
7	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 20%; padding: 2px;">If using:</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">PeG</td> <td style="padding: 2px;">Add 2 drops of PeG to each test tube</td> </tr> <tr> <td style="padding: 2px;">LISS</td> <td style="padding: 2px;">Add 2 drops of LISS to each test tube</td> </tr> </tbody> </table>	If using:	Then	PeG	Add 2 drops of PeG to each test tube	LISS	Add 2 drops of LISS to each test tube
If using:	Then						
PeG	Add 2 drops of PeG to each test tube						
LISS	Add 2 drops of LISS to each test tube						
8	Mix well and incubate at 37 \pm 1 for 10-30 minutes						
9	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 20%; padding: 2px;">If using:</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">PeG</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Examine for evidence of hemolysis Continue to next step </td> </tr> <tr> <td style="padding: 2px;">LISS</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Centrifuge according to calibration for LISS phase and observe for hemolysis Shake gently to resuspend the cell buttons and examine macroscopically for agglutination Read, grade and record reactions per SOP Grading Reactions </td> </tr> </tbody> </table>	If using:	Then	PeG	<ul style="list-style-type: none"> Examine for evidence of hemolysis Continue to next step 	LISS	<ul style="list-style-type: none"> Centrifuge according to calibration for LISS phase and observe for hemolysis Shake gently to resuspend the cell buttons and examine macroscopically for agglutination Read, grade and record reactions per SOP Grading Reactions
If using:	Then						
PeG	<ul style="list-style-type: none"> Examine for evidence of hemolysis Continue to next step 						
LISS	<ul style="list-style-type: none"> Centrifuge according to calibration for LISS phase and observe for hemolysis Shake gently to resuspend the cell buttons and examine macroscopically for agglutination Read, grade and record reactions per SOP Grading Reactions 						
10	Wash the tubes 3 times with saline						
11	Add 2 drops of Anti-IgG, mix well and centrifuge according to calibration for AHG phase testing						
12	Shake gently to resuspend the cell buttons, and examine macroscopically for agglutination						
13	<p>Read, grades, and record the reactions per SOP Grading Reactions</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 40%; padding: 2px;">If the reaction is</th> <th style="padding: 2px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 2px;">Negative</td> <td style="padding: 2px;"> <ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative antiglobulin test Mix gently and centrifuge according to calibration for AHG phase testing Shake gently to resuspend the cell button and examine for agglutination </td> </tr> <tr> <td style="padding: 2px;">Positive</td> <td style="padding: 2px;">Continue to the next step</td> </tr> </tbody> </table>	If the reaction is	Then	Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative antiglobulin test Mix gently and centrifuge according to calibration for AHG phase testing Shake gently to resuspend the cell button and examine for agglutination 	Positive	Continue to the next step
If the reaction is	Then						
Negative	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative antiglobulin test Mix gently and centrifuge according to calibration for AHG phase testing Shake gently to resuspend the cell button and examine for agglutination 						
Positive	Continue to the next step						
14	Interpret the AHG crossmatch correctly						

[Return to Table of Contents](#)

Antigen Typing		
1	Label tube(s) per SOP Labeling Tubes including saline/albumin controls and QC if required	
2	Prepare an approximate 3-4% cell suspension of patient/donor red cells to be tested	
3	Add the appropriate number of drops of antisera to labeled tubes as indicated by the manufacturer's package insert (or Instructions for Use)	
4	Add 1 drop of 3-4% cell suspension to labeled tubes	
5	Mix well and incubate at appropriate temperature and time for antisera as indicated by the manufacturer's package insert (or Instructions for Use)	
6	If required phase of testing is	Then
	Not AHG	Centrifuge using saline program
	AHG	<ul style="list-style-type: none"> Wash as indicated by the manufacturer's package insert (or Instructions for Use) Add 2 drops of Anti-IgG and centrifuge using AHG program
7	Shake gently to resuspend the cell button, and examine macroscopically for agglutination	
8	Read and grade per SOP Grading Reactions, and record on Phenotyping Worksheet and in Sunquest if applicable	
	If required phase of testing is	Then
	NOT AHG	No further action is required
	AHG	<ul style="list-style-type: none"> Add 1 drop of IgG coated control cells to each tube with a negative agglutination test Mix gently and centrifuge using the AHG program Shake gently to resuspend the cell button, and examine macroscopically for agglutination

[Return to Table of Contents](#)

Fetal Maternal hemorrhage RapidScreen	
1	Label tubes per SOP Labeling Tubes
2	Make an approximate 3-4% red cell suspension of well-mixed maternal blood in a properly labeled tube
3	Place one drop of the anti-D reagent from the kit in properly labeled tubes for patient and controls
4	Add one drop of the maternal red blood cell suspension to patient tube
5	Add one drop of each well mixed control to the appropriate tubes
6	Mix well and incubate for 5 minutes at room temperature (18-30°C)
7	Wash test tubes 4 times with blood bank saline
8	Add one drop of well-mixed Indicator Cells from the kit and mix gently
9	Centrifuge for time specified
10	Re-suspend red blood cell button completely and examine 5 fields on low-power (10X) microscopically for mixed field agglutination. This can be done in the tube or by transferring the contents of the test tubes onto a microscope slide. Counting should be performed from the slide only.
11	Read and grade the negative and positive controls with the specimen, per SOP <i>Grading Reactions</i>
12	Interpret FMH RapidScreen reactions correctly

SickleDex Hemoglobin S	
1	Label 12 x 75 mm tubes for each donor to be tested, one positive control, and one negative control per SOP Labeling Tubes
2	Place tubes in testing rack and fill each tube with working SickleDex Solubility Buffer to the red line on rack (approximately 2 mL)
3	Return working buffer to refrigerator immediately after use
4	Bring controls and buffer to room temperature (18°C to 30°C) for at least 10 minutes before use
5	Mix controls by holding vertically between hands and rolling the vials back and forth for 20-30 seconds followed by inverting end over-end 20 times
6	Add 1 drop of each control into the appropriately labeled control tube by inverting the control and holding it vertically directly over the test tube for accurate delivery. Wipe threads on each control, if necessary
7	Add 20µL of whole blood or 10µL of packed red blood cells to each tube
8	Mix contents of test tubes
9	Allow tubes to stand at room temperature for at least 6 minutes, no longer than 60 minutes
10	Read reaction macroscopically by looking through the test tubes at black lines on the back of the testing rack
11	Interpret Sickle Cell Screening test results correctly
12	Record results in Sunquest or manual form

[Return to Table of Contents](#)

Testing with DTT Treated Cells	
1	Label tubes for each reagent cell
2	Add DTT-treated reagent cells in Alsever's solution to tube
3	Centrifuge for 1 minute
4	Decant supernatant
5	Make a 3-4% cell suspension using saline
6	Perform daily reagent QC prior to patient testing

EGA Testing							
1	Label test tubes per SOP Labeling Tubes						
2	Wash patient red cells 3 times with blood bank saline, making a 3-4% cell suspension after final wash						
3	Place 30 drops of the 3-4% cell suspension in a labeled tube						
4	Centrifuge washed red cells to form a cell button and remove supernatant						
5	Prepare EDTA glycine by adding 16 drops of Solution 2 into a labeled tube and then 4 drops of Solution 1						
6	Add the combined EGA solutions to the packed washed cells and mix gently						
7	Incubate for 2 minutes at Room Temperature						
8	Add 4 drops of Solution 3, mix well, and centrifuge for 30 sec						
9	Remove and discard the supernatant, add saline to treated cells						
10	Wash red cells 3 times with blood bank saline						
11	Perform DAT on treated cells <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 40%; padding: 5px;">If DAT is</th> <th style="padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Positive</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Repeat treatment starting at step 3 At step 7 incubates for 1.5 minutes </td> </tr> <tr> <td style="padding: 5px;">Negative</td> <td style="padding: 5px;">Proceed to antigen typing</td> </tr> </tbody> </table>	If DAT is	Then	Positive	<ul style="list-style-type: none"> Repeat treatment starting at step 3 At step 7 incubates for 1.5 minutes 	Negative	Proceed to antigen typing
If DAT is	Then						
Positive	<ul style="list-style-type: none"> Repeat treatment starting at step 3 At step 7 incubates for 1.5 minutes 						
Negative	Proceed to antigen typing						

[Return to Table of Contents](#)

Eluate Testing	
1	Label tubes per SOP Labeling Tubes
2	Centrifuge patient specimen and remove most, if not all of the plasma, if necessary
3	Place an aliquot of patient red blood cells (a smaller volume is acceptable if specimen sample is insufficient) in a labeled tube and wash it once with saline.
4	Wash the aliquot four times with Working Wash Solution packing the fourth wash for at least two minutes
5	Remove supernatant and save the last wash in a labeled tube
6	Add 20 drops of washed red blood cells to a clean labeled tube (a smaller volume is acceptable if specimen sample is insufficient)
7	Add 20 drops of Eluting Solution to the washed red blood cells (use equal drops of Eluting Solution if the patient cells fall short of 20 drops) invert 4 times to mix.
8	Centrifuge promptly for 45-60 seconds
9	Transfer supernatant eluate into clean labeled tube, discard the treated red blood cells
10	Add Buffering Solution until it turns pale blue
11	Mix well and, if necessary, centrifuge for 30 seconds, transfer eluate to a clean labeled tube (eluate must be clear of any precipitate or cellular debris)

[Return to Table of Contents](#)