



University of Washington Medical Center 1959 NE Pacific Street. Seattle, WA 98195 Transfusion Services Laboratory Policies and Procedures Manual	Original Effective Date: 03-11-2016	Number: PC-0016.03
	Revision Effective Date: 08-15-2022	
TITLE: Receiving Blood Products into Inventory		

PURPOSE:

This procedure provides instructions for receiving routine blood products into the UWMC Transfusion Services Laboratory. Entry into the LIS and recording visual inspection are described. Donor segment retention and routing of donor samples for testing is also described.

PRINCIPLE & CLINICAL SIGNIFICANCE:

Receipt of blood products from a blood supplier is achieved through observation of packaging to maintain temperature, comparison of order quantities against quantities received, entry of the product into the LIS for tracking and a documented visual inspection of the blood product. When RBC containing products (whole blood, RBCs or granulocytes) are received, donor segment retention and ABO/Rh type confirmation are also required prior to making units available for allocation and issue.

POLICIES:

- Donor segments are retained for a minimum of two months
- Any shipments with questionable storage conditions must have the temperature verified and documented prior to accepting the shipment into inventory
- Donor units must be processed in a manner such that time out of controlled storage conditions is limited
- RBC containing products must be segregated from available inventory until the type confirmation is complete
- RBCs that will be stored in the AGNEG refrigerator will be labeled with a “STOP - NOT IRRADIATED” tag (see [Appendix A](#)) before being placed into the refrigerator.
- Platelets that are NOT Psoralen-Treated and arrive from the supplier NOT Irradiated will be labeled with a “STOP- NOT IRRADIATED” tag (see [Appendix A](#)) before being placed into the platelet incubator.

SPECIMEN REQUIREMENTS:

NA

REAGENTS/SUPPLIES/EQUIPMENT:

Reagents:	Supplies:	Equipment:
None	<ul style="list-style-type: none"> • Test tubes • Plastic bag • Retention date labels • Test tube rack • Scissors 	LIS with scanner

QUALITY CONTROL:

NA

INSTRUCTIONS:

TABLE of CONTENTS
[Accepting Delivery](#)
[Inspection of Blood Shipment](#)
[Blood Product Entry in LIS](#)

Accepting Delivery

STEP	ACTION						
1	Ensure shipment is delivered to the correct delivery location						
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 30%;">If delivery location is</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td style="border-bottom: none;">Correct</td> <td> <ul style="list-style-type: none"> Sign courier log if required </td> </tr> <tr> <td>Incorrect</td> <td> <ul style="list-style-type: none"> Inform courier and supplier of wrong location Do not sign for shipment or accept shipment </td> </tr> </tbody> </table>	If delivery location is	Then	Correct	<ul style="list-style-type: none"> Sign courier log if required 	Incorrect	<ul style="list-style-type: none"> Inform courier and supplier of wrong location Do not sign for shipment or accept shipment
	If delivery location is	Then					
Correct	<ul style="list-style-type: none"> Sign courier log if required 						
Incorrect	<ul style="list-style-type: none"> Inform courier and supplier of wrong location Do not sign for shipment or accept shipment 						
Communicate to courier any boxes and/or other items to be returned to blood supplier							

Inspection of Blood Shipment

STEP	ACTION												
1	Open the shipping container and time stamp or write the date and time of delivery on the packing slip as soon as possible upon opening the box												
2	Verify contents are packed appropriately and shipment appears undamaged												
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 40%;">If</th> <th style="width: 20%;">Packing condition</th> <th style="width: 40%;">*Temp Range</th> </tr> </thead> <tbody> <tr> <td style="border-bottom: none;">Red Blood Cells</td> <td>Wet ice is present</td> <td>1-10° C</td> </tr> <tr> <td style="border-bottom: none;">Platelets, Granulocytes</td> <td>Room temperature stabilizing packs</td> <td>20-24°C</td> </tr> <tr> <td>Fresh Frozen Plasma, Cryoprecipitate</td> <td>Dry Ice is present</td> <td>< -18°C</td> </tr> </tbody> </table>	If	Packing condition	*Temp Range	Red Blood Cells	Wet ice is present	1-10° C	Platelets, Granulocytes	Room temperature stabilizing packs	20-24°C	Fresh Frozen Plasma, Cryoprecipitate	Dry Ice is present	< -18°C
	If	Packing condition	*Temp Range										
	Red Blood Cells	Wet ice is present	1-10° C										
Platelets, Granulocytes	Room temperature stabilizing packs	20-24°C											
Fresh Frozen Plasma, Cryoprecipitate	Dry Ice is present	< -18°C											
* If temperature is in question, verify the product transport temperature to ensure the range has not been exceeded													
3	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 30%;">If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td style="border-bottom: none;">Shipment acceptable</td> <td> <ul style="list-style-type: none"> Go to next step </td> </tr> <tr> <td>Temperature not maintained, shipment leaking or otherwise damaged</td> <td> <ul style="list-style-type: none"> Notify shift lead or manager and complete QIM Report Notify blood supplier regarding the issue </td> </tr> </tbody> </table>	If	Then	Shipment acceptable	<ul style="list-style-type: none"> Go to next step 	Temperature not maintained, shipment leaking or otherwise damaged	<ul style="list-style-type: none"> Notify shift lead or manager and complete QIM Report Notify blood supplier regarding the issue 						
	If	Then											
	Shipment acceptable	<ul style="list-style-type: none"> Go to next step 											
Temperature not maintained, shipment leaking or otherwise damaged	<ul style="list-style-type: none"> Notify shift lead or manager and complete QIM Report Notify blood supplier regarding the issue 												
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 20%;">If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td style="border-bottom: none;">Temp not OK</td> <td> <ul style="list-style-type: none"> Use a NIST calibrated thermometer to verify the temperature by placing the thermometer between two components (if possible) or sandwich the single product and read temperature after 3-5 minutes Record shipment temperature or other shipment issue on the packing slip Quarantine all products if not immediately returned to supplier (refer to SOP <i>Quarantine and Final Disposition of Blood Products</i>) </td> </tr> </tbody> </table>	If	Then	Temp not OK	<ul style="list-style-type: none"> Use a NIST calibrated thermometer to verify the temperature by placing the thermometer between two components (if possible) or sandwich the single product and read temperature after 3-5 minutes Record shipment temperature or other shipment issue on the packing slip Quarantine all products if not immediately returned to supplier (refer to SOP <i>Quarantine and Final Disposition of Blood Products</i>) 									
If	Then												
Temp not OK	<ul style="list-style-type: none"> Use a NIST calibrated thermometer to verify the temperature by placing the thermometer between two components (if possible) or sandwich the single product and read temperature after 3-5 minutes Record shipment temperature or other shipment issue on the packing slip Quarantine all products if not immediately returned to supplier (refer to SOP <i>Quarantine and Final Disposition of Blood Products</i>) 												

STEP	ACTION							
		Shipment leaking <ul style="list-style-type: none"> Find source of the leak Record the condition of the box on the packing slip Quarantine all products if not immediately returned to supplier (refer to SOP <i>Quarantine and Final Disposition of Blood Products</i>) 						
4	Compare the components shipped with those listed on the packing slip and verify the following: <ul style="list-style-type: none"> Unit numbers match Components received match the order placed Notify shipper if any discrepancy is noted 							
5	Inspect each component according to the SOP <i>Visual Inspection of Blood Products</i> <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 30%; padding: 5px;">If visual inspection</th> <th style="width: 70%; padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Passes</td> <td style="padding: 5px;">Continue to next step</td> </tr> <tr> <td style="padding: 5px;">Doesn't pass</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Quarantine all products in the shipment until further investigation is complete Notify shift lead or manager and complete QIM Report </td> </tr> </tbody> </table>		If visual inspection	Then	Passes	Continue to next step	Doesn't pass	<ul style="list-style-type: none"> Quarantine all products in the shipment until further investigation is complete Notify shift lead or manager and complete QIM Report
If visual inspection	Then							
Passes	Continue to next step							
Doesn't pass	<ul style="list-style-type: none"> Quarantine all products in the shipment until further investigation is complete Notify shift lead or manager and complete QIM Report 							
6	For any blood products received that have antigen typing(s) on label or tag: <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 30%; padding: 5px;">If person performing entry is a</th> <th style="width: 70%; padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">MLS</td> <td style="padding: 5px;">Proceed to step 7</td> </tr> <tr> <td style="padding: 5px;">Not MLS</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Do not proceed with Blood Entry process for antigen typed unit(s) Give unit(s) to a MLS to perform Blood Product Entry and further testing Antigen types must be entered using appropriate Ab/Ag codes in SQ </td> </tr> </tbody> </table>		If person performing entry is a	Then	MLS	Proceed to step 7	Not MLS	<ul style="list-style-type: none"> Do not proceed with Blood Entry process for antigen typed unit(s) Give unit(s) to a MLS to perform Blood Product Entry and further testing Antigen types must be entered using appropriate Ab/Ag codes in SQ
If person performing entry is a	Then							
MLS	Proceed to step 7							
Not MLS	<ul style="list-style-type: none"> Do not proceed with Blood Entry process for antigen typed unit(s) Give unit(s) to a MLS to perform Blood Product Entry and further testing Antigen types must be entered using appropriate Ab/Ag codes in SQ 							
7	<table border="1" style="width: 100%; border-collapse: collapse; margin-top: 5px;"> <thead> <tr style="background-color: #d3d3d3;"> <th style="width: 30%; padding: 5px;">If</th> <th style="width: 70%; padding: 5px;">Then</th> </tr> </thead> <tbody> <tr> <td style="padding: 5px;">Receiving Platelets, Cryo, Plasma</td> <td style="padding: 5px;">Go to next section: Blood Product Entry in LIS</td> </tr> <tr> <td style="padding: 5px;">Receiving RBCs and Granulocytes</td> <td style="padding: 5px;"> <ul style="list-style-type: none"> Remove 2 segments from each unit Label one segment with a unit # sticker and place in a dated storage bag for retention Place one segment in a glass tube labelled with unit # and blood type for ABO/Rh confirmation. See SOP <i>Unit Type Confirmation</i> Go to next section: Blood Product Entry in LIS </td> </tr> </tbody> </table>		If	Then	Receiving Platelets, Cryo, Plasma	Go to next section: Blood Product Entry in LIS	Receiving RBCs and Granulocytes	<ul style="list-style-type: none"> Remove 2 segments from each unit Label one segment with a unit # sticker and place in a dated storage bag for retention Place one segment in a glass tube labelled with unit # and blood type for ABO/Rh confirmation. See SOP <i>Unit Type Confirmation</i> Go to next section: Blood Product Entry in LIS
If	Then							
Receiving Platelets, Cryo, Plasma	Go to next section: Blood Product Entry in LIS							
Receiving RBCs and Granulocytes	<ul style="list-style-type: none"> Remove 2 segments from each unit Label one segment with a unit # sticker and place in a dated storage bag for retention Place one segment in a glass tube labelled with unit # and blood type for ABO/Rh confirmation. See SOP <i>Unit Type Confirmation</i> Go to next section: Blood Product Entry in LIS 							

Blood Product Entry in LIS

STEP	ACTION	
1	Open the ' Blood Product Entry' (BPE) function in Sunquest (SQ)	

	Any component with antigen typing on the label or attached tag	<ul style="list-style-type: none"> Only MLS staff may receive these components in SQ Antigen types must be entered using appropriate Ab/Ag codes in SQ and only MLS staff may receive these products in SQ 						
	Autologous, Directed or HLA matched components	<ul style="list-style-type: none"> Segregate these from other components Receive individually in SQ and add the intended recipient as an assignee using the recipient's MRN <p>NOTE: If the patient does not have a MRN assigned in the LIS/HIS, quarantine the product and contact the ordering provider to obtain a pre-registration MRN for assignee purposes in Sunquest. The provide will need to register the patient if a MRN is not available.</p>						
	All other components	<ul style="list-style-type: none"> Go to next step 						
2	Open the 'Blood Product Entry' (BPE) in Sunquest (SQ)							
3	Scan the following barcodes from the component label <ul style="list-style-type: none"> Unit # <table border="1" style="width: 100%; border-collapse: collapse; margin-left: 20px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 30%;">If facility ID is</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Recognized</td> <td>Go to next step</td> </tr> <tr> <td>Not recognized</td> <td>Choose OTHER for supplier</td> </tr> </tbody> </table> <ul style="list-style-type: none"> Product Type ABO/Rh Expiration Date/Time CMV negative barcode, if applicable 		If facility ID is	Then	Recognized	Go to next step	Not recognized	Choose OTHER for supplier
If facility ID is	Then							
Recognized	Go to next step							
Not recognized	Choose OTHER for supplier							
4	Enter the component volume from the label for all products except single donation RBCs that should autofill with 350ml							
5	Add additional information when applicable <table border="1" style="width: 100%; border-collapse: collapse; margin-left: 20px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 30%;">If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Unit is low titer plasma</td> <td> <ul style="list-style-type: none"> Select the "Ag/Ab/Attribute" tab Type LTP in the Antigen/Antibody box Click <Add> </td> </tr> <tr> <td>Supplier was entered as "Other" in step 2</td> <td> <ul style="list-style-type: none"> Select the "Comments" tab Enter the collection facility name in the free text box – do not add a comment code Click <Add> </td> </tr> </tbody> </table>		If	Then	Unit is low titer plasma	<ul style="list-style-type: none"> Select the "Ag/Ab/Attribute" tab Type LTP in the Antigen/Antibody box Click <Add> 	Supplier was entered as "Other" in step 2	<ul style="list-style-type: none"> Select the "Comments" tab Enter the collection facility name in the free text box – do not add a comment code Click <Add>
If	Then							
Unit is low titer plasma	<ul style="list-style-type: none"> Select the "Ag/Ab/Attribute" tab Type LTP in the Antigen/Antibody box Click <Add> 							
Supplier was entered as "Other" in step 2	<ul style="list-style-type: none"> Select the "Comments" tab Enter the collection facility name in the free text box – do not add a comment code Click <Add> 							
6	Visually inspect the component according to the SOP <i>Visual Inspection of Blood Components</i> <table border="1" style="width: 100%; border-collapse: collapse; margin-left: 20px;"> <thead> <tr style="background-color: #e0e0e0;"> <th style="width: 30%;">If visual inspection</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Passes</td> <td> <ul style="list-style-type: none"> Select PASS from the drop-down menu Continue to next step </td> </tr> <tr> <td>Does not pass</td> <td> <ul style="list-style-type: none"> Do not receive component into inventory Contact the supplier and return or discard the component as instructed </td> </tr> </tbody> </table>		If visual inspection	Then	Passes	<ul style="list-style-type: none"> Select PASS from the drop-down menu Continue to next step 	Does not pass	<ul style="list-style-type: none"> Do not receive component into inventory Contact the supplier and return or discard the component as instructed
If visual inspection	Then							
Passes	<ul style="list-style-type: none"> Select PASS from the drop-down menu Continue to next step 							
Does not pass	<ul style="list-style-type: none"> Do not receive component into inventory Contact the supplier and return or discard the component as instructed 							

TITLE: Receiving Blood Products into Inventory	Number: PC-0016.03
---	-------------------------------

		NOTE: If the supplier does not want the component returned, receive the component in Sunquest, fail the visual inspection, and discard and dispose of appropriately.
7	Click the "Add" button closest to the "Unit Summary List" to allow batch entry of similar products	
8	Repeat steps 2-7 for each unit in the batch and click SAVE	
9	Compare the units listed on the screen with the packing list and ensure all information matches	
10	<ul style="list-style-type: none"> • Click <SAVE> • Click <EXIT> 	
11	Record the Worklist # for use in performing and documenting the ABO/Rh unit confirmation	
12	Place components in the appropriate storage area for the component type	
13	File the packing list in the appropriate file for verification of billing from the supplier	
14	If	Then
	Platelets, plasma or cryoprecipitate (non-cellular products)	No further action needed
	RBCs or Granulocytes	Route the segments to the testing area for ABO/Rh type confirmation (refer to SOP <i>Unit Type Confirmation</i>)

CALIBRATION:

NA

NOTES AND LIMITATIONS:

- Care must be taken to ensure that units are not out of monitored storage conditions for a prolonged time to prevent products from exceeding the acceptable temperature ranges.
- Apheresis RBC units must be tested manually and labelled with the product code or container # in addition to the unit #
- Some products must be entered into BPE separately. Examples include Autologous, Directed, and HLA-matched. These types of units should be sequestered and entered accordingly.

REFERENCES:

- Technical Manual. Bethesda, MD: AABB Press, current edition
- Standards for Blood Banks and Transfusion Services. Bethesda, MD: AABB Press, current edition

RELATED DOCUMENTS:

SOP Quarantine and Final Disposition of Blood Products
 SOP Visual Inspection of Blood Components
 SOP Unit Type Confirmation Using Tube Method
 SOP Specimen and Unit Segment Management

APPENDIX:

A - STOP- NOT IRRADIATED" tag

UWMC SOP Approval:		
UWMC CLIA Medical Director	Andrew Bryan, MD	Date
Transfusion Service Manager	Nina Sen	Date
QA Manger	Tayler Reeves	Date
Transfusion Service Medical Director	Monica B Pagano, MD	Date
UWMC Biennial Review:		
		Date
		Date

REVISION HISTORY:
04/22/2018: Updated to include changes due to Sunquest 8.1 upgrade. Most significant change is the visual inspection of each blood component can now be documented at the time the component is received into the LIS (Laboratory Information System). The Recording Visual Inspection Section of version PC-0016.01 was removed.

08/15/2022: Updated to include use of "STOP - NOT IRRADIATED" tag on antigen negative RBCs before placing units into the antigen negative refrigerator, and non-psoralen treated platelets received from the supplier without irradiation before placing units in the platelet incubator.

Appendix A

