



Veritor™ System

For Rapid Detection Flu A+B

CE 8087666(10)
2016-05

Laboratory kit configured for testing liquid nasopharyngeal wash, aspirate and swab in transport media samples.

Kit de laboratoire conçu pour l'analyse d'échantillons liquides de lavage, d'aspiration et d'écouvillonnage rhino-pharyngés en milieu de transport.

Labor-Testkit zum Testen von flüssigen Nasopharyngeal-Spülungen, -Aspiraten und -Abstrichen in Proben in Transportmedien.

Kit per laboratorio configurato per l'analisi di lavaggi, tamponi e aspirati nasofaringei nei campioni di terreni di trasporto.

Kit de laboratorio configurado para analizar muestras de lavados, aspirados y torundas nasofaríngeos en medio de transporte.

English: pages 2 – 19

Français: pages 19 – 36

Deutsch: Seiten 37 – 54

Italiano: pagine 54 – 72

Español: páginas 73 – 92

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyttä lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жепринікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

30

Determinations
Déterminations
Bestimmungen
Determinazioni
Determinaciones

BD Veritor™ System

For Rapid Detection of Flu A+B

Rx Only

For *in vitro* diagnostic use only.

INTENDED USE

The BD Veritor™ System for Rapid Detection of Flu A+B is a rapid chromatographic immunoassay for the direct and qualitative detection of influenza A and B viral nucleoprotein antigens from nasopharyngeal wash, aspirate and swab in transport media samples from symptomatic patients. The BD Veritor System for Rapid Detection of Flu A+B is a differentiated test, such that influenza A viral antigens can be distinguished from influenza B viral antigens from a single processed sample using a single device. The test is to be used as an aid in the diagnosis of influenza A and B viral infections. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay. Outside the U.S. a negative test is presumptive and it is recommended that these results be confirmed by viral culture or a molecular assay cleared for diagnostic use in the country of use. FDA has not cleared this device outside the U.S. Negative test results do not preclude influenza viral infection and should not be used as the sole basis for treatment or other patient management decisions. The test is not intended to detect influenza C antigens.

Performance characteristics for influenza A and B nasopharyngeal (NP) washes/aspirates were established during January through March of 2011 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2010-2011 Season, and Composition of the 2011-2012 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

Performance characteristics for influenza A and B NP swabs in transport media were established during January through April of 2012 when influenza viruses A/2009 H1N1, A/H3N2, B/Victoria lineage, and B/Yamagata lineage were the predominant influenza viruses in circulation according to the *Morbidity and Mortality Weekly Report* from the CDC entitled "Update: Influenza Activity—United States, 2011-2012 Season, and Composition of the 2012-2013 Influenza Vaccine." Performance characteristics may vary against other emerging influenza viruses.

If infection with a novel influenza virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Virus culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

SUMMARY AND EXPLANATION

Influenza illness classically presents with sudden onset of fever, chills, headache, myalgias, and a non-productive cough. Epidemics of influenza typically occur during winter months with estimated 114,000 hospitalizations¹ and 36,000 deaths² per year in the U.S. Influenza viruses can also cause pandemics, during which rates of illness and death from influenza-related complications can increase dramatically.

Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hours of onset of illness. It is important to rapidly distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Moreover, it is important to determine if influenza A or B is causing symptomatic disease in a particular institution (e.g., nursing home) or community, so that appropriate preventative intervention can be taken for susceptible individuals. It is therefore important to not only rapidly determine whether influenza is present, but also which type of influenza virus is present as severity and treatment can be different.³

Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction (PCR), serology, and viral culture.⁴⁻¹¹ Immunofluorescence assays entail staining of specimens immobilized on microscope slides using fluorescent-labeled antibodies for observation by fluorescence microscopy.^{6,12,13} Culture methods employ initial viral isolation in cell culture, followed by hemadsorption inhibition, immunofluorescence, or neutralization assays to confirm the presence of the influenza virus.¹³⁻¹⁵

The BD Veritor System for Rapid Detection of Flu A+B (also referred to as the BD Veritor System and BD Veritor System Flu A+B) is a digital immunoassay (DIA) to qualitatively detect influenza A or B nucleoprotein antigens from respiratory specimens of symptomatic patients with a time to result of 10 minutes. The speed and simplified workflow of the BD Veritor System for Rapid Detection of Flu A+B makes it applicable as a "STAT" influenza A and B antigen detection test providing relevant information to assist with the diagnosis of influenza. All BD Veritor System Flu A+B test devices are interpreted by a BD Veritor System Instrument, either a BD Veritor Reader or BD Veritor Plus Analyzer (the "Analyzer"). When using an Analyzer, procedures to evaluate test devices depend on the workflow configuration chosen. In **Analyze Now mode**, the instrument evaluates assay devices after manual timing of their development. In **Walk Away mode**, devices are inserted immediately after application of the specimen, and timing of assay development and analysis is automated. Additionally, connection of an Analyzer to a printer or IT system is possible if desired. Additional result documentation capabilities are possible with the integration of a BD Veritor InfoScan ("InfoScan") or BD Veritor InfoSync ("InfoSync") module. Please refer to the Analyzer *Instructions for Use* for details on how to implement these features. InfoSync is not available in all regions.

PRINCIPLES OF THE PROCEDURE

The **BD Veritor System** for Rapid Detection of Flu A+B is a qualitative digital immunoassay for the detection of influenza A and B viral antigens in samples processed from respiratory specimens. When specimens are processed and added to the test device, influenza A or B viral antigens bind to anti-influenza antibodies conjugated to detector particles in the A+B test strip. The antigen-conjugate complex migrates across the test strip to the reaction area and is captured by the line of antibody on the membrane. A positive result for influenza A is determined by the **BD Veritor System Instrument** when antigen-conjugate is deposited at the Test "A" position and the Control "C" position on the **BD Veritor System Flu A+B assay device**. A positive result for influenza B is determined by the **BD Veritor System Instrument** when antigen-conjugate is deposited at the Test "B" position and the Control "C" position in the **BD Veritor System Flu A+B assay device**. The instrument analyzes and corrects for non-specific binding and detects positives not recognized by the unaided eye to provide an objective digital result.

REAGENTS

The following components are included in the **BD Veritor System** for Rapid Detection of Flu A+B kit:

| | | |
|--|------------------------------|--|
| BD Veritor System Flu A+B Devices | 30 devices | Foil pouched device containing one reactive strip. Each strip has two test lines of monoclonal antibody specific to either Flu A or Flu B influenza viral antigen and murine monoclonal control line antibodies. |
| RV Reagent C | 30 tubes with 100 µL reagent | Detergent with < 0.1% sodium azide |
| 300 µL Pipette | 30 each | Transfer pipette |
| Control A+/B- Swab | 1 each | Flu A Positive and Flu B Negative Control Swab, influenza A antigen (inactive recombinant nucleoprotein) with < 0.1% sodium azide |
| Control B+/A- Swab | 1 each | Flu A Negative and Flu B Positive Control Swab, influenza B antigen (inactive recombinant nucleoprotein) with < 0.1% sodium azide |

Materials Required But Not Provided: **BD Veritor™ System Reader** (Cat. No. 256055) or **BD Veritor™ Plus Analyzer** (Cat. No. 256066), Timer, Tube Rack for specimen testing.

Optional Equipment: **BD Veritor™ InfoScan Module** (Cat. No. 256068), **BD Veritor™ InfoSync Module** (Cat. No. 256067), USB Printer Cable for **BD Veritor™ Analyzer** (Cat. No. 443907), Epson Printer model TM-T20 II.

Warnings and Precautions:

Warning



H302 Harmful if swallowed. **H402** Harmful to aquatic life. **H412** Harmful to aquatic life with long lasting effects.

P273 Avoid release to the environment. **P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

1. For *in vitro* Diagnostic Use.
2. Test results are not meant to be visually determined. **All test results must be determined using the BD Veritor System Instrument.**
3. If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to the state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.
4. Pathogenic microorganisms, including hepatitis viruses, Human Immunodeficiency Virus and novel influenza viruses, may be present in clinical specimens. "Standard Precautions"¹⁶⁻¹⁹ and institutional guidelines should be followed in handling, storing and disposing of all specimens and all items contaminated with blood and other body fluids.
5. Dispose of used **BD Veritor System** test devices as biohazardous waste in accordance with federal, state and local requirements.
6. Reagents contain sodium azide, which is harmful if inhaled, swallowed or exposed to skin. Contact with acids produces very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
7. Do not use kit components beyond the expiration date.
8. Do not reuse the **BD Veritor System** test device.
9. Do not use the kit if the Control A+/B- swab and Control B+/A- swab do not yield appropriate results.
10. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
11. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section.

12. FluMist® is made from attenuated live flu virus and although the concentration tested (1%) was non-interfering, it is possible when tested with higher concentrations that an influenza A and/or influenza B false positive may occur.
13. Specific training or guidance is recommended if operators are not experienced with specimen collection and handling procedures.

Storage and Handling: Kits may be stored at 2–30 °C. DO NOT FREEZE. Reagents and devices must be at room temperature (15–30 °C) when used for testing.

SPECIMEN COLLECTION AND HANDLING

Specimen Collection and Preparation: Acceptable specimens for testing with the BD Veritor System for Rapid Detection of Flu A+B include nasopharyngeal (NP) washes, aspirates and swab specimens in transport media. It is essential that correct specimen collection and preparation methods be followed. Specimens obtained early in the course of the illness will contain the highest viral titers.

Inadequate specimen collection, improper specimen handling and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to accurate test results.

Specimen Transport Media: The following transport media have been tested and found to be compatible using moderate positive samples with the BD Veritor System for Rapid Detection of Flu A+B:

- Modified Amies Medium (liquid), ESwab, Liquid Stuart Medium, Amies, Bartel ViraTrans™, BD Universal Transport, Hank's Balanced Salt solution, M4, M4-RT, M5, M6, Normal Saline, Phosphate Buffered Saline.

Samples in these transport media can be stored at 2–8 °C for up to 72 hours. Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Failure to do so may result in inconsistent flow.

Other transport media may be utilized if an appropriate validation exercise is performed.

Specimen Transport and Storage:

Freshly collected specimens should be processed within 1 hour. If necessary, specimens may be stored at 2–8 °C for up to 72 hours and then tested at room temperature. It is essential that correct specimen collection and preparation methods be followed. Do not centrifuge specimens prior to use, as the removal of cellular material may adversely affect test sensitivity.

Procedure for Nasopharyngeal Washes/Aspirates:

- For NP washes/aspirates, sample volumes of 1 to 3 mL are recommended. If transport medium is used, minimal dilution of specimens is recommended.
- Excessive wash volumes should be avoided as they may result in decreased test sensitivity.
- Process specimen as described in "Test Procedure."

Procedure for Nasopharyngeal Swabs in Transport Media:

- For NP swabs in transport media, a minimal volume of transport media (1 mL) is recommended to minimize dilution.
- Process specimen as described in "Test Procedure."

TEST PROCEDURE

NOTES: Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.

Prepare for testing


The following steps assume that users of a **BD Veritor Plus Analyzer** have chosen and set all configuration options, and that the Analyzer is ready to use. To choose or change these settings, see the **BD Veritor Plus Analyzer Instructions for Use**, section 4.7. A printer is not necessary to display results. However, if your facility has chosen to connect the **BD Veritor Plus Analyzer** to a printer, check that the printer is plugged into a power source, paper supply is adequate and any necessary network connections are enabled before testing.

For each patient specimen or control swab:

Step 1: Remove one **RV Reagent C** tube/tip and one **BD Veritor System Flu A+B** device from its foil pouch immediately before testing.

Step 2: Label one **BD Veritor System** device and one **RV Reagent C** tube for each specimen and control to be tested.

Step 3: Place the labeled **RV Reagent C** tube(s) in the designated area of the tube rack.




All reagents and specimens must be at room temperature prior to processing

Step 4: Process the Sample or Control:


a. For NP washes, aspirates and swab specimens in transport media:

1. Vortex or thoroughly mix specimen. Do not centrifuge.
2. Remove and discard the cap from the **RV Reagent C** tube corresponding to the sample to be tested.
3. Using the transfer pipette, transfer 300 µL of the specimen into the **RV Reagent C** tube. Discard pipette after use.




b. For kit swab controls:

1. Remove and discard the cap from the **RV Reagent C** tube corresponding to the sample to be tested.
2. Using the transfer pipette add 300 µL of distilled or deionized water to the **RV Reagent C** tube.
3. Insert the control swab into the tube and vigorously plunge the swab up and down in the fluid for a minimum of 15 seconds.
4. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.



Step 5:

- a. Press the attached tip firmly onto the **RV Reagent C** tube containing the processed specimen or control (threading/twisting not required).
- b. Vortex or mix thoroughly by swirling or flicking the bottom of the tube.




NOTE: Do not use tips from any other product, including other products from **BD** or other manufacturers.

| After step 5, choose from the model and workflow option below before continuing to step 6: | | | |
|--|---|--|---|
| | BD Veritor Reader or Analyzer in Analyze Now mode | BD Veritor Plus Analyzer in Walk Away mode | BD Veritor Plus Analyzer with either InfoScan or InfoSync modules In Analyze now mode---or---Walk Away mode |
| Instructions in section: | A | B | C |

A**Using a BD Veritor Reader or Analyzer in "Analyze Now" mode:****Step 6A: Adding the specimen**

- Invert the RV Reagent C tube and hold the tube vertically (approximately one inch above the labeled BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device.




NOTE: Squeezing the tube too close to the tip may cause leakage

Step 7A: Timing development

- After adding the sample, allow the test to run for 10 minutes before inserting into the BD Veritor Instrument.
- **NOTE:** If running test under laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8A: Using the BD Veritor Instrument:**

- During incubation time, turn the BD Veritor Instrument on by pressing the power button once.
- Insert the assay device when the 10 minute assay development time is complete. Follow the on screen prompts to complete the procedure.
- The status of the assay analysis process appears in the display window.



Do not touch the Instrument or remove the test device

Step 9A: Record the Result

- When analysis is complete, the test result appears in the display window.


ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

B**Using an Analyzer in "Walk Away" mode:
with no optional module installed**


To use Walk Away mode – connect the AC power adapter to the Analyzer and a power source

Step 6B: Starting Walk Away mode:

- Turn the Analyzer on by pressing the blue power button once.
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
 - Double-click the blue power button.

**Step 7B: Adding the specimen**

- When the display window reads "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY":
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device.



NOTE: Squeezing the tube close to the tip may cause leakage


Step 8B: Start the development and reading sequence

- Immediately insert the test device into the slot on the right side of the Analyzer.

The test device must remain horizontal to prevent spilling the specimen out of the sample well.

- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.

Do not touch the Analyzer or remove the test device during this time. Doing so will abort the assay analysis.


**Step 9B: Record the result**

- When analysis is complete, the test result appears in the display window.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (if the AC power adapter is connected).

C**Using an Analyzer in "Analyze Now" mode:
with InfoScan or InfoSync module installed****Step 6C: Adding the specimen**

- Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
- Gently squeeze the ridged body of the tube, dispensing three (3) drops of the processed specimen into the sample well of a labeled BD Veritor System Flu A+B device. **NOTE: Squeezing the tube close to the tip may cause leakage.**

**Step 7C: Timing development**



- Allow the test to develop for **10 minutes**.
- If running the test in a laminar flow hood or in an area with heavy ventilation, cover test device to avoid inconsistent flow.

**Step 8C: Using the Analyzer**


During incubation time, turn the BD Veritor Plus Analyzer on by pressing the blue button once.

The display window briefly shows "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.

- When assay development time is complete and the Analyzer display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE":
 - Insert the BD Veritor System Flu A+B device into the BD Veritor Plus Analyzer.

**Step 9C: Using the Bar Code scanner**

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER
 in accordance with site requirements and Analyzer settings.



- Prompts for each scanning step appear in the display window for only 10 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 8C. To restart this step, remove and reinsert the test device to initiate a new sequence.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

- After required scans are completed, the Analyzer displays a countdown timer and test analysis begins.
- Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.
- When analysis is complete a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed.

If a printer is not connected, record the result before removing the assay device.

ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 15 minutes (60 minutes if AC power adapter is connected).

Step 10C: Remove the test device


- Pull the device out. The display will show "INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE" to indicate the Analyzer is ready to perform another test.




- If an InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.
- In the event that the BD Veritor Plus Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

D**Using an Analyzer In "Walk Away" mode:
with InfoScan or InfoSync module installed****To use Walk Away mode - connect the AC power adapter to the Analyzer and a power source****Step 6D: Starting Walk Away mode**

- Turn the Analyzer on by pressing the blue power button.
The display window will briefly show "SCAN CONFIG BARCODE." This is an opportunity to change the configuration of the Analyzer. Please refer to the Analyzer *Instructions for Use* for configuration instructions. Ignore this message and postpone this process when an assay is awaiting analysis.
- When the display window reads: "INSERT TEST DEVICE OR DOUBLE-CLICK FOR WALK AWAY MODE"
 - Double-click the blue power button.

**Step 7D: Using the Bar Code scanner**

- Follow the prompts on the display window to complete any required barcode scans of:
 - OPERATOR ID
 - SPECIMEN ID and/or
 - KIT LOT NUMBER




According to site requirements and Analyzer settings.


- Prompts for each scanning step appear in the display window for only 10 seconds. Failure to complete scans during that time will cause the Analyzer to default to the beginning of step 6D. To restart this step, double-click the power button.
- Move the barcode slowly toward the window until a confirmation tone sounds. The scanned barcode value appears in the next display window.
- The Analyzer can record the kit Lot Number in the test record but does not restrict the use of expired or inappropriate reagents. Management of expired materials is the responsibility of the user. BD recommends against the use of expired materials.

Step 8D: Add the specimen to the test device

- When the display window reads: "ADD SPECIMEN TO TEST DEVICE AND INSERT IMMEDIATELY":
 - Invert the tube, holding it vertically (approximately one inch above the BD Veritor System Flu A+B device sample well).
 - Gently squeeze the ridged portion of the tube, allowing three (3) drops of the processed specimen to dispense into the sample well of a labeled BD Veritor System Flu A+B device. NOTE: Squeezing the tube close to the tip may cause leakage.


**Step 9D: Start the development and reading sequence**

- Immediately insert the test device into the slot on the right side of the Analyzer.
The test device must remain horizontal to prevent spilling the specimen out of the sample well.
- "DO NOT DISTURB TEST IN PROGRESS" appears in the display window. Automatic timing of the assay development, image processing and result analysis begins.
- A countdown timer in the display window shows the remaining analysis time.
Do not touch the Analyzer or remove the test device during this process. Doing so will abort the assay analysis.
- When analysis is complete, a result appears in the display window. If configured to display, the specimen ID barcode value also appears. If a printer is connected, specimen ID and result are automatically printed. **If a printer is not connected, record the result before removing the assay device.**



ATTENTION: TEST Results are NOT maintained in the display window when the device is removed or if the Analyzer is left unattended for more than 60 minutes (when AC power adapter is connected).

Step 10D: Remove the test device

- Pull the device out. The display will show INSERT TEST DEVICE OR DOUBLE-CLICK BUTTON FOR WALK AWAY MODE to indicate the Analyzer is ready to perform another test. Note that the Analyzer returns to Analyze Now mode at the conclusion of each read sequence.
-  If an InfoSync module is installed the ENVELOPE symbol will appear to indicate that results are transmitting.
- In the event that the BD Veritor Plus Analyzer does not detect adequate cellular network strength while the ENVELOPE symbol is still displayed, it will queue all results to be transmitted and continuously attempt to transmit them. If it is powered off during this time, it will attempt to transmit as soon as power is restored.

OPTIONAL TEST PROCEDURE: Use this procedure to test for both INFLUENZA A+B and RSV using a single NP wash, aspirate or swab specimen in transport medium.

Note: The BD Veritor™ System for Rapid Detection of RSV (Cat. # 256042) is required for this procedure in addition to the BD Veritor™ System for Rapid Detection of Flu A+B (Cat. # 256041).

IMPORTANT NOTE: This optional procedure allows for use of the remaining processed sample from Step 5 above to test additionally for RSV. The sample to be tested with the RSV kit must be from a patient less than 20 years of age as indicated in the BD Veritor RSV laboratory kit package insert. **THE PROCESSED SAMPLE SHOULD BE TESTED WITHIN 15 MINUTES.**

Allow specimens that have been stored at 2–8 °C to warm to room temperature prior to use. Thoroughly mix all specimens prior to removal of an aliquot for processing. Do not centrifuge specimens.

1. Collect specimen from the patient and follow Steps 1-5 of the test procedure above to prepare the sample for testing.
2. Using the sample from Step 5, continue the test procedure using the test device for RSV and the same workflow configuration used to obtain the Flu A+B result.
3. Refer to the product insert for **BD Veritor™ System for Rapid Detection of RSV (Cat. # 256042)** for the test procedure and full description of the **BD Veritor RSV** test. Follow the instructions in the insert and the Instrument on screen prompts to complete the test procedure and obtain results. Refer to the product insert for the **BD Veritor System RSV Kit** for result interpretation.

INTERPRETATION OF RESULTS

The **BD Veritor System Instrument** (purchased separately) must be used for interpretation of all test results. Operators should not attempt to interpret assay results directly from the test strip contained within the **BD Veritor System Flu A+B** assay device. With some specimens, up to four lines may be visible on the test device. The Instrument will appropriately interpret the result.

| Display | Interpretation |
|--------------------------|---|
| FLU A: + FLU B: - | Positive Test for Flu A (influenza A antigen present) |
| FLU A: - FLU B: + | Positive Test for Flu B (influenza B antigen present) |
| FLU A: - FLU B: - | Negative Test for Flu A and Flu B (no antigen detected) |
| RESULT INVALID | Result Invalid |
| POSITIVE CONTROL INVALID | Test Invalid. Repeat the test. |
| NEGATIVE CONTROL INVALID | Test Invalid. Repeat the test. |

Invalid Test – If the test is invalid, the **BD Veritor System Instrument** will display “RESULT INVALID” or “CONTROL INVALID” and the test or control must then be repeated. Because true dual positives are exceptionally rare, the **BD Veritor System Instrument** reports dual positive influenza A and influenza B results as “Result Invalid.” Specimens generating a “Result Invalid” should be retested. Upon retesting, if the specimen produces a “Result Invalid” the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

REPORTING OF RESULTS

Positive Test Positive for the presence of influenza A or influenza B antigen. A positive result may occur in the absence of viable virus.

Negative Test Negative for the presence of influenza A and influenza B antigen. Infection due to influenza cannot be ruled-out because the antigen present in the sample may be below the detection limit of the test. A negative test is presumptive and it is recommended that these results be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.

Invalid Test Test result is inconclusive. Do not report results. Repeat the test.

QUALITY CONTROL:

ATTENTION: To utilize the Analyzer's QC documentation capability, specimen barcode scanning must be enabled on an Analyzer equipped with either an InfoScan or InfoSync module. Please refer to the Analyzer *Instructions for Use*, section 4, to choose or change this configuration and for specific QC procedure steps.

Each BD Veritor System Flu A+B device contains both positive and negative internal/procedural controls:

1. The internal positive control validates the immunological integrity of the device, proper reagent function, and assures correct test procedure.
2. The membrane area surrounding test lines functions as a background check on the assay device.
The BD Veritor System Instrument evaluates the positive and negative internal/procedural controls after insertion of each BD Veritor System test device. The BD Veritor System Instrument prompts the operator if a quality issue occurs during assay analysis. Failure of the internal/procedural controls will generate an invalid test result. **NOTE:** The internal controls do not assess proper sample collection technique.

External Positive and Negative Controls:

Flu A positive/B negative and Flu B positive/A negative control swabs are supplied with each kit. These controls provide additional quality control material to assess that the test reagents and the BD Veritor System Instrument perform as expected. Prepare kit control swabs and test using the same procedure (either **Analyze Now** or **Walk Away** mode) as used for patient specimen swabs. When using the barcode scanning feature to document QC procedures, scan the barcode on the control swab packaging when prompted for a Specimen ID.

Your laboratory's standard Quality Control procedures and applicable local, state and/or federal regulations or accreditation requirements dictate the performance of external quality control procedures.

BD recommends external controls be run once for:

- each new kit lot,
- each new operator,
- each new shipment of test kits,
- as required by internal quality control procedures and in accordance with local, state and federal regulations or accreditation requirements.

LIMITATIONS OF THE PROCEDURE

- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- The contents of this kit are to be used for the qualitative detection of influenza type A and B antigens from NP wash, aspirate and swab in transport media specimens.
- The BD Veritor System for Rapid Detection of Flu A+B is capable of detecting both viable and non-viable influenza particles. The BD Veritor System for Rapid Detection of Flu A+B performance depends on antigen load and may not correlate with other diagnostic methods performed on the same specimen.
- Results from the BD Veritor System for Rapid Detection of Flu A+B test should be correlated with the clinical history, epidemiological data and other data available to the clinician evaluating the patient.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of influenza A or influenza B infection, and should be confirmed by viral culture or an FDA-cleared influenza A and B molecular assay.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not identify specific influenza A virus subtypes.
- Negative test results are not intended to rule out other non-influenza viral or bacterial infections.
- Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false positive results during periods of little/no influenza activity when disease prevalence is low. False negative test results are more likely during peak influenza activity when prevalence of disease is high.
- This device has been evaluated for use with human specimen material only.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, influenza A viruses that have undergone minor amino acid changes in the target epitope region.
- The analytical reactivity of this device has not been established for avian or swine origin influenza strains other than those included in the "Strain Reactivity" tables.
- The BD Veritor System Instrument reports dual positive influenza A and influenza B results as "Result Invalid." True dual positives are exceptionally rare. Specimens generating a "Result Invalid" should be retested. Upon retesting, if the specimen produces a "Result Invalid" the user may want to consider other methods to determine whether the sample is positive or negative for influenza virus.

EXPECTED VALUES

The rate of positivity observed in respiratory testing will vary depending on the method of specimen collection, handling/transport system employed, detection method utilized, the time of year, age of the patient, geographic location and most importantly, local disease prevalence.

The overall prevalence observed with an FDA-cleared influenza A and B molecular assay in the U.S. during the 2010–2011 clinical study was 23.9% for influenza A and 7.5% for influenza B. At the clinical site located in Hong Kong, the prevalence observed with the same FDA-cleared influenza A and B molecular assay was 7.2% for influenza A and 3.4% for influenza B. The overall prevalence observed with an FDA-cleared influenza A and B molecular assay in the U.S. during the 2011–2012 clinical study was 31.7% for influenza A and 4.5% for influenza B. At the clinical sites located in Japan, the prevalence observed with the same FDA-cleared influenza A and B molecular assay was 0% for influenza A and 89% for influenza B.

PERFORMANCE CHARACTERISTICS

Explanation of Terms:

PPA: Positive Percent Agreement = $a / (a+c) \times 100\%$

NPA: Negative Percent Agreement = $d / (b+d) \times 100\%$

P: Positive

N: Negative

C.I.: Confidence Interval

| New Test Method | Comparator Method | |
|-----------------|-------------------|-------|
| | P | N |
| P | a | b |
| N | c | d |
| Total | (a+c) | (b+d) |

Clinical Performance NP Washes/Aspirates 2010-2011:

Performance characteristics for the **BD Veritor System for Rapid Detection of Flu A+B** test were established using NP wash/aspirate specimens in multi-center clinical studies conducted at two U.S. trial sites and one Hong Kong trial site during the 2010-2011 respiratory season. A total of 1502 prospective specimens (1002 in the U.S and 500 in Hong Kong) were evaluated using the **BD Veritor System for Rapid Detection of Flu A+B** test and PCR. Five specimens were not evaluable because of data reconciliation issues, an additional 13 were excluded because of insufficient sample volume for reference method testing and 13 samples were excluded as "Result Invalid" (for an invalid rate of 0.9% [13/1484]).

The prospective specimens consisted of NP washes and aspirates from symptomatic patients. 49% of the samples were from females and 51% from males. 56.6% were from patients less than or equal to 5 years of age. 21.9% of the patients tested were in the 6-21 year age group, 5.7% were from 22-59 years of age and 15.8% were obtained from people greater than or equal to age 60 (the patient age was not provided for 0.1% of samples).

The performance of the **BD Veritor System for Rapid Detection of Flu A+B** test was compared to an FDA-cleared Influenza A and B molecular assay (PCR).

The performance is presented in Table 1 below.

Table 1: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for All NP Wash/Aspirate Specimens – All Sites

| Clinical kit: BD Flu A | Reference PCR | | |
|---|---------------|------|-------|
| | P | N | Total |
| P | 224 | 29 | 253 |
| N | 46 | 1172 | 1218 |
| Total | 270 | 1201 | 1471 |
| Reference Method: PCR PPA: 83.0% (95% C.I. 78.0%–87.0%) NPA: 97.6% (95% C.I. 96.6%–98.3%) | | | |

| Clinical kit: BD Flu B | Reference PCR | | |
|---|---------------|------|-------|
| | P | N | Total |
| P | 74 | 3 | 77 |
| N | 17 | 1377 | 1394 |
| Total | 91 | 1380 | 1471 |
| Reference Method: PCR PPA: 81.3% (95% C.I. 72.1%–88.0%) NPA: 99.8% (95% C.I. 99.4%–99.9%) | | | |

An additional 263 frozen retrospective specimens were evaluated with the **BD Veritor System for Rapid Detection of Flu A+B** test. Twelve samples were excluded because there was insufficient sample volume for reference method testing, one sample was excluded as a PCR "Unresolved" and one sample was excluded as "Result Invalid" (for an invalid rate of 0.4% [1/250]). The retrospective specimens consisted of NP washes and aspirates from symptomatic patients. 44.9% of the samples were from females and 55.1% from males. 87.5% were from patients less than or equal to 5 years of age.

The performance is presented in Table 2 below.

Table 2: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for Retrospective NP Wash/Aspirate Specimens

| Clinical kit: BD Flu A | Reference PCR | | | Clinical kit: BD Flu B | Reference PCR | | |
|---|---------------|-----|-------|---|---------------|-----|-------|
| | P | N | Total | | P | N | Total |
| P | 58 | 2 | 60 | P | 29 | 2 | 31 |
| N | 5 | 184 | 189 | N | 10 | 208 | 218 |
| Total | 63 | 186 | 249 | Total | 39 | 210 | 249 |
| Reference Method: PCR PPA: 92.1% (95% C.I. 82.7%–96.6%) NPA: 98.9% (95% C.I. 96.2%–99.7%) | | | | Reference Method: PCR PPA: 74.0% (95% C.I. 58.9%–85.4%) NPA: 99.0% (95% C.I. 96.6%–99.7%) | | | |

Clinical Performance NP Swabs in Transport Media 2011-2012; U.S. and Japan Combined

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at six clinical trial sites located in geographically diverse areas within the United States and five clinical sites in Japan using a total of 292 samples.

The combined results are presented in Table 3 below.

Table 3: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media - U.S. and Japan Combined

| Clinical kit: BD Flu A | Reference PCR | | | Clinical kit: BD Flu B | Reference PCR | | |
|---|---------------|-----|-------|---|---------------|-----|-------|
| | P | N | Total | | P | N | Total |
| P | 52 | 6 | 58 | P | 77 | 2 | 79 |
| N | 12 | 222 | 234 | N | 13 | 200 | 213 |
| Total | 64 | 228 | 292 | Total | 90 | 202 | 292 |
| Reference Method: PCR PPA: 81.3% (95% C.I. 70.0%–88.9%) NPA: 97.4% (95% C.I. 94.4%–98.8%) | | | | Reference Method: PCR PPA: 85.6% (95% C.I. 76.8%–91.4%) NPA: 99.0% (95% C.I. 96.5%–99.7%) | | | |

Clinical Performance NP Swabs in Transport Media 2011-2012; U.S

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at six clinical trial sites located in geographically diverse areas within the United States. A total of 217 prospective specimens were evaluated using the BD Veritor System for Rapid Detection of Flu A+B test and PCR. Two specimens could not be evaluated because of data reconciliation issues, one was eliminated because of an invalid control reading and 13 were excluded because the PCR results were unresolved.

The specimens consisted of NP swabs in transport media from symptomatic patients. 55.8% of the samples were from females and 44.2% from males. 16.1% were from patients less than or equal to 5 years of age, 25.3% were from patients 6-21 years of age, 47.5% were from patients 22-59 years of age and 11.1% were obtained from patients greater than or equal to 60 years of age.

The performance of the BD Veritor System for Rapid Detection of Flu A+B test was compared to an FDA-cleared Influenza A and B molecular assay (PCR).

The results are presented in Table 4 below.

Table 4: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media - U.S.

| Reference PCR | | | |
|---|----|-----|-------|
| Clinical kit: BD Flu A | P | N | Total |
| P | 52 | 6 | 58 |
| N | 12 | 131 | 143 |
| Total | 64 | 137 | 201 |
| Reference Method: PCR PPA: 81.3% (95% C.I. 70.0%–88.9%) NPA: 95.6% (95% C.I. 90.8%–98.0%) | | | |

| Reference PCR | | | |
|---|---|-----|-------|
| Clinical kit: BD Flu B | P | N | Total |
| P | 7 | 0 | 7 |
| N | 2 | 192 | 194 |
| Total | 9 | 192 | 201 |
| Reference Method: PCR PPA: 77.8% (95% C.I. 45.3%–93.7%) NPA: 100% (95% C.I. 98.0%–100%) | | | |

Clinical Performance NP Swabs in Transport Media 2011-2012; Japan

Performance characteristics for the BD Veritor System for Rapid Detection of Flu A+B test were established using NP swabs in transport media in multi-center studies conducted at five clinical trial sites in Japan. A total of 93 prospective specimens were evaluated using the BD Veritor System for Rapid Detection of Flu A+B test and PCR. Two specimens were excluded as the results were undetermined with the comparator assay.

The specimens consisted of NP swabs in transport media from symptomatic patients. 49.5% of the samples were from females and 50.5% from males. 31.2% were from patients less than or equal to 5 years of age, 63.4% were from patients 6-21 years of age and 5.4% were from patients 22-59 years of age (there were no specimens from patients greater than or equal to 60 years of age).

The results are presented in Table 5 below.

Table 5: Summary of the Performance of the BD Veritor System for Rapid Detection of Flu A+B Test Compared to PCR for NP Swabs in Transport Media - Japan

| Reference PCR | | | |
|---|---|----|-------|
| Clinical kit: BD Flu A | P | N | Total |
| P | 0 | 0 | 0 |
| N | 0 | 91 | 91 |
| Total | 0 | 91 | 91 |
| Reference Method: PCR No Data for PPA Calculation NPA: 100% (95% C.I. 95.9%–100%) | | | |

| Reference PCR | | | |
|---|----|----|-------|
| Clinical kit: BD Flu B | P | N | Total |
| P | 70 | 2 | 72 |
| N | 11 | 8 | 19 |
| Total | 81 | 10 | 91 |
| Reference Method: PCR PPA: 86.4% (95% C.I. 77.3%–92.2%) NPA: 80.0% (95% C.I. 49.0%–94.3%) | | | |

Reproducibility

The reproducibility of the BD Veritor System for Rapid Detection of Flu A+B test was evaluated at three clinical laboratory sites in 2010–2011. The reproducibility panel was composed of 30 simulated influenza A or B samples. These included moderate positive samples, low positive samples (near the assay limit of detection), high negative samples (i.e., containing very low concentrations of virus such that positive results occur ~5% of the time) and negative samples. The panel was tested by two operators at each site for five consecutive days. The results are summarized below.

| Reproducibility Results – Percent of Flu A Positives | | | | |
|--|---|---|---|---|
| Sample | Site 1 | Site 2 | Site 3 | Total |
| High negative H1N1 A | 3.3% (1/30) (95% C.I. 0.6%–16.7%) | 0.0% (0/30) (95% C.I. 0.0%–11.3%) | 0.0% (0/30) (95% C.I. 0.0%–11.3%) | 1.1% (1/90) (95% C.I. 0.2%–6.0%) |
| Low positive H1N1 A | 93.3% (28/30) (95% C.I. 78.7%–98.2%) | 86.7% (26/30) (95% C.I. 70.3%–94.7%) | 93.3% (28/30) (95% C.I. 78.7%–98.2%) | 91.1% (82/90) (95% C.I. 83.4%–95.4%) |
| Moderate positive H1N1 A | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 96.7% (29/30) (95% C.I. 83.3%–99.4%) | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 98.9% (89/90) (95% C.I. 94.0%–99.8%) |
| High negative H3N2 A | 16.7% (5/30) (95% C.I. 7.3%–33.6%) | 3.3% (1/30) (95% C.I. 0.6%–16.7%) | 0.0% (0/30) (95% C.I. 0.0%–11.3%) | 6.7% (6/90) (95% C.I. 3.1%–13.8%) |
| Low positive H3N2 A | 93.3% (28/30) (95% C.I. 78.7%–98.2%) | 86.7% (26/30) (95% C.I. 70.3%–94.7%) | 93.3% (28/30) (95% C.I. 78.7%–98.2%) | 91.1% (82/90) (95% C.I. 83.4%–95.4%) |
| Moderate positive H3N2 A | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 96.7% (29/30) (95% C.I. 83.3%–99.4%) | 98.9% (89/90) (95% C.I. 94.0%–99.8%) |
| Negatives | 0.8% (1/120) (95% C.I. 0.1%–4.6%) | 0.0% (0/120) (95% C.I. 0.0%–3.1%) | 0.0% (0/119) (95% C.I. 0.0%–3.1%) | 0.3% (1/359) (95% C.I. 0.0%–1.6%) |

| Reproducibility Results – Percent of Flu B Positives | | | | |
|--|---|---|---|---|
| Sample | Site 1 | Site 2 | Site 3 | Total |
| High negative B | 3.3% (1/30) (95% C.I. 0.6%–16.7%) | 0.0% (0/30) (95% C.I. 0.0%–11.3%) | 0.0% (0/30) (95% C.I. 0.0%–11.3%) | 1.1% (1/90) (95% C.I. 0.2%–6.0%) |
| Low positive B | 90.0% (27/30) (95% C.I. 74.4%–96.5%) | 63.3% (19/30) (95% C.I. 45.5%–78.1%) | 82.8% (24/29) (95% C.I. 65.5%–92.4%) | 78.7% (70/89) (95% C.I. 69.0%–85.9%) |
| Moderate positive B | 96.7% (29/30) (95% C.I. 83.3%–99.4%) | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 100.0% (30/30) (95% C.I. 88.6%–100.0%) | 98.9% (89/90) (95% C.I. 94.0%–99.8%) |
| Negatives | 0.0% (0/210) (95% C.I. 0.0%–1.8%) | 0.0% (0/210) (95% C.I. 0.0%–1.8%) | 0.0% (0/210) (95% C.I. 0.0%–1.8%) | 0.0% (0/630) (95% C.I. 0.0%–0.6%) |

Analytical Studies

Analytical Sensitivity (Limit of Detection)

The limit of detection (LOD) for the BD Veritor System for Rapid Detection of Flu A+B test was established for a total of 8 influenza strains: 5 influenza A and 3 influenza B. The LOD for each strain represents the lowest concentration producing a positivity rate of ≥95% based on testing 20 to 60 replicates.

| Type | Influenza Viral Strain | Calculated LOD (TCID ₅₀ /mL) | Calculated LOD (EID ₅₀ /mL) | No. Positive / Total | % Positive |
|------|--------------------------|---|--|----------------------|------------|
| A | A/Brisbane/10/2007 H3N2 | 7.27 x 10 ² | N/A | 57/60 | 95% |
| A | A/Brisbane/59/2007 H1N1 | 3.30 x 10 ² | N/A | 57/60 | 95% |
| A | A/California/7/2009 H1N1 | 5.00 x 10 ³ | N/A | 57/60 | 95% |
| A | A/Victoria/3/75 H3N2 | 3.11 x 10 ³ | N/A | 59/60 | 98.3% |
| A | A/Anhui/1/2013 H7N9 | N/A | 5.42 x 10 ⁶ | 59/60 | 98.3% |
| B | B/Brisbane/60/2008 | 7.42 x 10 ³ | N/A | 58/60 | 96.7% |
| B | B/Florida/4/2006 | 1.30 x 10 ³ | N/A | 58/60 | 96.7% |
| B | B/Lee/40 | 4.44 x 10 ⁴ | N/A | 20/20 | 100% |

TCID₅₀/mL = 50% Tissue Culture Infectious Dose

EID₅₀/mL = 50% Egg Infectious Dose

Strain Reactivity with Influenza A and B Viruses

The **BD Veritor** System for Rapid Detection of Flu A+B test was evaluated using a panel of influenza strains. Each strain was diluted and tested in triplicate until a point where not all of the replicates were positive. The dilution prior to that is provided in the table below as a minimal detected concentration. All influenza A strains showed positive Flu A test results and negative Flu B test results. Conversely, all of the influenza B strains showed positive Flu B test results and negative Flu A test results.

Although this test has been shown to detect novel avian influenza A (H7N9) and H3N2v cultured viruses the performance characteristics of this device with clinical specimens that are positive for novel avian influenza A (H7N9) and H3N2v influenza viruses have not been established. The **BD Veritor** System Flu A+B assay can distinguish between influenza A and B viruses, but it cannot differentiate influenza A subtypes.

| Strain | Subtype | Minimal Detected Concentration |
|--|---------|--|
| A/Brisbane/59/2007 | H1N1 | 3.3×10^2 TCID ₅₀ /mL* |
| A/California/7/2009 | H1N1 | 5.0×10^3 TCID ₅₀ /mL* |
| A/Denver/1/57 | H1N1 | 4.45×10^4 CEID ₅₀ /mL |
| A/FM/1/47 | H1N1 | 7.91×10^4 CEID ₅₀ /mL |
| A/Fujian-Gulou/1896/2009 | H1N1 | 4.5×10^5 CEID ₅₀ /mL |
| A/Mal/302/54 | H1N1 | 2.22×10^5 CEID ₅₀ /mL |
| A/New Caledonia/20/1999 | H1N1 | 2.5×10^3 TCID ₅₀ /mL |
| A/New Jersey/8/76 | H1N1 | 1.58×10^3 CEID ₅₀ /mL |
| A/NWS/33 | H1N1 | 1.58×10^4 CEID ₅₀ /mL |
| A/PR/8/34 | H1N1 | 6.31×10^2 TCID ₅₀ /mL |
| A/Solomon Island/03/2006 | H1N1 | 2.5×10^3 TCID ₅₀ /mL |
| A/Washington/24/2012 | H1N1 | 3.16×10^4 EID ₅₀ /mL |
| A>Weiss/43 | H1N1 | 7.03×10^6 CEID ₅₀ /mL |
| A/WWS/33 | H1N1 | 7.91×10^2 CEID ₅₀ /mL |
| A/Aichi/2/68 | H3N2 | 7.91×10^3 CEID ₅₀ /mL |
| A/Brisbane/10/2007 | H3N2 | 7.27×10^2 TCID ₅₀ /mL* |
| A/California/02/2014 | H3N2 | 1.45×10^2 TCID ₅₀ /mL |
| A/Hong Kong/8/68 | H3N2 | 8.89×10^4 CEID ₅₀ /mL |
| A/Moscow/10/99 | H3N2 | 5.8×10^6 TCID ₅₀ /mL |
| A/Perth/16/2009 | H3N2 | 1.0×10^6 TCID ₅₀ /mL |
| A/Port Chalmers/1/73 | H3N2 | 3.95×10^4 CEID ₅₀ /mL |
| A/Switzerland/9715293/2013 | H3N2 | 3.25×10^2 TCID ₅₀ /mL |
| A/Texas/50/2012 | H3N2 | 1.75×10^3 TCID ₅₀ /mL |
| A/Wisconsin/67/2005 | H3N2 | 2.5×10^5 TCID ₅₀ /mL |
| A/Victoria/3/75 | H3N2 | 3.11×10^3 TCID ₅₀ /mL* |
| A/Indiana/08/2011 | H3N2v | 1×10^4 TCID ₅₀ /mL |
| A/Indiana/10/2011 | H3N2v | 7.9×10^6 CEID ₅₀ /mL |
| A/Kansas/13/2009 | H3N2v | 1.0×10^3 TCID ₅₀ /mL |
| A/Minnesota/11/2010 | H3N2v | 7.9×10^5 CEID ₅₀ /mL |
| A/Pennsylvania/14/2010 | H3N2v | 1.26×10^6 CEID ₅₀ /mL |
| A/West Virginia/06/2011 | H3N2v | 7.9×10^3 TCID ₅₀ /mL |
| A/Anhui/01/2005 | H5N1 | 0.512 HA |
| A/Vietnam/1203/2004 | H5N1 | 0.512 HA |
| A/Northern Pintail/Washington/40964/2014 | H5N2 | 6.28×10^5 EID ₅₀ /mL |
| A/Pheasant/New Jersey/1355/1998 | H5N2 | 0.256 HA |
| A/Gyrfalcon/Washington/41088-6/2014 | H5N8 | 1.98×10^6 EID ₅₀ /mL |
| A/Mallard/Netherlands/12/2000 | H7N7 | 0.256 HA |
| A/Anhui/1/2013 | H7N9 | 5.42×10^6 CEID ₅₀ /mL* |
| A/Chicken/HongKong/G9/1997 | H9N2 | 1.024 HA |

*Values taken from preceding Analytical Limit of Detection Table

- a. EID₅₀ = 50% Egg Infectious Dose
- b. TCID₅₀ = 50% Tissue Culture Infectious Dose
- c. CEID₅₀ = 50% Chicken Embryo Infectious dose
- d. HA = Hemagglutination Assay

| Strain | Minimal Detected Concentration |
|--|--|
| B/Brazil/178/96 | 2.32 x 10 ⁴ TCID ₅₀ /mL |
| B/Brisbane/33/2008 (Victoria Lineage) | 2.45 x 10 ⁵ CEID ₅₀ /mL |
| B/Brisbane/60/2008 | 7.42 x 10 ³ TCID ₅₀ /mL* |
| B/Brisbane/72/97 | 1.00 x 10 ⁴ TCID ₅₀ /mL |
| B/Canada/548/99 | >0.64 HA |
| B/Egypt/393/99 | >1.28 HA |
| B/Florida/2/2006 | 1.08 x 10 ⁵ TCID ₅₀ /mL |
| B/Florida/4/2006 | 1.30 x 10 ³ TCID ₅₀ /mL* |
| B/Fujian/93/97 | 3.95 x 10 ⁵ TCID ₅₀ /mL |
| B/Fukushima/220/99 | 9.33 x 10 ² TCID ₅₀ /mL |
| B/Guangdong-Liwan/1133/2014 (Yamagata Lineage) | 9.0 x 10 ⁵ CEID ₅₀ /mL |
| B/Guangxi/547/98 | 2.32 x 10 ⁵ TCID ₅₀ /mL |
| B/Hawaii/01/97 | >6.4 HA |
| B/Hong Kong/5/72 | 1.11 x 10 ⁴ CEID ₅₀ /mL |
| B/Hong Kong/219/98 | >1 HA |
| B/Hong Kong/259/2010 (Victoria Lineage) | 1.35 x 10 ⁶ CEID ₅₀ /mL |
| B/Jiangsu/10/2003 | 1.16 x 10 ⁴ TCID ₅₀ /mL |
| B/Johannesburg/5/99 | 3.95 x 10 ⁴ TCID ₅₀ /mL |
| B/Lee/40 | 4.44 x 10 ⁴ CEID ₅₀ /mL* |
| B/Lisbon/03/96 | >0.08 HA |
| B/Malaysia/2506/2004 | 5.0 x 10 ⁴ TCID ₅₀ /mL |

| Strain | Minimal Detected Concentration |
|---|---|
| B/Maryland/1/59 | 3.51 x 10 ² CEID ₅₀ /mL |
| B/Massachusetts/2/2012 (Yamagata Lineage) | 1.25 x 10 ⁶ CEID ₅₀ /mL |
| B/Mass/3/66 | 1.58 x 10 ⁵ CEID ₅₀ /mL |
| B/Montana/5/2012 | 3.14 x 10 ⁵ EID ₅₀ /mL |
| B/Ohio/11/96 | >0.16 HA |
| B/Ohio/1/05 | 1.34 x 10 ⁵ TCID ₅₀ /mL |
| B/Phuket/3073/2013 | 6.08 x 10 ³ TCID ₅₀ /mL |
| B/Puerto Mont/10427/98 | 0.02 HA |
| B/Russia/69 | 3.9 x 10 ² TCID ₅₀ /mL |
| B/Shandong/7/97 | 1.58 x 10 ⁶ CEID ₅₀ /mL |
| B/Shanghai/04/97 | 1.58 x 10 ⁵ TCID ₅₀ /mL |
| B/Shenzhen/135/97 | 3.16 x 10 ⁴ TCID ₅₀ /mL |
| B/Sichuan/116/96 | 0.016 HA |
| B/Taiwan/2/62 | 2.81 x 10 ² CEID ₅₀ /mL |
| B/Texas/06/2011 (Yamagata Lineage) | 6.2 x 10 ⁵ CEID ₅₀ /mL |
| B/Texas/02/2013 (Victoria Lineage) | 2.75 x 10 ⁴ CEID ₅₀ /mL |
| B/Utah/09/2014 (Yamagata Lineage) | 6.3 x 10 ³ CEID ₅₀ /mL |
| B/Victoria/504/00 | 4.64 x 10 ⁴ TCID ₅₀ /mL |
| B/Wisconsin/01/2010 (Yamagata Lineage) | 7.0 x 10 ² CEID ₅₀ /mL |
| B/Yamagata/16/88 | 9.75 x 10 ³ TCID ₅₀ /mL |
| B/Yamanashi/166/98 | 4.88 x 10 ⁴ TCID ₅₀ /mL |

*Values taken from preceding Analytical Limit of Detection Table

- a. EID₅₀ = 50% Egg Infectious Dose
- b. TCID₅₀ = 50% Tissue Culture Infectious Dose
- c. CEID₅₀ = 50% Chicken Embryo Infectious dose
- d. HA = Hemagglutination Assay

Analytical Specificity (Cross Reactivity)

The BD Veritor System for Rapid Detection of Flu A+B test was evaluated with a total of 51 microorganisms. The 37 bacteria and yeast were tested at a target concentration of approximately 10⁷ CFU/mL (CFU – Colony Forming Units) with the exception of *Staphylococcus aureus*, which was tested at a final concentration of 10⁶ CFU/mL. The 14 viruses were evaluated at concentrations of 10³ to 10¹⁰ TCID₅₀/mL. Of the 51 microorganisms tested, none showed cross-reactivity in either the Flu A or Flu B tests.

| |
|------------------------------------|
| <i>Bacteriodes fragilis</i> |
| <i>Bordetella pertussis</i> |
| <i>Candida albicans</i> |
| <i>Chlamydia pneumoniae</i> |
| <i>Corynebacterium diphtherium</i> |
| <i>Escherichia coli</i> |
| <i>Fusobacterium nucleatum</i> |
| <i>Haemophilus influenzae</i> |
| <i>Haemophilus parainfluenzae</i> |
| <i>Kingella kingae</i> |
| <i>Klebsiella pneumoniae</i> |
| <i>Lactobacillus</i> sp. |
| <i>Legionella</i> sp. |
| <i>Moraxella catarrhalis</i> |
| <i>Mycobacterium tuberculosis</i> |
| <i>Mycoplasma pneumoniae</i> |
| <i>Neisseria gonorrhoeae</i> |

| |
|--|
| <i>Neisseria meningitides</i> |
| <i>Neisseria mucosa</i> |
| <i>Neisseria</i> sp. (<i>Neisseria perflaus</i>) |
| <i>Neisseria subflava</i> |
| <i>Peptostreptococcus anaerobius</i> |
| <i>Porphyromonas asaccharolyticus</i> |
| <i>Prevotella oralis</i> |
| <i>Propionibacterium acnes</i> |
| <i>Proteus mirabilis</i> |
| <i>Pseudomonas aeruginosa</i> |
| <i>Serratia marcescens</i> |
| <i>Staphylococcus aureus</i> |
| <i>Staphylococcus epidermidis</i> |
| <i>Streptococcus mutans</i> |
| <i>Streptococcus pneumoniae</i> |
| <i>Streptococcus pyogenes</i> |
| <i>Streptococcus</i> sp. Group C |

| |
|------------------------------------|
| <i>Streptococcus</i> sp. Group G |
| <i>Streptococcus salivarius</i> |
| <i>Veillonella parvula</i> |
| Adenovirus, type 1 |
| Adenovirus, type 7 |
| Cytomegalovirus |
| Enterovirus |
| Epstein Barr Virus |
| HSV Type 1 |
| Human Coronavirus OC43 |
| Human Coronavirus 2229E |
| Human metapneumovirus (HMPV-27 A2) |
| Human Parainfluenza |
| Measles virus |
| Mumps virus |
| Respiratory syncytial virus |
| Rhinovirus |

Interfering Substances

Various substances were evaluated with the BD Veritor System for Rapid Detection of Flu A+B test. These substances included whole blood (2%) and various medications. No interference was noted with this assay for any of the substances tested.

| Substance | Concentration | Substance | Concentration |
|--------------------------|---------------|------------------------------|---------------|
| 4-Acetamidophenol | 10 mg/mL | Homeopathic Allergy Medicine | 10 mg/mL |
| Acetylsalicylic acid | 20 mg/mL | Ibuprofen | 10 mg/mL |
| Albuterol | 0.083 mg/mL | Loratidine | 100 ng/mL |
| Amantadine Hydrochloride | 500 ng/mL | Menthol Throat Lozenges | 10 mg/mL |
| Ayr Saline Nasal Gel | 10 mg/mL | Mometasone | 500 ng/mL |
| Beclomethasone | 500 ng/mL | Mupirocin | 500 ng/mL |
| Budesonide | 500 ng/mL | Oseltamivir | 500 ng/mL |
| Chlorpheniramine maleate | 5 mg/mL | Oxymetazoline | 0.05 mg/mL |
| Dexamethasone | 10 mg/mL | Phenylephrine | 1 mg/mL |
| Dextromethorphan | 10 mg/mL | Pseudoephedrine HCl | 20 mg/mL |
| Diphenhydramine HCl | 5 mg/mL | Purified Mucin Protein | 1 mg/mL |
| Fexofenadine | 500 ng/mL | Ribavirin | 500 ng/mL |
| FluMist | 1% | Rimantadine | 500 ng/mL |
| Flunisolide | 500 ng/mL | Three OTC mouthwashes | 5 % |
| Fluticasone | 500 ng/mL | Tobramycin | 500 ng/mL |
| Four OTC nasal sprays | 10 % | Triamcinolone | 500 ng/mL |
| Four OTC throat drops | 25 % | Whole Blood | 2% |
| Guaiaicol Glyceryl Ether | 20 mg/mL | Zanamivir | 1 mg/mL |

Of the 44 substances tested in this study, none exhibited interfering reactions when tested with influenza A and influenza B positive samples. Based on the data, the substances tested at the indicated concentration levels did not interfere with the BD Veritor System for Rapid Detection of Flu A+B test.

AVAILABILITY

| Cat. No. | Description |
|----------|--|
| 256041 | BD Veritor™ System for Rapid Detection of Flu A+B, 30 tests |
| 256042 | BD Veritor™ System for Rapid Detection of RSV, 30 tests |
| 256051 | BD Veritor™ System Flu A+B Control Swab Set, 10 pairs of swabs |
| 256055 | BD Veritor™ System Reader |
| 256066 | BD Veritor™ Plus Analyzer |
| 256067 | BD Veritor™ InfoSync Module |
| 256068 | BD Veritor™ InfoScan Module |
| 443907 | USB Printer Cable for BD Veritor™ Analyzer |

REFERENCES

1. Simonsen L., Fukuda K, Schonberger LB, Cox NJ. Impact of influenza epidemics on hospitalizations. *J. Infect. Dis.* 2000; *181*:831-7.
2. Thompson WW, Shay DK, Weintraub E, et al. Mortality associated with influenza and respiratory syncytial virus in the United States. *JAMA* 2003; *289*: 179-86.
3. Treanor, J.J., Hayden, F.G., Vrooman, P.S., et al. 2000. Efficacy and safety of the oral neuraminidase inhibitor oseltamivir in treating acute influenza: a randomized controlled trial. *JAMA.* 283:1016-1024.
4. Kaiser, L., Couch, R.B., Galasso, G.J., Glezen, W.P., Webster, R.G., Wright, P.F., and Hayden, F.G. 1999. First international symposium on influenza and other respiratory viruses: summary and overview Kapalua, Maui, Hawaii, December 4-6, 1998. *Antiviral Res.*, *42*:149-176.
5. Cox, N.J., and Bender, C.A. 1995. The molecular epidemiology of influenza viruses. *Virology*, *6*:359-370.
6. Todd, S.J., Minnich, L., and Waner, J.L. 1995. Comparison of rapid immunofluorescence procedure with TestPack RSV and Directigen Flu A for diagnosis of respiratory syncytial virus and influenza A virus. *J. Clin. Microbiol.* *33*:1650-1651.
7. Harris, P.O. 1989. Clinical relevance and efficient detection of seven major respiratory viruses. *ACL.* p. 15-19.
8. McElhaney, J.E., Gravenstein, S., Krause, P., Hooton, J.W., Upshaw, C.M., and Drinka, P. 1998. Assessment of markers of the cell-mediated immune response after influenza virus infection in frail older adults. *Clin. Diag. Lab. Immunol.* *5*:840-844.
9. Fan, J., Henrickson, K.J., and Savatski, L.L. 1998. Rapid simultaneous diagnosis of infections with respiratory syncytial viruses A and B, influenza viruses A and B, and human parainfluenza virus types 1, 2, and 3 by multiplex quantitative reverse transcription-polymerase chain reaction-hybridization assay (hexaplex). *Clin. Infect. Disease* *26*:1397-1402.

10. Wright, K.E., Wilson, G.A.R., Novosad, D., Dimock, C., Tan, D., and Weber, J.M. 1995. Typing and subtyping of influenza viruses in clinical samples by PCR. *J. Clin. Microbiol.*33:1180-1184.
11. Kendal, A.P. 1985. Influenza Viruses. p. 341-357. *Laboratory Diagnosis of Viral Infections*, In H. Lennette, (ed.) Marcel Dekker, Inc., New York.
12. McQuillen, J., Madeley, C.R., and Kendal, A.P. 1985. Monoclonal antibodies for the rapid diagnosis of influenza A and B virus infections by immunofluorescence. *Lancet*.ii: 911-914.
13. Guenther, S.H., and Linnemann, C.C., Jr. 1988. Indirect immunofluorescence assay for rapid diagnosis of influenza virus. *Laboratory Medicine*.19: 581-583.
14. Minnick, L.L., and Ray, C.G. 1986. Early testing of cell cultures for detection of hemadsorbing viruses. *J. Clin. Microbiol.*25: 421-422.
15. Schmidt, N.J., Ota, M., Gallo, D., and Fox, V.L. 1982. Monoclonal antibodies for rapid, strain specific identification of influenza virus isolates. *J. Clin. Microbiol.*16: 763-765.
16. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
17. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.*17:53-80.
18. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S Government Printing Office, Washington, D.C.
19. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). *Office Journal L262*, 17/10/2000, p.021-0045.

Technical Information: In the United States, contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

BD Veritor System

Français

For Rapid Detection of Flu A+B

Pour le diagnostic *in vitro* uniquement.

APPLICATION

Le test BD Veritor System for Rapid Detection of Flu A+B (système pour la détection rapide de Flu A+B) est un dosage immunologique chromatographique rapide conçu pour la détection qualitative et directe des antigènes de nucléoprotéine viraux d'influenza A et B à partir d'échantillons rhino-pharyngés prélevés par lavage, aspiration et écouvillonnage chez des patients présentant des symptômes. Le test BD Veritor System for Rapid Detection of Flu A+B est un test de différenciation, tel que les antigènes viraux d'influenza A peuvent être distingués de ceux d'influenza B à partir d'un même échantillon analysé sur un même dispositif. Le test devrait servir comme aide au diagnostic des infections par les virus influenza A ou B. Un test négatif est présomptif et il est recommandé de confirmer ces résultats par culture sur cellules virales ou un essai moléculaire pour l'influenza A et B agréé par la FDA. En dehors des États-Unis, un test négatif est présomptif et il est recommandé de confirmer ces résultats par culture sur cellules virales ou un essai moléculaire à des fins diagnostiques dans le pays d'utilisation. La FDA n'a pas donné son agrément pour l'utilisation de ce dispositif en dehors des États-Unis. Les résultats négatifs à l'influenza n'écartent pas la possibilité d'une infection virale et ne doivent pas servir de seule base à un traitement ou à d'autres décisions thérapeutiques. Le test n'est pas conçu pour la détection des antigènes de l'influenza de type C.

Les caractéristiques de performances des échantillons rhino-pharyngés prélevés par lavage/aspiration (LA) pour l'influenza A et B ont été établies de janvier à mars 2011 lorsque les virus d'influenza A/2009 H1N1, A/H3N2, lignée B/Victoria et lignée B/Yamagata étaient les principaux virus d'influenza en circulation d'après le rapport *Morbidity and Mortality Weekly Report* du CDC intitulé « Update: Influenza Activity—United States, 2010-2011 Season, and Composition of the 2011-2012 Influenza Vaccine ». Les caractéristiques de performances peuvent varier par rapport à d'autres nouveaux virus d'influenza.

Les caractéristiques de performances des échantillons LA placés dans un milieu de transport pour l'influenza A et B ont été établies de janvier à avril 2012 lorsque les virus d'influenza A/2009 H1N1, A/H3N2, lignée B/Victoria et lignée B/Yamagata étaient les principaux virus d'influenza en circulation d'après le rapport *Morbidity and Mortality Weekly Report* du CDC intitulé « Update: Influenza Activity—United States, 2011-2012 Season, and Composition of the 2012-2013 Influenza Vaccine ». Les caractéristiques de performances peuvent varier par rapport à d'autres nouveaux virus d'influenza.

Si, sur la base des critères cliniques et de dépistage épidémiologique actuellement recommandés par les autorités de santé publique, on soupçonne une infection avec un nouveau virus influenza, il convient de prélever des échantillons en prenant les précautions de contrôle de l'infection appropriées pour des nouveaux virus influenza virulents et de les envoyer aux laboratoires de dépistage nationaux ou locaux pour des tests de confirmation. Dans de tels cas, une culture virale ne doit pas être tentée à moins de disposer d'une unité BSL 3+ pour recevoir et mettre en culture les échantillons.