PROCEDURE

Title: Retyping for History					
Procedure #: 2015BLOODBANK82					
Institution: Highlands Regional Medical Center					
Address: 3600 Highlands Avenue, Sebring Florida	a 33870				
Prepared by: Anita Smith	Date: 6/12/2015				
Title: Laboratory Administrative Director	,				
Accepted by: Date: 6/12/15					
Title: Laboratory Medical Director					
Date Patient Testing Implemented: 1/1/2009					
Review of procedure every two years					
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Discontinued testing date:	·				



Policy Name: Retyping/Redrawing Blood Bank Patients Department: Blood Bank - Lab

Departmental Review: Policy #:

INITIATE DATE DATE REVIEWED/REVISED PAGE 1 of 2

1/2009 02/2013

PURPOSE:

Mistransfusion occurs from misidentification of the intended recipient at the time of collection, labeling error of the pretransfusion testing sample, during laboratory testing or preparation of units to be issued. Misidentification at sample collection occurs approximately once in every 1000 samples and one in every 12,000 transfusions the recipient receives the wrong unit. To reduce risk of misidentification and consequently, of a hemolytic transfusion reaction historic record checks will be done on each blood bank sample received for testing.

POLICY:

Patients without historic data for blood bank testing will have a second specimen drawn to confirm the ABO/Rh of the patient before blood products may be issued for transfusion. Repeat antibody screen is not required.

PROCEDURE:

- 1. All patients requiring blood products must be properly identified and banded with the blood band identification system. This unique number will appear on the patient's armband and specimen as well as any documentation in the laboratory.
- 2. Check Patient's blood bank records for historical data.
- 3. If historical data is not available, the patient must be retyped with a second specimen drawn at a separate phlebotomy.
- 4. Results from ABO type of the second phlebotomy specimen are recorded on the daily blood bank log and in the LIS.
- 5. Any discrepancies will result in the withholding of blood products until resolved.

REFERENCES:

AABB Technical Manual

CAP Requirements: TRM.30575, TRM.40300, TRM.40670, TRM.40820



Policy Name: Rety	ping/Redrawing Blood Bank	Realients	Department:	Blood Bank – Lab	
Departmental Revie	w:	Policy #:			
INITIATE DATE 1/2009		DATE REVIEWED/REVISED PAGE 2 of 2 02/2013			
Reviewed by	Reviewed Date 8.6.14 5.27.15 SAIT		ed by	Reviewed Date	
Initial Implementation	on Date:				
Taken out of Service	ce:				
Reason:					
Reviewed by: Asome Date: 6/10/13					
Reviewed by: Department Supervisor Date: 6/4/13 Department Adm. Director					
Reviewed by:	epartment Chief Technolog	Da ist	te:	·	
Reviewed and Approved by:					