

PROCEDURE

Title: Soft Bank Down Time Procedure
Procedure #: 2015BLOODBANK66


Institution: Highlands Regional Medical Center

Address: 3600 Highlands Avenue, Sebring Florida 33870

Prepared by: Anita Smith

Date: 6/12/2015

Title: Laboratory Administrative Director

Accepted by:  Date: 6/12/15

Title: Laboratory Medical Director

Date Patient Testing Implemented: 5/20/2014

Review of procedure every two years

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Reviewed by: _____ Date: _____

Discontinued testing date: _____



Policy Name: SoftBank Downtime Procedure

Department: Blood Bank

Departmental Review:

Policy #:

INITIATE DATE
05/20/14

DATE REVIEWED/REVISED

PAGE 1 of 2

PURPOSE:

The purpose of this procedure is to manage the blood bank workflow when laboratory information (LIS) is down and to provide a system of recovery for once the LIS is back on line.

PROCEDURE:

1. Print downtime forms
 - a. Computer Downtime Worksheet
 - b. Computer Downtime - Phenotyping
 - c. Computer Downtime - Product Issue Record
 - d. Computer Downtime - Daily QC
2. Access SoftBank backup file to look up patient blood bank history.
3. Print patient order labels.
4. Print and fill out blank transfusion slip.
5. Record results on downtime forms.
6. When the LIS system is on line, gather all forms and log results into SoftBank.
7. Ensure that all forms are reviewed and signed.
8. Save all manual form in the downtime logbook.

REFERENCES:

- SoftBank super user Manual, SCC, FL

Computer Downtime Worksheet

PATIENT NAME:		MEDICAL RECORD NO:		SAMPLE DATE:		TIME:	
BIRTH DATE:		ARM BAND NO:		Data Entered into the computer:			
DOWNTIME BARCODE LABEL		ORDER LABEL		Tech:		Date:	
				Reviewed by:		Date:	

1). Card checked: _____ YES _____ NO _____ 2). Historical ABO/Rh: _____ NO _____ YES _____ ABO/Rh _____ Date _____

3). Significant antibody/-ies: _____ 4). Significant adverse reaction: _____

5. Special instructions: _____

ABO/Rh Forward-Reverse Typing					Antibody Screen																							
anti-A	anti-B	anti-D	D cont.	Al cells	B cells	TYPE/Rh	TECH	IS			37°C			AHG				INTERP	TECH									
								I	II	III	AC	I	II	III	AC	I	CC			II	CC	III	CC	AC	CC			
RETYPE																												

Compatibility Testing										Direct Antiglobulin Testing						
Date	Donor No.	Exp.	Type	IS	37°C	AHG	INTERP	TECH	ANTIGEN	TECH	ANTISERA	REACTION	CC	INTERP	TECH	
											Poly					
											Anti-IgG					
											Anti-C3					

Computer Downtime Worksheet – Phenotyping

Patient Name: _____

ABO/Rh: _____

Medical Record #: _____

Date of Sample: _____

Date Tested: _____

Tech: _____

Antisera	C	E	c	e	C ^w	K	k	Fy ^a	FY ^b	Jk ^a	Jk ^b	S	s	M	N	Le ^a	Le ^b	P ₁															
Patient																																	
Positive Control Cell# / result																																	
Negative Control Cell#/result																																	
Antisera																																	
Manufacturer: Immucor / Ortho																																	
Lot #																																	
Expiration Date																																	

Other Manufacturer: _____

Positive and Negative Control Cells Manufacturer / Lot No: _____

DATE:

Date	Poly IgG	CC IgG	CC IgG	CC C3b,- C3d	CCC	RESULT

Reviewed by/date: _____

Computer Downtime Worksheet – Product Issue Record

PATIENT LABEL	UNIT NO. PROD.	ISSUED TO DATE/TIME LOCATION	VISUAL CHECK AND CONSENT SIGNED Y/N	TECH	ENTERED INTO COMPUTER BY/DATE

Reviewed by/date: _____

Computer Downtime Worksheet – Daily Quality Control

TECH: _____

DATE: _____

Reagent	Lot No.	Exp.	Reagent	Lot No.	Exp.
Anti-A			MTS2		
Anti-B			MTS2P		
Anti-D			MTS A/B/D/R		
Rh Control			MTS Anti-IgG		
A1 Cells			MTS Buffer		
B Cells			0.8%Afirmagen A1/B cells		
N-Hance			Surgiscreen I, II, III		
PeG			Confidence QC		
Saline					
Panoscreen I,II,III					
Check Cells					
Comp. Check Cells					
Anti-IgG					
Anti-IgG,-C3d					
Anti-C3b,-C3d					
corQC Red Cells					
corQC Antisera					

Reviewed by: _____ Date: _____

Computer Downtime Worksheet – Daily Quality Control

TECH: _____

DATE: _____

Reagent	Control	QCTest	A/P	QCInt P/F	Reagent	Control	QCTest	A/P	QCInt P/F
Anti-A	CORC				A-GEL	CONC1			
Anti-B	CORC				B-GEL	CONC1			
Anti-D	CORC				D-GEL	CONC1			
RhCon	CORC				DC-GEL	CONC1			
A1C	CORA				A1C-G	MTS2P			
BC	CORA				BC-G	MTS2P			
Anti-D	A1C/BC				A-GEL	CONC2			
SC1	SALIN				B-GEL	CONC2			
SC2	SALIN				D-GEL	CONC2			
SC3	SALIN				DC-GEL	CONC2			
SC1N	CORA				A1C-G	CONAB			
SC2N	CORA				BC-G	CONAB			
SC3N	CORA				SUR1	MTS2			
SC1P	CORA				SUR2	MTS2			
SC2P	CORA				SUR3	MTS2			
SC3P	CORA				SUR1	CONAB			
IGG	CC				SUR2	CONAB			
POLY	CC				SUR3	CONAB			
C3	C3CC				BUFRG	MTS2			
					BUFRG	CONAB			

Reviewed by: _____ Date: _____



Policy Name: SoftBank Downtime Procedure

Department: Blood Bank

Departmental Review:

Policy #:

INITIATE DATE
05/20/14

DATE REVIEWED/REVISED

PAGE 2 of 2

Reviewed by	Reviewed Date	Reviewed by	Reviewed Date
<i>Jaw</i>	5-25-15		
<i>CA</i>	5-29-15		

Initial Implementation Date: 05/20/14

Reviewed by: *Christel Pondel* Date: 5-19-14
Department Supervisor

Reviewed by: *Regula Lauster* Date: 5/19/14
Department Adm. Director

Reviewed by: NA Date: _____
Department Chief Technologist

Reviewed and Approved by: _____ Date: _____
Department Medical Director