PROCEDURE

Title: Reference Specimens to One Blood					
Procedure #: 2015BLOODBANK91					
Institution: Highlands Regional Medical Center					
Address: 3600 Highlands Avenue, Sebring Florid	a 33870				
Prepared by: Anita Smith	Date: 6/12/2015				
Title: Laboratory Administrative Director					
Accepted by:	M Date: 6/12/15				
Title: Laboratory Medical Director					
Date Patient Testing Implemented: 4/5/2015					
Review of procedu	re every two years				
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:					
Reviewed by:					
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:					
Reviewed by:	Date:				
Reviewed by:	Date:				
Reviewed by:	Date:				
Discontinued testing date:					



Policy Name: Referral and Reference Study Requests

Department: Blood Bank

Departmental Review: Policy #:

INITIATE DATE

DATE REVIEWED/REVISED

PAGE 1 of 3

04/05/2015

PRINCIPLE:

To aid in the resolution of difficult crossmatching, antibody/ies identification problems or to provide specialized services to patients and physicians, which are beyond the scope of our Blood Bank department; specimens are referred to One Blood Inc. One Blood Inc. Reference Laboratory is accredited as an AABB Immunohematology Reference Lab. They provide consultation services, assist in resolving complex antibody problems, offer transfusion recommendations and locate rare or unusual blood units to meet the transfusion needs of the hospital's patients. The Reference Laboratory functions as a channel to the Rare Donor files should units need to be acquired from outside the immediate area. The HLA Laboratory performs routine HLA typing for patients requiring HLA matched products.

SPECIMEN:

I. Reference Lab:

• TSXM, Antibody work-up: 2 to 3 large EDTA blood sample

1 serum sample (if available)

• HDN work-up: EDTA cord blood sample, baby and mother's

EDTA samples

• Transfusion Reaction work-up: Pre and post transfusion reaction samples

Send the bag and tubing if available or labeled

Segments from the transfused units

II. HLA Lab:

• All patient samples: A,B typing:2 yellow (10 ml) ACD tubes

B27 typing: 1 yellow ACD tube

HLA ab screen: 1 plain red top tube (1 ml serum)

Bone Marrow Donor: 2 yellow (10ml) ACD tubes

For HLA testing, keep all tubes at room temperature after collection and send to the HLA lab.

PROCEDURE:

- 1. Fill out the "Reference Laboratory Consultation Request Form" or "HLA Consultation Request Form" completely. The required identification on the tubes must match the request form.
- 2. Call the OBI Reference Lab at (407)226-3824 and fax a copy of the request form and blood bank laboratory data to (407)226-3830 or (407)641-9961.



Policy Name: Referral and Reference Study Requests

Department: Blood Bank

Departmental Review: Policy #:

INITIATE DATE

DATE REVIEWED/REVISED

PAGE 2 of 3

04/05/2015

3. Transport specimen in a leak proof container and send along with request form to OneBlood Reference Lab at 32 West Gore St. Orlando, FL.

4. If possible send the patient sample with the next available OneBlood courier run. For STAT request contact our courier for shipping.

REPORTING RESULTS:

Results of the testing requested will be phoned and faxed to the hospital Blood Bank as soon as they are available. A hard copy of the written report will follow by courier or mail. One copy of the report is to be sent to the patient's physician, another to the patient's chart or medical record and the original filed in the Reference Reports Book in the Blood Bank. In addition, any special notes such as identification of the antibody present are to be noted on the front of the Patient History Card.

PROCEDURE NOTES:

- 1. OneBlood technical staff is available from 7 am to 11 pm Monday through Friday. If assistance is required at other times, a staff member is on call for emergency transfusion request and/or consultation services.
- 2. Samples arriving after 9 pm Monday through Friday and during weekends and holidays will be processed if they are emergencies and cannot be postponed until the next staffed workday.
- 3. Notify the nurse in charge of the delay in blood processing and availability of blood product requested. Document communication with RN.

REFERENCES:

OneBlood Inc.

Reference Laboratory Consultation Request

Facility Name:	Phone:	Fax:
Ordering Physician:	Submitted By:	Sample Draw Date:
Date/ Time Sample Sent:	. 1	Priority: □ STAT □ ASAP □ Routine
Patient Information Last Name:	First Name:	Middle:
		Gender: □Male □Female
		ABO/Rh:
Race:	□Native American □Caucasian □Oth	ner:HGB/HCT:
Medication(s):		
	s: □Yes □No Blood Products:	Date(s):
		Date(s):
		· · · · · · · · · · · · · · · · · · ·
		RhIG: □Yes, Date: □No
Reason for Consult		IVIG: □Yes, Date: □No patible Crossmatch □DAT/Elution
		yping DOther:
Blood Products Blood Product Requested:	Qi.	antity: Crossmatch: □Yes □No
□Leukoreduced □CMV □Other:	Negative - DHgB S Negati	ve □Irradiated □Washed
□Phenotype Matched □Antigen N	Vegative for:	
Submitting Facility Checklist ☐Please attach copies of your Block		
□Notify Laboratory by phone, fax	or electronically of specimen sh	ipment.
□Send at least 2 full 7ml EDTA ar	nd (if available) Serum sample	S.
□Samples must be labeled with fu	ll name, permanent identifier, da	ite drawn and collector's identifier.
□Transport in secure container.		
Comments:		



Document

FORM-26

Version

2



$\begin{array}{c} \textbf{HIGHLANDS REGIONAL MEDICAL CENTER} \\ \textbf{Sebring, FL} \\ \textbf{\textit{Laboratory}} \end{array}$

DOCUMENT CHANGE RECORD

Document Name: Re-	formal of beforence Shidy Reguests.	
Document Section:	ferral + leference Shidy Reguests. Blood Bank	
Author: Maritel	Porfula	
Please circle one of the	following: NEW REVISIO	N ARCHIVE
Effective		9
Description of docume	ent, changes, and rationale: wydate on OB	I I come information
	<u> </u>	<i>V</i>
Formal Training of sta	iff required: None	
Attach email sent to st	aff about new procedures or changes to proced	lure if applicable
Method Validation req	uired (attach documents): Whe	
List any changes to the	E Lab Information system:	
		i
Review and Approval	Signature	Date
Author	Maristel Ponymela	4.5-13
Chief Technologist		
Admin. Lab Director		
Laboratory Director		

Implementation occurs after signature by Laboratory Director.



Policy Name: Referral and Reference Study Requests Department: Blood Bank

Departmental Review:

Policy #:

INITIATE DATE

DATE REVIEWED/REVISED 04/05/2015

PAGE 3 of 3

Reviewed by	Reviewed Date	Reviewed by	Reviewed Date
Troviow od og			
Initial Implementation Da	ate:		
Reviewed hv		Date:	
Reviewed by: Department Supe	ervisor		
Reviewed by:		Date:	
Department Adm			
Reviewed by:		Date:	
Department Chie	f Technologist		
D 1 1 1 4	1	Data	
Reviewed and Approved Department Med		Date.	