

6- month Soft ID Training Checklist

Name: *Dee Dee Contreras*

Standard Operation	Date	Initials	Trainer
	Performed		
Turn on MC75 via red button on upper right corner	5/9/16	DC	✓
Turn on QLN 220 printer via green button on lower left of printer	5/9/16	DC	✓
Tap upper left corner to activate Soft ID software	5/9/16	DC	✓
Enter Log in ID or scan badge	5/9/16	DC	✓
Enter password	5/9/16	DC	✓
Select Assignment, Select "Okay"	5/9/16	DC	✓
-Verify via verbal interview or wristband that this is the correct patient	5/9/16	DC	✓
-Select "Verify"	5/9/16	DC	✓
Review Orders	5/9/16	DC	✓
Select "Print" key	5/9/16	DC	✓
-Select "Print" again	5/9/16	DC	✓
-Labels will print	5/9/16	DC	✓
Draw blood and/or collect specimen	5/9/16	DC	✓
Attach specimen labels to the correct tube or container designated by the	5/9/16	DC	✓
container type on lower right of label	5/9/16	DC	✓
Scan each specimen label	5/9/16	DC	✓
Select "Submit"	5/9/16	DC	✓
Scan patient wristband to Submit	5/9/16	DC	✓
Change Assignment	5/9/16	DC	✓
Select on "Session" at lower left of screen	5/9/16	DC	✓
Select "Change Assignment"	5/9/16	DC	✓
Choose an Assignment	5/9/16	DC	✓
Reschedule Function (used when a draw is unsuccessful)	5/9/16	DC	✓
Follow scenario above through label printing	5/9/16	DC	✓
Scan only labels for which you have successful draws	5/9/16	DC	✓
Hold stylus or cursor over test(s) to be rescheduled	5/9/16	DC	✓
Select reason for rescheduling	5/9/16	DC	✓
Select "1-OK"	5/9/16	DC	✓
No Print Function	5/9/16	DC	✓
Follow standard orders through Select "print" key	5/9/16	DC	✓
Tap on check mark on left side of order to cancel label print of that test.	5/9/16	DC	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/9/16	DC	✓
Hold stylus or cursor over test(s) with draw type to change	5/9/16	DC	✓
Select draw type	5/9/16	DC	✓
Select "1-OK"	5/9/16	DC	✓
Changing Specimen Type for Micro Specimens	5/9/16	DC	✓
When reviewing Micro specimens, Click on "Specimen"	5/9/16	DC	✓
Follow prompts to change Source and/or Site	5/9/16	DC	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5/9/16	DC	DC
Select test to highlight	5/9/16	DC	DC
Select "Action" in lower left corner	5/9/16	DC	DC
Select "Draw History"	5/9/16	DC	DC
Date, Time, Order#, and tube types previously collected appear	5/9/16	DC	DC
View Comments	5/9/16	DC	DC
Scan patient wristband	5/9/16	DC	DC
Select "7-Instr."	5/9/16	DC	DC
Patient comments, stay comments, and order comments appear	5/9/16	DC	DC
View patient information	5/9/16	DC	DC
Scan patient wristband	5/9/16	DC	DC
Select "... " in upper right corner	5/9/16	DC	DC
Name, MR#, Sex, DOB, and Financial number will appear	5/9/16	DC	DC
Add Extra Tubes	5/9/16	DC	DC
Scan patient wristband	5/9/16	DC	DC
Select ""5-Add"	5/9/16	DC	DC
Scroll list of extra tubes and select those you wish to draw	5/9/16	DC	DC
Select "Append" and continue to select tubes until done	5/9/16	DC	DC
Select "1-Add" at bottom of screen	5/9/16	DC	DC
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/9/16	DC	DC
Select test to highlight	5/9/16	DC	DC
Select "Action" in lower left corner	5/9/16	DC	DC
Scan badge or type alternate phlebotomist name	5/9/16	DC	DC
Select "OK"	5/9/16	DC	DC

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Wanda D Davis

Standard Operation	Date	Initials	Trainer
	Performed		
Turn on MC75 via red button on upper right corner	5-10-16	WD	✓
Turn on QLN 220 printer via green button on lower left of printer	5-10-16	WD	✓
Tap upper left corner to activate Soft ID software	5-10-16	WD	✓
Enter Log in ID or scan badge	5-10-16	WD	✓
Enter password	5-10-16	WD	✓
Select Assignment, Select "Okay"	5-10-16	WD	✓
-Verify via verbal interview or wristband that this is the correct patient	5-10-16	WD	✓
-Select "Verify"	5-10-16	WD	✓
Review Orders	5-10-16	WD	✓
Select "Print" key	5-10-16	WD	✓
-Select "Print" again	5-10-16	WD	✓
-Labels will print	5-10-16	WD	✓
Draw blood and/or collect specimen	5-10-16	WD	✓
Attach specimen labels to the correct tube or container designated by the	5-10-16	WD	✓
container type on lower right of label	5-10-16	WD	✓
Scan each specimen label	5-10-16	WD	✓
Select "Submit"	5-10-16	WD	✓
Scan patient wristband to Submit	5-10-16	WD	✓
Change Assignment	5-10-16	WD	✓
Select on "Session" at lower left of screen	5-10-16	WD	✓
Select "Change Assignment"	5-10-16	WD	✓
Choose an Assignment	5-10-16	WD	✓
Reschedule Function (used when a draw is unsuccessful)	5-10-16	WD	✓
Follow scenario above through label printing	5-10-16	WD	✓
Scan only labels for which you have successful draws	5-10-16	WD	✓
Hold stylus or cursor over test(s) to be rescheduled	5-10-16	WD	✓
Select reason for rescheduling	5-10-16	WD	✓
Select "1-OK"	5-10-16	WD	✓
No Print Function	5-10-16	WD	✓
Follow standard orders through Select "print" key	5-10-16	WD	✓
Tap on check mark on left side of order to cancel label print of that test.	5-10-16	WD	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-10-16	WD	✓
Hold stylus or cursor over test(s) with draw type to change	5-10-16	WD	✓
Select draw type	5-10-16	WD	✓
Select "1-OK"	5-10-16	WD	✓
Changing Specimen Type for Micro Specimens	5-10-16	WD	✓
When reviewing Micro specimens, Click on "Specimen"	5-10-16	WD	✓
Follow prompts to change Source and/or Site	5-10-16	WD	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-10-16	WD	✓
Select test to highlight	5-10-16	WD	✓
Select "Action" in lower left corner	5-10-16	WD	✓
Select "Draw History"	5-10-16	WD	✓
Date, Time, Order#, and tube types previously collected appear	5-10-16	WD	✓
View Comments	5-10-16	WD	✓
Scan patient wristband	5-10-16	WD	✓
Select "7-Instr."	5-10-16	WD	✓
Patient comments, stay comments, and order comments appear	5-10-16	WD	✓
View patient information	5-10-16	WD	✓
Scan patient wristband	5-10-16	WD	✓
Select "... " in upper right corner	5-10-16	WD	✓
Name, MR#, Sex, DOB, and Financial number will appear	5-10-16	WD	✓
Add Extra Tubes	5-10-16	WD	✓
Scan patient wristband	5-10-16	WD	✓
Select ""5-Add"	5-10-16	WD	✓
Scroll list of extra tubes and select those you wish to draw	5-10-16	WD	✓
Select "Append" and continue to select tubes until done	5-10-16	WD	✓
Select "1-Add" at bottom of screen	5-10-16	WD	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-10-16	WD	✓
Select test to highlight	5-10-16	WD	✓
Select "Action" in lower left corner	5-10-16	WD	✓
Scan badge or type alternate phlebotomist name	5-10-16	WD	✓
Select "OK"	5-10-16	WD	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Imke Foulakis*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/25/16	IF	JMW
Turn on QLN 220 printer via green button on lower left of printer	5/25/16	IF	JMW
Tap upper left corner to activate Soft ID software	5/25/16	IF	JMW
Enter Log in ID or scan badge	5/25/16	IF	JMW
Enter password	5/25/16	IF	JMW
Select Assignment, Select "Okay"	5/25/16	IF	JMW
-Verify via verbal interview or wristband that this is the correct patient	5/25/16	IF	JMW
-Select "Verify"	5/25/16	IF	JMW
Review Orders	5/25/16	IF	JMW
Select "Print" key	5/25/16	IF	JMW
-Select "Print" again	5/25/16	IF	JMW
-Labels will print	5/25/16	IF	JMW
Draw blood and/or collect specimen	5/25/16	IF	JMW
Attach specimen labels to the correct tube or container designated by the	5/25/16	IF	JMW
container type on lower right of label	5/25/16	IF	JMW
Scan each specimen label	5/25/16	IF	JMW
Select "Submit"	5/25/16	IF	JMW
Scan patient wristband to Submit	5/25/16	IF	JMW
Change Assignment	5/25/16	IF	JMW
Select on "Session" at lower left of screen	5/25/16	IF	JMW
Select "Change Assignment"	5/25/16	IF	JMW
Choose an Assignment	5/25/16	IF	JMW
Reschedule Function (used when a draw is unsuccessful)	5/25/16	IF	JMW
Follow scenario above through label printing	5/25/16	IF	JMW
Scan only labels for which you have successful draws	5/25/16	IF	JMW
Hold stylus or cursor over test(s) to be rescheduled	5/25/16	IF	JMW
Select reason for rescheduling	5/25/16	IF	JMW
Select "1-OK"	5/25/16	IF	JMW
No Print Function	5/25/16	IF	JMW
Follow standard orders through Select "print" key	5/25/16	IF	JMW
Tap on check mark on left side of order to cancel label print of that test.	5/25/16	IF	JMW
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/25/16	IF	JMW
Hold stylus or cursor over test(s) with draw type to change	5/25/16	IF	JMW
Select draw type	5/25/16	IF	JMW
Select "1-OK"	5/25/16	IF	JMW
Changing Specimen Type for Micro Specimens	5/25/16	IF	JMW
When reviewing Micro specimens, Click on "Specimen"	5/25/16	IF	JMW
Follow prompts to change Source and/or Site	5/25/16	IF	JMW

6- month Soft ID Training Checklist

Standard Operation	Date		
	Performed	Initials	Trainer
View Patient Draw History	5/26/16	JH	JMW
Select test to highlight	5/26/16	JH	JMW
Select "Action" in lower left corner	5/26/16	JH	JMW
Select "Draw History"	5/26/16	JH	JMW
Date, Time, Order#, and tube types previously collected appear	5/26/16	JH	JMW
View Comments	5/26/16	JH	JMW
Scan patient wristband	5/26/16	JH	JMW
Select "7-Instr."	5/26/16	JH	JMW
Patient comments, stay comments, and order comments appear	5/26/16	JH	JMW
View patient information		JH	
Scan patient wristband	5/26/16	JH	JMW
Select "... " in upper right corner	5/26/16	JH	JMW
Name, MR#, Sex, DOB, and Financial number will appear	5/26/16	JH	JMW
Add Extra Tubes			
Scan patient wristband	5/26/16	JH	JMW
Select ""5-Add"	5/26/16	JH	JMW
Scroll list of extra tubes and select those you wish to draw	5/26/16	JH	JMW
Select "Append" and continue to select tubes until done	5/26/16	JH	JMW
Select "1-Add" at bottom of screen	5/26/16	JH	JMW
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/26/16	JH	JMW
Select test to highlight	5/26/16	JH	JMW
Select "Action" in lower left corner	5/26/16	JH	JMW
Scan badge or type alternate phlebotomist name	5/26/16	JH	JMW
Select "OK"	5/26/16	JH	JMW

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Sue Forke*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/6/16	SKF	JMW
Turn on QLN 220 printer via green button on lower left of printer	5/6/16	SKF	JMW
Tap upper left corner to activate Soft ID software	5/6/16	SKF	JMW
Enter Log in ID or scan badge	5/6/16	SKF	JMW
Enter password	5/6/16	SKF	JMW
Select Assignment, Select "Okay"	5/6/16	SKF	JMW
-Verify via verbal interview or wristband that this is the correct patient	5/6/16	SKF	JMW
-Select "Verify"	5/6/16	SKF	JMW
Review Orders	5/6/16	SKF	JMW
Select "Print" key	5/6/16	SKF	JMW
-Select "Print" again	5/6/16	SKF	JMW
-Labels will print	5/6/16	SKF	JMW
Draw blood and/or collect specimen	5/6/16	SKF	JMW
Attach specimen labels to the correct tube or container designated by the	5/6/16	SKF	JMW
container type on lower right of label	5/6/16	SKF	JMW
Scan each specimen label	5/6/16	SKF	JMW
Select "Submit"	5/6/16	SKF	JMW
Scan patient wristband to Submit	5/6/16	SKF	JMW
Change Assignment	5/6/16	SKF	JMW
Select on "Session" at lower left of screen	5/6/16	SKF	JMW
Select "Change Assignment"	5/6/16	SKF	JMW
Choose an Assignment	5/6/16	SKF	JMW
Reschedule Function (used when a draw is unsuccessful)	5/6/16	SKF	JMW
Follow scenario above through label printing	5/6/16	SKF	JMW
Scan only labels for which you have successful draws	5/6/16	SKF	JMW
Hold stylus or cursor over test(s) to be rescheduled	5/6/16	SKF	JMW
Select reason for rescheduling	5/6/16	SKF	JMW
Select "1-OK"	5/6/16	SKF	JMW
No Print Function	5/6/16	SKF	JMW
Follow standard orders through Select "print" key	5/6/16	SKF	JMW
Tap on check mark on left side of order to cancel label print of that test.	5/6/16	SKF	JMW
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/6/16	SKF	JMW
Hold stylus or cursor over test(s) with draw type to change	5/6/16	SKF	JMW
Select draw type	5/6/16	SKF	JMW
Select "1-OK"	5/6/16	SKF	JMW
Changing Specimen Type for Micro Specimens	5/6/16	SKF	JMW
When reviewing Micro specimens, Click on "Specimen"	5/6/16	SKF	JMW
Follow prompts to change Source and/or Site	5/6/16	SKF	JMW

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5/6/16	SKF	JMW
Select test to highlight	5/6/16	SKF	JMW
Select "Action" in lower left corner	5/6/16	SKF	JMW
Select "Draw History"	5/6/16	SKF	JMW
Date, Time, Order#, and tube types previously collected appear	5/6/16	SKF	JMW
View Comments	5/6/16	SKF	JMW
Scan patient wristband	5/6/16	SKF	JMW
Select "7-Instr."	5/6/16	SKF	JMW
Patient comments, stay comments, and order comments appear	5/6/16	SKF	JMW
View patient information	5/6/16	SKF	JMW
Scan patient wristband	5/6/16	SKF	JMW
Select "... " in upper right corner	5/6/16	SKF	JMW
Name, MR#, Sex, DOB, and Financial number will appear	5/6/16	SKF	JMW
Add Extra Tubes	5/6/16	SKF	JMW
Scan patient wristband	5/6/16	SKF	JMW
Select ""5-Add"	5/6/16	SKF	JMW
Scroll list of extra tubes and select those you wish to draw	5/6/16	SKF	JMW
Select "Append" and continue to select tubes until done	5/6/16	SKF	JMW
Select "1-Add" at bottom of screen	5/6/16	SKF	JMW
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/6/16	SKF	JMW
Select test to highlight	5/6/16	SKF	JMW
Select "Action" in lower left corner	5/6/16	SKF	JMW
Scan badge or type alternate phlebotomist name	5/6/16	SKF	JMW
Select "OK"	5/6/16	SKF	JMW

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Julie Frackowiak*

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-16-16	JF	✓
Turn on QLN 220 printer via green button on lower left of printer	5-16-16	JF	✓
Tap upper left corner to activate Soft ID software	5-16-16	JF	✓
Enter Log in ID or scan badge	5-16-16	JF	✓
Enter password	5-16-16	JF	✓
Select Assignment, Select "Okay"	5-16-16	JF	✓
-Verify via verbal interview or wristband that this is the correct patient	5-16-16	JF	✓
-Select "Verify"	5-16-16	JF	✓
Review Orders	5-16-16	JF	✓
Select "Print" key	5-16-16	JF	✓
-Select "Print" again	5-16-16	JF	✓
-Labels will print	5-16-16	JF	✓
Draw blood and/or collect specimen	5-16-16	JF	✓
Attach specimen labels to the correct tube or container designated by the	5-16-16	JF	✓
container type on lower right of label	5-16-16	JF	✓
Scan each specimen label	5-16-16	JF	✓
Select "Submit"	5-16-16	JF	✓
Scan patient wristband to Submit	5-16-16	JF	✓
Change Assignment	5-16-16	JF	✓
Select on "Session" at lower left of screen	5-16-16	JF	✓
Select "Change Assignment"	5-16-16	JF	✓
Choose an Assignment	5-16-16	JF	✓
Reschedule Function (used when a draw is unsuccessful)	5-16-16	JF	✓
Follow scenario above through label printing	5-16-16	JF	✓
Scan only labels for which you have successful draws	5-16-16	JF	✓
Hold stylus or cursor over test(s) to be rescheduled	5-16-16	JF	✓
Select reason for rescheduling	5-16-16	JF	✓
Select "1-OK"	5-16-16	JF	✓
No Print Function	5-16-16	JF	✓
Follow standard orders through Select "print" key	5-16-16	JF	✓
Tap on check mark on left side of order to cancel label print of that test.	5-16-16	JF	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-16-16	JF	✓
Hold stylus or cursor over test(s) with draw type to change	5-16-16	JF	✓
Select draw type	5-16-16	JF	✓
Select "1-OK"	5-16-16	JF	✓
Changing Specimen Type for Micro Specimens			✓
When reviewing Micro specimens, Click on "Specimen"	5-16-16	JF	✓
Follow prompts to change Source and/or Site	5-16-16	JF	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-16-16	JR	C
Select test to highlight	5-16-16	JR	C
Select "Action" in lower left corner	5-16-16	JR	C
Select "Draw History"	5-16-16	JR	C
Date, Time, Order#, and tube types previously collected appear	5-16-16	JR	C
View Comments	5-16-16	JR	C
Scan patient wristband	5-16-16	JR	C
Select "7-Instr."	5-16-16	JR	C
Patient comments, stay comments, and order comments appear	5-16-16	JR	C
View patient information	5-16-16	JR	C
Scan patient wristband	5-16-16	JR	C
Select "..." in upper right corner	5-16-16	JR	C
Name, MR#, Sex, DOB, and Financial number will appear	5-16-16	JR	C
Add Extra Tubes	5-16-16	JR	C
Scan patient wristband	5-16-16	JR	C
Select "5-Add"	5-16-16	JR	C
Scroll list of extra tubes and select those you wish to draw	5-16-16	JR	C
Select "Append" and continue to select tubes until done	5-16-16	JR	C
Select "1-Add" at bottom of screen	5-16-16	JR	C
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-16-16	JR	C
Select test to highlight	5-16-16	JR	C
Select "Action" in lower left corner	5-16-16	JR	C
Scan badge or type alternate phlebotomist name	5-16-16	JR	C
Select "OK"	5-16-16	JR	C

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Sybil Her gas*

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/10/16	ly	✓
Turn on QLN 220 printer via green button on lower left of printer	5/10/16	ly	✓
Tap upper left corner to activate Soft ID software	5/10/16	ly	✓
Enter Log in ID or scan badge	5/10/16	ly	✓
Enter password	5/10/16	ly	✓
Select Assignment, Select "Okay"	5/10/16	ly	✓
-Verify via verbal interview or wristband that this is the correct patient	5/10/16	ly	✓
-Select "Verify"	5/10/16	ly	✓
Review Orders	5/10/16	ly	✓
Select "Print" key	5/10/16	ly	✓
-Select "Print" again	5/10/16	ly	✓
-Labels will print	5/10/16	ly	✓
Draw blood and/or collect specimen	5/10/16	ly	✓
Attach specimen labels to the correct tube or container designated by the	5/10/16	ly	✓
container type on lower right of label	5/10/16	ly	✓
Scan each specimen label	5/10/16	ly	✓
Select "Submit"	5/10/16	ly	✓
Scan patient wristband to Submit	5/10/16	ly	✓
Change Assignment	5/10/16	ly	✓
Select on "Session" at lower left of screen	5/10/16	ly	✓
Select "Change Assignment"	5/10/16	ly	✓
Choose an Assignment	5/10/16	ly	✓
Reschedule Function (used when a draw is unsuccessful)	NA	-	✓
Follow scenario above through label printing	NA	-	✓
Scan only labels for which you have successful draws	NA	-	✓
Hold stylus or cursor over test(s) to be rescheduled	NA	-	✓
Select reason for rescheduling	NA	-	✓
Select "1-OK"	NA	-	✓
No Print Function			✓
Follow standard orders through Select "print" key	5/10/16	ly	✓
Tap on check mark on left side of order to cancel label print of that test.	5/10/16	ly	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	NA	-	✓
Hold stylus or cursor over test(s) with draw type to change	NA	-	✓
Select draw type	NA	-	✓
Select "1-OK"	NA	-	✓
Changing Specimen Type for Micro Specimens	NA	-	✓
When reviewing Micro specimens, Click on "Specimen"	I	-	✓
Follow prompts to change Source and/or Site	I	-	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5/10/16	ls	cr
Select test to highlight	5/10/16	ls	cr
Select "Action" in lower left corner	5/10/16	ls	cr
Select "Draw History"	5/10/16	ls	cr
Date, Time, Order#, and tube types previously collected appear	5/10/16	ls	cr
View Comments	5/10/16	ls	cr
Scan patient wristband	5/10/16	ls	cr
Select "7-Instr."	5/10/16	ls	cr
Patient comments, stay comments, and order comments appear	5/10/16	ls	cr
View patient information	5/10/16	ls	cr
Scan patient wristband	5/10/16	ls	cr
Select "... " in upper right corner	5/10/16	ls	cr
Name, MR#, Sex, DOB, and Financial number will appear	5/10/16	ls	cr
Add Extra Tubes	NA	-	cr
Scan patient wristband	NA	-	cr
Select "'5-Add"	NA	-	cr
Scroll list of extra tubes and select those you wish to draw	NA	-	cr
Select "Append" and continue to select tubes until done	NA	-	cr
Select "1-Add" at bottom of screen	NA	-	cr
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/10/16	ls	cr
Select test to highlight	5/10/16	ls	cr
Select "Action" in lower left corner	5/10/16	ls	cr
Scan badge or type alternate phlebotomist name	5/10/16	ls	cr
Select "OK"	5/10/16	ls	cr

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Marla Harris

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/9/16	MAH	SKF
Turn on QLN 220 printer via green button on lower left of printer	5/9/16	MAH	SKF
Tap upper left corner to activate Soft ID software	5/9/16	MAH	SKF
Enter Log in ID or scan badge	5/9/16	MAH	SKF
Enter password	5/9/16	MAH	SKF
Select Assignment, Select "Okay"	5/9/16	MAH	SKF
-Verify via verbal interview or wristband that this is the correct patient	5/9/16	MAH	SKF
-Select "Verify"	5/9/16	MAH	SKF
Review Orders	5/9/16	MAH	SKF
Select "Print" key	5/9/16	MAH	SKF
-Select "Print" again	5/9/16	MAH	SKF
-Labels will print	5/9/16	MAH	SKF
Draw blood and/or collect specimen	5/9/16	MAH	SKF
Attach specimen labels to the correct tube or container designated by the	5/9/16	MAH	SKF
container type on lower right of label	5/9/16	MAH	SKF
Scan each specimen label	5/9/16	MAH	SKF
Select "Submit"	5/9/16	MAH	SKF
Scan patient wristband to Submit	5/9/16	MAH	SKF
Change Assignment	5/9/16	MAH	SKF
Select on "Session" at lower left of screen	5/9/16	MAH	SKF
Select "Change Assignment"	5/9/16	MAH	SKF
Choose an Assignment	5/9/16	MAH	SKF
Reschedule Function (used when a draw is unsuccessful)	5/9/16	MAH	SKF
Follow scenario above through label printing	5/9/16	MAH	SKF
Scan only labels for which you have successful draws	5/9/16	MAH	SKF
Hold stylus or cursor over test(s) to be rescheduled	5/9/16	MAH	SKF
Select reason for rescheduling	5/9/16	MAH	SKF
Select "1-OK"	5/9/16	MAH	SKF
No Print Function	5/9/16	MAH	SKF
Follow standard orders through Select "print" key	5/9/16	MAH	SKF
Tap on check mark on left side of order to cancel label print of that test.	5/9/16	MAH	SKF
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/9/16	MAH	SKF
Hold stylus or cursor over test(s) with draw type to change	5/9/16	MAH	SKF
Select draw type	5/9/16	MAH	SKF
Select "1-OK"	5/9/16	MAH	SKF
Changing Specimen Type for Micro Specimens	5/9/16	MAH	SKF
When reviewing Micro specimens, Click on "Specimen"	5/9/16	MAH	SKF
Follow prompts to change Source and/or Site	5/9/16	MAH	SKF

Marla Harris

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5/9/16	MMH	SKF
Select test to highlight	5/9/16	MMH	SKF
Select "Action" in lower left corner	5/9/16	MMH	SKF
Select "Draw History"	5/9/16	MMH	SKF
Date, Time, Order#, and tube types previously collected appear	5/9/16	MMH	SKF
View Comments	5/9/16	MMH	SKF
Scan patient wristband	5-9-16	MMH	SKF
Select "7-Instr."	5-9-16	MMH	SKF
Patient comments, stay comments, and order comments appear	5-9-16	MMH	SKF
View patient information	5-9-16	MMH	SKF
Scan patient wristband	5-9-16	MMH	SKF
Select "... " in upper right corner	5-9-16	MMH	SKF
Name, MR#, Sex, DOB, and Financial number will appear	5-9-16	MMH	SKF
Add Extra Tubes	5-9-16	MMH	SKF
Scan patient wristband	5-9-16	MMH	SKF
Select ""5-Add"	5-9-16	MMH	SKF
Scroll list of extra tubes and select those you wish to draw	5-9-16	MMH	SKF
Select "Append" and continue to select tubes until done	5-9-16	MMH	SKF
Select "1-Add" at bottom of screen	5-9-16	MMH	SKF
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-9-16	MMH	SKF
Select test to highlight	5-9-16	MMH	SKF
Select "Action" in lower left corner	5-9-16	MMH	SKF
Scan badge or type alternate phlebotomist name	5-9-16	MMH	SKF
Select "OK"	5-9-16	MMH	SKF

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Nancy KERNER*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-10-16	nl	✓
Turn on QLN 220 printer via green button on lower left of printer	5-10-16	nl	✓
Tap upper left corner to activate Soft ID software	5-10-16	nl	✓
Enter Log in ID or scan badge	5-10-16	nl	✓
Enter password	5-10-16	nl	✓
Select Assignment, Select "Okay"	5-10-16	nl	✓
-Verify via verbal interview or wristband that this is the correct patient	5-10-16	nl	✓
-Select "Verify"	5-10-16	nl	✓
Review Orders	5-10-16	nl	✓
Select "Print" key	5-10-16	nl	✓
-Select "Print" again	5-10-16	nl	✓
-Labels will print	5-10-16	nl	✓
Draw blood and/or collect specimen	5-10-16	nl	✓
Attach specimen labels to the correct tube or container designated by the	5-10-16	nl	✓
container type on lower right of label	5-10-16	nl	✓
Scan each specimen label	5-10-16	nl	✓
Select "Submit"	5-10-16	nl	✓
Scan patient wristband to Submit	5-10-16	nl	✓
Change Assignment	5-10-16	nl	✓
Select on "Session" at lower left of screen	5-10-16	nl	✓
Select "Change Assignment"	5-10-16	nl	✓
Choose an Assignment	5-10-16	nl	✓
Reschedule Function (used when a draw is unsuccessful)			✓
Follow scenario above through label printing	5-10-16	nl	✓
Scan only labels for which you have successful draws	5-10-16	nl	✓
Hold stylus or cursor over test(s) to be rescheduled	5-10-16	nl	✓
Select reason for rescheduling	5-10-16	nl	✓
Select "1-OK"	5-10-16	nl	✓
No Print Function			✓
Follow standard orders through Select "print" key	5-10-16	nl	✓
Tap on check mark on left side of order to cancel label print of that test.	5-10-16	nl	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-10-16	nl	✓
Hold stylus or cursor over test(s) with draw type to change	5-10-16	nl	✓
Select draw type	5-10-16	nl	✓
Select "1-OK"	5-10-16	nl	✓
Changing Specimen Type for Micro Specimens	5-10-16	nl	✓
When reviewing Micro specimens, Click on "Specimen"	NA		✓
Follow prompts to change Source and/or Site	NA		✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History			<i>C</i>
Select test to highlight	5-10-16	<i>ml</i>	<i>C</i>
Select "Action" in lower left corner	5-10-16	<i>ml</i>	<i>C</i>
Select "Draw History"	5-10-16	<i>ml</i>	<i>C</i>
Date, Time, Order#, and tube types previously collected appear	5-10-16	<i>ml</i>	<i>C</i>
View Comments			<i>C</i>
Scan patient wristband	5-10-16	<i>ml</i>	<i>C</i>
Select "7-Instr."	5-10-16	<i>ml</i>	<i>C</i>
Patient comments, stay comments, and order comments appear	5-10-16	<i>ml</i>	<i>C</i>
View patient information			<i>C</i>
Scan patient wristband	5-10-16	<i>ml</i>	<i>C</i>
Select "... " in upper right corner	5-10-16	<i>ml</i>	<i>C</i>
Name, MR#, Sex, DOB, and Financial number will appear	5-10-16	<i>ml</i>	<i>C</i>
Add Extra Tubes			<i>C</i>
Scan patient wristband	5-10-16	<i>ml</i>	<i>C</i>
Select ""5-Add"	5-10-16	<i>ml</i>	<i>C</i>
Scroll list of extra tubes and select those you wish to draw	5-10-16	<i>ml</i>	<i>C</i>
Select "Append" and continue to select tubes until done	5-10-16	<i>ml</i>	<i>C</i>
Select "1-Add" at bottom of screen	5-10-16	<i>ml</i>	<i>C</i>
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)			<i>C</i>
Select test to highlight	5-10-16	<i>ml</i>	<i>C</i>
Select "Action" in lower left corner	5-10-16	<i>ml</i>	<i>C</i>
Scan badge or type alternate phlebotomist name	5-10-16	<i>ml</i>	<i>C</i>
Select "OK"	5-10-16	<i>ml</i>	<i>C</i>

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

Euse

6- month Soft ID Training Checklist

Name:

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-10-16	OME	✓
Turn on QLN 220 printer via green button on lower left of printer	5-10-16	OME	✓
Tap upper left corner to activate Soft ID software	5-10-16	OME	✓
Enter Log in ID or scan badge	5-10-16	OME	✓
Enter password	5-10-16	OME	✓
Select Assignment, Select "Okay"	5-10-16	OME	✓
-Verify via verbal interview or wristband that this is the correct patient	5-10-16	OME	✓
-Select "Verify"	5-10-16	OME	✓
Review Orders	5-10-16	OME	✓
Select "Print" key	5-10-16	OME	✓
-Select "Print" again	5-10-16	OME	✓
-Labels will print	5-10-16	OME	✓
Draw blood and/or collect specimen	5-10-16	OME	✓
Attach specimen labels to the correct tube or container designated by the	5-10-16	OME	✓
container type on lower right of label	5-10-16	OME	✓
Scan each specimen label	5-10-16	OME	✓
Select "Submit"	5-10-16	OME	✓
Scan patient wristband to Submit	5-10-16	OME	✓
Change Assignment	5-10-16	OME	✓
Select on "Session" at lower left of screen	5-10-16	OME	✓
Select "Change Assignment"	5-10-16	OME	✓
Choose an Assignment	5-10-16	OME	✓
Reschedule Function (used when a draw is unsuccessful)	5-10-16	OME	✓
Follow scenario above through label printing	5-10-16	OME	✓
Scan only labels for which you have successful draws	5-10-16	OME	✓
Hold stylus or cursor over test(s) to be rescheduled	5-10-16	OME	✓
Select reason for rescheduling	5-10-16	OME	✓
Select "1-OK"	5-10-16	OME	✓
No Print Function	5-10-16	OME	✓
Follow standard orders through Select "print" key	5-10-16	OME	✓
Tap on check mark on left side of order to cancel label print of that test.	5-10-16	OME	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-10-16	OME	✓
Hold stylus or cursor over test(s) with draw type to change	5-10-16	OME	✓
Select draw type	5-10-16	OME	✓
Select "1-OK"	5-10-16	OME	✓
Changing Specimen Type for Micro Specimens	5-10-16	OME	✓
When reviewing Micro specimens, Click on "Specimen"	5-10-16	OME	✓
Follow prompts to change Source and/or Site	5-10-16	OME	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-10-16	DM	✓
Select test to highlight	5-10-16	DM	✓
Select "Action" in lower left corner	5-10-16	DM	✓
Select "Draw History"	5-10-16	DM	✓
Date, Time, Order#, and tube types previously collected appear	5-10-16	DM	✓
View Comments	5-10-16	DM	✓
Scan patient wristband	5-10-16	DM	✓
Select "7-Instr."	5-10-16	DM	✓
Patient comments, stay comments, and order comments appear	5-10-16	DM	✓
View patient information	5-10-16	DM	✓
Scan patient wristband	5-10-16	DM	✓
Select "..." in upper right corner	5-10-16	DM	✓
Name, MR#, Sex, DOB, and Financial number will appear	5-10-16	DM	✓
Add Extra Tubes	5-10-16	DM	✓
Scan patient wristband	5-10-16	DM	✓
Select ""5-Add"	5-10-16	DM	✓
Scroll list of extra tubes and select those you wish to draw	5-10-16	DM	✓
Select "Append" and continue to select tubes until done	5-10-16	DM	✓
Select "1-Add" at bottom of screen	5-10-16	DM	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-10-16	DM	✓
Select test to highlight	5-10-16	DM	✓
Select "Action" in lower left corner	5-10-16	DM	✓
Scan badge or type alternate phlebotomist name	5-10-16	DM	✓
Select "OK"	5-10-16	DM	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Ched Masfin*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-7-16	<i>u</i>	KR
Turn on QLN 220 printer via green button on lower left of printer	5-7	<i>u</i>	KR
Tap upper left corner to activate Soft ID software	5-7	<i>u</i>	KR
Enter Log in ID or scan badge	5-7	<i>u</i>	KR
Enter password	5-7	<i>u</i>	KR
Select Assignment, Select "Okay"	5-7	<i>u</i>	KR
-Verify via verbal interview or wristband that this is the correct patient	5-7	<i>u</i>	KR
-Select "Verify"	5-7	<i>u</i>	KR
Review Orders	5-7	<i>u</i>	KR
Select "Print" key	5-7	<i>u</i>	KR
-Select "Print" again	5-7	<i>u</i>	KR
-Labels will print	5-7	<i>u</i>	KR
Draw blood and/or collect specimen	5-7	<i>u</i>	KR
Attach specimen labels to the correct tube or container designated by the	5-7	<i>u</i>	KR
container type on lower right of label	5-7	<i>u</i>	KR
Scan each specimen label	5-7	<i>u</i>	KR
Select "Submit"	5-7	<i>u</i>	KR
Scan patient wristband to Submit	5-7	<i>u</i>	KR
Change Assignment	5-7	<i>u</i>	KR
Select on "Session" at lower left of screen	5-7	<i>u</i>	KR
Select "Change Assignment"	5-7	<i>u</i>	KR
Choose an Assignment	5-7	<i>u</i>	KR
Reschedule Function (used when a draw is unsuccessful)	5-7	<i>u</i>	KR
Follow scenario above through label printing	5-7	<i>u</i>	KR
Scan only labels for which you have successful draws	5-7	<i>u</i>	KR
Hold stylus or cursor over test(s) to be rescheduled	5-7	<i>u</i>	KR
Select reason for rescheduling	5-7	<i>u</i>	KR
Select "1-OK"	5-7	<i>u</i>	KR
No Print Function	5-7	<i>u</i>	KR
Follow standard orders through Select "print" key	5-7	<i>u</i>	KR
Tap on check mark on left side of order to cancel label print of that test.	5-7	<i>u</i>	KR
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-7	<i>u</i>	KR
Hold stylus or cursor over test(s) with draw type to change	5-7	<i>u</i>	KR
Select draw type	5-7	<i>u</i>	KR
Select "1-OK"	5-7	<i>u</i>	KR
Changing Specimen Type for Micro Specimens	5-7	<i>u</i>	KR
When reviewing Micro specimens, Click on "Specimen"	5-7	<i>u</i>	KR
Follow prompts to change Source and/or Site	5-7	<i>u</i>	KR

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-7	✓	KR
Select test to highlight	5-7	✓	KR
Select "Action" in lower left corner	5-7	✓	KR
Select "Draw History"	5-7	✓	KR
Date, Time, Order#, and tube types previously collected appear	5-7	✓	KR
View Comments	5-7	✓	KR
Scan patient wristband	5-7	✓	KR
Select "7-Instr."	5-7	✓	KR
Patient comments, stay comments, and order comments appear	5-7	✓	KR
View patient information	5-7	✓	KR
Scan patient wristband	5-7	✓	KR
Select "... " in upper right corner	5-7	✓	KR
Name, MR#, Sex, DOB, and Financial number will appear	5-7	✓	KR
Add Extra Tubes	5-7	✓	KR
Scan patient wristband	5-7	✓	KR
Select ""5-Add"	5-7	✓	KR
Scroll list of extra tubes and select those you wish to draw	5-7	✓	KR
Select "Append" and continue to select tubes until done	5-7	✓	KR
Select "1-Add" at bottom of screen	5-7	✓	KR
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-7	✓	KR
Select test to highlight	5-7	✓	KR
Select "Action" in lower left corner	5-7	✓	KR
Scan badge or type alternate phlebotomist name	5-7	✓	KR
Select "OK"	5-7	✓	KR

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Deh Powell*

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/9/16	DGB	<input checked="" type="checkbox"/>
Turn on QLN 220 printer via green button on lower left of printer	5/9/16	DGB	<input checked="" type="checkbox"/>
Tap upper left corner to activate Soft ID software	5/9/16	DGB	<input checked="" type="checkbox"/>
Enter Log in ID or scan badge	5/9/16	DGB	<input checked="" type="checkbox"/>
Enter password	5/9/16	DGB	<input checked="" type="checkbox"/>
Select Assignment, Select "Okay"	5/9/16	DGB	<input checked="" type="checkbox"/>
-Verify via verbal interview or wristband that this is the correct patient	5/9/16	DGB	<input checked="" type="checkbox"/>
-Select "Verify"	5/9/16	DGB	<input checked="" type="checkbox"/>
Review Orders	5/9/16	DGB	<input checked="" type="checkbox"/>
Select "Print" key	5/9/16	DGB	<input checked="" type="checkbox"/>
-Select "Print" again	5/9/16	DGB	<input checked="" type="checkbox"/>
-Labels will print	5/9/16	DGB	<input checked="" type="checkbox"/>
Draw blood and/or collect specimen	5/9/16	DGB	<input checked="" type="checkbox"/>
Attach specimen labels to the correct tube or container designated by the	5/9/16	DGB	<input checked="" type="checkbox"/>
container type on lower right of label	5/9/16	DGB	<input checked="" type="checkbox"/>
Scan each specimen label	5/9/16	DGB	<input checked="" type="checkbox"/>
Select "Submit"	5/9/16	DGB	<input checked="" type="checkbox"/>
Scan patient wristband to Submit	5/9/16	DGB	<input checked="" type="checkbox"/>
Change Assignment	5/9/16	DGB	<input checked="" type="checkbox"/>
Select on "Session" at lower left of screen	5/9/16	DGB	<input checked="" type="checkbox"/>
Select "Change Assignment"	5/9/16	DGB	<input checked="" type="checkbox"/>
Choose an Assignment	5/9/16	DGB	<input checked="" type="checkbox"/>
Reschedule Function (used when a draw is unsuccessful)	5/9/16	DGB	<input checked="" type="checkbox"/>
Follow scenario above through label printing	5/9/16	DGB	<input checked="" type="checkbox"/>
Scan only labels for which you have successful draws	5/9/16	DGB	<input checked="" type="checkbox"/>
Hold stylus or cursor over test(s) to be rescheduled	5/9/16	DGB	<input checked="" type="checkbox"/>
Select reason for rescheduling	5/9/16	DGB	<input checked="" type="checkbox"/>
Select "1-OK"	5/9/16	DGB	<input checked="" type="checkbox"/>
No Print Function	5/9/16	DGB	<input checked="" type="checkbox"/>
Follow standard orders through Select "print" key	5/9/16	DGB	<input checked="" type="checkbox"/>
Tap on check mark on left side of order to cancel label print of that test.	5/9/16	DGB	<input checked="" type="checkbox"/>
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/9/16	DGB	<input checked="" type="checkbox"/>
Hold stylus or cursor over test(s) with draw type to change	5/9/16	DGB	<input checked="" type="checkbox"/>
Select draw type	5/9/16	DGB	<input checked="" type="checkbox"/>
Select "1-OK"	5/9/16	DGB	<input checked="" type="checkbox"/>
Changing Specimen Type for Micro Specimens	5/9/16	DGB	<input checked="" type="checkbox"/>
When reviewing Micro specimens, Click on "Specimen"	5/9/16	DGB	<input checked="" type="checkbox"/>
Follow prompts to change Source and/or Site	5/9/16	DGB	<input checked="" type="checkbox"/>

6- month Soft ID Training Checklist

Standard Operation	Date		
	Performed	Initials	Trainer
View Patient Draw History	5/19/2016	DGP	✓
Select test to highlight	5/19/16	DGP	✓
Select "Action" in lower left corner	5/19/16	DGP	✓
Select "Draw History"	5/19/16	DGP	✓
Date, Time, Order#, and tube types previously collected appear	5/19/16	DGP	✓
View Comments	5/19/16	DGP	✓
Scan patient wristband	5/19/16	DGP	✓
Select "7-Instr."	5/19/16	DGP	✓
Patient comments, stay comments, and order comments appear	5/19/16	DGP	✓
View patient information	5/19/16	DGP	✓
Scan patient wristband	5/19/16	DGP	✓
Select "... " in upper right corner	5/19/16	DGP	✓
Name, MR#, Sex, DOB, and Financial number will appear	5/19/16	DGP	✓
Add Extra Tubes	5/19/16	DGP	✓
Scan patient wristband	5/19/16	DGP	✓
Select ""5-Add"	5/19/16	DGP	✓
Scroll list of extra tubes and select those you wish to draw	5/19/16	DGP	✓
Select "Append" and continue to select tubes until done	5/19/16	DGP	✓
Select "1-Add" at bottom of screen	5/19/16	DGP	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/19/16	DGP	✓
Select test to highlight	5/19/16	DGP	✓
Select "Action" in lower left corner	5/19/16	DGP	✓
Scan badge or type alternate phlebotomist name	5/19/16	DGP	✓
Select "OK"	5/19/16	DGP	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Becky R.*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-18-16	<i>BR</i>	<i>C</i>
Turn on QLN 220 printer via green button on lower left of printer	5-18-16	<i>BR</i>	<i>C</i>
Tap upper left corner to activate Soft ID software	5-18-16	<i>BR</i>	<i>C</i>
Enter Log in ID or scan badge	5-18-16	<i>BR</i>	<i>C</i>
Enter password	5-18-16	<i>BR</i>	<i>C</i>
Select Assignment, Select "Okay"	5-18-16	<i>BR</i>	<i>C</i>
-Verify via verbal interview or wristband that this is the correct patient	5-18-16	<i>BR</i>	<i>C</i>
-Select "Verify"	5-18-16	<i>BR</i>	<i>C</i>
Review Orders	5-18-16	<i>BR</i>	<i>C</i>
Select "Print" key	5-18-16	<i>BR</i>	<i>C</i>
-Select "Print" again	5-18-16	<i>BR</i>	<i>C</i>
-Labels will print	5-18-16	<i>BR</i>	<i>C</i>
Draw blood and/or collect specimen	5-18-16	<i>BR</i>	<i>C</i>
Attach specimen labels to the correct tube or container designated by the	5-18-16	<i>BR</i>	<i>C</i>
container type on lower right of label	5-18-16	<i>BR</i>	<i>C</i>
Scan each specimen label	5-18-16	<i>BR</i>	<i>C</i>
Select "Submit"	5-18-16	<i>BR</i>	<i>C</i>
Scan patient wristband to Submit	5-18-16	<i>BR</i>	<i>C</i>
Change Assignment	5-18-16	<i>BR</i>	<i>C</i>
Select on "Session" at lower left of screen	5-18-16	<i>BR</i>	<i>C</i>
Select "Change Assignment"	5-18-16	<i>BR</i>	<i>C</i>
Choose an Assignment	5-18-16	<i>BR</i>	<i>C</i>
Reschedule Function (used when a draw is unsuccessful)	5-18-16	<i>BR</i>	<i>C</i>
Follow scenario above through label printing	5-18-16	<i>BR</i>	<i>C</i>
Scan only labels for which you have successful draws	5-18-16	<i>BR</i>	<i>C</i>
Hold stylus or cursor over test(s) to be rescheduled	5-18-16	<i>BR</i>	<i>C</i>
Select reason for rescheduling	5-18-16	<i>BR</i>	<i>C</i>
Select "1-OK"	5-18-16	<i>BR</i>	<i>C</i>
No Print Function	5-18-16	<i>BR</i>	<i>C</i>
Follow standard orders through Select "print" key	5-18-16	<i>BR</i>	<i>C</i>
Tap on check mark on left side of order to cancel label print of that test.	5-18-16	<i>BR</i>	<i>C</i>
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-18-16	<i>BR</i>	<i>C</i>
Hold stylus or cursor over test(s) with draw type to change	5-18-16	<i>BR</i>	<i>C</i>
Select draw type	5-18-16	<i>BR</i>	<i>C</i>
Select "1-OK"	5-18-16	<i>BR</i>	<i>C</i>
Changing Specimen Type for Micro Specimens	5-18-16	<i>BR</i>	<i>C</i>
When reviewing Micro specimens, Click on "Specimen"	5-18-16	<i>BR</i>	<i>C</i>
Follow prompts to change Source and/or Site	5-18-16	<i>BR</i>	<i>C</i>

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-18-16	pr	ca
Select test to highlight	5-18-16	pr	ca
Select "Action" in lower left corner	5-18-16	pr	ca
Select "Draw History"	5-18-16	pr	ca
Date, Time, Order#, and tube types previously collected appear	5-18-16	pr	ca
View Comments	5-18-16	pr	ca
Scan patient wristband	5-18-16	pr	ca
Select "7-Instr."	5-18-16	pr	ca
Patient comments, stay comments, and order comments appear	5-18-16	pr	ca
View patient information	5-18-16	pr	ca
Scan patient wristband	5-18-16	pr	ca
Select "..." in upper right corner	5-18-16	pr	ca
Name, MR#, Sex, DOB, and Financial number will appear	5-18-16	pr	ca
Add Extra Tubes	5-18-16	pr	ca
Scan patient wristband	5-18-16	pr	ca
Select ""5-Add"	5-18-16	pr	ca
Scroll list of extra tubes and select those you wish to draw	5-18-16	pr	ca
Select "Append" and continue to select tubes until done	5-18-16	pr	ca
Select "1-Add" at bottom of screen	5-18-16	pr	ca
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-18-16	pr	ca
Select test to highlight	5-18-16	pr	ca
Select "Action" in lower left corner	5-18-16	pr	ca
Scan badge or type alternate phlebotomist name	5-18-16	pr	ca
Select "OK"	5-18-16	pr	ca

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Kim Runca

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-7	KR	✓
Turn on QLN 220 printer via green button on lower left of printer	5-7	KR	✓
Tap upper left corner to activate Soft ID software	5-7	KR	✓
Enter Log in ID or scan badge	5-7	KR	✓
Enter password	5-7	KR	✓
Select Assignment, Select "Okay"	5-7	KR	✓
-Verify via verbal interview or wristband that this is the correct patient	5-7	KR	✓
-Select "Verify"	5-7	KR	✓
Review Orders	5-7	KR	✓
Select "Print" key	5-7	KR	✓
-Select "Print" again	5-7	KR	✓
-Labels will print	5-7	KR	✓
Draw blood and/or collect specimen	5-7	KR	✓
Attach specimen labels to the correct tube or container designated by the	5-7	KR	✓
container type on lower right of label	5-7	KR	✓
Scan each specimen label	5-7	KR	✓
Select "Submit"	5-7	KR	✓
Scan patient wristband to Submit	5-7	KR	✓
Change Assignment	5-7	KR	✓
Select on "Session" at lower left of screen	5-7	KR	✓
Select "Change Assignment"	5-7	KR	✓
Choose an Assignment	5-7	KR	✓
Reschedule Function (used when a draw is unsuccessful)	5-7	KR	✓
Follow scenario above through label printing	5-7	KR	✓
Scan only labels for which you have successful draws	5-7	KR	✓
Hold stylus or cursor over test(s) to be rescheduled	5-7	KR	✓
Select reason for rescheduling	5-7	KR	✓
Select "1-OK"	5-7	KR	✓
No Print Function	5-7	KR	✓
Follow standard orders through Select "print" key	5-7	KR	✓
Tap on check mark on left side of order to cancel label print of that test.	5-7	KR	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-7	KR	✓
Hold stylus or cursor over test(s) with draw type to change	5-7	KR	✓
Select draw type	5-7	KR	✓
Select "1-OK"	5-7	KR	✓
Changing Specimen Type for Micro Specimens	5-7	KR	✓
When reviewing Micro specimens, Click on "Specimen"	5-7	KR	✓
Follow prompts to change Source and/or Site	5-7	KR	✓

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-7	KR	C
Select test to highlight	5-7	KR	C
Select "Action" in lower left corner	5-7	KR	C
Select "Draw History"	5-7	KR	C
Date, Time, Order#, and tube types previously collected appear	5-7	KR	C
View Comments	5-7	KR	C
Scan patient wristband	5-7	KR	C
Select "7-Instr."	5-7	KR	C
Patient comments, stay comments, and order comments appear	5-7	KR	C
View patient information	5-7	KR	C
Scan patient wristband	5-7	KR	C
Select "... " in upper right corner	5-7	KR	C
Name, MR#, Sex, DOB, and Financial number will appear	5-7	KR	C
Add Extra Tubes	5-7	KR	C
Scan patient wristband	5-7	KR	C
Select ""5-Add"	5-7	KR	C
Scroll list of extra tubes and select those you wish to draw	5-7	KR	C
Select "Append" and continue to select tubes until done	5-7	KR	C
Select "1-Add" at bottom of screen	5-7	KR	C
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-7	KR	C
Select test to highlight	5-7	KR	C
Select "Action" in lower left corner	5-7	KR	C
Scan badge or type alternate phlebotomist name	5-7	KR	C
Select "OK"	5-7	KR	C

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name:

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-9-14	MPS	✓
Turn on QLN 220 printer via green button on lower left of printer	5-9-14	MPS	✓
Tap upper left corner to activate Soft ID software	5-9-14	MPS	✓
Enter Log in ID or scan badge	5-9-14	MPS	✓
Enter password	5-9-14	MPS	✓
Select Assignment, Select "Okay"	5-9-14	MPS	✓
-Verify via verbal interview or wristband that this is the correct patient	5-9-14	MPS	✓
-Select "Verify"	5-9-14	MPS	✓
Review Orders	5-9-14	MPS	✓
Select "Print" key	5-9-14	MPS	✓
-Select "Print" again	5-9-14	MPS	✓
-Labels will print	5-9-14	MPS	✓
Draw blood and/or collect specimen	5-9-14	MPS	✓
Attach specimen labels to the correct tube or container designated by the	5-9-14	MPS	✓
container type on lower right of label	5-9-14	MPS	✓
Scan each specimen label	5-9-14	MPS	✓
Select "Submit"	5-9-14	MPS	✓
Scan patient wristband to Submit	5-9-14	MPS	✓
Change Assignment	5-9-14	MPS	✓
Select on "Session" at lower left of screen	5-9-14	MPS	✓
Select "Change Assignment"	5-9-14	MPS	✓
Choose an Assignment	5-9-14	MPS	✓
Reschedule Function (used when a draw is unsuccessful)	5-9-14	MPS	✓
Follow scenario above through label printing	5-9-14	MPS	✓
Scan only labels for which you have successful draws	5-9-14	MPS	✓
Hold stylus or cursor over test(s) to be rescheduled	5-9-14	MPS	✓
Select reason for rescheduling	5-9-14	MPS	✓
Select "1-OK"	5-9-14	MPS	✓
No Print Function	5-9-14	MPS	✓
Follow standard orders through Select "print" key	5-9-14	MPS	✓
Tap on check mark on left side of order to cancel label print of that test.	5-9-14	MPS	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-9-14	MPS	✓
Hold stylus or cursor over test(s) with draw type to change	5-9-14	MPS	✓
Select draw type	5-9-14	MPS	✓
Select "1-OK"	5-9-14	MPS	✓
Changing Specimen Type for Micro Specimens	5-9-14	MPS	✓
When reviewing Micro specimens, Click on "Specimen"	5-9-14	MPS	✓
Follow prompts to change Source and/or Site	5-9-14	MPS	✓

Melissa
Staranowicz

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-9-14	MRS	✓
Select test to highlight	5-9-14	MRS	✓
Select "Action" in lower left corner	5-9-14	MRS	✓
Select "Draw History"	5-9-14	MRS	✓
Date, Time, Order#, and tube types previously collected appear	5-9-14	MRS	✓
View Comments	5-9-14	MRS	✓
Scan patient wristband	5-9-14	MRS	✓
Select "7-Instr."	5-9-14	MRS	✓
Patient comments, stay comments, and order comments appear	5-9-14	MRS	✓
View patient information	5-9-14	MRS	✓
Scan patient wristband	5-9-14	MRS	✓
Select "... " in upper right corner	5-9-14	MRS	✓
Name, MR#, Sex, DOB, and Financial number will appear	5-9-14	MRS	✓
Add Extra Tubes	5-9-14	MRS	✓
Scan patient wristband	5-9-14	MRS	✓
Select ""5-Add"	5-9-14	MRS	✓
Scroll list of extra tubes and select those you wish to draw	5-9-14	MRS	✓
Select "Append" and continue to select tubes until done	5-9-14	MRS	✓
Select "1-Add" at bottom of screen	5-9-14	MRS	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5-9-14	MRS	✓
Select test to highlight	5-9-14	MRS	✓
Select "Action" in lower left corner	5-9-14	MRS	✓
Scan badge or type alternate phlebotomist name	5-9-14	MRS	✓
Select "OK"	5-9-14	MRS	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Crystal Southcombe 5/23/14

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/23	CJS	✓
Turn on QLN 220 printer via green button on lower left of printer	5/23	CJS	✓
Tap upper left corner to activate Soft ID software	5/23	CJS	✓
Enter Log in ID or scan badge	5/23	CJS	✓
Enter password	5/23	CJS	✓
Select Assignment, Select "Okay"	5/23	CJS	✓
-Verify via verbal interview or wristband that this is the correct patient	5/23	CJS	✓
-Select "Verify"	5/23	CJS	✓
Review Orders	5/23	CJS	✓
Select "Print" key	5/23	CJS	✓
-Select "Print" again	5/23	CJS	✓
-Labels will print	5/23	CJS	✓
Draw blood and/or collect specimen	5/23	CJS	✓
Attach specimen labels to the correct tube or container designated by the	5/23	CJS	✓
container type on lower right of label	5/23	CJS	✓
Scan each specimen label	5/23	CJS	✓
Select "Submit"	5/23	CJS	✓
Scan patient wristband to Submit	5/23	CJS	✓
Change Assignment	5/23	CJS	✓
Select on "Session" at lower left of screen	5/23	CJS	✓
Select "Change Assignment"	5/23	CJS	✓
Choose an Assignment	5/23	CJS	✓
Reschedule Function (used when a draw is unsuccessful)	5/23	CJS	✓
Follow scenario above through label printing	5/23	CJS	✓
Scan only labels for which you have successful draws	5/23	CJS	✓
Hold stylus or cursor over test(s) to be rescheduled	5/23	CJS	✓
Select reason for rescheduling	5/23	CJS	✓
Select "1-OK"	5/23	CJS	✓
No Print Function	5/23	CJS	✓
Follow standard orders through Select "print" key	5/23	CJS	✓
Tap on check mark on left side of order to cancel label print of that test.	5/23	CJS	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/23	CJS	✓
Hold stylus or cursor over test(s) with draw type to change	5/23	CJS	✓
Select draw type	5/23	CJS	✓
Select "1-OK"	5/23	CJS	✓
Changing Specimen Type for Micro Specimens	5/23	CJS	✓
When reviewing Micro specimens, Click on "Specimen"	5/23	CJS	✓
Follow prompts to change Source and/or Site	5/23	CJS	✓

6- month Soft ID Training Checklist

Standard Operation	Date		
	Performed	Initials	Trainer
View Patient Draw History	5/23	CJS	✓
Select test to highlight	5/23	CJS	✓
Select "Action" in lower left corner	5/23	CJS	✓
Select "Draw History"	6/23	CJS	✓
Date, Time, Order#, and tube types previously collected appear	5/23	CJS	✓
View Comments	5/23	CJS	✓
Scan patient wristband	5/23	CJS	✓
Select "7-Instr."	6/23	CJS	✓
Patient comments, stay comments, and order comments appear	5/23	CJS	✓
View patient information	6/23	CJS	✓
Scan patient wristband	5/23	CJS	✓
Select "..." in upper right corner	5/23	CJS	✓
Name, MR#, Sex, DOB, and Financial number will appear	5/23	CJS	✓
Add Extra Tubes	6/23	CJS	✓
Scan patient wristband	5/23	CJS	✓
Select ""5-Add"	5/23	CJS	✓
Scroll list of extra tubes and select those you wish to draw	5/23	CJS	✓
Select "Append" and continue to select tubes until done	5/23	CJS	✓
Select "1-Add" at bottom of screen	5/23	CJS	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/23	CJS	✓
Select test to highlight	5/23	CJS	✓
Select "Action" in lower left corner	5/23	CJS	✓
Scan badge or type alternate phlebotomist name	5/23	CJS	✓
Select "OK"	5/23	CJS	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Kim Thompson

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-16	[initials]	[initials]
Turn on QLN 220 printer via green button on lower left of printer	5-16	[initials]	[initials]
Tap upper left corner to activate Soft ID software	5-16	[initials]	[initials]
Enter Log in ID or scan badge	5-16	[initials]	[initials]
Enter password	5-16	[initials]	[initials]
Select Assignment, Select "Okay"	5-16	[initials]	[initials]
-Verify via verbal interview or wristband that this is the correct patient	5-16	[initials]	[initials]
-Select "Verify"	5-16	[initials]	[initials]
Review Orders	5-16	[initials]	[initials]
Select "Print" key	5-16	[initials]	[initials]
-Select "Print" again	5-16	[initials]	[initials]
-Labels will print	5-16	[initials]	[initials]
Draw blood and/or collect specimen	5-16	[initials]	[initials]
Attach specimen labels to the correct tube or container designated by the	5-16	[initials]	[initials]
container type on lower right of label	5-16	[initials]	[initials]
Scan each specimen label	5-16	[initials]	[initials]
Select "Submit"	5-16	[initials]	[initials]
Scan patient wristband to Submit	5-16	[initials]	[initials]
Change Assignment	5-16	[initials]	[initials]
Select on "Session" at lower left of screen	5-16	[initials]	[initials]
Select "Change Assignment"	5-16	[initials]	[initials]
Choose an Assignment	5-16	[initials]	[initials]
Reschedule Function (used when a draw is unsuccessful)	5-16	[initials]	[initials]
Follow scenario above through label printing	5-16	[initials]	[initials]
Scan only labels for which you have successful draws	5-16	[initials]	[initials]
Hold stylus or cursor over test(s) to be rescheduled	5-16	[initials]	[initials]
Select reason for rescheduling	5-16	[initials]	[initials]
Select "1-OK"	5-16	[initials]	[initials]
No Print Function	5-16	[initials]	[initials]
Follow standard orders through Select "print" key	5-16	[initials]	[initials]
Tap on check mark on left side of order to cancel label print of that test.	5-16	[initials]	[initials]
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-16	[initials]	[initials]
Hold stylus or cursor over test(s) with draw type to change	5-16	[initials]	[initials]
Select draw type	5-16	[initials]	[initials]
Select "1-OK"	5-16	[initials]	[initials]
Changing Specimen Type for Micro Specimens	5-16	[initials]	[initials]
When reviewing Micro specimens, Click on "Specimen"	5-16	[initials]	[initials]
Follow prompts to change Source and/or Site	5-16	[initials]	[initials]

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5-16	[Signature]	✓
Select test to highlight	5/16	[Signature]	✓
Select "Action" in lower left corner	5/16	[Signature]	✓
Select "Draw History"	5/16	[Signature]	✓
Date, Time, Order#, and tube types previously collected appear	5/16	[Signature]	✓
View Comments	5/16	[Signature]	✓
Scan patient wristband	5/16	[Signature]	✓
Select "7-Instr."	5/16	[Signature]	✓
Patient comments, stay comments, and order comments appear	5/16	[Signature]	✓
View patient information	5/16	[Signature]	✓
Scan patient wristband	5/16	[Signature]	✓
Select "..." in upper right corner	5/16	[Signature]	✓
Name, MR#, Sex, DOB, and Financial number will appear	5/16	[Signature]	✓
Add Extra Tubes	5/16	[Signature]	✓
Scan patient wristband	5/16	[Signature]	✓
Select ""5-Add"	5/16	[Signature]	✓
Scroll list of extra tubes and select those you wish to draw	5/16	[Signature]	✓
Select "Append" and continue to select tubes until done	5/16	[Signature]	✓
Select "1-Add" at bottom of screen	5/16	[Signature]	✓
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/16	[Signature]	✓
Select test to highlight	5/16	[Signature]	✓
Select "Action" in lower left corner	5/16	[Signature]	✓
Scan badge or type alternate phlebotomist name	5/16	[Signature]	✓
Select "OK"	5/16	[Signature]	✓

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: Matthew Wells

Standard Operation	Date		
	Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5-18-16	mw	✓
Turn on QLN 220 printer via green button on lower left of printer	5-18-16	mw	✓
Tap upper left corner to activate Soft ID software	5-18-16	mw	✓
Enter Log in ID or scan badge	5-18-16	mw	✓
Enter password	5-18-16	mw	✓
Select Assignment, Select "Okay"	5-18-16	mw	✓
-Verify via verbal interview or wristband that this is the correct patient	5-18-16	mw	✓
-Select "Verify"	5-18-16	mw	✓
Review Orders	5-18-16	mw	✓
Select "Print" key	5-18-16	mw	✓
-Select "Print" again	5-18-16	mw	✓
-Labels will print	5-18-16	mw	✓
Draw blood and/or collect specimen	5-18-16	mw	✓
Attach specimen labels to the correct tube or container designated by the	5-18-16	mw	✓
container type on lower right of label	5-18-16	mw	✓
Scan each specimen label	5-18-16	mw	✓
Select "Submit"	5-18-16	mw	✓
Scan patient wristband to Submit	5-18-16	mw	✓
Change Assignment	5-18-16	mw	✓
Select on "Session" at lower left of screen	5-18-16	mw	✓
Select "Change Assignment"	5-18-16	mw	✓
Choose an Assignment	5-18-16	mw	✓
Reschedule Function (used when a draw is unsuccessful)	5-18-16	mw	✓
Follow scenario above through label printing	5-18-16	mw	✓
Scan only labels for which you have successful draws	5-18-16	mw	✓
Hold stylus or cursor over test(s) to be rescheduled	5-18-16	mw	✓
Select reason for rescheduling	5-18-16	mw	✓
Select "1-OK"	5-18-16	mw	✓
No Print Function	5-18-16	mw	✓
Follow standard orders through Select "print" key	5-18-16	mw	✓
Tap on check mark on left side of order to cancel label print of that test.	5-18-16	mw	✓
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5-18-16	mw	✓
Hold stylus or cursor over test(s) with draw type to change	5-18-16	mw	✓
Select draw type	5-18-16	mw	✓
Select "1-OK"	5-18-16	mw	✓
Changing Specimen Type for Micro Specimens	5-18-16	mw	✓
When reviewing Micro specimens, Click on "Specimen"	5-18-16	mw	✓
Follow prompts to change Source and/or Site	5-18-16	mw	✓

6- month Soft ID Training Checklist

Standard Operation	Date		
	Performed	Initials	Trainer
View Patient Draw History	5.18.14	ll	ll
Select test to highlight	5.18.14	ll	ll
Select "Action" in lower left corner	5.18.14	ll	ll
Select "Draw History"	5.18.14	ll	ll
Date, Time, Order#, and tube types previously collected appear	5.18.14	ll	ll
View Comments	5.18.14	ll	ll
Scan patient wristband	5.18.14	ll	ll
Select "7-Instr."	5.18.14	ll	ll
Patient comments, stay comments, and order comments appear	5.18.14	ll	ll
View patient information	5.18.14	ll	ll
Scan patient wristband	5.18.14	ll	ll
Select "... " in upper right corner	5.18.14	ll	ll
Name, MR#, Sex, DOB, and Financial number will appear	5.18.14	ll	ll
Add Extra Tubes	5.18.14	ll	ll
Scan patient wristband	5.18.14	ll	ll
Select ""5-Add"	5.18.14	ll	ll
Scroll list of extra tubes and select those you wish to draw	5.18.14	ll	ll
Select "Append" and continue to select tubes until done	5.18.14	ll	ll
Select "1-Add" at bottom of screen	5.18.14	ll	ll
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5.18.14	ll	ll
Select test to highlight	5.18.14	ll	ll
Select "Action" in lower left corner	5.18.14	ll	ll
Scan badge or type alternate phlebotomist name	5.18.14	ll	ll
Select "OK"	5.18.14	ll	ll

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training

6- month Soft ID Training Checklist

Name: *Jody Wilmer*

Standard Operation	Date Performed	Initials	Trainer
Turn on MC75 via red button on upper right corner	5/6/16	JMW	SKF
Turn on QLN 220 printer via green button on lower left of printer	5/6/16	JMW	SKF
Tap upper left corner to activate Soft ID software	5/6/16	JMW	SKF
Enter Log in ID or scan badge	5/6/16	JMW	SKF
Enter password	5/6/16	JMW	SKF
Select Assignment, Select "Okay"	5/6/16	JMW	SKF
-Verify via verbal interview or wristband that this is the correct patient	5/6/16	JMW	SKF
-Select "Verify"	5/6/16	JMW	SKF
Review Orders	5/6/16	JMW	SKF
Select "Print" key	5/6/16	JMW	SKF
-Select "Print" again	5/6/16	JMW	SKF
-Labels will print	5/6/16	JMW	SKF
Draw blood and/or collect specimen	5/6/16	JMW	SKF
Attach specimen labels to the correct tube or container designated by the	5/6/16	JMW	SKF
container type on lower right of label	5/6/16	JMW	SKF
Scan each specimen label	5/6/16	JMW	SKF
Select "Submit"	5/6/16	JMW	SKF
Scan patient wristband to Submit	5/6/16	JMW	SKF
Change Assignment	5/6/16	JMW	SKF
Select on "Session" at lower left of screen	5/6/16	JMW	SKF
Select "Change Assignment"	5/6/16	JMW	SKF
Choose an Assignment	5/6/16	JMW	SKF
Reschedule Function (used when a draw is unsuccessful)	5/6/16	JMW	SKF
Follow scenario above through label printing	5/6/16	JMW	SKF
Scan only labels for which you have successful draws	5/6/16	JMW	SKF
Hold stylus or cursor over test(s) to be rescheduled	5/6/16	JMW	SKF
Select reason for rescheduling	5/6/16	JMW	SKF
Select "1-OK"	5/6/16	JMW	SKF
No Print Function	5/6/16	JMW	SKF
Follow standard orders through Select "print" key	5/6/16	JMW	SKF
Tap on check mark on left side of order to cancel label print of that test.	5/6/16	JMW	SKF
Change draw type (Used to change from line to venipuncture or from venipuncture to line draw)	5/6/16	JMW	SKF
Hold stylus or cursor over test(s) with draw type to change	5/6/16	JMW	SKF
Select draw type	5/6/16	JMW	SKF
Select "1-OK"	5/6/16	JMW	SKF
Changing Specimen Type for Micro Specimens	5/6/16	JMW	SKF
When reviewing Micro specimens, Click on "Specimen"	5/6/16	JMW	SKF
Follow prompts to change Source and/or Site	5/6/16	JMW	SKF

6- month Soft ID Training Checklist

Standard Operation	Date Performed	Initials	Trainer
View Patient Draw History	5/6/16	JMW	SKF
Select test to highlight	5/6/16	JMW	SKF
Select "Action" in lower left corner	5/6/16	JMW	SKF
Select "Draw History"	5/6/16	JMW	SKF
Date, Time, Order#, and tube types previously collected appear	5/6/16	JMW	SKF
View Comments	5/6/16	JMW	SKF
Scan patient wristband	5/6/16	JMW	SKF
Select "7-Instr."	5/6/16	JMW	SKF
Patient comments, stay comments, and order comments appear	5/6/16	JMW	SKF
View patient information	5/6/16	JMW	SKF
Scan patient wristband	5/6/16	JMW	SKF
Select "..." in upper right corner	5/6/16	JMW	SKF
Name, MR#, Sex, DOB, and Financial number will appear	5/6/16	JMW	SKF
Add Extra Tubes	5/6/16	JMW	SKF
Scan patient wristband	5/6/16	JMW	SKF
Select ""5-Add"	5/6/16	JMW	SKF
Scroll list of extra tubes and select those you wish to draw	5/6/16	JMW	SKF
Select "Append" and continue to select tubes until done	5/6/16	JMW	SKF
Select "1-Add" at bottom of screen	5/6/16	JMW	SKF
Select Alternate Phlebotomist (Used when another is collecting and you are labeling as in CSF, or bone marrow collection, etc.)	5/6/16	JMW	SKF
Select test to highlight	5/6/16	JMW	SKF
Select "Action" in lower left corner	5/6/16	JMW	SKF
Scan badge or type alternate phlebotomist name	5/6/16	JMW	SKF
Select "OK"	5/6/16	JMW	SKF

S: Soft 6 month checklists\SOFT 6 month checklist for Soft ID Training