

HOWARD UNIVERSITY HOSPITAL
DEPARTMENT OF PATHOLOGY AND LABORATORY MEDICINE
STANDARD OPERATING POLICY AND PROCEDURE MANUAL

Special Hematology/102

Leuko EZ- Fecal Lactoferrin

PURPOSE

The *LEUKO EZ VUE*® test is an immunochromatographic test for the qualitative detection of elevated levels of fecal lactoferrin, a marker for fecal leukocytes and an indicator of intestinal inflammation. The *LEUKO EZ VUE*® test detects lactoferrin in liquid, semi-solid, and solid fecal specimens. A positive test result indicates an increased level of fecal lactoferrin and warrants additional testing.

The *LEUKO EZ VUE*® overcomes the problems of microscopy by utilizing immunochromatography technology and provides results in 10 minutes. The assay detects elevated levels of lactoferrin in fecal samples. Lactoferrin is very stable and is not degraded during infections by the toxins of pathogens such as *C. difficile* (6).

PRINCIPLE OF THE TEST

The *LEUKO EZ VUE*® test utilizes rabbit anti-lactoferrin antibodies that are conjugated directly to gold particles. The Membrane Cassette contains two stripes of immobilized antibodies. One stripe contains anti-lactoferrin antibodies. The other, representing a control stripe, contains anti-IgG antibodies. The diluted sample and gold conjugate migrate by capillary action when the sample is added to the well. If elevated lactoferrin is present in the sample, gold conjugate-lactoferrin complexes form and are captured by the immobilized anti-lactoferrin antibodies in the stripe. The lactoferrin-conjugate-antibody complexes appear as a single red line in the test portion of the Results Window. In the control stripe, conjugate binds to the immobilized anti-IgG antibodies, demonstrating correct migration of the sample and conjugate along the membrane. The conjugate-anti-IgG antibodies appear as a single red line in the control portion of the Results Window.

HOWARD UNIVERSITY HOSPITAL

DEPARTMENT OF PATHOLOGY AND LABORATORY MEDICINE

STANDARD OPERATING POLICY AND PROCEDURE MANUAL

Special Hematology/102

Leuko EZ- Fecal Lactoferrin

assay, and that the sample migrated properly through the Membrane Cassette. A clear background in the result area is considered an internal negative control. If the test had been performed correctly and reagents are working properly, the background will be clear to give a discernible result.

External: The reactivity of the LEUKO EZ VUE® test should be verified on receipt using the Positive Control and negative control (Diluent). The Positive Control is supplied with the kit (red-capped bottle). The Positive Control confirms the reactivity of the other reagents associated with the assay. Diluent is used for the negative control.

- The test should not be used if control tests do not produce the correct results.
- Positive and Negative controls run on each day of patient testing.
- Failure of the internal and external controls to produce the expected results suggests the test was not performed correctly (i.e., incorrect volume of reagents added, incorrect incubation temperature or times used, or that reagents were not brought to room temperature prior to testing).
- Repeat the control tests as the first step in determining the cause of the failure

VISUAL INTERPRETATION OF RESULTS

Positive result



Negative result



Invalid test



Invalid test



HOWARD UNIVERSITY HOSPITAL
DEPARTMENT OF PATHOLOGY AND LABORATORY MEDICINE
STANDARD OPERATING POLICY AND PROCEDURE MANUAL

Special Hematology/102

Leuko EZ- Fecal Lactoferrin

SHELF LIFE AND STORAGE

- The expiration date of the kit is given on the outside of the box.
- Expiration dates for each component are listed on the individual labels.
- The kit containing the reagents should be stored between 2° and 30°C (refrigerated or room temperature).
- Membrane Cassettes should be kept in the sealed pouches until used.

LIMITATIONS OF THE PROCEDURE

1. The LEUKO EZ VUE® test detects elevated levels of lactoferrin released from fecal leukocytes as a marker of intestinal inflammation. The test may not be appropriate in immunocompromised persons.
2. The 1:50 dilution of fecal specimen has been evaluated in clinical trials and found to be optimal for fecal dilutions. The use of lower dilutions may result in positive reactions due to the presence of normal lactoferrin levels. Therefore, only the dilution recommended in the brochure should be used.
3. The LEUKO EZ VUE® test has not been clinically evaluated for detecting leukocytes in other types of clinical specimens.
4. The intensity of a positive sample test line does not indicate the amount of lactoferrin or severity of disease.
5. Fecal samples from breast fed infants should not be used with this assay.

Refer to the product insert for performance characteristics

