



Fecal Occult Blood

Identifier:	Fecal Occult Blood Procedure	Version #:	2
Folder:	MML POCT WVD	Type:	Procedure
Subfolder:		Effective on:	08/14/2019

CLINICAL SIGNIFICANCE:

The Fecal Occult Blood test is recommended for use as a diagnostic aide during routine physical examination to monitor for bleeding in patients with iron deficiency anemia or recuperating from surgery, peptic ulcer, ulcerative colitis and other conditions, and in screening programs for colorectal cancer.

PRINCIPLE:

Seracult Slides contain a specially prepared guaiac-impregnated paper. When the stool specimen containing occult blood is applied to the test paper, the hemoglobin comes in contact with the guaiac. Application of Seracult Developer (a stabilized hydrogen peroxide and denatured alcohol solution) creates a guaiac-peroxidase-like reaction, which turns the test paper blue if occult blood is present. This reaction requires that the blood cells be hemolyzed. When blood is present, hemolysis is made possible by substances in the stool – primarily water and salts.

EQUIPMENT:

- Seracult Slides: Seracult slides and developer are to be stored at room temperature, (15-30C) and are stable until the expiration date printed on the box or bottle. Slides should be protected from heat, sunlight, florescent light and ultra-violet radiation.
- Seracult Developer: Aqueous solution of hydrogen peroxide and ethanol should be kept away from heat and light and the bottle tightly capped when not in use.

Do not use Seracult slides or developer after expiration date!

- Gloves
- Wood Application sticks

SPECIMEN:

Fresh stool specimen

Written by:	DIANE HAGERMAN	Approved by (sign.):	
Reviewed by:	MARGARET MORAN	Approved by (name):	BRADLEY ROPP
Reviewed on:	05/07/2019	Approved on:	08/13/2019
Renewed by:		Late date:	
Renewed on:			

Patient Preparation

If possible, the patient should be placed on a red meat-free, high residue diet to avoid the possibility of a false-positive reaction and to help uncover silent lesions, which may bleed only intermittently. This special diet should be started two days before testing and continued through the testing period. The patient should inform the physician if any of the dietary recommendations cause discomfort. The patient should always consult the physician before interrupting or discontinuing any prescription medication.

An alternative to this procedure is to omit the special diet for initial tests and to impose it on patients whose stools yield positive results and are to be retested.

Special Diet: This diet, if employed, should be started two days prior to and continued through the test period.

Consume:

- Generous amounts of cooked and uncooked vegetables such as lettuce, corn and spinach.
- Plenty of fruits, such as plums, grapes and apples.
- Moderate amounts of bran cereal, peanuts and popcorn (high fiber foods).
- Well-cooked pork, poultry and fish.

Do not consume:

- Rare and lightly cooked red meats.
- Cauliflower, horseradish, red radishes, turnips, broccoli and cantaloupe.
- Aspirin and medications which may cause gastrointestinal irritation.
- Vitamin C in excess of 250 mg/day and iron-rich supplements.
- Iron rich supplements.

CALIBRATION:

None required

QUALITY CONTROL:

Performance Control – Internal Control

Each Seracult slide has a **Performance Control Area** to be tested after the patient specimen. Apply one drop of Seracult developer to this area. A blue color is a verification of correct product performance. A lack of blue color indicates that the slide test is not performing to product specifications and the test result may be invalid. The testing should be repeated. Before recollection, check expiration date of the Seracult slide and developer. If recollected test and control area perform in the same manner, contact www.proppermfg.com or call 1-800-832-4300.

PROCEDURE:**1. Preparation of slide**

- a. Collect a small specimen of stool on one end of applicator stick. Apply a very thin smear to one window.
- b. From another part of the stool, collect another small specimen using the same applicator as before and apply to the other window.
- c. Close front flap of slide. Slides may be prepared and developed immediately or once prepared; they may be stored for up to 8 days before being developed. Keep prepared slides away from heat and light.

2. Development of slide

- a. Open perforated cover on back of slide.
- b. Apply two drops of developer solution to each smear in the Specimen Test Area.
- c. Read results within 30-60 seconds.
- d. **Any trace of blue color is positive for occult blood. Color begins to fade after 2-4 minutes.**
- e. Develop the performance control ONLY AFTER SPECIMEN TESTS HAVE BEEN COMPLETED AND INTERPRETED.
- f. Apply only one drop of developer solution to the Performance Control Area.
- g. The performance control test yields a blue color within 30 seconds when the Seracult reagents and developer are reacting according to product specifications.

RESULTS:

Results must be obtained visually **30 – 60 seconds** after application of the developer solution.

POSITIVE
NEGATIVE

ANY trace of blue color in the Specimen Test Area
No presence of blue color in the Specimen Test Area

Since any developed blue color may fade after 2-4 minutes, it is imperative that the developed slide is read within the recommended time period.

Note: Because this test is visually read and requires color differentiation it should not be interpreted by color impaired personnel.

Record the Quality control and patient test results in the patient chart or EMR (Electronic Medical Record).

LIMITATIONS OF THE PROCEDURE:

Interfering Substances:

- Vitamin C in excess of 250 mg/day (False negative).
- Oral iron preparations (False positive).
- Certain Oral medications may cause gastrointestinal irritations and bleeding. Aspirin, indomethacin, phenylbutazone, and corticosteroid and reserpine should, with physician consent, be discontinued for two days prior to and during the testing period.

Patients experiencing hemorrhoidal bleeding, having a menstrual period, or bleeding from the nose, gums, etc. should delay testing for at least 48 hours from the time that all such bleeding as stopped.

To increase the chances of detecting intermittent gastrointestinal bleeding, it is recommended that stool samples be collected from three consecutive bowel movements and that two smears are made from two different areas of each bowel movement.

REFERENCES:

1. Package Insert Seracult Propper Mfg. Co., Inc. Long Island City, N.Y. 11101, January 2010

6/17 Removed initials/date results
3/19 QC in patient EMR