

Piccolo Xpress chemistry analyzer Operator's Manual

for In-Vitro Diagnostic Use Only

**Customer and Technical Support:
1-800-822-2947**

10/2020



Abaxis, Inc.
3240 Whipple Road
Union City, CA 94587
USA



ABAXIS Europe GmbH
Bunsenstr. 9-11
64347 Griesheim
Germany

Customer and Technical Support:
USA/Canada: 1-800-822-2947 or AB-PiccoloTechSupport@zoetis.com
International: +49 (6155) 780 21 300 or int.piccolosupport@zoetis.com

PN: 1100-7108-1 Manual Text, Rev.E

© 2020, Abaxis, Inc.

Piccolo is a registered trademark of Abaxis, Inc.




10.5.2 Running HyperTerminal

1. Prepare HyperTerminal to capture results by loading a saved session file or by creating a new session file as described in “[Capturing Results with HyperTerminal](#)” on page 10-9.
2. Click **Transfer > Capture Text** on the menu. Select or type in a path and file name, then click **Start**.



Any results transmitted to the computer from this point forward will be captured in the designated text file.

3. Transmit results to the computer following the instructions in “[Transmitting All Results](#)” on page 5-9. The results will appear in the HyperTerminal window, as shown in this example. In addition, the results can be manually copied and pasted into other applications by highlighting the results and using the **Edit** menu.



Note: *If the results do not appear in the HyperTerminal window, an incorrect port may have been selected in setup. Repeat the setup and select a different port as described in “[Capturing Results with HyperTerminal](#)” on page 10-9.*

4. After the results are transmitted, select **Transfer > Capture Text > Stop** in the menu. The captured text file can be used by most word processing, spreadsheet, and database applications.
5. Close HyperTerminal by selecting **File > Exit** in the menu.