

 <input checked="" type="checkbox"/> <b>POLICY</b> <input checked="" type="checkbox"/> <b>PROCEDURE</b> <input type="checkbox"/> <b>STANDARD OF CARE</b> <input type="checkbox"/> <b>STANDARDIZED PROCEDURE</b> <input type="checkbox"/> <b>GUIDELINE</b> <input type="checkbox"/> <b>OTHER</b>	<b>APPROVAL DATE</b>  <b>4/3/2019</b>	<b>MANUAL:</b> <b>Laboratory Administrative</b>
		<b>TRACKING #</b> <b>AD038</b>
	<b>TITLE:</b> <b>Safe Handling of Dry Ice</b>	
	<b>PERFORMED BY:</b> <b>Laboratory Staff</b>	

## 1. PURPOSE

Dry ice can be hazardous if not handled properly. This procedure describes procedures for safely handling, storage, and disposal of dry ice.

## 2. PROPERTIES AND HAZARDS

- 2.1. Dry ice is the solid form of Carbon Dioxide and is available in flakes, pellets or block form. It has a temperature of -109.3 degrees<sup>o</sup> F (-78.5<sup>o</sup>C) and can cause burns or frostbite to skin in short periods of time.
- 2.2. Dry ice will sublimate at a temperature of -109.3<sup>o</sup>F / -78.5<sup>o</sup> C or higher at approximately 5 to 10 pounds every 24 hours (blocks last longer) in a typical storage cooler.
- 2.3. Use of dry ice in poorly ventilated areas can result in the depletion of the oxygen level resulting in asphyxiation.
- 2.4. Placing dry ice into a tightly sealed container can produce sufficient gas build up to cause an explosion.

## 3. SAFE HANDLING, STORAGE, AND DISPOSAL

- 3.1. Because dry ice sublimates into large quantities of carbon dioxide (CO<sub>2</sub>) gas, it must be stored and used in well-ventilated areas.
- 3.2. Store dry ice in an ultra-low freezer (-70<sup>o</sup>C or below). Do not store dry ice in confined areas such as a walk-in refrigerator.

- 3.3. Small quantities of dry ice for benchtop use may be stored in a styrofoam container with a lid loosely in place at all times.
- 3.4. The following Personal Protective Equipment must be worn when handling dry ice:
  - 3.4.1. Eyewear/safety goggles or faceshield.
  - 3.4.2. Insulated gloves. Nitrile gloves will not provide enough protection. Never handle dry ice with bare hands. Using bare hands can result in burns/frostbite to the skin in a short period of time. Use a scoop or tongs to handle dry ice when possible.
  - 3.4.3. Lab coat or apron (arms and legs must be covered).
- 3.5. To dispose of unneeded dry ice, allow the unused portion sublimate in a well-ventilated area. Sublimation may take several hours or days depending on the volume.
  - 3.5.1. Never dispose of dry ice in a sink, toilet, drain, or trash container.
  - 3.5.2. Never leave surplus dry ice in the hallways or outside of laboratory work areas.

#### **4. TRAINING**

- 4.1. Personnel that handle dry ice are trained on the hazards and safe handling, storage, and disposal.
- 4.2. Personnel responsible for shipping packages containing dry ice must be properly trained in the Department of Transportation (DOT) and/or International Air Transportation Association shipping requirements.

#### **5. REFERENCES**

OSHA Quick Facts Laboratory Safety - Dry Ice and Cryogens.

SDS: Dry Ice (Proxair SDS P-4575).