



Ortho VISION Maintenance Record

Instructions for Use: The date of each activity must be documented. Add the initials of the operator in the box to document completion of the task(s).

Site (Select one)	<input type="radio"/> RHBB <input type="radio"/> MHBB <input type="radio"/> UHBB	Model (Select one)	<input type="radio"/> VISION Max <input type="radio"/> VISION Swift		Month																													
		Name/Serial Number							Year																									
DAILY MAINTENANCE		DATE	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Daily Probe Maintenance																																		
Perform Liquid System Monitoring (empty waste and re-fill saline)																																		
Rotate Diluent and Reagent RBC (Every 24 hrs.) <i>Visually verify the outdate of diluent and diluent is clear without clouding or particulate matter.</i>																																		
Verify Health Check Report (Status OK)																																		
Run Daily QC Profiles: ABORH, ABO3CELL, DONOR, ABORH2, 3CELL, Weak D, Weak D2, Buffered Gel Card, and Diluent 2/XM QC																																		
No Testing Performed on the analyzer – Initial here when no testing is performed or analyzer is out of service																																		
Daily Quality Control Peer Review																																		

WEEKLY MAINTENANCE	WEEK 1		WEEK 2		WEEK 3		WEEK 4		WEEK 5	
	Initials	Date								
Liquid System Decontamination and Pump Test										
Perform Archive Backup										
Weekly Instrument Cleaning										

MONTHLY MAINTENANCE	Operator's Initials	Date
Monthly Instrument Cleaning		

Reviewed by: _____ Date: _____