Title	Units	Normal
WBC	x10^9/L	9.9
PLT	x10^9/L	260
НВ	g/L	131
RBC	x10^12/L	4.27
MCV	fL	86
НСТ	L/L	0.367
МСН	pg	30.7
МСНС	g/L	357
RDW	%	13.6

Presented are the parameters of a NORMAL FBE.

HB, RBC, MCV and RDW are measured directly on the analyser.

The **HCT** is calculated from RBC x MCV / 1000.

The **MCH** is calculated from HB / RBC.

The **MCHC** calculation is expressed in a number of ways

- a) HB / (RBC x MCV) x 1000
- b) HB / (HCT)
- c) MCH / MCV x 1000