






# Assay Technical Training: Xpert<sup>®</sup> CT/NG

For CE-IVD and US-IVD Use



 In Vitro Diagnostic Medical Device

  In Vitro Diagnostic Medical Device



# Training Agenda

- **Xpert® CT/NG Training**

- Intended use
- Reagents
- Sample collection
- Kit storage and handling
- Preparing the cartridge
- Quality Controls
- Results Analysis

- **Discussion**





# Training Objectives

- **At the end of the training, users will be able to:**
  - Store and handle the Xpert® CT/NG cartridge kit and Sample collection kits
  - Follow proper laboratory safety precautions
  - Collect appropriate specimen types and transport specimens
  - Perform the cartridge set up and run the assay
  - Report the various software-generated results
  - Understand the assay control strategy



# The Cepheid Solution



- Simultaneous detection of CT and a dual target NG (NG2/NG4)
- On-board controls for each individual sample
  - Probe Check Control (PCC)
  - Specimen Processing Control (SPC)
  - Sample Adequacy Control (SAC)
- Results in approximately 90 minutes
- Closed cartridge system minimizes risk of contamination
- On-demand results
- Random access



# Intended Use

The Xpert<sup>®</sup> CT/NG Assay is an automated in vitro diagnostic test for qualitative detection and differentiation of DNA from *Chlamydia trachomatis* (CT) and *Neisseria gonorrhoeae* (NG).

The assay may be used to test the following specimens from asymptomatic and symptomatic individuals:

- **Specimen**

- Urine Specimen (male and female)
- Pharyngeal Swab Specimen (male and female)
- Rectal Specimen (male and female)
- Patient-collected Vaginal Specimen (collected in a clinical setting)
- Endocervical Specimen

- **Detects:**

- CT target sequence
- NG (NG2/NG4) target sequences



# Targets

5 targets are detected:

- CT1 located on genomic DNA (also present in the genome of the swedish variant strains of *Chlamydia trachomatis*)
- NG2 independent and unique target of *Neisseria gonorrhoeae*
- NG4 independent and unique target of *Neisseria gonorrhoeae*
- SPC
- SAC





# Assay Requirements

## GeneXpert® Systems

- GeneXpert® Dx Software **v4.3** or higher/ Xpertise™ Software **v6.0** or higher

## Test Kits

- GXCT/NG-10 and GXCT/NG-120 (US-IVD)
- GXCT/NGX-CE-10 and GXCT/NGX-CE-120 (CE-IVD)

## Sample Collection

- SWAB/A-50
- SWAB/G-50
- URINE/A-50

## Other materials

- Personal Protective Equipment (PPE)
- 1:10 dilution of household bleach
- 70% ethanol or denatured ethanol

## Optional

- Uninterruptible Power Supply /Surge Protector
- Printer



# Good Laboratory Practice

## Personal Protective Equipment (PPE)

- Wear clean lab coats, safety glasses and gloves
- Change gloves between processing samples

## Lab Bench area

- Clean work surfaces routinely with:
  - ✓ 1:10 dilution of household bleach\*
  - ✓ 70% Ethanol Solution
- After cleaning, ensure work surfaces are dry

*\*Final Active Chlorine concentration should be 0.5% regardless of the household bleach concentration in your country*

## Specimens, Samples, and Kits Storage

- Store specimens and sample away from kit to prevent contamination

## Equipment

- Use filtered tips when recommended
- Follow the manufacturer's requirements for calibration and maintenance of equipment(s)



# Kit Handling





# Xpert<sup>®</sup> CT/NG Assay Kit Contents

Catalog Number	GXCT/NG-10 or GXCT/NGX-CE-10 GXCT/NG-120 or GXCT/NGX-CE-120
Tests Per Kit	10 or 120
Kit CD	Assay Definition File (ADF)
	Assay Import Instructions
	Package Insert (PDF)
Transfer Pipettes	10 or 120
Storage	2-28 °C

*Note: Sample Reagent contains guanidinium thiocyanate, which is harmful if swallowed (H303) and Irritating to eyes and skin (H315, H319).*

*Cartridges contain chemically hazardous substances-please see Package Insert and Safety Data Sheet for more detailed information.*





# Warnings and Precautions

- Store the Xpert® CT/NG Assay cartridges and reagents at 2–28°C
- Follow your institution's safety procedures for working with chemicals and handling biological samples
- Do not use Collection Reagent tubes that have not been validated by Cepheid
- Open the Assay cartridge lid only when adding the Sample, close the lid and proceed with processing
- With the GeneXpert System, start the test within 30 minutes after adding the sample to the cartridge
- With the Infinity System, place the cartridge on the conveyor within 30 minutes of adding the sample.





# Warnings and Precautions

- Do not shake the cartridge
- Do not use a cartridge that...:
  - appears wet, has leaked or if the lid seal appears to have been broken
  - appears damaged
  - has been dropped after removing it from packaging
  - has been dropped or shaken after adding the sample to it
  - has a damaged reaction tube
  - has been used: each cartridge is single-use to process one test
  - is expired
  - do not reuse disposable pipettes

*Dispose Xpert® CT/NG Assay cartridges and reagents according to your institution's and country's guidelines for disposal of hazardous materials*

# Warnings and Precautions

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions.
- Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures.
- If national or regional regulations do not provide clear direction on proper disposal, biological specimens and used cartridges should be disposed per WHO [World Health Organization] medical waste handling and disposal guidelines.



# Xpert<sup>®</sup> CT/NG Assay Limitations

- The Xpert CT/NG test has been validated with the following specimen types, collected with the Xpert Swab Specimen Collection Kit, Xpert CT/NG Urine Specimen Collection Kit or Xpert Urine Specimen Collection Kit.
  - Endocervical swabs
  - Patient-collected vaginal swabs
  - Male and female pharyngeal swabs
  - Male and female rectal swabs
  - Male and female urine
- Presence of vaginal discharge, tampons, douching, and Cepheid non-validated specimen collection have not been determined.
- Collection and testing of urine specimens with the Xpert<sup>®</sup> CT/NG Assay is not intended to replace a cervical exam and endocervical sampling for diagnosis of urogenital infections. Other genitourinary tract infections can be caused by other infectious agents.
- The performance of Xpert<sup>®</sup> CT/NG has not been evaluated for patients **less than 14 years** of age or **in patients with a history of hysterectomy**.



# Specimen Collection, Storage and Transport





# Urine Specimen Collection

- **Urine**
  - Refer to urine specimen collection package insert
  - 20 to 50 ml of first-catch urine should be collected in a sterile urine collection cup with no preservative (not provided by Cepheid), from which 7 mL is transferred to the Urine Sample tube containing the preservative
  - Use only Xpert Urine Sample Collection kit for processing male and female urine prior to testing in the GeneXpert CT/NG assay
  - Urine specimen must be collected and tested before the expiration date of the Xpert Urine Sample Collection kit


URINE/A-50 is designed to preserve and transport CT and NG DNA in first-catch male and female urine specimens



*\*Cepheid catalog # URINE/A-50*




# Urine Specimen Collection (First Catch)



Direct patient to provide first catch urine (20-50 mL) into a urine collection cup.

Note: The patient should not have urinated for at least 1 hour before. Patient should not cleanse the genital area prior to collecting specimen.




Transfer approximately 7 mL of urine into the transport tube, using the disposable transfer pipette. The correct volume is marked by the black dashed line on the pipette.




The Xpert® Urine Specimen Collection kit contains


- Ⓐ Large transfer pipette
- Ⓑ Urine Transport Reagent tube




Open the package of disposable transfer pipette.




Remove the yellow cap from the transportation tube.



Replace the yellow cap on the transport tube and tighten securely.



Invert the transport tube 3-4 times to ensure that the specimen and reagent are well mixed.



Return the tube as instructed by your doctor, nurse or health care provider.







*Under or over dispensing of urine into the Xpert Urine Transport Reagent tubes may affect assay performance*



# Specimen Collection, Transport and Storage

- Unprocessed Urine samples**

Pictures : greencrossvet.com





		Temperature (°C)	Storage Time
	Female urine 	Room temperature	24 hours
	Male urine 		3 days
	Male or Female Urine  		8 days



# Specimen Collection, Transport and Storage

- **Urine samples transferred to Xpert® CT/NG Urine transport tube**

*Cepheid catalog # URINE/A-50*

		Temperature (°C)	Storage Time
	Female urine sample ♀	+2  +15 °C	45 days
	Female urine sample ♀	+2  +30 °C	3 days
	Male urine sample ♂	+2  +30 °C	45 days



# Vaginal and Endocervical Collection

- **Swab**

- Pharyngeal, rectal, vaginal and endocervical specimens are collected from patients using flocked swabs included in the kit
- Vaginal samples are collected by the patient. Conversely, endocervical samples are collected by a clinician
- Swabs are broken off into the transport reagent tubes to elute organisms and stabilize DNA
- Swab specimens are then transported to the laboratory for testing on the GeneXpert® Instrument

*Volume of transport medium: 2,4 mL*




*\*Cepheid catalog # SWAB/A-50 for vaginal & endocervical samples*

*or Cepheid catalog # SWAB/G-50 kits for pharyngeal, rectal, vaginal & endocervical sample*


SWAB/A-50 and SWAB/G-50 kits are designed to collect, preserve and transport endocervical and vaginal specimens from symptomatic and asymptomatic individuals to the lab prior to analysis with the Xpert® CT/NG Assay




# Pharyngeal Collection




Open the individual rectal specimen collection package that contains the pink-capped swab transport tube and individually wrapped collection swab. Discard the larger swab.




Open the collection swab wrapper by peeling the top of the wrapper




Hold the swab in your hand, placing your thumb and forefinger in the middle of the swab shaft across the scoreline.




Instruct the patient to open mouth widely. Position the tongue toward the bottom of the mouth. Swab areas of the pharynx (tonsil, posterior wall, uvula, posterior wall).



While holding the swab in the same hand, unscrew the cap from the Xpert Swab Transport Reagent tube.



Carefully break the swab shaft against the side of the tube at the scoreline




Re-cap the transport tube and tighten the cap securely.


Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming. Label the transport tube with the sample ID, including date of the collection, as required.

**Avoid splashing contents of the transport tube on the skin. Wash with soap and water if exposed.**


# Rectal Collection




Open the individual rectal specimen collection package that contains the pink-capped swab transport tube and individually wrapped collection swab. Discard the larger swab.




Open the collection swab wrapper by peeling the top of the wrapper




Hold the swab in your hand, placing your thumb and forefinger in the middle of the swab shaft across the scoreline.




Carefully insert the swab approximately 1 cm beyond the anal sphincter, so that the fiber tips are no longer visible and rotate the swab.



While holding the swab in the same hand, unscrew the cap from the Xpert Swab Transport Reagent tube.



Carefully break the swab shaft against the side of the tube at the scoreline





Re-cap the transport tube and tighten the cap securely.

Invert or gently shake the tube 3-4 times to elute material from the swab. Avoid foaming. Label the transport tube with the sample ID, including date of the collection, as required.


**Avoid splashing contents of the transport tube on the skin.**  
**Wash with soap and water if exposed.**

# Assessing a correct sample

**Figure 1. Examples of Acceptable Rectal Swabs**




**Figure 2. Examples of Unacceptable Rectal Swabs**







# Vaginal Collection (Patient Collected)




Open the individual vaginal/endocervical specimen collection package **A** that contains the pink-capped swab transport tube and individually wrapped collection swab. Set the tube aside. Discard the larger **B** swab




Open the collection swab by peeling open the top of the wrapper. Remove the swab, taking care not to touch the tip or lay it down.




Hold the swab in your hand, placing your thumb and forefinger in the middle of the swab shaft across the scoreline.




Carefully insert the swab into your vagina about 2 inches/5cm inside the opening of the vagina




Gently rotate the swab for 10-30 seconds. Ensure the swab touches the walls of the vagina so that the moisture is absorbed by the swab. Withdraw the swab and continue to hold it in your hand.



Unscrew the cap from the transport tube. Immediately place the collection swab into the transport tube.



Identifying the scoreline, break the swab shaft against the side of the tube. Discard the top portion of the swab shaft. Avoid splashing contents on the skin. Wash with soap and water if exposed.



Re-cap the transport tube and tighten the cap securely. Return the tube as instructed by your doctor, nurse or health care provider.




# Endocervical Collection (Clinician-collected)




The Vaginal /Endocervical Specimen Collection kit contains


- A** Individual Collection swab
- B** Cleaning Swab




Remove excess mucus from the cervix and surrounding area using the large individually wrapped cleaning swab. Discard the swab. **B**




Open the package that contains the pink-capped Xpert Swab Transport tube and individually wrapped collection swab. Open the collection swab wrapper by peeling the top of the wrapper




Hold the swab in your hand, placing your thumb and forefinger in the middle of the swab shaft.




Insert the collection swab into the endocervical canal. Rotate the swab for 30 seconds in the endocervical canal. Withdraw the swab carefully.



Unscrew the cap from the transport tube. Immediately place the collection swab into the transport tube.



Align the small groove against the edge (rim) of the tube and break it off. If needed, gently rotate the shaft to complete the breakage. Discard the top part of the swab shaft.



Re-cap the transport tube and tighten the cap securely. Label the transport tube with the sample ID and date of collection, as required.








# Specimen Collection, Transport and Storage

- Swab samples transferred to Xpert<sup>®</sup> CT/NG swab transport tube**

*Cepheid catalog # SWAB/A-50*

*or Cepheid catalog # SWAB/G-50 kits*

	Swab Samples	Validated Collection Tool	Temperature (°C)	Storage Time
	Endocervical collection Swab	SWAB/A-50 SWAB/G-50 kits		60 days
	Vaginal collection Swab	SWAB/A-50 SWAB/G-50 kits		60 days
	Pharyngeal collection Swab	SWAB/G-50 kits		60 days
	Rectal collection Swab	SWAB/G-50 kits		60 days





# Cartridge Preparation






# Cartridge Preparation – Urine or Swab

1




Obtain appropriately collected urine or swab specimen in Xpert Sample Collection Kit

2




Take one Xpert CT/NG cartridge and the provided transfer pipette

3



Label the side of the cartridge with the same ID as the collection tube


4



Open the cartridge lid


5

x3-4



Gently mix by inverting the transport tube 3-4 times

6



Pipette at least 1 mL of the sample using the provided pipette\*

7




Slowly empty the pipette into the sample chamber of the cartridge

8



Close the lid firmly. Start the test within the time frame specified in the package insert.



\*CT/NG Assay  
Transfer pipette



# Run a Test


## 1 Create Test

GeneXpert



Start the test within **30 minutes** after adding the sample to the cartridge

## 2 Scan cartridge barcode message




*By default, do not click on Manual Entry or Cancel*

## 3 Scan the cartridge



GeneXpert  
Infinity



Place the cartridge on the conveyor within **30 minutes** of adding the sample.


Order Test - Assay

[Scan Cartridge Barcode](#)

Cartridge barcode is successfully scanned when you hear the beep.



Patient ID	patientid
Sample ID	sampleid
Last Name	patient
Priority	Normal
First Name	id



*"For complete details on how to run a test, refer to the Package Insert and the GeneXpert® Dx or Xpertise™ Operator Manuals.*



# Create a Test on GeneXpert<sup>®</sup> Dx Software

4 Complete the fields as required

5 Select the appropriate Assay Protocol from the drop-down list: you will get only the assay related result

6 The module is selected automatically

7 Click on Start Test

8 A green light will flash on the module  
Load the cartridge into module and close the door

The screenshot shows the 'Create Test' window with the following fields and options:

- Patient ID: [Text Field]
- Sample ID: [Text Field]
- Patient ID 2: [Text Field]
- Last Name: [Text Field]
- First Name: [Text Field]
- Select Assay: [Drop-down menu showing 'Xpert CT']
- Select Module: [Drop-down menu showing 'Xpert NG', 'Xpert CT\_NG', 'Xpert CT']
- Reagent Lot ID\*: [Text Field]
- Test Type: [Drop-down menu showing 'Specimen']
- Sample Type: [Drop-down menu showing 'Other']
- Other Sample Type: [Text Field]
- Notes: [Text Area]
- Start Test: [Button]
- Scan Cartridge Barcode: [Button]





# Create a Test on Xpertise™ Software – Assay Selection

- 4 Select the appropriate Assay Protocol from the drop-down list: you will get only the assay related result

- 5 Click on Continue

**Order Test - Assay Selection**

Assay	Version
Xpert NG	3
Xpert CT_NG	3
Xpert CT	3

**Patient ID**  
patientid

**Sample ID**  
sampleid

**Last Name**  
patient

**First Name**  
id

**Priority**  
Normal





# Create a Test on Xpertise™ Software – Test Information

6 Review and complete the test information →


**Order Test - Test Information**

<b>Patient ID</b>	
<input type="text" value="patientid"/>	
<b>Sample ID</b>	
<input type="text" value="sampleid"/>	
<b>Last Name</b>	<b>First Name</b>
<input type="text" value="patient"/>	<input type="text" value="id"/>
<b>Assay*</b>	
<input type="text" value="Xpert CT_NG"/>	
<b>Version*</b>	
<input type="text" value="3"/>	
<b>Reagent Lot ID*</b>	<b>Cartridge S/N*</b>
<input type="text" value="12102"/>	<input type="text" value="282769448"/>
<b>Expiration Date*</b>	<b>Priority</b>
<input type="text" value="2018/11/04"/>	<input type="text" value="Normal"/>
<b>Test Type</b>	
<input type="text" value="Specimen"/>	
<b>Sample Type</b>	<b>Other Sample Type</b>
<input type="text" value="Other"/>	<input type="text"/>
<b>Notes</b>	
<input type="text"/>	

7 Click on SUBMIT →



8 Place the cartridge on the conveyor belt







# Automated Xpert<sup>®</sup> Protocol

1

Sample is added to the cartridge

2

The cartridge is loaded into the instrument

3

Nucleic acids are purified

Purified nucleic acids mix with the PCR reagents

4

Simultaneous amplification and detection occurs

5

Results are ready to view

6



# Quality Controls





# Cepheid Control Strategy

- **Instrument System Control – Check Status**

- System control checks the optics, temperature of the module and mechanical integrity of each cartridge.
- If the system controls fail, an ERROR test result will be reported.

- **Assay Quality Controls**

- Each Xpert cartridge is a self-contained test device
- Cepheid designed specific molecular methods to include internal controls that enable the system to detect specific failure modes within each cartridge
  - Sample Adequacy Control (SAC)
  - Specimen Processing Control (SPC)
  - Probe Check Controls (PCC)



# Internal Quality Controls

- **Sample Adequacy Control (SAC)** HMBS (Hydroxymethylbilane synthase)
  - Verifies that human cells are present in the sample
- **Probe Check Controls (PCC)**
  - Before the PCR step, fluorescence signal is measured on all probes and compared with pre-established factory settings to monitor
    - bead rehydration
    - reaction tube filling
    - probe integrity
    - dye stability
- **Sample Processing Controls (SPC)**
  - Genomic DNA of *Bacillus globigii* in each cartridge
    - Verifies adequate sample processing
    - Verifies lysis, presence of the organism and detects PCR inhibition
    - Must be positive in a negative sample
    - Can be positive or negative in a positive sample



# Commercially Available External Controls

Part Number	Description	Configuration	Storage
NATCT(434)-6MC	CT positive control	1 mL x 6 vials	2-8°C
NATNG-6MC	NG positive control	1 mL x 6 vials	2-8°C
NATCT/NGNEG-6MC	CT and NG Negative controls	1 mL x 6 vials	2-8°C

<http://www.zeptometrix.com>

1. *Take 1 vial of the control material*
2. *Vortex for 10 seconds*
3. *Using the kit pipette, aspirate control mix to the mark (1mL) and add it to the sample chamber of the cartridge*
4. *Close the lid and launch the test on the GeneXpert®*

- *External controls should be used in accordance with local, state accrediting organizations, as applicable*
- *NATtrol™ products are Research Use Only and not for in-vitro diagnostic use.*




# Result Interpretation







# Result Interpretation Algorithm






# Xpert<sup>®</sup> CT/NG – All possible results

Result displayed	CT1	NG2	NG4	SPC	SAC
<b>CT DETECTED</b>	+	+	+	+/-	+/-
<b>NG DETECTED</b>					
<b>CT DETECTED</b>	+	+	-	+/-	+/-
<b>NG NOT DETECTED</b>					
<b>CT DETECTED</b>	+	-	+	+/-	+/-
<b>NG NOT DETECTED</b>					
<b>CT NOT DETECTED</b>	-	+	+	+/-	+/-
<b>NG DETECTED</b>					
<b>CT NOT DETECTED</b>	-	-	+	+/-	+/-
<b>NG NOT DETECTED</b>					
<b>CT NOT DETECTED</b>	-	-	-	+	+
<b>NG NOT DETECTED</b>					
<b>INVALID</b>	-	-	-	-	+/-
<b>INVALID</b>	-	-	-	+/-	-



# CT DETECTED NG DETECTED

CT DETECTED;  
NG DETECTED




- The targets CT1, NG2 and NG4 are detected and the Ct values are within the valid range
- SAC: NA (not applicable)
  - SAC is ignored because a target amplification occurred
- SPC: NA (not applicable)
  - SPC is ignored because a target amplification occurred
- Probe Check: PASS



# CT DETECTED; NG NOT DETECTED

CT DETECTED;  
NG NOT DETECTED




- The target CT1 is detected and the Ct value is within the valid range
- None of the NG targets are detected
- SAC: NA (not applicable)
  - SAC is ignored because an amplification occurred
- SPC: NA (not applicable)
  - SPC is ignored because the CT1 target amplification occurred
- Probe Check: PASS



# CT NOT DETECTED; NG DETECTED

CT NOT DETECTED;  
NG DETECTED




- The targets NG2 and NG4 are detected and the Ct values are within the valid range
- The target CT1 is not detected
- SAC: NA (not applicable)
  - SAC is ignored because the NG target amplification occurred
- SPC: NA (not applicable)
  - SPC is ignored because the NG target amplification occurred
- Probe Check: PASS



# CT NOT DETECTED; NG NOT DETECTED

CT NOT DETECTED;  
NG NOT DETECTED



- The targets CT1, NG2 and NG4 are NOT detected
- SAC: PASS
  - SAC has a Ct value within the valid range
- SPC: PASS
  - SPC has a Ct value within the valid range
- Probe Check: PASS



# Troubleshooting






# Factors That Negatively Affect Results

- Improper specimen collection
  - The number of organisms in the specimen is below the detection limit of the test
  - Performance with other collection devices and specimen types has not been assessed
- Improper transport or storage of collected specimen
  - Storage and transport conditions are specimen specific
  - Refer to the Package Insert for the appropriate handling instructions
- Improper testing procedure
  - Modification to the testing procedures may alter the performance of the test.
  - Careful compliance with the package insert is necessary to avoid erroneous results



# INVALID

INVALID




- Presence or absence of the CT1 and NG2/4 targets can not be determined
- SAC: PASS
  - SPC has a Ct value within the valid range
- SPC: FAIL
  - SAC Ct value is not within the valid range
- Probe Check: PASS



# INVALID

INVALID




- Presence or absence of the CT1 and NG2/4 targets can not be determined
- SAC: FAIL
  - SAC Ct value is not within the valid range
- SPC: PASS
  - SPC has a Ct value within the valid range
- Probe Check: PASS



# INVALID

INVALID



- Presence or absence of the CT1 and NG2/4 targets can not be determined
- SAC: FAIL
  - SPC Ct value is not within the valid range
- SPC: FAIL
  - SAC Ct value is not within the valid range
- Probe Check: PASS



# INVALID

- INVALID result with failing **SPC** and/or **SAC**

## Origin(s)

- PCR was inhibited due to interfering substances
- Inadequate sample was used
- Improper specimen storage/collection/preparation
- Improper kit storage conditions

## Solution(s)

- Use the correct **specimen type**
- Check the **sample quality** (Blood, Mucin, topical medication...)
- Follow recommended **instructions** on sample collection, preparation and storage
- Check **kit storage** conditions and shelf life
- Collect a new sample when necessary and retest





# INVALID

- INVALID result with failing SAC only

- Origin(s)

- Inadequate sample was used
    - Improper specimen collection
    - Improper sample storage or preparation
    - Improper kit storage conditions

- Solution(s)

- Use the correct specimen type
    - Check the collection: Urine first catch must be collected to ensure a proper epithelial cell concentration – Proper swabbing must be performed (according to illustrated collection instructions)
    - Follow recommended instructions on sample collection, preparation and storage
    - Check kit storage conditions and shelf life
    - Collect a new sample in the appropriate conditions, when necessary, and retest



# ERROR

ERROR

Test Result	Analyte Result	Detail	Errors	History	Support
Assay Name	Xpert CT_NG		Version	3	
Test Result	ERROR				
For In Vitro Diagnostic Use Only.					

- The Test result tab displays “ERROR”
- The error code and description can be found in the “Errors” Tab
- The test must be re-run, after corrective actions



# NO RESULT

NO RESULT

The screenshot shows a software window with a title bar. Inside, there's a section labeled 'Test Result' with a dropdown menu currently set to 'NO RESULT'. Below this, a large text area contains the phrase 'For In Vitro Diagnostic Use Only.' The window has a standard scrollbar at the bottom.

## NO RESULT

- Test could not be completed and insufficient data was collected

## ORIGIN(S)

- Power failure during test
- “Stop Test” function was used.
- Computer freeze or crash during test



# Re-test Procedure

1

Discard used cartridge

Follow your institution's safety guidelines for disposal of cartridges

2



Obtain the residual treated sample from either

- Swab Transport Reagent or
- Urine Transport Reagent tube

If the leftover sample volume is insufficient, or the retest continues to return an INVALID, ERROR, or NO RESULT, collect a new sample

3



Obtain a new cartridge


Label appropriately as retest on the new cartridge

Process the sample per the package insert

4



Run the test on the GeneXpert® System





# Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
  - Product name
  - Lot number
  - Serial number of the System
  - Error messages (if any)
  - Software version and, if applicable, Computer Service Tag number
- Log your complaint using the following link <http://www.cepheid.com/us/support>

Region	Telephone	Technical Support Email
US	+ 1 888 838 3222	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>
Australia and New Zealand	+ 1800 130 821 + 0800 001 028	<a href="mailto:techsupportANZ@cepheid.com">techsupportANZ@cepheid.com</a>
Brazil and Latin America	+ 55 11 3524 8373	<a href="mailto:latamsupport@cepheid.com">latamsupport@cepheid.com</a>
China	+ 86 400 821 0728	<a href="mailto:techsupportchina@cepheid.com">techsupportchina@cepheid.com</a>
France	+ 33 563 825 319	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Germany	+ 49 69 710 480 480	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
India, Bangladesh, Bhutan, Nepal and Sri Lanka	+ 91 11 48353010	<a href="mailto:techsupportindia@cepheid.com">techsupportindia@cepheid.com</a>
Italy	+ 39 800 902 567	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
South Africa	+ 27 861 22 76 35	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
United Kingdom	+ 44 3303 332 533	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Belgium, Luxembourg and Netherlands	+33 563 825 3319	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other European, Middle East, and African countries	+ 33 563 825 319 + 971 4 253 3218	<a href="mailto:support@cepheideurope.com">support@cepheideurope.com</a>
Other countries not listed	+1 408 400 8495	<a href="mailto:techsupport@cepheid.com">techsupport@cepheid.com</a>



Thank You.



[www.Cepheid.com](http://www.Cepheid.com)

