



# KAISER PERMANENTE®

<b>DOCUMENT NUMBER:</b> SCPMG-PPP-0113
<b>DOCUMENT TITLE:</b> Instructions For Use: Incubator Carrier Tote
<b>DOCUMENT NOTES:</b>

<b>LOCATION:</b> SCPMG-rel	<b>VERSION:</b> 03
<b>DOC TYPE:</b> Microbiology	<b>STATUS:</b> Release

<b>EFFECTIVE DATE:</b> 09 Mar 2017	<b>NEXT REVIEW DATE:</b> 09 Mar 2019
<b>RELEASE DATE:</b> 09 Mar 2017	<b>EXPIRATION DATE:</b>

<b>AUTHOR:</b> H099970	<b>PREVIOUS NUMBER:</b>
<b>OWNER:</b> RRL Dir Lab Svcs, AP	<b>CHANGE NUMBER:</b> SCPMG-CR-0308

## Instructions for Use: Incubator Carrier Tote

**Purpose**

The purpose of this procedure is to provide instructions for use of the portable incubator tote and the temperature data logger which are used to provide a controlled temperature environment [30- 37°C] for the transportation of plated samples for culture.

**Specimen Transport**

- Agar plates streaked with sample for culture are placed in a CO<sub>2</sub> biobag
- The CO<sub>2</sub> biobag is placed in an incubator tote for transport to the Regional Lab


**Note:** The incubator tote is not for transport of blood culture bottles or raw patient samples.

**Maintaining the tote when not in use**

**Charging the battery**

**NOTE:** Unit must be off when charging the battery


Step	Action
1	Insert the plug on the green battery charger into the plug on the front of the incubator [labeled CHARGING PORT] Set the charge current switch to 1.8



*Continued on next page*

## Instructions for Use: Incubator Carrier Tote, Continued

**Maintaining the tote when not in use, continued**      **Charging the battery, continued**


Step	Action
	<p>Plug the GREEN BATTERY CHARGER into a wall socket</p> <ul style="list-style-type: none"> <li>• The light on the charger will be red when charging</li> <li>• The light will turn green when cycle is complete</li> <li>• <i>Once cycle is complete the charger will not restart unless it is unplug and plugged in again</i></li> <li>• <i>This should be done daily to keep the batteries fully charged</i></li> </ul> <p><b>Note:</b> A completely drained battery will require 20-25 hours to charge to capacity.</p>
3	<p>Press the black button on the Battery Monitor on the front of the incubator to check the charge.</p> <ul style="list-style-type: none"> <li>• A MINIMUM of 4 lights must be illuminated otherwise the battery is not fully charged.</li> <li>• The Battery Monitor is separate from the power port. It will only display battery status when the unit is off.</li> </ul> <div style="text-align: center;">  </div>

*Continued on next page*

## Instructions for Use: Incubator Carrier Tote, Continued

### Preparing the Tote for use-

### Turning on the Heating Unit

Step	Action
1	Insert the end of the plug from the black power supply into the power port on the left end of the incubator <ul style="list-style-type: none"> <li>• The temperature display will light up.</li> <li>• The motor will start up</li> <li>• The power switch will automatically turn “on”</li> </ul> 
2	Allow the incubator carrier unit to come up to temperature and stabilize before placing samples inside [approximately 1 hour].  <b>Note: The unit must be turned on, fully charged, and allowed to come up to temperature prior to use.</b>
3	Load samples for transport.

*Continued on next page*

## Instructions for Use: Incubator Carrier Tote, Continued

### Preparing the tote for transport-

Step	Action
1	Unplug the power supply from the power port
2	Push power switch in to turn on the unit. The temperature display should light up
3	Verify that there is a temperature logger inside the incubator
4	Record the temperature from the display on the lower left of the incubator.
5	Send via courier to Sherman Way lab.

### Receipt in testing department-

Step	Action
1	Record the temperature from the lower left digital display
2	Push power switch in to turn off the unit.
3	Remove the samples
4	Return to Medical Center

### Quality Control

- The temperature of the unit will be recorded prior to shipment to SWL
- The temperature and time of receipt of the unit will be recorded upon receipt in Bacteriology
- The incubator contains a calibrated data logger device which will record the internal temperature of the incubator carrier at 30 minute intervals
  - If outside of the functional temperature range [34-37°C], inform a manager
  - The incubator carrier will be taken out of service until the data from the logger can be downloaded and reviewed.
- The unit will be pulled from service once a year for calibration

---

## Instructions for Use: Incubator Carrier Tote, Continued

---

### Downloading The Data from the data logger

Perform this procedure according to the following steps.

Note: Requires DeltaTrak FlashPDF software

Step	Action
1	Retrieve data logger from incubator carrier
2	Plug into USB port of computer. A pop up will appear displaying the progress of the download.
3	Once complete, the PDF will automatically display <ul style="list-style-type: none"> <li>• Save a copy of the PDF on the shared drive</li> <li>• Print a copy of the PDF for review and signature. Record the IR number of the incubator carrier unit and the date the data was retrieved.</li> </ul>
4	Once the data has been retrieved, right click on the DeltaTrak icon on the taskbar and select “settings” <ul style="list-style-type: none"> <li>• Advance to screen that contains the alarm parameters and verify the alarm parameters</li> <li>• Return to initial screen and click “reactivate”. A pop up will appear that asks “reconfigure token?”- select “yes”.</li> <li>• Upon completion another message will pop up and ask if another token is to be configured, reply “no”.</li> </ul>
5	Remove token from the USB port and return to the incubator carrier
6	Place a sticker on the incubator carrier with the date the data was downloaded and the due date for the next data retrieval.
7	Return the carrier to service.

### Non-controlled Documents:

The following non-controlled documents support this instruction:

- Nicholson, Ron. Incubator Carrier Instructions. October 09, 2015.
- DeltaTrak FlashPDF instructions. Copyright 2012.

*End*

**Signature Manifest****Document Number:** SCPMG-PPP-0113**Revision:** 03**Title:** Instructions For Use: Incubator Carrier Tote

All dates and times are in Pacific Standard Time.

**New Lab Director - Micro****Change Request**

Name/Signature	Title	Date	Meaning/Reason
Paulette Medina (K088673)	ASST DIR REGL LAB		
Vahe Khanlian (O532803)	RRL DIR OF LAB SVCS, MIC		
Matthew Jones (F754627)	Systems Consultant	06 Mar 2017, 01:30:40 PM	Approved

**Collaboration**

Name/Signature	Title	Date	Meaning/Reason
Matthew Jones (F754627)	Systems Consultant	06 Mar 2017, 01:34:33 PM	Complete

**Initial Approval**

Name/Signature	Title	Date	Meaning/Reason
Charles Park (K239415)	Director of Operations	06 Mar 2017, 02:17:50 PM	Approved
Ken Van Horn (K660731)	Asst. Technical Director Bact	07 Mar 2017, 04:22:46 PM	Approved

**Final Approval**

Name/Signature	Title	Date	Meaning/Reason
David Quam (P092597)	Rgnl Mg Admn-Pmg Executive	09 Mar 2017, 11:41:57 AM	Approved

**Set Effective Date**

Name/Signature	Title	Date	Meaning/Reason
Paulette Medina (K088673)	ASST DIR REGL LAB		
Matthew Jones (F754627)	Systems Consultant	09 Mar 2017, 01:55:25 PM	Approved



**Notify Users**

Name/Signature	Title	Date	Meaning/Reason
Paulette Medina (K088673)	ASST DIR REGL LAB	09 Mar 2017, 01:55:25 PM	Email Sent

## Instructions for Use: Incubator Carrier Tote

---

Reviewed and approved by (for Medical Center Area Approval Only):

SIGNATURE	DATE
	5.10.18
Name: <u>Janice M. Wolf</u> Operations Director, Area Laboratory	
	5/10/18
Name: <u>Sony Wirio, MD</u> CLIA Laboratory Director	

---