

## Blood Parasite Identification

### BP-06 Introduction

Thick and thin Giemsa-stained smears were obtained from 51-year-old male from Wisconsin with fever and chills. The specimen contained *Babesia* sp. A response of "*Babesia* sp." and "*Plasmodium* sp./*Babesia* sp. seen, referred for identification" were considered satisfactory.


Parasite Identification	Referees (32)		Participants (447)	
	No.	%	No.	%
<i>Babesia</i> sp.	32	100.0	401	89.7
<i>Plasmodium falciparum</i>	-	-	40	8.9

Parasite Screen	Referees (20)		Participants (770)	
	No.	%	No.	%
<i>Plasmodium</i> sp./ <i>Babesia</i> sp. seen, referred for identification	20	100.0	758	98.4

BP-06

If you have identified a *Plasmodium* sp. or *Babesia* sp.: Report percent of infected RBCs seen (number of infected RBCs per 100 red blood cells). (Ungraded)



### Summary of Key Features for Identification:

1. *Babesia* spp. infect red blood cells of all sizes.
2. Only the ring stages are identified on blood smears.
3. Differentiating *Babesia* spp. from *P. falciparum* is sometimes difficult. The distinguishing features of *Babesia* spp. include pleomorphic ring forms, the "Maltese cross pattern" and extracellular forms.

Upon initial presentation of the patient, the parasite may present in such low numbers in blood that they cannot be seen in thin films. The diagnosis then depends on finding parasites in the thick films. If only ring forms are present on the thick films, it may be impossible to distinguish between *P. falciparum* and *Babesia* species. The clinical history along with collection of multiple blood samples then are of extreme importance. Molecular methods (PCR) testing or malarial antigen testing may also be useful in these cases. Confirmed and suspected cases should be reported to the department of health.

---

## Discussion

### Causal Agents

Babesiosis is caused by apicomplexan parasites of the genus *Babesia*. Species most-commonly attributed to human disease include *B. microti* (endemic to northeastern and northern midwestern United States; introduced elsewhere), *B. duncani* (northwestern North America), *B. divergens* (Old World), and *Babesia* sp. MO-1 (midwestern United States).

### Biology and Life Cycle

*Babesia* spp. have a two-host life cycle. The definitive hosts are ticks in the genus *Ixodes*, while the intermediate hosts are mammals (usually rodents in nature). Infective sporozoites are released into the mammalian hosts when an infected *Ixodes* takes a blood meal. Sporozoites enter erythrocytes and undergo asexual reproduction by budding. Early ring-form trophozoites give rise to merozoites (including the classic 'tetrad form') that either produce more trophozoites or gametes. The latter is a dead-end stage in the mammalian host, and gametes must be ingested by a tick in order for sexual reproduction to occur. Humans can become infected either by being fed upon by an infected tick or through blood transfusions. Rare congenital cases have also been reported.

### Diagnosis

Human babesiosis is typically diagnosed initially by the finding of trophozoites and merozoites on stained blood films. Merozoites displayed in a 'Maltese-cross' formation (tetrads) are typically diagnostic for *Babesia* spp. Trophozoites may be ring-shaped, pyriform (pear-shaped), or pleomorphic and vacuolated, and arranged singly or in short chains. Multiply-infected RBCs are not uncommon. Ring-form trophozoites need to be distinguished from *Plasmodium* spp., especially *P. falciparum*. *Plasmodium* spp. produce pigment, which is never produced by *Babesia* spp. In cases of babesiosis, extracellular ring-forms may be observed, singly or in clusters.

Species-level identification of *Babesia* cannot be reliably accomplished by microscopy alone. Patient travel history can be useful information for determining the species of *Babesia* present, but confirmatory diagnosis should be performed by molecular (PCR, or PCR in combination with sequencing analysis) or serologic methods. Serology can also play a very important role in screening potential donors in transfusion-acquired cases.

### Clinical significance

The clinical outcome can vary from asymptomatic infection to death depending on infective organism, presence of other tick-borne infections such as borreliosis, age, host immune status, and other underlying factors (e.g., splenectomy or a history of blood transfusions). Presenting symptoms include malaise, chills, myalgia, fatigue, anemia, and high-grade fever not unlike acute malaria. More chronic symptoms described in cases are nausea, vomiting, night sweats, loss of weight, and bloody or dark urine.

### Treatment

The majority of babesiosis cases self-resolve without need for drug therapy. The standard of care for more severe cases is clindamycin combined with quinine. An alternative regimen is atovaquone combined with azithromycin. In severe disease, or in splenectomized patients, exchange transfusion may be needed in addition to antimicrobial therapy.

---

### References

1. Ash LP, Orihel TC. *Atlas of Human Parasitology*, 5<sup>th</sup> ed. Chicago, IL: ASCP Press: 2007.
2. Centers for Disease Control and Prevention, Division of Parasitic Diseases (DPD). *Laboratory Identification of Parasites of Public Concern*. <http://www.cdc.gov/dpdx>
3. Mandell GL, Bennett JE, Dolin R. *Mandell, Douglas, and Bennett's Principles and Practice of Infectious Diseases*, 7th ed. Philadelphia, PA: Churchill Livingstone Elsevier, 2009

## BP-07 Introduction

Thick and thin Giemsa-stained smears were obtained from a Moroccan Peace Corps worker with fever. A response of "No parasite(s) seen" or "Specimen screened for blood parasites, no organisms seen" was considered satisfactory.

BP-07	<b>Parasite Identification</b>	<b>Referees (32)</b>	<b>Participants (428)</b>
		<b>No. %</b>	<b>No. %</b>
	No parasite(s) seen	34 100.0	425 99.3
	<b>Parasite Screen</b>	<b>Referees (18)</b>	<b>Participants (790)</b>
		<b>No. %</b>	<b>No. %</b>
	Specimen screened for blood parasites, no organisms seen	18 100.0	778 98.5

## Discussion

### Identification

Careful examination of multiple thin and thick blood films is imperative to exclude the diagnosis of blood parasites, particularly for patients living in endemic areas. For thin films: (1) all blood components (erythrocytes, white blood cells, and platelets) should be intact, (2) the background should be clean and free from debris, (3) erythrocytes should stain a pale grayish-pink, and neutrophilic leukocytes should have deep purple nuclei and well defined granules, and (4) erythrocytes at the terminal, feathered end of the film should be adjacent, but not overlap (one layer thick). For thick films: (1) the background should be clean, free from debris, with a pale mottled-gray color derived from lysed erythrocytes, (2) leukocytes should stain deep purple with pale purple cytoplasm, and (3) eosinophilic granules should stain a bright purple-red and neutrophilic granules should stain deep pink-purple.

Thick films are most useful for screening since they provide a larger quantity of blood for examination. Thin films, on the other hand, are most useful for speciation since they provide the best red blood cell (RBC) and parasite morphology. All requests for peripheral blood smear examination to detect *Plasmodium* spp. should be performed without delay. Both thick and thin films should first be fully screened at low power (ie, using the 10x objective) to detect microfilaria which may be present in low numbers anywhere on the slides and which may not be detected in the standard 300 field slide review at higher magnification.

Due to the severe implications of a misdiagnosis, laboratory personnel should then examine at least 300 oil immersion fields (using the 100X oil immersion objective) for each thick and thin blood film. In addition, one set of blood films is not sufficient to exclude the diagnosis of malaria and the laboratory should recommend collection of multiple blood specimens approximately at 6-8 hour intervals to definitively exclude the presence of blood parasitemia. This comment should accompany the final report "No blood parasites seen."

## References

1. Garcia LS. 2007. *Diagnostic Medical Parasitology*, 5<sup>th</sup> ed., Washington, DC. ASM Press.
2. CDC: Malaria surveillance - United States 2008. *MMWR* 2010;59(SS07);1-15.
3. Greenwood BM, Bojang K, Whitty CJ, Targett GA. Malaria. *Lancet*. 2005;365:1487-1498.

**BP-08 Introduction**

Thick and thin Giemsa-stained smears were obtained from a 35-year old missionary worker in sub-Saharan Africa with fever and diarrhea. The specimen contained *Plasmodium falciparum*. A response of "*Plasmodium falciparum*," or "*Plasmodium sp.*, NOS would refer or request another specimen, or perform additional molecular testing" and "*Plasmodium sp./Babesia sp.* seen, referred for identification" were considered satisfactory.

Parasite Identification	Referees (33)		Participants (450)	
	No.	%	No.	%
<i>Plasmodium falciparum</i>	27	81.8	343	76.7
<i>Plasmodium sp.</i> , NOS would refer or request another specimen, or perform additional molecular testing	3	9.1	56	12.5
<i>Plasmodium vivax</i>	1	3.0	12	2.7
<i>Plasmodium sp.</i> , not <i>P. falciparum</i> , referred for identification	2	6.1	14	3.1


  

Parasite Screen	Referees (19)		Participants (775)	
	No.	%	No.	%
<i>Plasmodium sp./Babesia sp.</i> seen, referred for identification	17	89.5	733	94.6

BP-08

BP-08

If you have identified a *Plasmodium sp.* or *Babesia sp.*: Report percent of infected RBCs seen (number of infected RBCs per 100 red blood cells). (Ungraded)



Key morphologic features on thin blood film that suggests a diagnosis of *P. falciparum* may include:

- Normal size and shape of infected erythrocytes
- Smaller, more delicate ring stages (approximately 1/5 the diameter of the erythrocyte) frequently with two chromatin dots (so-called "head phone" forms)
- Erythrocytes infected with multiple ring forms
- Presence of ring forms at the edge of the erythrocyte (appliqué forms)
- Absence of mature trophozoites and schizonts in the peripheral blood film (may be seen if there is a delay in processing the blood specimen)
- Presence of crescent-shaped gametocytes (not always seen)
- Absence of Schüffner's stippling. Larger, comma-shaped dots (Maurer's clefts) may be seen, especially when the stain buffer is at a pH of 7.2

Distinguishing *P. falciparum* from *Babesia* spp. can be challenging, given that both have a predominance of small ring forms, infect RBCs of all ages, and can produce multiply infected RBCs. The presence of Maurer's clefts, malarial pigment, and banana-shaped gametocytes eliminate *Babesia* infection from consideration. Furthermore, *Babesia* parasites are usually more pleomorphic with spindled, elliptical and ovoid forms. Extracellular forms are also more common in babesiosis. Finally, identification of the classic "Maltese cross" or tetrad form of *Babesia* sp. is diagnostic for babesiosis. Molecular or antigen-detection methods, in addition to clinical/travel history may be useful adjuncts for distinguishing between these two similar appearing parasites.

---

## Discussion

### Causal Agents:

There are four species of *Plasmodium* that cause human malaria: *P. falciparum*, *P. malariae*, *P. ovale*, and *P. vivax*. In addition, there are at least six species of simian *Plasmodium* that have been documented to cause zoonotic infections in humans, the most notable being *P. knowlesi*, infections of which appear to be increasing on the Malaysian peninsula.

*Plasmodium falciparum* occurs nearly worldwide in the tropics and subtropics, particularly in Africa and Southeast Asia. *Plasmodium malariae* also occurs nearly worldwide in the tropics and subtropics, but has a more patchy distribution than *P. falciparum*; most common in tropical Africa, Indian subcontinent, and Southeast Asia. *Plasmodium ovale* occurs primarily in tropical western Africa, but also New Guinea and Southeast Asia; *P. ovale* has not yet been documented from the New World. *Plasmodium vivax* occurs nearly worldwide in the tropics, subtropics, and some regions of northern and eastern Africa, the Middle East, the Indian subcontinent, Southeast Asia, and the Americas.

### Biology and Life Cycle:

*Plasmodium* spp. are transmitted by mosquitoes in the genus *Anopheles*. Infected female mosquitoes inject sporozoites when taking a blood meal. Sporozoites are carried via blood to the liver where they invade hepatocytes and form schizonts. The liver schizonts rupture, releasing large numbers of merozoites that then invade erythrocytes starting the erythrocytic cycle. Early ring forms develop into mature trophozoites and take one of two pathways: 1) they develop into schizonts (which rupture and continue the erythrocytic cycle) or 2) develop into gametocytes. Gametocytes are a dead-end stage in the human host but are required for sexual reproduction in the mosquito. In the mosquito host, microgametocytes (=males) exflagellate and fertilize macrogametocytes (=females), resulting in an ookinete. Ookinetes further develop into oocysts, which when mature rupture and release the infective sporozoites. In *P. falciparum*, late trophozoites and schizonts express a protein on the surface of the erythrocytic membrane causing the infected erythrocyte to adhere to the endothelial lining of capillaries in internal organs. Thus, only ring forms and gametocytes are usually seen in well-prepared peripheral blood smears.

### Diagnosis:

The ideal specimen for laboratory identification of malaria is fresh capillary blood from a finger or heel stick, with immediate preparation of thick and thin blood films. Since this is not feasible in most settings, however, venous blood collected in EDTA anticoagulant is also acceptable. It is important to transport the blood as quickly as possible to the laboratory for examination since prolonged exposure to EDTA may result in distortion of the malaria parasites and compromised morphology.

Malaria is primarily diagnosed by the identification of *Plasmodium* parasites on thin and thick blood films stained with Giemsa, Wright, or Wright-Giemsa stain. Molecular methods such as PCR may be employed when an identification cannot be made morphologically or there is morphologic evidence of a mixed infection. Rapid diagnostic tests are also commonly used to distinguish *P. falciparum* from other malaria infections. Serology is not used for routine diagnosis but may be helpful during transfusion investigations.

Typically, thick films are used for the recognition of *Plasmodium*, with a species-level identification performed on the thin film. Thin films should be read at 1000x magnification with oil for at least 100 microscopic fields. Immunologically naïve patients (eg, returning travelers born in non-endemic areas) may present with stronger clinical manifestations at a lower parasitemia. Severe malaria is currently defined as a parasitemia  $\geq 2\%$  in immunologically naïve patients and  $\geq 5\%$  in non-naïve patients.

### Morphologic Identification:

Two important questions regarding morphologic identification of malaria are:

- 1) Is it malaria?
- 2) Is it *Plasmodium falciparum*?

Recognition of *Plasmodium* is based on observing stages of the parasite inside infected red blood cells. In a well-prepared specimen, the cytoplasm will stain blue and the chromatin red. Pigment (which is absent in *Babesia* infections) will present as golden-brown to black flecks. In some species, structures such as Schüffner's stippling or Maurer's clefts may be present when stained at an appropriate pH. Identification of *Plasmodium* to the species level is extremely important for patient management, as different species may be treated differently (for example, it is important to target the liver stages of *P. ovale* and *P. vivax* to prevent relapse of the disease).

The following table compares the morphologic features of the four stages of human *Plasmodium* spp.

Morphologic Criteria	<i>Plasmodium falciparum</i>	<i>Plasmodium malariae</i>	<i>Plasmodium ovale</i>	<i>Plasmodium vivax</i>
Size of infected RBC	Normal	Normal to smaller	Enlarged	Enlarged
Rings (early trophozoites)	Common; usually with thin, delicate cytoplasm and double chromatin dots; often multiple rings per infected RBC; appliqué forms common	Cytoplasm sturdy, usually with single, large chromatin dot; occasional 'birds-eye' forms	Cytoplasm sturdy, with 1-2 large chromatin dots	Large, sturdy cytoplasm, usually with large, single chromatin dot
Developing trophozoites	Rare, but may be seen if a delay in processing; form compact, pigment usually evident	Variable; may be compact to elongate (band-form) or pleomorphic and vacuolated (basket-form); pigment coarse	Compact to slightly amoeboid with dark pigment; elongation and fimbriation may be observed	Pleomorphic to grossly amoeboid; pigment diffuse and golden-brown to nearly black
Schizonts	Rare in peripheral blood; 8-24 small merozoites when mature; pigment dark, compact	6-12 merozoites when mature, often in a rosette pattern around central mass of pigment	6-14 merozoites when mature; pigment dark brown to black and discrete when mature; elongation and fimbriation may be present	12-24 merozoites when mature; may fill entire RBC; often noticeable enlarged
Gametocytes	Crescent-shaped; chromatin discrete (macrogametocyte) or diffuse (microgametocyte); Laveran's bib may be present	Small, round, compact; pigment coarse and diffuse	Round to oval, compact; if elongated and fimbriated may not fill entire infected RBC; pigment coarse, dark	Large and round to pleomorphic (may 'hug' surrounding RBCs); may fill most of infected RBC; pigment golden-brown to nearly black
Other Features	Maurer's clefts may be present; ring-form trophozoites usually predominate	Generally smaller; pigment coarse; Ziemann's stippling may be present; all stages seen	Schüffner's stippling may be present at appropriate pH; elongation and fimbriation may be observed; all stages seen	Schüffner's stippling may be present at appropriate pH; enlargement of infected RBCs usually pronounced; all stages seen



### Calculating Percent Parasitemia:

The percent parasitemia is very important to calculate for prognostic purposes and also to evaluate response to antimalarial therapy.

Parasitemia can be calculated on a thin blood film as follows:

1. Count the number of infected RBCs per 100 RBCs in different oil immersion fields.
2. Apply the formula:

$$\frac{\text{\# of infected RBCs}}{\text{total \# of RBCs counted}} \times 100 = \% \text{ parasitemia}$$

Notes:

- 1) At least 500 RBC's should be counted, with counting 2000 or more RBCs providing the most accurate estimation of parasitemia
- 2) An infected RBC containing multiple parasites is calculated only once
- 3) Fields devoid of parasites should be included, if encountered
- 4) Gametocytes should not be included in the count. Justification is because: 1) many antimalarial drugs are not gametocidal and the presence of gametocytes post-treatment is not indicative of the effectiveness of the treatment and 2) gametocytes are a dead-end stage in the human host.

### Clinical Significance:

In 2014, ninety seven countries and territories had ongoing malaria transmission. Over half a million people die from malaria each year. Most malaria cases and deaths occur in children in sub-Saharan Africa. In 2011, a 40-year high of 1,925 cases of malaria were reported to the CDC in the United States, almost all in recent travelers and immigrants. Although the *Anopheles* mosquito is endemic in parts of North America, malaria transmission was largely eradicated in the 1940s through public health efforts.

Malaria infection can be classified as either uncomplicated or severe (complicated). In uncomplicated infections, patients present with nonspecific symptoms including fever, chills, sweats, headaches, nausea/vomiting, body aches and malaise. Symptoms classically (but infrequently observed) recur either in a two-day cycle (*P. falciparum*, *P. vivax* and *P. ovale*) or in three-day cycle (*P. malariae*). In severe infections, organ failure and/or metabolic abnormalities occur including severe anemia, acute respiratory distress syndrome, acute kidney failure, metabolic acidosis, cerebral infection, and coagulation abnormalities. Severe infection is considered a medical emergency requiring urgent treatment. *Plasmodium falciparum* can cause severe illness and death whereas *P. vivax*, *P. malariae* and *P. ovale* tend to cause less severe illness. The hypnozoite form of *P. vivax* and *P. ovale* can remain dormant in a patient's liver and cause relapsing infection.

### Treatment:

Treatment of malaria should ideally wait until a laboratory diagnosis has been made. Treating "presumptively" should occur only when no other option exists. Therapy is guided by the infecting species of *Plasmodium*, the clinical status of the patient and the drug susceptibility of the infecting parasites (dependent on geographic area and previous anti-malarial treatment). Because of the rapid progression of *P. falciparum* infections and a high risk of fatality, urgent treatment is essential.

If the infection is uncomplicated, oral anti-malarial medication can provide effective treatment. However, severe infections necessitate parenteral therapy. *Plasmodium falciparum* and *P. vivax* have different drug resistance patterns in different geographic regions. Although not readily available in North America, the WHO recommends artemisinin-based combination therapy as first-line treatment in uncomplicated *P. falciparum* malaria (oral administration), severe malarial infections (intravenous administration) and *P. vivax* infections in areas of known chloroquine resistance. Other, non-artemisinin based combination treatments include sulfadoxine-pyrimethamine

plus either chloroquine, amodiaquine, or atovaquone-proguanil. In recent years, resistance to artemisinins has been detected in Cambodia, Laos, Myanmar, Thailand and Vietnam. For confirmed *P. vivax* and *P. ovale* infections, radical cure can be achieved with treatment using primaquine and in order to prevent relapse due to the hypnozoite form. In high-transmission settings re-infection with *P. vivax* is likely. Mixed-species malarial infections are not common but may be underestimated by routine microscopy.

---

## References

1. Centers for Disease Control and Prevention. Treatment of Malaria: Guideline for Clinicians. Available at: [http://www.cdc.gov/malaria/diagnosis\\_treatment/clinicians3.html](http://www.cdc.gov/malaria/diagnosis_treatment/clinicians3.html). Accessed June 23, 2015.
2. World Health Organization. Guidelines for the Treatment of Malaria. 3<sup>rd</sup> ed. Geneva, 2015.
3. World Health Organization. Malaria: Fact Sheet #94. Updated December 2014. Available at: <http://www.who.int/mediacentre/factsheets/fs094/en/> Accessed online October 26, 2015.
4. Garcia LS. *Diagnostic Medical Parasitology*. 5<sup>th</sup> ed. Washington, DC. ASM Press; 2007.
5. CDC: Malaria surveillance – United States 2005. *MMWR* 2007;56(SS06);23-38.
6. Greenwood BM, Bojang K, Whitty CJ, Targett GA. Malaria. *Lancet*.2005;365:1487-1498.

## BP-09 Introduction

Thick and thin Giemsa-stained smears were obtained from 22-year-old Peace Corps worker in India with peripheral eosinophils and fever. The specimen contained *Microfilaria-Brugia* sp. A response of "Microfilaria-*Brugia* sp", "Microfilaria, NOS, referred for identification" and "Blood or tissue parasite, not *Plasmodium* sp. or *Babesia* sp., referred for identification" was considered satisfactory.

BP-09	Parasite Identification	Referees (34)		Participants (384)	
		No.	%	No.	%
	<i>Microfilaria-Brugia</i> sp	30	88.2	316	82.5
	<i>Microfilaria-Wuchereria bancrofti</i>	2	5.9	47	12.3
	<i>Microfilaria-Loa loa</i>	2	5.9	15	3.9
BP-09	Parasite Screen	Referees (18)		Participants (840)	
		No.	%	No.	%
	Microfilaria, NOS, referred for identification	18	100.0	745	88.7
	Blood or tissue parasite, not <i>Plasmodium</i> sp. or <i>Babesia</i> sp., referred for identification	-	-	68	8.1

## Discussion

### Causal Agent:

Lymphatic filariasis is caused by the filarid nematodes *Brugia malayi*, *B. timori*, and *Wuchereria bancrofti*. *Wuchereria bancrofti* is distributed nearly worldwide in the tropics. *Brugia malayi* is distributed in Southeast Asia and the Indian subcontinent, while *B. timori* is endemic to the Lesser Sunda Islands of the Indonesian archipelago.

### Biology and Life Cycle:

All three species have a similar life cycle, and adults of all three species reside in the lymphatic vessels of the human definitive host. Gravid females release sheathed microfilariae which circulate in the blood at night, exhibiting nocturnal periodicity (except for some populations of *W. bancrofti* in Southeast Asia which do not express specific periodicity). An appropriate mosquito intermediate host becomes infected while ingesting microfilariae during the course of a blood meal. Microfilariae migrate from the midgut of the mosquito to the flight muscles where they develop into infective L3 larvae in approximately two weeks. L3 larvae migrate through the hemocoel of the mosquito to the head and mouthparts. Humans become infected when a mosquito deposits L3 larvae onto the skin while taking a blood meal. Larvae migrate to the lymphatics system where it takes several months to develop into sexually mature adults.

### Diagnosis:

The diagnosis of all three species is based primarily on the identification of microfilariae in thick and thin blood films stained with Giemsa, Wright stain, or hematoxylin. Concentration procedures, such the Knott's method, may increase sensitivity. Because all three species exhibit nocturnal periodicity, the optimal time to collect blood specimens from a patient is between 10 PM and 2 AM.

All three species may possess a sheath, although the sheath may be absent in stained blood smears so the absence of a sheath should not in itself rule-out any of these species. The most important features for identifying

these nematodes to the genus level are the nuclear arrangements in the head and tail. The following table summarizes the important morphologic features.

Species	Size (in stained blood films)	Sheath (color, when properly stained with Giemsa)*	Head Space (distance between anterior end of nuclear column and tip of worm)	Tail nuclei
<i>Wuchereria bancrofti</i>	244-296 $\mu\text{m}$	Colorless*	Short	Tail anucleate
<i>Brugia malayi</i>	177-230 $\mu\text{m}$	Bright pink*	Long	Terminal and subterminal nuclei present, with gaps in between
<i>Brugia timori</i>	310 $\mu\text{m}$ avg. length	Colorless*	Long	Terminal and subterminal nuclei present, with gaps in between

\*Sheath color is pH-dependent, and at times the sheath of *B. malayi* may not stain bright pink. Likewise, on rare occasions, the sheath of *W. bancrofti* has been known to stain bright pink.

There are no routine molecular or rapid tests available for lymphatic filariasis in the United States. A rapid format immunochromatographic test is available outside the U.S. however. An EIA is available for detecting circulating antibodies in blood. Unlike with microscopy, blood does not need to be collected at night to perform the EIA. This test is reliable for *W. bancrofti* and *B. malayi*, but has not been properly validated for *B. timori*. There is also some cross-reactivity with *Onchocerca volvulus* and *Loa loa*.

Clinical Significance:

Most microfilarial infections are asymptomatic with subclinical tortuosity and dilation of lymphatics. The spectrum of disease for those with symptoms includes lymphedema, hydrocele, acute attacks of febrile lymphangitis and, less frequently, pulmonary tropical eosinophilia syndrome or chyluria. The range of clinical presentations varies slightly with species and geography. For example, involvement of the genital lymphatics occurs almost exclusively with *W. bancrofti* infection. Acute symptoms are often more intense in patients from non-endemic areas. With low worm burden and a good immune response, long-term sequelae in these patients are rare. In contrast, for those who live in endemic areas and sustain repeated bites by infected mosquitos, worm burdens are higher and lymphatics are more likely to become obstructed leading to chronic lymphedema. Lymphedema occurs more commonly in the lower extremities but can also involve the upper extremities, breasts in females and scrotum in males. Subsequent skin thickening and fissuring invites recurrent bacterial infection. With time, the lymphedema and skin changes can progress to elephantiasis.

Treatment:

The treatment of choice for active lymphatic filariasis is diethylcarbamazine (DEC) because it is both microfilaricidal and active against the adult worm. Adult worms must be killed in order to prevent relapse. However, DEC is contraindicated in patients with onchocerciasis co-infection and should be used with extreme caution in those with *Loa loa* infections. There is also some evidence that treatment targeting *Wolbachia*, the rickettsial endosymbiont bacteria that lives inside *Wuchereria* and *Brugia* spp., may stop microfilarial production. Due to low prevalence of the disease, DEC is no longer FDA-approved in the United States but can be obtained through the Centers for Disease Control and Prevention. Other therapeutic options include ivermectin (kills only microfilariae), and albendazole (has some macrofilarial activity). If lymphedema is already established, antifilarial medication has not been shown to be of benefit. Instead, management of symptoms includes exercise, elevation and local skin care.

---

## References

1. Ash LP, Orihel TC. *Atlas of Human Parasitology*, 5th ed. Chicago, IL: ASCP Press: 2007.
2. Centers for Disease Control and Prevention, Division of Parasitic Diseases (DPD). *Laboratory Identification of Parasites of Public Concern*. <http://www.cdc.gov/dpdx>
3. Garcia LS. 2007. *Diagnostic Medical Parasitology*, 5<sup>th</sup> ed., Washington, DC. ASM Press.
4. Chatterjee S, Nutman TB. "Filarial Nematodes." In *Manual of Clinical Microbiology*, ed. Jorgensen JH et al., 2461-2470. Washington, DC. ASM Press, 2015.


## BP-10 Introduction

Thick and thin Giemsa-stained smears were obtained from a 44-year-old Nigerian male with relapsing fevers, chills, and night sweats. The specimen contained *Plasmodium ovale*. A response of "*Plasmodium ovale*", "*Plasmodium vivax/ovale* NOS", "*Plasmodium* sp., not *P. falciparum* referred for identification", "*Plasmodium* sp., NOS would refer or request another specimen, or perform additional molecular testing", and "*Plasmodium* sp./*Babesia* sp. seen referred for identification" would have been considered satisfactory. However, consensus was not met by participants or referees for Parasite Identification.

Parasite Identification*	Referees (33)		Participants (445)	
	No.	%	No.	%
<i>Plasmodium ovale</i>	1	3.0	39	8.8
<i>Plasmodium vivax/ovale</i> , NOS	5	15.2	69	15.5
<i>Plasmodium</i> sp., not <i>P. falciparum</i> , referred for identification	16	48.5	137	30.9
<i>Plasmodium</i> sp., NOS would refer or request another specimen, or perform additional molecular testing	2	6.1	20	4.5
<i>Plasmodium malariae</i>	9	27.3	114	25.7
<i>Plasmodium vivax</i>	-	-	62	14.0
Parasite Screen	Referees (18)		Participants (778)	
<i>Plasmodium</i> sp./ <i>Babesia</i> sp. seen, referred for identification	18	94.7	743	95.5

\* Parasite identification was not graded due to lack of participant and referee consensus

If you have identified a *Plasmodium* sp. or *Babesia* sp.: Report percent of infected RBCs seen (number of infected RBCs per 100 red blood cells). (Ungraded)



Both thick and thin blood films should be prepared when the diagnosis of malaria is suspected. Examination of the thick blood film is considered the gold standard for identification since a larger blood volume can be examined than with the thin film, thus enabling the detection of low levels of parasitemia. In contrast, thin blood films provide the best morphology for species differentiation. In this Survey, the challenge was to recognize the presence of *Plasmodium ovale*. Consensus was not reached due to a number of laboratories identifying the species as *Plasmodium malariae* or *Plasmodium vivax*.

Key morphologic features on thin and thick blood films that characterize *P. ovale* are:

1. Enlarged size of the infected red blood cells (seen with both *P. ovale* and *P. vivax*).
2. Infected red blood cells with fimbriations and an ovoid shape (seen in up to 1/3 of infected cells)
3. Small red-pink granules in the infected red blood cell known as Schüffner's dots (seen in both *P. ovale* and *P. vivax*). Note: These may not be seen in early ring forms in *P. ovale*
4. Mature schizonts with 8-12 merozoites
5. Compact ring compared to the more amoeboid trophozoite of *P. vivax*.

Distinguishing *P. ovale* from *P. vivax* may be challenging when the trophozoites do not characteristic features and fimbriations are not abundant. The absence of schizont forms also makes the differentiation between the two species more difficult. In these cases a diagnosis of *Plasmodium vivax/ovale* may be acceptable. If available, molecular testing can be used confirm the species when morphology is not definitive.

**Note:** *The ideal specimen for laboratory identification of malaria is fresh capillary blood from a finger or heel stick, with immediate preparation of thick and thin blood films. Since this is not feasible in most settings, however, venous blood collected in EDTA anticoagulant is also acceptable. It is important to transport the blood as quickly as possible to the laboratory for examination since prolonged exposure to EDTA may result in distortion of the malaria parasites and compromised morphology.*

---

## Discussion

See BP-08 page 7 for general Malaria discussion.



## Actions Laboratories Should Take when a PT Result is Not Graded

The College uses Exception Reason Codes that signify the proficiency testing (PT) for an analyte has not been graded. The Exception Reason Code is located on the evaluation report in brackets to the right of the result. Your laboratory must identify all of the analytes with an Exception Reason Code and investigate the acceptability of performance with the same rigor as if it were an unacceptable performance. The actions accredited laboratories should take include but are not limited to:

<b>Code</b>	<b>Exception Reason Code Description</b>	<b>Action Required</b>
11	Unable to analyze.	Document why the specimens were not analyzed (eg, instrument not functioning or reagents not available). Perform and document alternative assessment (ie, split samples) for the period that commercial PT was not tested to the same level and extent that would have been tested.
20	No appropriate target/response; cannot be graded.	Document that the laboratory performed a self-evaluation using the data presented in the Participant Summary and compared its results to a similar method, all method, or all participant statistics if provided. If comparison is not available, perform and document alternative assessment (ie, split samples) for the period that commercial PT was not tested to the same level and extent that would have been tested.
21	Specimen problem.	Document that the laboratory has reviewed the proper statistics supplied in the Participant Summary. Perform and document alternative assessment for the period that commercial PT was not tested to the same level and extent that would have been tested. Credit is not awarded in these cases.
22	Result is outside the method/ instrument reportable range.	Document the comparison of results to the proper statistics supplied in the Participant Summary. Verify detection limits.
24	Incorrect response due to failure to provide a valid response code.	Document the laboratory's self-evaluation against the proper statistics and evaluation criteria supplied in the Participant Summary. Perform and document the corrective action of any unacceptable results. Document corrective action to prevent future failures.
25	Inappropriate use of antimicrobial.	Document the investigation of the result as if they were unacceptable and review the proper reference documents to gain knowledge of the reason your response is not appropriate.
26	Educational challenge.	Response to the CAP is not required. Laboratory should document its review.
27,31	Lack of participant or referee consensus.	Document that the laboratory performed a self-evaluation and compared its results to the intended response when provided in the Participant Summary. If comparison is not available, perform and document alternative assessment (ie, split samples) for the period that commercial PT reached non-consensus to the same level and extent that would have been tested.
28	Response qualified with a greater than or less than sign; unable to quantitate.	Document that the laboratory performed a self-evaluation and compared its results to the proper statistics supplied in the Participant Summary. Verify detection limits.
30	Scientific Committee decision.	Document that the laboratory has reviewed the proper statistics supplied in the Participant Summary.
33	Specimen determined to be unsatisfactory after contacting the CAP.	Document that the laboratory has contacted the CAP and no replacements specimens were available. Perform and document alternative assessment (ie, split samples) for the period that commercial PT was not tested to the same level and extent that would have been tested.
40	Results for this kit were not received.	Document why results were not received, corrective action to prevent recurrence and the laboratory's self-evaluation of the results by comparing results to the proper statistics and evaluation criteria supplied in the Participant Summary. If PT specimens were not analyzed, perform and document alternative assessment (ie, split samples) for the period that commercial PT was not tested to the same level and extent that would have been tested.
41	Results for this kit were received past the evaluation cut-off date.	Document why results were not received, corrective action to prevent recurrence and the laboratory's self-evaluation of the results by comparing results to the proper statistics and evaluation criteria supplied in the Participant Summary. If PT specimens were not analyzed, perform and document alternative assessment (ie, split samples) for the period that commercial PT was not tested to the same level and extent that would have been tested.
42	No credit assigned due to absence of response.	The Participant Summary indicates which tests are graded (see evaluation criteria) and which tests are Not Evaluated/Educational. Updates to grading will also be noted. If a test is educational, the laboratory is not penalized for leaving a result(s) blank. The code 42 that appears on the evaluation is <b>not</b> a penalty. However, if a test is graded (regulated and non-regulated analytes) and your laboratory performs that test, results cannot be left blank. The laboratory is required to submit results for <b>all</b> challenges within that test or use an appropriate exception code or indicate test not performed/not applicable/not indicated. Exceptions may be noted in the Kit Instructions and/or the Result Form. Document corrective actions to prevent future failures.
44	This drug is not included in our test menu. Use of this code counts as a correct response.	Verify that the drug is not tested on patient samples and document to ensure proper future reporting.
45	Antimicrobial agent is likely ineffective for this organism or site of infection.	Document that the laboratory performed a self-evaluation of written protocols and practices for routine reporting of antimicrobial susceptibility reports to patient medical records. Document that routine reporting of this result to clinicians for patient care is compliant with specific recommendations of relevant Medical Staff and Committees (eg, infectious Diseases, Pharmacy and Therapeutics, Infection Control). Response to the CAP is not required.
77	Improper use of the exception code for this mailing.	Document the identification of the correct code to use for future mailings.
91	There was an insufficient number of contributing challenges to establish a composite grade.	Document the investigation of the result as if it were an unacceptable result. Perform and document the corrective action if required.
35, 43, 88, 92	Various codes.	No action required.



## Attestation of Participation for Self-Reported Training\*

We the participants below have completed the review of the CAP BP-B 2016 Participant  
Product Mailing, Year

Summary/Final Critique report, and can self report the recommended 0.5 hours towards  
Education Hours

fulfilling education and certification of maintenance requirements.

Participant	Date	Participant	Date
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____
_____	_____	_____	_____

\_\_\_\_\_  
**Director (or Designee) Signature** - I have verified that the individuals listed      Date  
above have successfully participated in this activity.

**Retain this page for record-keeping and auditing purposes.**

Individuals can also track their participation of educational activities through the CAP Learning Management System (LMS).

1. Log in to [www.cap.org](http://www.cap.org), using your User ID and Password. If you don't have an online account, you will need to create one.
2. Click **Learning**, select **Learning Transcript**
3. Click **'Add My Own Activity'**
4. Enter the required information, and click **Save** when complete

For assistance, call our Customer Contact Center at 800-323-4040 or 847-832-7000 option 1.

*\*CAP Self-Reported Training activities do not offer CE credit, but can be used towards fulfilling requirements for certification of maintenance by agencies such as the American Society of Clinical Pathology (ASCP). Please verify with your certifying agency to determine your education requirements.*