

## MC 6.013 *Acinetobacter* spp. Susceptibility Reporting

---

**Purpose** This procedure provides instruction and guidance for routine testing and selective and cascade susceptibility reporting on *Acinetobacter* spp.

---

**Principal and Clinical Significance** The decisions for the most appropriate antimicrobial agents to test and report are made with input from Pharmacy, Infectious Disease and the Clinical Laboratory. The goal is to provide clinically relevant information that will decrease the chance of developing antibiotic resistance, harmful effects of inappropriate antimicrobial use and avoid reporting results that could adversely affect patient care using selective and cascade reporting rules.

---

**Policy Statements** This procedure applies to Microbiologists who perform susceptibility testing.

---

**Special Safety Precautions** Microbiologists are subject to occupational risks associated with specimen handling.

- [Biohazard Containment](#)
- [Biohazardous Spills](#)
- [Safety in the Microbiology Laboratory](#)

---

**Procedure**

1. Antibiotics appropriate for routine testing and reporting for *Acinetobacter* spp. include:
  - Ampicillin-Sulbactam
  - Ceftazidime
  - Cefepime
  - Ciprofloxacin
  - Levofloxacin
  - Gentamicin
  - Tobramycin
  - Imipenem
  - Meropenem
  - Amikacin-KBS with day of use QC
  - Piperacillin-Tazobactam
  - Trimethoprim-Sulfamethoxazole
  - Minocycline
  - Cefiderocol
  - Doxycycline
  - Cefotaxime
  - Ceftriaxone
2. Perform susceptibilities using Vitek cards N806 and XN30, MicroScan NUC101 or Kirby Bauer Method.
3. Not all antibiotics are available on every panel.
  - Vitek will be the primary method of testing. Both Vitek cards N806 and XN30 will need to be tested for *Acinetobacter* spp. isolates.
  - MicroScan and Kirby Bauer are back up if testing fails.
  - Do not perform testing on multiple methods to cover all the antibiotics. It is acceptable if cascaded antibiotics are not reported.

---

- Cefiderocol is not available in house and would need to be sent out if requested.
4. Vitek results will be accepted under the Online Tab. Modifications will be under the **VITMIC** keyboard under the Susceptibility tab.
  5. MicroScan and Kirby Bauer results will be entered manually under the **MMIC** and **KB** keyboards respectively under the Susceptibility tab.

**Selective Reporting**

1. Antibiotics will be reported in a specific order indicating the first and subsequent preferences of the Antimicrobial Stewardship Committee and Infectious Disease physicians
2. Antibiotics will be reported for all sources.
3. There will be exceptions based on the method used for testing.
4. There will be 6 antibiotics routinely reported for all sources.
  - Ampicillin/Sulbactam
  - Cefepime
  - Ciprofloxacin
  - Levofloxacin
  - Tobramycin
  - Gentamicin

Online Instrument Data Organism #2 - ACINETOBACTER BAUMANNII COMPLEX			
#2.1	ID-GN95 (VK2)	2.1	SUS-GN95 (VK2)
	ACINETOBACTER BAUMANNII	SS	AS(<=2), CEFE(1), CP(0.5), LEVO(<=0.12), TO(<=1), GM(<=1)
#2.2	ID-XN08 (VK2)	HIDE	CAZ(4-SS), IMP(<=0.25-SS), PIPT(<=4-SS), TS(<=20-SS)
	ACINETOBACTER BAUMANNII COMPLEX		


Figure 1 Selective reporting with 6 antibiotics reported

5. If any carbapenem is resistant, confirm result by alternate method. Report in Sunquest if confirmed resistant and submit isolate to MDH. Resistant carbapenems will be reported regardless of selective or cascade reporting. Refer to [MCVI 4.0 MDRO IP Notification Flow Chart](#) for further instructions.


**Cascade Reporting-Vitek Method**

If resistance is encountered, additional antibiotics will automatically be released.

1. For Vitek method, hidden antibiotics will be released when Accept Selected is clicked.



The released results will be visible in Results Summary.



2. Results can be edited under the **VITMIC** keyboard but is not needed to report the cascaded antibiotics.

**Method: Vitek**

- If Cefepime is R → report Meropenem
    - If Meropenem is R → report Imipenem
      - If Gentamicin is R → report Amikacin\*\*
      - If Tobramycin is R → report Amikacin\*\*
3. \*\*Perform Amikacin testing using the KB method and perform day of use QC.



Organism #2	- ACINETOBACTER BAUMANNII COMPLEX	
- VITMIC -		
	SS	AS(<=2),CP(0.5),LEVO(<=0.12),TO(<=1),GM(<=1),MERO(1)
	R	CEFE(32)
	HIDE	DOXY(1-SS),IMP(2-SS),MINO(1-SS),TS(<=20-SS),CAZ(4-SS),CAX(4-SS)

Figure 2 Meropenem released due to resistant Cefepime

- All other antibiotics on the panel that are not part of the cascade are hidden. Do not release and report. Hidden antibiotics include:
  - Trimethoprim/Sulfa
  - Minocycline
  - Ceftriaxone
  - Ceftazidime
  - Doxycycline
  - Piperacillin/Tazobactam
- With a provider request, antibiotics may be released.

**Cascade Reporting-MicroScan Method**

If resistance is encountered, additional antibiotics will automatically be released.


- For MicroScan method, hidden antibiotics will be released when File  is clicked. The released Results will be visible in Results Summary. 
- Enter **all** results manually under the **MMIC** keyboard. Antibiotics will be released following the cascade rules.

**Method: MicroScan**

- If Cefepime is R → report Meropenem

Org #2.	- ACINETOBACTER BAUMANNII COMPLEX	
- MMIC -		
	SS	AS(<=2),CP(0.5),LEVO(<=0.12),TO(<=1),GM(<=4)
	R	CEFE(64)
	HIDE	MERO(1-SS),MINO(1-SS),TS(<=20-SS),CAZ(4-SS)

Figure 3 Meropenem hidden

- It will appear as the results are hidden but Sunquest will release from HIDE automatically after results are filed. 

Organism #2	- ACINETOBACTER BAUMANNII COMPLEX	
- MMIC -		
	SS	AS(<=2),CP(0.5),LEVO(<=0.12),TO(<=1),GM(<=4),MERO(1)
	R	CEFE(64)
	HIDE	MINO(1-SS),TS(<=20-SS),CAZ(4-SS)

Figure 4 Meropenem released due to resistant cefepime

- All other antibiotics on the panel that are not part of the cascade are hidden. Do not release and report. Hidden antibiotics include:
  - Trimethoprim Sulfa
  - Minocycline

- Ceftazidime
- Piperacillin Tazobactam

5. With a provider request, antibiotics may be released.

**Cascade  
Reporting-Kirby  
Bauer Method**

If resistance is encountered, additional antibiotics will automatically be released.

1. Enter **all** results manually under the **KB** keyboard. Antibiotics will be released following the cascade rules.

**Method: Kirby Bauer**

- If Cefepime is R → report Meropenem
  - If Meropenem is R → report Imipenem

Org #3.	2+ ACINETOBACTER BAUMANNII COMPLEX	
- KB -	SS	AS(20), CP(22), LEVO(18), TO(16), GM(16)
	R	CEFE(14)
	HIDE	CAX(15-I), CAZ(16-I), IMP(23-SS), MERO(14-R), TS(18-SS)


Figure 5 Cefepime and Meropenem resistant

Organism #3	- 2+ ACINETOBACTER BAUMANNII COMPLEX	
- KB -	SS	AS(20), CP(22), LEVO(18), TO(16), GM(16), IMP(23)
	R	CEFE(14), MERO(14)
	HIDE	CAX(15-I), CAZ(16-I), TS(18-SS)

Figure 6 Meropenem and Imipenem released


2. All other antibiotics on the panel that are not part of the cascade are hidden. Do not release and report. Hidden antibiotics include:
  - Trimethoprim Sulfa
  - Ceftazidime
  - Ceftriaxone
  - Piperacillin Tazobactam
3. With a provider request, antibiotics may be released.

**Appendix**



Alternative methods:  
when reported drugs fail, have a product limitation or a physician wants to add a drug\*\*

The only drugs that can be reported on Acinetobacter per CLSI			NF KBS		
Disk (Abbrev.)	Sunquest code	Drug Name	Disk (Abbrev.)	Sunquest code	Drug Name
SAM	AS	Ampicillin/Sulbactam	ATM	AZTR	Aztreonam
AK	AK	Amikacin-KBS with QC	CAZ	CAZ	Ceftazidime
FEP	CEFE	Cefepime	CRO	CAX	Ceftriaxone
CIP	CFT	Cefotaxime	CIP	CP	Ciprofloxacin
CIP	CP	Ciprofloxacin	FEP	CEFE	Cefepime
CAZ	CAZ	Ceftazidime	GM	GM	Gentamicin
CRO	CAX	Ceftriaxone	IMP	IMP	Imipenem
	DOXY	Doxycycline	LUX	LEVO	Levofloxacin
GM	GM	Gentamicin	MEM	MERO	Meropenem
IMP	IMP	Imipenem	TZP	PIPT	Pip/Taz (ZOSYN)
LUX	LEVO	Levofloxacin	SXT	TS	Trimethaprim/Sulfa
MEM	MERO	Meropenem	NN	TO	Tobramycin
	Mino	Minocycline			
TZP	PIPT	Pip/Taz (zosyn)			
NN	TO	Tobramycin			
SXT	TS	Trimethaprim/Sulfa			
		Cefiderocol-Need to send to UM, if requested			



EBAC

Disk (Abbrev.)	Sunquest code	Drug Name
SAM	AS	Ampicillin/Sulbactam
CRO	CAX	Ceftriaxone
CIP	CP	Ciprofloxacin
FEP	CEFE	Cefepime
GM	GM	Gentamicin
IMP	IMP	Imipenem
LUX	LEVO	Levofloxacin
MEM	MERO	Meropenem
TZP	PIPT	Pip/Taz (ZOSYN)
SXT	TS	Trimethaprim/Sulfa
NN	TO	Tobramycin

MicroScan MIC Neg Urine Combo (NUC) 101

MSCN (Abbrev.)	Sunquest code	Drug Name
A/S	AS	AMP/SULB (UNASYN)
CAZ	CAZ	Ceftazidime
Cax	CAX	Ceftriaxone
Cpe	CEFE	Cefepime
Cp	CP	Ciprofloxacin
GM	GM	Gentamicin
Lux	LEVO	Levofloxacin
Mer	MERO	Meropenem
Min	Mino	Minocycline
P/T	PIPT	Pip/Taz (zosyn)
To	TO	Tobramycin
SXT	TS	Trimethaprim/Sulfa

\*\*If a physician wants to add-on a drug that we do not have, send to UM Fairview Medical Center\*\*

**References**

1. bioMerieux Vitek 2 AST-N806 Gram Negative Susceptibility Card 424709 2023-07
2. bioMerieux Vitek 2 AST-XN30 Gram Negative Susceptibility Card 424639 20235-04
3. Beckman Coulter Diagnostics. 250 South Kraemer Boulevard. Brea, CA 92821-6232 USA, MicroScan® Dried Gram Negative (8/2022).
4. CLSI M100 edition 34 Performance Standards for Antimicrobial Susceptibility Testing 2024

**Training Plan/  
Competency  
Assessment**

Training Plan	Initial Competency Assessment
<ol style="list-style-type: none"> <li>1. Employee must read the procedure.</li> <li>2. Employee will observe trainer performing the procedure.</li> <li>3. Employee will demonstrate the ability to perform procedure, record results and document corrective action after instruction by the trainer.</li> </ol>	<ol style="list-style-type: none"> <li>1. Direct observation.</li> </ol>

**Historical  
Record**

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	Susan DeMeyere	8/27/2024	Initial version