

MC 6.015 Enterobacterales with concern for AmpC Susceptibility Reporting

Purpose This procedure provides instruction and guidance for routine testing and selective and cascade susceptibility reporting on Enterobacterales, with concern for derepression of AmpC β -lactamase.

Principal and Clinical Significance The decisions for the most appropriate antimicrobial agents to test and report are made with input from Pharmacy, Infectious Disease and the Clinical Laboratory. The goal is to provide clinically relevant information that will decrease the chance of developing antibiotic resistance, harmful effects of inappropriate antimicrobial use and avoid reporting results that could adversely affect patient care using selective and cascade reporting rules.

Some Enterobacterales may develop resistance during therapy with third generation cephalosporins as a result of derepression on AmpC β -lactamase. This derepression is most commonly seen with *Citrobacter freundii* complex, *Enterobacter cloacae* complex and *Klebsiella aerogenes*. Isolates that are initially susceptible may become resistant within a few days after initiation of therapy. The decision to restrict reporting of Ampicillin, Ampicillin Sulbactam, Cephalosporins I, II, III and Piperacillin Tazobactam was made by the Antimicrobial Stewardship Committee for the organisms listed below.

Policy Statements This procedure applies to Microbiologists who perform susceptibility testing.

Special Safety Precautions Microbiologists are subject to occupational risks associated with specimen handling.

- [Biohazard Containment](#)
- [Biohazardous Spills](#)
- [Safety in the Microbiology Laboratory](#)

Procedure

1. Enterobacterales with concern for derepression of AmpC β -lactamase include:
 - *Citrobacter braakii*
 - *Citrobacter freundii*
 - *Citrobacter freundii* complex
 - *Citrobacter* species not listed in MC 6.014 without concern for AmpC
 - *Citrobacter werkmanii*
 - *Citrobacter youngae*
 - *Enterobacter* species
 - *Hafnia* species
 - *Klebsiella aerogenes*
 - *Kluyvera intermedia*
 - *Lelliottia* species
 - *Pantoea agglomerans*
 - *Pluralibacter* species
2. Antibiotics appropriate for routine testing and reporting for Enterobacterales with concern for derepression of AmpC β -lactamase include:
 - Cefepime
 - Gentamicin
 - Tobramycin
 - Ciprofloxacin
 - Levofloxacin
 - Trimethoprim Sulfamethoxazole

- Ertapenem
 - Imipenem
 - Meropenem
 - Cefiderocol
 - Ceftazidime-avibactam
 - Ceftolozane-tazobactam
 - Imipenem-relebactam
 - Meropenem-vaborbactam
 - Aztreonam
 - Amikacin
 - Nitrofurantoin-Urine only
3. Perform susceptibilities using Vitek cards N806 and XN30, MicroScan NUC101 or Kirby Bauer Method.
 4. Not all antibiotics are available on every panel.
 - Vitek will be the primary method of testing.
 - In most cases, only the N806 cards needs to be performed.
 - **Perform XN30 for respiratory cultures/CF patients for Tobramycin.**
 - MicroScan and Kirby Bauer are back up if testing fails.
 - Do not perform testing on multiple methods to cover all the antibiotics. It is acceptable if cascaded antibiotics are not reported.
 - Cefiderocol is not available in house and would need to be sent out if requested.
 5. Vitek results will be accepted under the Online Results tab. Modifications will be under the **VITMIC** keyboard under the Susceptibility tab.
 6. MicroScan and Kirby Bauer results will be entered manually under the **MMIC** and **KB** keyboards respectively under the Susceptibility tab.

Selective Reporting

1. Antibiotics will be reported in a specific order indicating the first and subsequent preferences of the Antimicrobial Stewardship Committee and Infectious Disease physicians.
2. Antibiotics will be reported depending on the source, if the sample is a urine or non-urine source.
3. There will be exceptions based on the method used for testing, type of culture and diagnosis of the patient.
4. For **urine sources**, these 5 antibiotics will be reported routinely, in this order.
 1. Cefepime
 2. Ciprofloxacin
 3. Gentamicin
 4. Trimethoprim-Sulfa
 5. Nitrofurantoin

Figure 1 -Only CEFE, CP, GM, TS, FD are reported. All other antibiotics are in HIDE.

Org #2.	>100000 COL/ML CITROBACTER FREUNDII
- VITMIC -	
SS	CEFE(2),CP(0.25),GM(2),TS(40),FD(32)
HIDE	ETP(0.5-SS),MERO(1-SS),LEVO(0.5-SS)

5. For **non-urine sources**, these 4-5 antibiotics will be reported routinely, in this order.
 1. Cefepime
 2. Ciprofloxacin
 3. Gentamicin
 4. Trimethoprim-Sulfa
 5. Tobramycin-respiratory cultures/CF patients

Figure 2-Only CEFE, CP, GENT, TS are reported. All other antibiotics are in HIDE.


Org #2.	3+ ENTEROBACTER CLOACAE COMPLEX
- VITMIC -	
SS	CEFE (2) ,CP (0.25) ,GM (2) ,TS (40)
HIDE	ETP (0.5-SS) ,MERO (1-SS) ,LEVO (0.5-SS) ,TO (2-SS)

6. If any carbapenem is resistant, confirm result by alternate method. Report in Sunquest if confirmed resistant and submit isolate to MDH. See Carbapenemase Detection section of procedure for further instructions. Resistant carbapenems will be reported regardless of selective or cascade reporting.


Cascade Reporting- Vitek Method

If resistance is encountered, additional antibiotics will automatically be released.

1. For Vitek method, hidden antibiotics will be released when Accept Selected is clicked.



The released results will be visible in Results Summary.



2. Results can be edited under the **VITMIC** keyboard but is not needed to report the cascaded antibiotics.

Method: Vitek -urine source

- If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin

Figure 3-Ertapenem, Meropenem and Levofloxacin released from HIDE due to resistance.

Organism #2	>100000 COL/ML CITROBACTER FREUNDII
- VITMIC -	
SS	GM (2) ,TS (40) ,FD (32) ,MERO (1) ,LEVO (0.5)
R	CEFE (16) ,CP (2) ,ETP (2)

Method: Vitek: -non-urine sources

- If culture type is respiratory → report Tobramycin
- If CF patient → report Tobramycin
 - If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin



Figure 4-Ertapenem and Levofloxacin released from HIDE due to resistance. Tobramycin released due to CF patient.

Organism #2		3+ ENTEROBACTER CLOACAE COMPLEX
- VITMIC - (ZZ01)		
SS		GM(2), TS(40), ETP(0.5), LEVO(0.5), TO(2)
R		CEFE(16), CP(1)
HIDE		MERO(1-SS), AZTR(0.5-SS), CZA(8-SS), CZT(2-SS), DOXY(1-SS), IMP(1-SS), IMIR(1-SS), MINO(4-SS), MVM(4-SS)

- All other antibiotics on the XN30 panel that are not part of the cascade are hidden. Do not release and report. Hidden antibiotics include:
 - Aztreonam
 - Imipenem
 - Ceftazidime-avibactam
 - Ceftolozane-tazobactam
 - Imipenem-relebactam
 - Meropenem-vaborbactam
 - Doxycycline
 - Minocycline
- Additional antibiotics may be released with a provider request.

Cascade Reporting-MicroScan Method


If resistance is encountered, additional antibiotics will automatically be released.

- Enter **all** results manually under the **MMIC** keyboard. Antibiotics will be released following the cascade rules.
- For MicroScan method, hidden antibiotics will be released when File  is clicked. Results will be visible in Results Summary. 

Method: MicroScan – urine source

- If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin

Figure 5- Ertapenem, Meropenem and Levofloxacin still in HIDE

Org #4.		>100000 COL/ML ENTEROBACTER CLOACAE COMPLEX (STRAIN 2)
- MMIC -		
SS		GM(2), TS(40)
R		CEFE(16), CP(1)
HIDE		ETP(2-R), MERO(1-SS), LEVO(0.5-SS) 


- It will appear as the results are hidden but Sunquest will release from HIDE automatically after results are filed. 

Figure 6- Ertapenem, Meropenem and Levofloxacin released from HIDE due to resistance

Organism #4 - >100000 COL/ML ENTEROBACTER CLOACAE COMPLEX (STRAIN 2)	
- MMIC -	
SS	GM(2), TS(40), MERO(1), LEVO(0.5)
R	CEFE(16), CP(1), ETP(2)

Method: MicroScan – non-urine sources

- If culture type is respiratory → report Tobramycin
- If CF patient → report Tobramycin
 - If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin



Figure 7 Example of CF patient with resistance.

Organism #2 - 4+ ENTEROBACTER CLOACAE COMPLEX	
- MMIC -	
SS	GM(2), TS(40), MERO(1), LEVO(0.5), TO(0.5)
R	CEFE(16), CP(1), ETP(2)
HIDE	AZTR(1-SS), CZA(4-SS), MVM(4-SS), MINO(4-SS)

4. All other antibiotics on the panel that are not part of the cascade are hidden. Do not release and report. Hidden antibiotics include:
 - Aztreonam
 - Minocycline
 - Ceftazidime-avibactam
 - Meropenem-vaborbactam
5. Additional antibiotics may be released with a provider request.

Cascade Reporting- Kirby Bauer Method

If resistance is encountered, additional antibiotics will automatically be released.

1. Enter **all** results manually under the **KB** keyboard. Antibiotics will be released following the cascade rules.
2. For Kirby Bauer method, hidden antibiotics will be released when File  is clicked. The released Results will be visible in Results Summary. 

Method: Kirby Bauer – urine source

- If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin

Figure 8 -Cefepime resistance

Org #3. >100000 COL/ML ENTEROBACTER CLOACAE COMPLEX	
- KB -	
SS	CP(26), GM(18), TS(16)
R	CEFE(18)
HIDE	LEVO(22-SS), TO(18-SS), MERO(24-SS), ETP(23-SS)


3. It will appear as the results are hidden but Sunquest will release from HIDE automatically after results are filed. 

Figure 9-Ertapenem released from HIDE

- KB -	
SS	CP (26) , GM (18) , TS (16) , ETP (23)
R	CEFE (18)
HIDE	LEVO (22-SS) , TO (18-SS) , MERO (24-SS)

Method: Kirby Bauer – non-urine source

- If culture type is respiratory → report Tobramycin
- If CF patient → report Tobramycin
 - If Cefepime is R → report Ertapenem
 - If Ertapenem is R → report Meropenem
 - If Ciprofloxacin is R → report Levofloxacin

Figure 10 -Cefepime and Ciprofloxacin resistance and CF patient

Org #2. 4+ ENTEROBACTER CLOACAE COMPLEX	
- KB -	
SS	GM (20) , TS (18)
R	CEFE (8) , CP (20)
HIDE	ETP (22-SS) , MERO (23-SS) , LEVO (22-SS) , TO (20-SS)

Figure 11 -Ertapenem and Levofloxacin released, Tobramycin released due to CF patient

- KB -	
SS	GM (20) , TS (18) , ETP (22) , LEVO (22) , TO (20)
R	CEFE (8) , CP (20)
HIDE	MERO (23-SS)

4. All other antibiotics on the panel that are not part of the cascade are hidden. Do not release and report.
5. Requested antibiotics may be released with a provider request.

Method Performance Specifications

1. Usually, only the Vitek card AST-N806 will be required for testing Enterobacterales with concern for AmpC.
2. Report Tobramycin on respiratory cultures and CF patients
 - Set up the XN30 to obtain result for Tobramycin and release only the Tobramycin. Refer to MC 6.02 AST-XN30 for further reporting instructions.
3. Mucoid isolate susceptibility testing should be performed with the Kirby Bauer method.

Carbapenemase Detection

1. If a carbapenem exhibits resistance, confirm the resistant result with an alternate method before reporting results.
 - Confirmation is not needed if patient has known resistance.
2. If resistance is not confirmed, report the MicroScan or KB result.
3. With known or confirmed resistance:
 - Report MIC results as tested. Do not change interpretations.
 - Label as **MDRO**
 - Notify provider
 - Send isolate to MDH Project 1380

- Add MDHADD comment
 - Freeze isolate
4. Notify Infection Control with positive mCIM result
 5. Report MDH results with NCPO with negative mCIM result or the appropriate Carbapenemase Producer code. e.g. VIM, KPC, NDM, etc.

Carbapenemase Result Reporting

1. Result Enterobacterales with resistance to a carbapenem and negative mCIM result as: example-ENCL-MDRO-NCPO
2. Result Enterobacterales with resistance to a carbapenem and positive mCIM result as: example-ENCL-MDRO-OXA
3. Add SCAND when MDH reports are scanned.
4. Add PRAC with known or previous resistance.

Vitek Product Limitations

Results for an antibiotic/organism combination may have limitations and may be suppressed from reporting. Refer to table below for specific limitations.

AST-N806 card for Enterobacterales with concern for AmpC

Antibiotic	Product Limitations
Cefepime	Perform KB on <i>Hafnia alvei</i>
Ertapenem	Perform alternate method with MIC of 0.25-0.5 Perform alternate method on <i>Hafnia alvei</i>
Gentamicin	If resistance is observed on <i>Klebsiella aerogenes</i> , <i>Enterobacter cloacae</i> . Send isolate to MDH.
Meropenem	Perform alternate method on <i>Enterobacter cloacae</i> with MIC \geq 16

AST-XN30 card for Enterobacterales with concern for AmpC


Antibiotic	Product Limitations
Amikacin	Not validated. Perform KB and day of use QC
Imipenem	Perform alternate method on <i>Klebsiella aerogenes</i>

Appendix


Additional Tables and Flowcharts

Antibiotics that can be reported for Enterobacterales with concern for AmpC	
N806 Card	Cefepime
	Ciprofloxacin
	Ertapenem
	Gentamicin
	Levofloxacin
	Meropenem
	Trimethoprim Sulfamethoxazole
	Nitrofurantoin-Urine only
XN30 Card	Amikacin
	Aztreonam
	Ceftazidime-avibactam
	Ceftolozane-tazobactam
	Imipenem
	Imipenem-relebactam
	Meropenem-vaborbactam
	Tobramycin
	Cefpodoxime-urine only, if cefazolin R and ESBL negative on EC, KLPN, PRMI
	Cefiderocol- UM sendout


Enterobacterales with concern for AmpC
Citrobacter species- unless listed in the without AmpC procedure
Citrobacter braakii
Citrobacter freundii
Citrobacter freundii complex
Citrobacter werkmanii
Citrobacter youngae
Enterobacter species
Hafnia species
Klebsiella aerogenes
Kluyvera intermedia
Lelliottia species
Pantoea agglomerans
Pluralibacter species



**Microscan-
Urine and Non-urine**
Results will appear hidden until they are accepted in Sunquest



**Kirby-Bauer (KB)-
EBAC stamper**
Urine and Non-Urine
Results will appear hidden until they are accepted in Sunquest



References

1. bioMerieux Vitek 2 AST-N806 Gram Negative Susceptibility Card 424709 2023-07
2. bioMerieux Vitek 2 AST-XN30 Gram Negative Susceptibility Card 424639 20235-04
3. Beckman Coulter Diagnostics. 250 South Kraemer Boulevard. Brea, CA 92821-6232 USA, MicroScan® Dried Gram Negative (8/2022).
4. CLSI M100 edition 34 Performance Standards for Antimicrobial Susceptibility Testing 2024

**Training Plan/
 Competency
 Assessment**

Training Plan	Initial Competency Assessment
<ol style="list-style-type: none"> 1. Employee must read the procedure. 2. Employee will observe trainer performing the procedure. 3. Employee will demonstrate the ability to perform procedure, record results and document corrective action after instruction by the trainer. 	<ol style="list-style-type: none"> 1. Direct observation.

**Historical
 Record**

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	Susan DeMeyere	11/12/2024	Initial version
2			