

Agents of Bioterrorism Competency Exam

Instructions:

- Review power point
- Reference Guides (both are by the mail boxes):
 1. “Bioterrorism Guide for Clinical Labs” from the CDC
 2. Microbiology: Bioterrorism Procedure Manual (gray binder)
- Answer all questions on exam (**on Medtraining.org**)
- **DUE DATE: December 5, 2014**

Laboratory Response Network (LRN)

Classification of Laboratories

Sentinel Labs: Clinical labs that operate under BSL2 levels, some may have a BSL3 lab that is used for Cat A infectious substances and also for possible bioterrorism agents. Rule out or forward isolates. Hospital based, clinical or commercial diagnostic labs

Reference Labs: State Public Health Labs, BSL2 with BSL3 practice. Definitive confirmation of isolates

National Labs : BSL4 practice. Highest level of organism characterization and containment. Includes CDC, Military Research and Infectious Disease institutes

(taken from the CDC-Bioterrorism Response Guide for Clinical Labs)

Brucella spp.

- Endemic to Mediterranean basin, Middle East, western Asia, Africa and Latin America
- Bioterrorism CDC list
- **Best Sources:**
 - Blood (most sensitive)
 - Bone marrow
 - Tissue (lymph node)
- Symptoms:
 - profound weakness
 - diurnal fevers
 - generalized muscle aches
- Disease:
 - Brucellosis
 - Undulant fever
 - Malta fever




Brucella sp.

- Gram stain: tiny gram neg coccobacilli (GNCB) that stain faintly
- Growth: BAP, Choc and MTM, no growth on MAC (specialized media)
- Some species require 5-10% CO₂
- Slow grower: 48-72 hours, hold culture 7-10 days
- Colonies appear pinpoint and are not hemolytic
- Suspected: work up must be done in **BSL3 lab** with appropriate PPE (gloves, lab coat, N95 mask)

Brucella sp

Tiny Gram negative
coccobacilli



Small gamma, gray
colonies at 48-72 hrs



Key Characteristics of Brucella spp. (Rule in/out process)

If organism is:

- Slow Grower: 48-72 hrs or longer in CO₂
- Small gray, non-hemolytic
- Gram negative coccobacillus, stains faintly
- Catalase: Pos
- Oxidase: Pos
- Urea: Pos (strong < 2 hrs)
- Non-motile
- Presumptive ID-→Send to MDCH for confirmation


Francisella tularensis

- Isolated from 100's of animals:
- mammals (wild rabbits, sheep, beaver, muskrat)
 - avian spp.
 - arthropods (ticks, deerflies, mosquitos)
 - domestic animals



Franciscella tularensis

Faint staining tiny
GNCB



Glucose Cystine Blood
agar



Key Characteristics of *Francisella tularensis*

(Rule in/out process)

If organism is:

- Growth on: BAP with cystine, CHOC, MTM, BCYE, no growth on MAC
- 5% sheep blood: no growth or weak
- 35-37C , no CO2 needed
- Slow Grower: 48-72 hrs or longer, hold 7 days
- Small gray, alpha-hemolytic
- Very tiny Gram negative coccobaccilli, stains faintly
- Catalase: Neg or weak pos
- Oxidase: Neg
- Urea: Neg
- Non-motile
- Work in **BSL3 lab**
- Presumptive ID-→Send to MDCH for confirmation


Yersinia pestis

- Gram stain : Gram negative rods (bipolar/safety pin)
- Growth: best at 4C or room temperature(21-25C)
- Selective media = CIN(cefsulodin-irgasan-novobiocin)
- Mac and Hek: slow-lactose fermenter
- CIN morphology: small, red center, transparent edge,(bulls eye)
- Suspect--→ **work in BSL3 lab**

Yersinia pestis


Gram stain

(may or may not see safety pin)



Wright- Giemsa stain

(safety pin)



Yersinia pestis

MAC at 24-48hours




BAP at 48-72 hrs



Yersinia pestis

- CIN agar at RT



Key Characteristics of *Yersinia pestis*

(Rule in/out process)

If organism is:

- Gram negative bacillus(may resemble “safety pin”)
- Grows at 35-37C, faster at room temperature
- Catalase : positive
- Non-motile at 37C and RT (other species of *Yersinia* are motile at RT)
- Oxidase: negative
- Identification by enteric biochemicals

Presumptive ID-→ send to MDCH


Bacillus anthracis

- Non-motile
- Non-hemolytic
- Have irregular edges to colonies
- Tenacious
- “Medusa head”
- Samples:
 - ✓ Blood
 - ✓ Vesicular fluid
- Suspect on culture: **BSL3 lab**
- Send to MDCH to rule out anthracis



Gram stain of Bacillus species

No spores seen



Spores seen (stressed)




Bacillus anthracis on BAP---tenacious



B. anthracis vs. B. species on BAP

**Bacillus anthracis
(medusa head)**



*Courtesy of Larry Stauffer,
Oregon State Public Health Laboratory*

Bacillus sp. not anthracis



ASM MicrobeLibrary.org © Buxton

Botulinum Toxin

- Lab must be notified if suspected
- MDCH must be notified immediately
- **Responsibility of the Level A** lab is limited to advising medical staff on specimen selection, packing, shipping and notification of recipient lab about suspected case.
- Specimens: Enema fluid, Food sample, Nasal swab, Serum, Stool, Enviromental (if intentionally released)
- **Sentinel labs should not attempt to culture this organism!!**

Clostridium botulinum

ANA BAP/spready colony



Gram stain/terminal or subterminal spores



Agents of Bioterrorism Competency Exam

Complete exam attached with the
powerpoint on MTS

This is an open book exam

DUE DATE: 12/5/2014