

Neisseria meningitidis

Workup

Process Improvement Project
By Milad Latif

LEARNING OBJECTIVES

After completing this *N. meningitidis* workup training session, the user will be able to:

- 1- Describe the culture growth conditions for *N. meningitidis*
- 2- Identify the proper steps required for a suspected *N. meningitidis* culture.


NATURE OF THE PROBLEM

Exposure to *Neisseria meningitidis* is a serious hazard in the microbiology laboratory. Medical technologists may become exposed to this pathogen without knowing while working up patient specimens.

The goal of this project is to inform the staff about the proper procedure for handling specimens suspicious of containing *N. meningitidis*.

Ensuring the staff safety was the main reason behind choosing this issue to work on.

The project will address the following points:

- Description of *Neisseria meningitidis* pathogenicity
 - Comparison to other morphologically similar organisms
 - Breakdown for the procedure
- 
- A series of three parallel white diagonal lines in the bottom right corner of the slide, extending from the middle of the right edge towards the bottom left.

NEISSERIA MENINGITIDIS

- ▶ *Neisseria meningitidis* is an opportunistic pathogen which can colonize the mucous membranes of humans and may cause significant infection.
- ▶ Medical laboratory technologists are at increased risk of acquiring infections when dealing with *N. meningitidis* isolates.
- ▶ Requires 35-37C, with 5-10% CO₂ for growth in the laboratory
- ▶ Growth on Chocolate, Modified Thayer Martin (MTM) and Blood agars

NEISSERIA MENINGITIDIS

- ▶ *Neisseria meningitidis* is a leading cause of bacterial meningitis and sepsis in the United States. It can also cause focal disease, such as pneumonia and arthritis.



Meningococcal Disease (very serious and can be fatal)

Symptoms include sudden fever, headache, and stiff neck. Other symptoms can include nausea, vomiting, increased sensitivity to light, and confusion.

MORPHOLOGY



- ▶ Chocolate agar:
1-2 mm creamy and gray colonies
- ▶ Blood agar:
Grey and non pigmented and appear
round, smooth, moist, glistening, and convex

N. gonorrhoeae may have similar morphology on Chocolate agar, but shouldn't grow on blood agar except in rare cases.



GRAM STAIN

Gram-negative diplococcus
May be seen intracellularly in PMN's



WORKUP

Suspected colony on Chocolate agar + Growth on blood agar



Move to BSC
hood



Perform Oxidase test on the suspected organism




Positive



Perform Gram stain on the suspected morphology

Workup - continued



- If MALDI is not available, perform RapID NH.

WORKUP - CONTINUED

ID is *N. gonorrhoeae*



Report as presumptive *N. gonorrhoeae* and send to MDHHS

ID is *N. meningitidis*



Report result according to appropriate procedure. Send to MDHHS only if from sterile site

ID is Other organism



Workup as appropriate

Identification Flowcharts **Attachment 10:** **Gram Negative Diplococci** **and Gram Negative** **Coccobacilli**



Start Here

Suspect colony growing on Chocolate

(+) Growth on BAP

(-) No Growth on BAP

Workup in the RED box must be done in the BSC.^c

Oxidase = Positive

Gram Stain = Gram Negative Diplococci

Gram Stain = Gram Negative Coccobacilli

Gram Stain = Other

Oxidase = Negative
Work up as appropriate.

POSSIBLE *N. MENINGITIDIS*!
Perform MALDI^a.

Possible *Haemophilus*
Perform MALDI.
Report as appropriate.

Work up as appropriate.

ID = *N. meningitidis*

ID = *N. gonorrhoeae*

ID = Other organism:
Work up as appropriate.

Report result according to appropriate procedure.
Send to MDHHS.^b

Report as Presumptive *N. gonorrhoeae*.
Send to MDHHS^b.

Oxidase = Positive

Gram Stain = Gram Negative Diplococci

Gram Stain = Gram Negative Coccobacilli

Oxidase = Negative

Gram Stain = Gram Negative Coccobacilli

Gram Stain = Other

Possible *N. gonorrhoeae*
Perform MALDI.

Possible *Haemophilus*
Perform MALDI.
Report as appropriate.

Work up as appropriate.

ID = *N. gonorrhoeae*

ID = Other organism:
Work up as appropriate.

Genital Source:
Report as *N. gonorrhoeae*.
Send to MDHHS^b.

Non-Genital Source:
Report as Presumptive *N. gonorrhoeae*.
Send to MDHHS^b.

KEY

a: Until *N. meningitidis* has been ruled out, the MALDI target must be inoculated with test organism within the BSC hood. Once organism is dry on target, reagents can be added on the bench. All preparation steps of the extractions must be performed in the hood.

b: MDHHS (Michigan Department of Health and Human Services). *N. meningitidis* from sterile sources and all isolates of *N. gonorrhoeae* must be forwarded to MDHHS for confirmatory testing.

c: BSC (Biological Safety Cabinet) must be utilized until *N. meningitidis* is successfully ruled out.