

Mycobacteriology Service Implementation

Module 3
Version 2019

*Sheryl R. Stuckey, MLS (ASCP)
Manager, Microbiology Laboratory
Holy Cross Hospital
Silver Spring, Maryland*

SPECIMEN PROCESSING



Processing of CSF and Sterile Body Fluid

- ▶ During specimen quality assessment check for companion routine bacterial cultures
- ▶ Note the volume and appearance of the fluid
- ▶ Check the gram stain/culture results
- ▶ If the fluid is of sufficient volume, is cloudy or organisms are seen in the gram stain, **process the specimen as non-sterile**
- ▶ If the fluid is clear and no organisms are seen in the gram stain, **centrifuge the specimen and use the sediment to directly inoculate media and microscope slides**

Processing of Nonsterile Specimens

This is a 3 part process

1. Digestion – uses mucolytic agents to liquefy the specimen and release the AFB
2. Decontamination – uses agents toxic to normal flora contaminants to prevent their overgrowth and allow mycobacteria to flourish
3. Centrifugation – used to concentrate AFB in a specimen sediment for inoculation to solid agar and broth culture media

Processing of Nonsterile Specimens

- ▶ Digestion/Decontamination reagents are toxic to mycobacteria
- ▶ Processing must be precisely timed to limit decreasing the yield of mycobacteria
- ▶ Centrifugation creates high heat that also can decrease the yield of mycobacteria
- ▶ Must monitor culture positivity to detect if process is too harsh or centrifugation is poor
- ▶ Must monitor contamination rates to detect if process reagents are not strong enough

Processing of Nonsterile Specimens

- ▶ Basic Processing Steps include:
 - ❖ Disinfecting the workspace
 - ❖ Lining the workspace with a disinfectant soaked absorbent pad
 - ❖ Placing specimen tubes in a rack with a space between each tube
 - ❖ Process no more than a centrifuge load of specimens in each processing batch
 - ❖ Open and handle only one specimen container at a time

Digestion / Decontamination Methods

- ▶ N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) – Principle method
- ▶ NaOH
- ▶ Oxalic Acid – *Pseudomonas species*; Cystic Fibrosis patients
- ▶ Sulfuric Acid – Urines, Gastric Lavage, Watery fluids
- ▶ Bleach – May be used for direct smears

General Specimen Processing Preparation

- ▶ Perform specimen quality assessment, noting quantity and quality of the specimen
- ▶ Contact the primary caregiver prior to processing suboptimal specimens (i.e. QNS may require recollection)
- ▶ Let all refrigerated items warm to room temperature
- ▶ Make sure that you have sufficient digestion, decontamination and buffer reagents
- ▶ Make sure you have sufficient PPE supplies (i.e. gloves for changing if needed)
- ▶ Make sure you have adequate labels and usable pen/markers for labeling

General Specimen Processing Preparation

- ▶ Prepare the BSC workspace
- ▶ Label microscope slides, media, individual reagent vials for each specimen
- ▶ Label centrifuge tubes as needed for the specimen (submitted in non-50mL conical centrifuge container), digestion aliquot

General Digestion/ Decontamination Procedure

- ▶ Add digestion/decontamination reagent equal to the volume of the specimen to the tube
- ▶ Vortex to mix and liquefy the sample (10 – 30 seconds)
- ▶ Let stand for 15 – 20 minutes, vortexing intermittently to mix
- ▶ Q.S. to 50 mL volume with buffer reagent.
- ▶ Centrifuge for 15 minutes at 3000 x g
- ▶ Decant the supernatant in to a second 50 mL conical tube; discard
- ▶ Re-suspend the sediment with buffer
- ▶ Prepare a smear with 1 – 2 drops (~100 uL) and heat fix.
- ▶ Inoculate the broth culture media with 0.5 mL of the sediment
- ▶ Inoculate solid agar media with 1 – 3 drops of sediment
- ▶ Incubate media appropriately for the specimen source

General Processing

- ▶ <https://www.youtube.com/watch?v=YIFFSbAoCdM>
- ▶ <https://www.youtube.com/watch?v=SG9WbK9rZVI>

SPECIMEN PROCESSING METHODS

- » Goal: Increase test sensitivity and decrease probability of false positives

N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) Decontamination Method

- ▶ NALC = Mucolytic Agent – dissolved in diluent
- ▶ Very susceptible to air; must be used within 24 hours
- ▶ Additional NALC can be added to very mucoid specimens
- ▶ NaOH acts as decontamination reagent; mixed with Na Citrate that binds heavy metal ions that can inactivate NALC

NaOH Decontamination Method

- ▶ Normal decontamination method uses a 2% – 3% NaOH method
- ▶ NaOH at a FINAL concentration of $\geq 2.0\%$ can be lethal to mycobacteria (may see decrease culture sensitivity in smear negative specimens)
- ▶ 4% NaOH method is used for specimens with contaminating flora resistant to the routine NALC–NaOH method
- ▶ Na Citrate buffer may be included as part of commercial solutions
- ▶ Sterile physiological saline used in lieu of chemical buffers (M15 Phosphate, N67, etc.)

Oxalic Acid

- ▶ 5% Oxalic Acid is used in lieu of the normal reagents
- ▶ This method is used with specimens that are contaminated with *Pseudomonas* species
- ▶ Sterile physiological saline is used in lieu of chemical buffers
- ▶ Phenol red used to check pH as NaOH is used to neutralize and re-suspend the pellet

Sulfuric Acid

- ▶ 4% Sulfuric Acid is used in lieu of the normal reagents
- ▶ This method is used with urines and watery body fluid specimens that remain contaminated using alkaline methods
- ▶ Sterile water is used in lieu of chemical buffers
- ▶ Phenol red used to check pH as NaOH is used to neutralize and re-suspend the pellet

Bleach Decontamination

- ▶ Used only for the preparation of direct specimen smears in an emergency situation to ensure non-viability of AFB
- ▶ Aliquot of sample should be retained for further studies prior to decontamination procedure
- ▶ Strict timing required to avoid destruction of all AFB
- ▶ Sterile water is used in lieu of chemical buffers to wash the sediment

SPECIMEN PROCESSING QUALITY CONTROL

- »» Goal: Improve the ability to detect false test results and increase the accuracy of the test result

Specimen Processing Quality Control

- ▶ Leave 1 – 2 spaces between each tube in the batch (this includes split sample tubes for the same sample)
- ▶ Never batch process more than a full centrifuge load (8 – 12 specimens)
- ▶ Include a Negative Processing Control (10 ml sterile water or buffer) with each batch of specimens processed (*This sample goes through the entire process and media inoculation as the patient samples*)

Avoiding Cross Contamination

- ▶ Open the specimen tubes very gently to avoid aerosol generation.
- ▶ When adding reagents to the tube, open one tube at a time. Do not keep all the tubes open at the same time.
- ▶ Place distance between tubes in the rack
- ▶ Change gloves often as needed
- ▶ Never set up proficiency test samples in the same batch or just before setting up patient samples
- ▶ Allow at least 15 – 20 minutes of airflow exchanges between batches



to be continued