

Leucine Aminopeptidase (LAP) Test

I. Principle

The LAP Test is a rapid assay for the detection of the enzyme leucine aminopeptidase in bacteria cultured on laboratory media. It is used as one of the tests for the presumptive identification of catalase-negative gram positive cocci. Leucine- β -naphthylamide impregnated discs serve as a substrate for the detection of leucine aminopeptidase. Following hydrolysis of the substrate by the enzyme, the resulting β -naphthylamine produces a red color upon the addition of cinnamaldehyde reagent.

II. Reagents and Materials

- A. LAP discs: reagent discs impregnated with a solution of leucine- β -naphthylamide (0.7%) then dried (Murex). Store at 2-8°C until stated expiration date.
- B. Color developer: 0.01% p-dimethylamino-cinnamaldehyde
- C. Sterile deionized water
- D. Wooden applicator stick or inoculating loop
- E. Forceps

III. Quality Control

- A. Quality control is performed upon receipt of each new lot or shipment. The control organisms and their expected reactions are as follows:
 - 1. *Enterococcus faecalis* ATCC 29212 = positive (pink/red color development)
 - 2. *Aerococcus viridans* ATCC 11563 = negative (no color/slight yellow color development)
- B. If controls do not display expected results, quality control must be repeated. Notify the supervisor.

IV. Procedure

- A. Using forceps, place a LAP disc on a glass slide or in a petri dish.
- B. Moisten the LAP disc with 10 μ L of sterile deionized water.
- C. Select 5-10 pure or isolated colonies of catalase-negative, gram positive cocci with a loop or stick and rub onto the disc.
- D. Incubate at room temperature for 5 min.
- E. Add one drop of the color developer, and examine for up to 1 min for pink to red color development.

V. Interpretation

- A. Pink/red = positive
- B. No color change/slight yellow = negative

VI. Limitations

- A. The LAP Test is only part of the overall scheme for identifying catalase-negative, gram positive cocci. Further biochemical characterization and serological grouping may be necessary for specific identification. False negatives may result from using too small an inoculum.
- B. When Gram stains are prepared from agar growth, some *Leuconostoc* and *Streptococcus* strains may appear coccobacillary or even rod-shaped and may be confused with members of the genus *Lactobacillus*. The most consistent Gram stains are prepared from broth media.

VII. References

- A. Coleman, G., and L.C. Ball. 1984. Identification of Streptococci in the Medical Laboratory, J. Appl. Microbiol. 57:1-14.
- B. Facklam, R.R. and J.A. Washington. 1990. Streptococci and Related Catalase Negative Gram-Positive Cocci, Manual of Clinical Microbiology, 5th Ed., ASM, Washington, DC.
- C. Lennette, E.H., et al. 1985. Manual of Clinical Microbiology, 4th Ed., ASM, Washington, DC.
- D. Murex Diagnostics LAP package insert. Sept 1993.

Effective 01/24/2006

Reviewed by Microbiology Director, Dr. Ann Robinson: 01/24/2006

Reviewed by Medical Director, Dr. Joseph Schappert: 03/10/2010

Reviewed by Microbiology Supervisor, Jerry Claridge: 01/24/2006, 11/2006, 10/2007, 05/2008, 05/2009, 04/01/2011, 03/2013, Jason Ammons 05/2015

Updates and Revisions: