

# MYCOLOGY CULTURE EVALUATION

## I. Principle

Mycology plates are interpreted frequently to optimize the recognition and identification of yeasts and molds.

## II. Specimen

Inoculated mycology culture plates

## III. Materials

- A. Glass slides
- B. 22 X 22 mm coverslips
- C. Scotch tape
- D. Light Microscope with 10X and 40X objectives
- E. LPCB (lactophenol cotton blue) stain, store at room temperature until the expiration date.
- F. Sterile water, store at room temperature until the expiration date.
- G. Wooden applicator sticks

## IV. Schedule

- A. Blood Isolators
  - 1. Examine the SAB and BHIA with gent and chloro Monday, Wednesday, and Friday.
  - 2. Refer to the "**Fungal Blood Cultures: DuPont Isolator System Method**" procedure for the workup of positive cultures.
- B. Routine Mycology cultures
  - 1. Additional Friday duties
    - a. Interpret the mycology cultures (day 1 through week 4).
    - b. Perform all of the workups on new and old isolates.
  - 2. Additional Monday duties
    - a. Interpret all of the mycology cultures (day 1 through week 4).
    - b. Perform all of the workups on new and old isolates.
    - c. Move the cultures ahead 1 week to the next spot in the incubator.
    - d. Finalize the reports on the 4 week old cultures under MNG (Micro No Growth) if negative, and discard the media into appropriate receptacles.
    - e. Review the Incomplete log.
  - 3. Additional monthly duties
    - a. Subculture all of the stock yeast and mold isolates on the last Friday of every month.
- C. Cultures referred from Bacteriology and Mycobacteriology
  - 1. Perform identifications on these isolates to the genus and species level.
  - 2. Report in Mysis, and charge appropriately.

## V. Procedure

- A. Remove the positive plates from the 30° incubator, and place them on the Mycology cart.
- B. Review each culture, and bring positive cultures up on rounds.
- C. Plates with yeast-like or bacteria-like growth: perform a wet prep examination on all of the different colony types.
  - 1. Wet prep procedure

- a. Place a drop of water on an appropriately labeled glass slide.
  - b. Touch a suspected colony with an applicator stick, and inoculate it into the drop of water. Apply a cover slip.
  - c. Examine the smear under the 10X and 40X objectives for bacteria/yeast.
  - d. If yeast is observed, refer to the Yeast Identification Procedure.
- D. Yeasts for antifungal susceptibility testing
1. Refer to the *Candida* Disk Diffusion Procedure for routine testing for fluconazole and voriconazole.
  2. If other antifungal agents are requested, the antibiotics should be specified by the requesting physician.
  3. Subculture the yeast onto 2 SAB slants.
  4. Give 1 SAB slant and the test request to PAML or SHMC sendouts to be sent to ARUP for MIC testing.
  5. Store the other SAB slant at 2-8°C in the Microbiology.
- E. *Nocardia* for Susceptibility Testing
1. This testing is only performed on physician request. The physician must request what antibiotics should be tested.
  2. Subculture the organism to 2 SAB slants.
  3. Submit 1 SAB slant with the request to the PAML or SHMC sendout department.
  4. *Nocardia* isolates are sent for susceptibility testing to ARUP.
  5. Store the second SAB slant in the mycology lab until the results are final.
- F. Plates with mold:
1. Refer to Figure 1 for mold isolate workup.
  2. Determine the specimen type.
    - a. Sterile body fluids and nonsterile sites except for hair/skin/nails with one colony of mold or one colony of several molds: bring it up on rounds before performing the lactophenol cotton blue preparation.
    - b. Sterile body fluids with numerous colonies of a similar mold on one or both plates: perform the lactophenol cotton blue scotch tape preparation (refer to "**Lactophenol Cotton Blue Preparations**" procedure), and bring it up on rounds.
    - c. Nonsterile body fluids with numerous colonies of a similar mold on one or both plates: perform a lactophenol cotton blue scotch tape preparation, and bring it up on rounds.
    - d. Hair/skin/nails: perform a lactophenol cotton blue preparation on all mold colonies.
      - i. If the mold shows sterile hyphae on the LPCB prep, subculture the isolate to potato flake agar.
  3. Bring up all significant mold isolates on rounds.
  4. Mold isolates should be identified, if necessary, using Figures 2-7.
  5. Record all of the isolate's characteristics in the Misys computer.
    - a. Indicate the color of both the surface and the reverse of the colony.
    - b. Describe any typical microscopic morphologic characteristics observed.
- G. Notify Rounds of any significant bacterial isolates on sterile body sites.
- H. Follow the appropriate procedures for the identification of yeasts and molds.

## VI. Result Reporting

- A. Report "No fungus isolated" on all final mycology cultures that exhibit no fungal growth after the appropriate length of incubation.

- B. Send a preliminary report of "*Candida albicans* or yeast, not *Candida albicans*" on cultures demonstrating yeast.
- C. Report preliminary results on any mold isolated.
- D. Send out a preliminary report with an organism identification if the culture work-up is completed. Send out final reports when the plate incubation is complete. Report the genus and species of the fungus isolated.
  - 1. Example:     *Candida albicans*  
                  *Aspergillus fumigatus*
- E. Critical values should be called to the patient care unit and/or physician of record. (Refer to "**Critical Values - Division of Microbiology**").

## VII. Reference

- A. Isenberg, HD, et al. (ed.). 1992. Clinical Microbiology Procedure Handbook, Vol 1. p 6.6.1-6.8.1. American Society for Microbiology, Washington, DC.



Figure 1

Mold Isolate Workup







Figure 3

ZYGOMYCETES

Woolly, White/gray, Rapidly Growing Colony with Broad Aseptate or Sparsely Septate Hyphae



OTHER POTENTIAL PATHOGENS: *Conidiobolus coronatus*, *Basidiobolus sp.*, *Syncephalastrum racemosum*, *Cokeromyces recurvatus*, *Saksenaea vasiformis*, *Apophysomyces elegans*

Figure 4

DEMATIACEOUS MOLDS  
Rapidly Growing  
(mature  $\leq$  7 days)






Figure 4A

**RAPIDLY GROWING  
DEMATIACEOUS MOLDS  
WITH SINGLE CELLED CONIDIA**

**Single Celled Conidia**





Figure 4B

**SLOW GROWING  
DEMATIACEOUS MOLDS**

**Slow Growing**



**Figure 5**  
**RAPIDLY GROWING HYALINE MOLDS**




Figure 5A

ASPERGILLUS SPECIES




Figure 5B

**RAPIDLY GROWING HYALINE MOLDS  
WITH CONIDIA BORNE SINGLY OR  
IN CLUSTERS**





Figure 6

DERMATOPHYTES  
SITE: HAIR, SKIN, NAILS

Produces macro and/or microconidia




\*Other hair, skin, nail potential pathogens include: *Bipolaris* (skin), *Exserohium* (skin), *Scopulariopsis* (nail), *Acremonium* (nail), *Fusarium* (skin, nail), *Chrysosporium* (nail)

Figure 7

Dimorphic Fungi

Yeast in Tissue (35°C)  
Mold at 30°C



Document Control

Effective 5/10/2000

Medical Director Approval: Reviewed by Dr. Schappert 3/10/2010.

Microbiology Director Approval: Dr. Ann Robinson 05/10/2000

Microbiology Supervisor Reviews: Jerry Claridge 05/01/2000, 06/14/2001, 04/01/2002, 03/2003, 04/2004, 06/2005, 06/2006, 06/2007, 05/2008, 07/2009, 04/01/2011, 03/2013, Jason Ammons 07/2015

Revisions & Updates: updated procedure to delete use of corn meal for subculturing molds. Deleted references for staining procedures no longer in use. Changed reference laboratory for send out testing to ARUP.