

STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 1 of 6	
Approved by: C. Case, Manager of Diagnostic Services	Signed by: <i>Cheryl Case</i>	

PURPOSE:

This system is used for the identification of non-fastidious, non-enteric Gram negative rods. The strip consists of 20 microtubes containing dehydrated substrates. The reactions are read and inputted into the database for identification.

SAMPLE INFORMATION:

Storage Requirements	Store at 2-8°C
-----------------------------	----------------

REAGENTS and/or MEDIA:

- API 20NE strip and incubation box (bioMerieux Inc, REF 20 050)
- Blood Agar Plate (BAP)
- 5 mL of 0.85% NaCl
- 5 mL API AUX medium (supplied)
- ~5 mL sterile water
- Ferric Chloride Reagent
- James or Kovacs Reagent
- NIT1 and KOH Reagent
- Zn Powder
- Mineral Oil
- 29°C Incubator

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50215API-NEPRO.doc	PRINT DATE: 6 March 2014

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 2 of 6	

SUPPLIES:

- Sterile Pipette
- Densichek

SPECIAL SAFETY PRECAUTIONS:

- Standard precautions should always be followed.

QUALITY CONTROL:

Performed on each shipment/lot # received:

1. Aeromonas hydrophilia ATCC35654
2. Alcaligenes faecalis ATCC35655

QC	Org	NO3	TRP	GLU	ADH	URE	ESC	GEL	PNPG	GLUJ	JARAJ	IMNEI	IMANI	INAGI	IMALI	IGNTI	ICAP	JADI	IMLTI	ICITI	IPACI	OX
1		+	+	+	+	-	+	+	+	+	+	+	+	+	+	+	+	-	+	*	-	+
2		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	+	+	+	+

Generate TQC order via TQC Order Entry – result QC results in TQC.

PROCEDURE INSTRUCTIONS:


Step	Action
Setting Up an API-20E	
1	In your plate log – Order ^NE.

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50215API-NEPRO.doc	PRINT DATE: 6 March 2014

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 3 of 6	

2	Prepare the incubation tray by adding 5 mL of sterile water – Filling up the honeycomb wells.
3	Remove the strip from its packaging and place it in the incubation tray.
4	Write the specimen number on the flap attached to the tray and date.
5	Order and perform an oxidase (^OX) and record the result on the result sheet.
6	Aliquot approximately 3 mLs of 0.85% saline into a plastic test tube.
7	Prepare a 0.5 McFarland suspension of the organism – the culture should be pure and 18-24 hours old.
8	Tilt the API strip and, using a sterile pipette, slowly distribute the bacterial suspension into the tubes. <ul style="list-style-type: none"> • Hold the pipette tip against the top side of the well to minimize bubble formation.
9	Inoculate tests NO3 → PNPG filling only the tubes.
10	Open the API AUX Medium ampule and add approximately 200 µL of the 0.5 McFarland – mix well avoiding the formation of bubbles.
11	Fill the tubes AND cupules from <u>IGLUI</u> → <u>IPACI</u> with the suspension ensuring a slightly convex meniscus.
12	Overlay <u>GLU</u> , <u>ADH</u> , and <u>URE</u> with mineral oil until a convex meniscus is formed.
13	Using one drop of the suspension – streak out a BA Purity Plate.
14	Incubate at 29°C +/- 2°C for 24 hours.
15	Inspect the Purity Plate – if not pure, repeat API using pure culture.
16	The following wells require the addition of reagent: <ol style="list-style-type: none"> 1. NO3: Add 1 drop NIT1 and 1 drop 40% KOH – wait 5 minutes <ul style="list-style-type: none"> • Red colour develops: Positive • No colour develops: add 2-3 mg of Zn to the cupule – wait 5 minutes <ul style="list-style-type: none"> ○ After the addition of Zn – Pink colour develops – test is negative as the Zn reacts with the nitrates left in the tube ○ After the addition of Zn – No Pink develops – test is positive as all of the NO3 has been converted to nitrogen gas and has dissipated into the air

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 4 of 6	

	2. TRP: Add 1 drop of James reagent	
17	Refer to the Reference Table for colour reactions or the supplied package insert. Assimilation Tests: Observe for bacterial growth. An opaque cupule indicates a positive reaction.	
18	Log in to the apiweb: https://apiweb.biomerieux.com	
19	Login name: NSTANTONTERRITORIALHOSPITAL Password: YKNIFE Hit Go .	
20	Select the appropriate API item (ie. APINE).	
21	Input reactions and hit " CONFIRM ".	
22	Print out ID sheet and evaluate the outcome.	
23	Re incubation for an additional 24 hours is necessary if: <ul style="list-style-type: none"> 1. Low discrimination 2. Unacceptable/doubtful profile 3. "Identification not valid before 48 hour incubation" 	
24	If re-incubation is required – remove the NIT1, KOH and James Reagents by suction with a pipette.	
25	Re-incubate at 29°C for an additional 24 hours.	
26	All tests from <u>ADH</u> → <u>IPACI</u> can be re-read at 48 hours.	

EXPECTED RESULTS:

Consult the Identification Table at the end of the package insert for a range of expected results.

REFERENCES:

- bioMerieux. (2006, 02). api 20E.

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50215API-NEPRO.doc	PRINT DATE: 6 March 2014

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 5 of 6	

REFERENCE TABLE:

TESTS	REACTIONS	RESULTS	
		NEGATIVE	POSITIVE
NO ₃	Reduction of nitrates to nitrites	NIT 1 + KOH ~5 minutes	
		Colorless*	Pink-red
	Reduction of nitrates to nitrogen	*Addition of Zn ~5 mins	
		Pink	colorless
TRP	Indole Production	JAMES - IMMEDIATELY	
		Colorless Pale green/Yellow	Pink
GLU	Fermentation	Blue to Green	Yellow
ADH	Arginine DiHydrolase	Yellow	Orange/pink/red
URE	UREase	Yellow	Orange/pink/red
ESC	Hydrolysis(B-glucosidase)	Yellow	Grey/brown/black
GEL	Hydrolysis (protease)	No pigment diffusion	Diffusion of black pigment
PNPG	B-galactosidase	Colorless	Yellow
GLU	Assimilation (GLUcose)	Transparent	Opaque
ARA	Assimilation (ARAbinose)	Transparent	Opaque
MNE	Assimilation (MANNosE)	Transparent	Opaque
MAN	Assimilation (MANnitrol)	Transparent	Opaque
NAG	Assimilation (N-Acetyl-)	Transparent	Opaque
MAL	Assimilation (MALtose)	Transparent	Opaque
GNT	Assimilation(potassium GlucoNate)	Transparent	Opaque
CAP	Assimilation (CAPric Acid)	Transparent	Opaque
ADI	Assimilation (Adipic Acid)	Transparent	Opaque
MLT	Assimilation (MaLaTe)	Transparent	Opaque
CIT	Assimilation (trisodium CITrate)	Transparent	Opaque
PAC	Assimilation (PhenylAcetic acid)	Transparent	Opaque
OX	Cytochrome oxidase	Transparent	Opaque

TITLE: API-NE	Revision Date: 06-March-2016	Issue Date: 06-March-2014
Document Number: MIC50215	Status: Approved	
Distribution: Microbiology Test Manual	Page: 6 of 6	

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Jul13	Initial Release	A. Darrach
1.1	06Mar14	Document control number changed from MTE10215 to MIC50215	C. Russell