



# STANTON TERRITORIAL HEALTH AUTHORITY

## Yellowknife, Northwest Territories

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 1 of 6</b>	
<b>Approved by:</b> Cheryl Case, Manager of Diagnostic Services	<b>Signed by:</b> <i>Cheryl Case</i>	

### PURPOSE:

The system is used for identification of *Neisseria spp*, *Haemophilus spp* and *Moraxella catarrhalis*. It utilizes 10 microtubes containing dehydrated substrates to test for 12 separate enzymatic reactions or sugar fermentations, as well as the presence of a penicillinase. Only organisms belonging to the above genera should be tested.

### SAMPLE INFORMATION:

<b>Storage Requirements</b>	Kit (strip + saline): 2-8°C James/ZymB reagent: 2-8°C in the dark, expire 1 month after opening
-----------------------------	--

### REAGENTS and/or MEDIA:

- apiNH strips (bioMerieux Inc, REF 10 400)
- 0.85% NaCl (provided with kit)
- Incubation box
- ~ 5 mls sterile water
- James reagent
- Zym B reagent
- Mineral oil
- Densichek
- Chocolate Agar Plate

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC50230APINHPRO.doc</b>	<b>PRINT DATE: 6 March 2014</b>

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 2 of 6</b>	

- 35-37C Incubator

**SUPPLIES:**

- Sterile pipette

**SPECIAL SAFETY PRECAUTIONS:**

Standard precautions should always be followed

**QUALITY CONTROL:**

Performed on each shipment/lot# received using the following strain:

- *Neisseria gonorrhoeae* ATCC31426

<u>PEN</u>	<u>GLU</u>	<u>FRU</u>	<u>MAL</u>	<u>SAC</u>	<u>ODC</u>	<u>URE</u>	<u>LIP</u>	<u>PAL</u>	<u>βGAL</u>	<u>ProA</u>	<u>GGT</u>	<u>IND</u>
+	+	-	-	-	-	-	-	-	-	+	-	-

Generate TQC order via TQC Order Entry – result QC results in TQC

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Setting Up and Interpreting apiNH Strips</b>	
<b>1</b>	In the plate log – Order ^NH
<b>2</b>	Prepare the incubation tray by adding 5mL of sterile water – filling up the honeycomb wells
<b>3</b>	Remove the strip from its packaging and place it in the incubation tray
<b>4</b>	Write the Specimen# and date on the flap attached to the tray
<b>5</b>	Using the provided 0.85% NaCl prepare a suspension of the organism in question to the equivalent of 4.0 McFarland - the culture should be pure and 18-24hrs old

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC50230APINHPRO.doc</b>	<b>PRINT DATE: 6 March 2014</b>

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 3 of 6</b>	

<b>6</b>	Tilt the API strip and using a sterile pipette slowly distribute the bacterial suspension into the tubes <ul style="list-style-type: none"> <li>• Hold the pipette tip against the top side of the well to minimize bubble formation</li> </ul>
<b>7</b>	Fill the tube of the first 7 microtubules <u>PEN</u> → <u>URE</u>
<b>8</b>	Fill the tube and cupules of the last 3 microtubules   <u>LIP/ProA</u>  ,   <u>PAL/GGT</u>  ,   <u>BGAL/IND</u>   Avoid the formation of a convex meniscus
<b>9</b>	Cover all underlined tests with mineral oil ( <u>PEN</u> → <u>URE</u> )
<b>10</b>	Cover and Incubate the tray for 2-2.5hrs at 36C in O2
<b>11</b>	Using 1 drop of your suspension, inoculate a Chocolate Agar Plate for purity – Incubate in CO2
<b>12</b>	After the incubation period, read the reactions and record the results on the resulting sheet
<b>13</b>	To the following wells add the reagents:   <u>ProA</u>  ,   <u>GGT</u>   - 1 drop of ZymB   <u>IND</u>   - 1 drop of James Wait 3 minutes then read the reactions See Appendix for reagent reconstitution instructions
<b>14</b>	**If after 2 hrs reactions are doubtful, re-incubate for an additional 2 hours
<b>15</b>	Refer to the attached “ <b>Reference Table</b> ” for color reactions or the supplied package insert
<b>16</b>	Log into the apiweb: <a href="https://apiweb.biomerieux.com">https://apiweb.biomerieux.com</a>
<b>17</b>	Login name: <b>NSTANTONTERRITORIALHOSPITAL</b> Password: <b>YKNIFE</b> Hit <b>Go</b>
<b>18</b>	Select the appropriate API item (ie. APINH)
<b>19</b>	Input reactions and hit “ <b>CONFIRM</b> ”
<b>20</b>	Print out ID sheet and evaluate the outcome
<b>21</b>	Limitations of test: 1. If the ID is <i>N.gonorrhoeae</i> , ensure that the   <u>ProA</u>   is positive – the ID must be confirmed at Dynalife

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC50230APINHPRO.doc</b>	<b>PRINT DATE: 6 March 2014</b>

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 4 of 6</b>	

	2. Certain species of <i>Moraxella</i> may be wrongly identified as <i>N. meningitidis</i> or <i>N. gonorrhoeae</i> – confirm with serological testing at Dynalife
--	--

**INTERPRETATION OF RESULTS:**

TESTS	REACTIONS	RESULTS	
		NEGATIVE	POSITIVE
1) <u>PEN</u>	PENicillinase	Blue (penicillinase absent)	Yellow Yellow-green Yellow-blue (penicillinase present)
2) <u>GLU</u> 3) <u>FRU</u> 4) <u>MAL</u> 5) <u>SAC</u>	GLUcose (Acidification) FRUctose (Acidification) MALtose (Acidification) Saccharose/Sucrose (Acidification)	red red-orange	yellow orange
6) <u>ODC</u>	Ornithine DeCarboxylase	yellow-green grey-green	blue
7) <u>URE</u>	UREase	yellow	pink-violet
8a) <u>LIP</u>	LIPase	colorless pale grey	blue (+ precipitate)
9a) <u>PAL</u>	ALKaline Phosphatase	colorless pale yellow	yellow
10a) <u>BGAL</u>	Beta GALactosidase	colorless	yellow
8b) <u>ProA</u>	Proline Arylamidase if LIP is +	ZYM B / 3 min	

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC50230APINHPRO.doc</b>	<b>PRINT DATE: 6 March 2014</b>

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 5 of 6</b>	

	ProA is always -	yellow pale orange (brown if LIP +)	orange
9b) <u>[GGT]</u>	Gamma Glutamyl Transferase	ZYM B / 3 min	
		yellow pale orange (yellow-orange if PAL +)	orange
10b) <u>[IND]</u>	INDole	JAMES / 3 min	
		Colorless	Pink

### **APPENDIX:**

Type	ZYM B Reagent	James Reagent
<b>Storage Requirements</b>	2-8°C in the dark	
<b>Stability</b>	<ul style="list-style-type: none"> <li>Stable until the expiration date indicated on the packaging.</li> <li>The reagents may be kept for up to 1 month after the ampoules have been opened and the reagents transferred into the dropper bottles.</li> </ul>	
<b>Reagent Preparation</b>	<ol style="list-style-type: none"> <li>Break the ampoule cap.</li> <li>Turn the ampoule upside down and maintain it in a vertical position.</li> <li>Squeeze gently on the cap to release a drop to transfer all the reagent into the dropper bottle.</li> </ol>	<ol style="list-style-type: none"> <li>Break open the ampoule of solvent associated with the James reagent (R1).</li> <li>Take up the contents of the ampoule using a completely dry pipette and transfer this solvent into the dropper bottle (R2).</li> <li>Fit the dropper to the bottle, carefully close it and shake the bottle containing the dehydrated active ingredient.</li> <li>Wait 10 minutes until it is completely dissolved.</li> </ol>

<b>TITLE: API-NH</b>	<b>Revision Date:</b> 06-March-2016	<b>Issue Date:</b> 06-March-2014
<b>Document Number: MIC50230</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology</b>	<b>Page: 6 of 6</b>	

<b>Notes</b>	These reagents are VERY sensitive to light: check the appearance of the reagents before transferring them into the dropper bottles and before each use. After transferring the contents of the ampoules into the dropper-bottles, put the bottles into a used Vitek card foil package (or other aluminum foil).	
	The Zym B reagent is normally yellow to amber in colour. Dispose of the reagent if any tint of pink is observed (sign of deterioration). Exposure to laboratory lighting for a short period of time (approx. 1 hour) will damage the reagent.	The James reagent may only be used if it is pale yellow. If a pink colour appears when the reagent is reconstituted with the solvent, wait until this pink colour has completely disappeared before using the reagent.

**REFERENCES:**

- bioMerieux. (2010, 02). apiNH.

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Dec2013	Initial Release	A.Darrach