

# STANTON TERRITORIAL HEALTH AUTHORITY

## Yellowknife, Northwest Territories

<b>TITLE: Oxidase</b>	<b>Revision Date:</b> 11-March-2016	<b>Issue Date:</b> 11-March-2014
<b>Document Number: MIC51400</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 1 of 4</b>	
<b>Approved by:</b> Cheryl Case, Manager of Diagnostic Services	<b>Signed by:</b> <i>Cheryl Case</i>	

### INTRODUCTION:

This test is useful in differentiating organisms, especially member of genera *Neisseria*, *Pseudomonas*, *Moraxella*, *Vibrio* and *Aeromonas* which are oxidase positive.

### PRINCIPLE:

This test detects the presence of intracellular oxidase enzymes (cytochromes) which play a part in the electron transport system of respiration in some aerobic and facultative bacteria. The reagent is based upon the oxidation of tetramethyl-p-phenylenediamine by bacterial cytochromes in the presence of atmospheric oxygen to form a **purple** coloured compound (Wurster's **blue**).

### SAMPLE INFORMATION:

<b>Type</b>	One well isolated colony
-------------	--------------------------

### REAGENTS and/or MEDIA:

<b>Type</b>	Pro-Lab Test Oxidase Reagent, Cat#PL.390
<b>Storage Requirements</b>	<ul style="list-style-type: none"><li>• Store at controlled room temperature (15°C – 30°C) in the original container.</li><li>• Do not freeze or overheat.</li><li>• Protect from light.</li></ul>

**NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.**

**FILENAME:** MIC51400OxidasePRO.doc

**PRINT DATE:** 11 March 2014

<b>TITLE: Oxidase</b>	<b>Revision Date:</b> 11-March-2016	<b>Issue Date:</b> 11-March-2014
<b>Document Number: MIC51400</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 2 of 4</b>	

<b>Stability</b>	<ul style="list-style-type: none"> <li>Keep the screw cap tightly closed.</li> </ul>
	Product stored under the above conditions will be stable until expiry date shown on the label.

**SUPPLIES:**

- Inoculating loop or sticks
- Filter paper strips or pads
- Incubator
- QC organisms

**SPECIAL SAFETY PRECAUTIONS:**

Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

Quality Control is set up each day the test is performed using the following control organisms:

**Positive:**     *Pseudomonas aeruginosa* ATCC # 27853

**Negative:**    *Klebsiella pneumoniae*     ATCC # 13883

- A TQC order is automatically generated to record the QC results

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME:</b> MIC51400OxidasePRO.doc	<b>PRINT DATE:</b> 11 March 2014

<b>TITLE: Oxidase</b>	<b>Revision Date:</b> 11-March-2016	<b>Issue Date:</b> 11-March-2014
<b>Document Number: MIC51400</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 3 of 4</b>	

### **PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Direct Colony Method</b>	
<b>1</b>	In the plate log – Order <b>^OX</b>
<b>2</b>	Add one drop of reagent to a well-isolated colony on the surface of recommended agar medium.
<b>3</b>	Observe the colony for a colour change within 30 seconds. (If test isolate produces excessively mucoid or slimy colonies, allow up to 1 minute for colour development.)

Step	Action
<b>Filter Paper Method</b>	
<b>1</b>	In the plate log – Order <b>^OX</b>
<b>2</b>	Add 1 to 2 drops of reagent to any convenient size of filter paper. Wait 1 to 2 minutes for proper reagent redistribution.
<b>3</b>	Using a wooden mixing stick or disposable inoculating loop (nichrome wire loops are not recommended), remove a medium size colony from the surface of the recommended agar medium and rub the inoculum onto the reagent-saturated area of the filter paper.
<b>4</b>	Observe the filter paper for colour change within 30 seconds.


### **INTERPRETATION OF RESULTS:**

#### **Positive test:**

The production of a distinct **blue** or **purple** colour.

#### **Negative test:**

The absence of a distinct **blue** or **purple** colour.



<b>TITLE: Oxidase</b>	<b>Revision Date:</b> 11-March-2016	<b>Issue Date:</b> 11-March-2014
<b>Document Number: MIC51400</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 4 of 4</b>	

### **NOTES AND PRECAUTIONS:**

1. Avoid contact with skin, eyes and clothing. Flammable.
2. The reagent should appear colourless, cloudy or very **light tan**. Do not use if the reagent is **purple**.
3. Do not take organisms off of media with dyes or indicators such as MacConkey since it will interfere with the colour reaction.

### **REFERENCES:**

- Test Oxidase Reagent package insert, 2009
- Clinical Microbiology Procedures Handbook Henry D. Isenberg – Editor in Chief 2004

### **REVISION HISTORY:**

<b>REVISION</b>	<b>DATE</b>	<b>Description of Change</b>	<b>REQUESTED BY</b>
1.0	25/01/11	Initial Release	J. Whitson
1.1	31Jul13	Addition of Computer Steps and Illustration	A. Darrach
1.2	11Mar14	Changed from Document control number MTE11400 to MIC51400	C. Russell

**NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.**

**FILENAME:** MIC51400OxidasePRO.doc

**PRINT DATE:** 11 March 2014