

STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

TITLE: Mycobacteria Media Inoculation & Direct Smear Preparation	Revision Date: 07-April-2017	Issue Date: 07-April-2015
Document Number: MIC80800	Status: Approved	
Distribution: Mycobacteria Manual	Page: 1 of 5	
Approved by: Gloria Badari, Director, Corporate Services and Chief Financial Officer	Signed by: (Original Signed Copy in Microbiology)	

PURPOSE:

To standardize the inoculation of Mycobacteria culture media and the preparation of Direct Smears.

SPECIAL SAFETY PRECAUTIONS:

- Handle all patient samples and testing reagent using “Routine Practices”
- Please refer to the Northwest Territories Infection Prevention and Control Manual, March 2012
- Prior to testing all patient are to be identified as per I-0500 Use of Two Patient Identifiers.

SUPPLIES:

- 50 mL conical FALCON tubes
- 5 ml unsterile Pasteur pipettes
- 1 mL graduated Transfer pipettes, individually wrapped, sterile
- Yellow Waste bucket (1/3rd full with fresh Accel TB)
- Vortex mixer
- Alcohol wipes
- Clean frosted glass slides
- 5% phenol alcohol
- Pencil
- Paper towels
- 4x4 gauze
- 50 mL beaker





NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAMEMIC80800 - Mycobacteria Media Inoculation & Direct Smear PreparationPRO.doc

PRINT DATE: 9 March 2015

TITLE: Mycobacteria Media Inoculation & Direct Smear	Revision Date: 07-April-2017	Issue Date: 07-April-2015
Document Number: MIC80800	Status: Approved	
Distribution: Mycobacteria Manual	Page: 2 of 5	

MEDIA & EQUIPMENT:

Criteria	Information
List of supplies, reagents, media and equipment required for Mycobacteria processing:	
LJ slants 	<ul style="list-style-type: none"> • Manufacturer: Remel • Vendor: Unipath/Oxoid • Storage: 2- 8 Celsius
MGIT 	<ul style="list-style-type: none"> • Storage: In-use box kept by Allegra with the lid closed to protect from light (light degrades the fluorescent compound inside tube). • Refer to PANTA document for additional information
MGIT 960 (analyzer) 	<ul style="list-style-type: none"> • Manufacturer: Becton Dickinson • Refer to Equipment document for use and maintenance
BSC (CL-2 protection) 	<ul style="list-style-type: none"> • Manufacturer: LABCONCO • Refer to Equipment document for use and maintenance

MEDIA INOCULATION & DIRECT SMEAR PROCEDURE:

Step	Action
Inoculate in the following order: MGIT tube → LJ → Slide	
Inoculate samples and slide one sample at a time (only one sample is open).	
1	Prepare work area: Bring labeled MGIT tube and LJ slants into BSC. Labeled slides should be sitting in the metal tray in the BSC (prepared before sample processing). Ensure the Accel TB work-surface is still moist with disinfectant from the processing steps. Place empty spare rack on work-surface and use to hold the tubes steady as they are being inoculated. Work on patient samples one-by-one on the moistened work-surface.
2	Place sample concentrate conical tube in 50mL beaker on the work surface.

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC80800 - Mycobacteria Media Inoculation & Direct Smear PreparationPRO.doc

PRINT DATE: 9 March 2015

TITLE: Mycobacteria Media Inoculation & Direct Smear Preparation	Revision Date: 07-April-2017	Issue Date: 07-April-2015
Document Number: MIC80800	Status: Approved	
Distribution: Mycobacteria Manual	Page: 3 of 5	

Obtain prepared MGIT and LJ and place into rack on the work-surface.



Figure 1 (above): Work-surface area; specimen concentrate ready for culture

3	Unwrap a sterile Transfer pipette. The same pipette can be used on MGIT and LJ cultures per one patient sample.
4	Unscrew MGIT tube. Using your pinkie, keep cap in hand so it doesn't touch the surface of the BSC.
5	Using the graduated marks on the pipette withdraw 0.5 mL of patient sample → inoculate MGIT tube. → Place pipette into patient sample tube → Close MGIT cap with free hand.
6	Invert several times to mix. Return MGIT to the plastic rack.
7	Unscrew LJ cap. Using your pinkie, keep cap in hand so it doesn't touch the surface of the BSC.
8	With the same pipette, inoculate LJ slant by slanting the tube and dispensing two drops of patient sample onto slant → close cap LOOSELY. Return LJ slant to the plastic rack.
9	Discard pipette in yellow discard bucket with Accel TB

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC80800 - Mycobacteria Media Inoculation & Direct Smear PreparationPRO.doc

PRINT DATE: 9 March 2015

TITLE: Mycobacteria Media Inoculation & Direct Smear Preparation	Revision Date: 07-April-2017	Issue Date: 07-April-2015
Document Number: MIC80800	Status: Approved	
Distribution: Mycobacteria Manual	Page: 4 of 5	

10	<p>Locate slide in metal tray (leave in metal tray but know its location).</p> <p>With a yellow 10 uL inoculating loop prepare a smear on the pre-labeled patient slide: Dip loop inside conical tube and touch the loop on the center of the labeled slide to dispense the drop on the slide.</p> <p>Thin the drop out with the loop to cover a dime sized area.</p> <p><i>Note: A slide for direct smear does not require a drop of albumin since there is plenty of protein in the patient sample.</i></p>
11	<p>Discard loop in yellow discard bucket with Accel TB.</p> <p>Close cap on the conical centrifuge tube.</p> <p>Move on to next patient sample.</p>
12	Repeat steps 1 → 9 for each sample.
13	<p>Allow smears to completely air dry before fixing.</p> <p>Leave BSC sash open while drying (and fixing) for proper protective airflow.</p> <p>Refer to Mycobacteria Smear Fixing.</p>
14	<p>Wipe or spray sample concentrate tube with Accel TB before removal from BSC.</p> <p>Wipe or spray down inoculated MGIT and LJ's before removal from BSC.</p> <ul style="list-style-type: none"> • Refer to MGIT 960 procedure for insertion into MGIT analyzer • Refer to LJ Culture Procedure for LJ slant incubation <p>Wipe or spray down the "TB Concentrates" bag before removal from BSC.</p> <p>Save the concentrates in the media fridge on bottom shelf in a Labeled Biohazard Bag "TB Concentrates" or similar label for 48 hours.</p> <p>Direct Concentrates may be sent for PCR if Direct Smear positive, at the request of the physician.</p>

RELATED DOCUMENTS:

- MIC 80810 Mycobacteria Smear Fixing
- MIC 81500 Auroamine-Rhodamine Stain
- MIC 81600 Kinyoun Stain

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC80800 - Mycobacteria Media Inoculation & Direct Smear PreparationPRO.doc	PRINT DATE: 9 March 2015

TITLE: Mycobacteria Media Inoculation & Direct Smear Preparation	Revision Date: 07-April-2017	Issue Date: 07-April-2015
Document Number: MIC80800	Status: Approved	
Distribution: Mycobacteria Manual	Page: 5 of 5	

REFERENCES:

- Central Public Health Laboratory. (2003). *Mycobacteriology Bench Manual*. Ottawa.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31-Jan-15	Initial Release	L. Driedger
	03Feb2015	Review	S. Webber

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC80800 - Mycobacteria Media Inoculation & Direct Smear PreparationPRO.doc	PRINT DATE: 9 March 2015