


STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 1 of 7	
Approved by: S. Asmussen, Manager of Diagnostic Services	Signed by: 	

PURPOSE:

The system is used for identification of *Neisseria spp.*, *Haemophilus spp.* and *Moraxella catarrhalis*. It utilizes 10 microtubes containing dehydrated substrates to test for 12 separate enzymatic reactions or sugar fermentations, as well as the presence of a penicillinase. Only organisms belonging to the above genera should be tested.

SAMPLE INFORMATION:

Storage Requirements	Kit (strip + saline): 2-8°C James/ZymB reagent: 2-8°C in the dark, expire 1 month after opening
-----------------------------	--

REAGENTS and/or MEDIA:

- apiNH strips (bioMerieux Inc, REF 10 400)
- 0.85% NaCl (provided with kit)
- Incubation box
- ~ 5 mLs sterile water
- James reagent
- Zym B reagent
- Mineral oil
- Densichek
- Chocolate Agar Plate
- 35-37C Incubator

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC50230API-NHPRO.doc

PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 2 of 7	

SUPPLIES:

- Sterile pipette

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

QUALITY CONTROL:

Performed on each shipment/lot# received using the following strain:

- *Neisseria gonorrhoeae* ATCC31426

<u>PEN</u>	<u>GLU</u>	<u>FRU</u>	<u>MAL</u>	<u>SAC</u>	<u>ODC</u>	<u>URE</u>	<u>LIP</u>	<u>PAL</u>	<u>βGAL</u>	<u>ProA</u>	<u>GGT</u>	<u>IND</u>
+	+	-	-	-	-	-	-	-	-	+	-	-

Generate TQC order via TQC Order Entry – result QC results in TQC

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50230API-NHPRO.doc	PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 3 of 7	

PROCEDURE INSTRUCTIONS:

Step	Action
Setting Up and Interpreting apiNH Strips	
1	In the plate log – Order ^NH
2	Prepare the incubation tray by adding 5mL of sterile water – filling up the honeycomb wells
3	Remove the strip from its packaging and place it in the incubation tray
4	Write the Specimen# and date on the flap attached to the tray
5	Using the provided 0.85% NaCl prepare a suspension of the organism in question to the equivalent of 4.0 McFarland - the culture should be pure and 18-24hrs old
6	Tilt the API strip and using a sterile pipette slowly distribute the bacterial suspension into the tubes <ul style="list-style-type: none"> • Hold the pipette tip against the top side of the well to minimize bubble formation
7	Fill the tube of the first 7 microtubules <u>PEN</u> → <u>URE</u>
8	Fill the tube and cupules of the last 3 microtubules <u>LIP/ProA</u> , <u>PAL/GGT</u> , <u>BGAL/IND</u> Avoid the formation of a convex meniscus
9	Cover all underlined tests with mineral oil (<u>PEN</u> → <u>URE</u>)
10	Cover and Incubate the tray for 2-2.5hrs at 36C in O2
11	Using 1 drop of your suspension, inoculate a Chocolate Agar Plate for purity – Incubate in CO2
12	After the incubation period, read the reactions and record the results on the resulting sheet
13	To the following wells add the reagents: <u>ProA</u> , <u>GGT</u> - 1 drop of ZymB <u>IND</u> - 1 drop of James Wait 3 minutes then read the reactions See Appendix for reagent reconstitution instructions
14	**If after 2 hrs reactions are doubtful, re-incubate for an additional 2 hours
15	Refer to the attached “ Reference Table ” for color reactions or the supplied package

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50230API-NHPRO.doc	PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 4 of 7	

	insert
16	Log into the apiweb: https://apiweb.biomerieux.com
17	Login name: NSTANTONTERRITORIALHOSPITAL Password: YKNIFE Hit Go
18	Select the appropriate API item (ie. APINH)
19	Input reactions and hit " CONFIRM "
20	Print out ID sheet and evaluate the outcome
21	Limitations of test: <ol style="list-style-type: none"> 1. If the ID is <i>N.gonorrhoeae</i>, ensure that the <u>ProA</u> is positive – the ID must be confirmed at Dynalife 2. Certain species of <i>Moraxella</i> may be wrongly identified at <i>N. meningitis</i> or <i>N. gonorrhoeae</i> – confirm with serological testing at Dynalife

INTERPRETATION OF RESULTS:

TESTS	REACTIONS	RESULTS	
		NEGATIVE	POSITIVE
1) <u>PEN</u>	PENicillinase	Blue (penicillinase absent)	Yellow Yellow-green Yellow-blue (penicillinase present)
2) <u>GLU</u>	GLUcose (Acidification)	red red-orange	yellow orange
3) <u>FRU</u>	FRUctose (Acidification)		
4) <u>MAL</u>	MALtose (Acidification)		
5) <u>SAC</u>	MALtose (Acidification)		

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50230API-NHPRO.doc	PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 5 of 7	

	Saccharose/Sucrose (Acidification)		
6) <u>ODC</u>	Ornithine DeCarboxylase	yellow-green grey-green	blue
7) <u>URE</u>	UREase	yellow	pink-violet
8a) <u>LIP</u>	LIPase	colourless pale grey	blue (+ precipitate)
9a) <u>PAL</u>	ALKaline Phosphatase	colourless pale yellow	yellow
10a) <u>BGAL</u>	Beta GALactosidase	colourless	yellow
8b) <u>ProA</u>	Proline Arylamidase if LIP is + ProA is always -	ZYM B / 3 min	
		yellow pale orange (brown if LIP +)	orange
9b) <u>GGT</u>	Gamma Glutamyl Transferase	ZYM B / 3 min	
		yellow pale orange (yellow-orange if PAL +)	orange
10b) <u>IND</u>	INDole	JAMES / 3 min	
		Colourless	Pink

APPENDIX:

Type	ZYM B Reagent	James Reagent
Storage Requirements	2-8°C in the dark	
Stability	<ul style="list-style-type: none"> Stable until the expiration date indicated on the packaging. The reagents may be kept for up to 1 month after the ampoules have been opened and the reagents transferred into the dropper bottles. 	
Reagent Preparation	<ol style="list-style-type: none"> Break the ampoule cap. Turn the ampoule upside 	<ol style="list-style-type: none"> Break open the ampoule of solvent associated with the

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50230API-NHPRO.doc	PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 6 of 7	

	<p>down and maintain it in a vertical position.</p> <p>3. Squeeze gently on the cap to release a drop to transfer the entire reagent into the dropper bottle.</p>	<p>James reagent (R1).</p> <p>2. Take up the contents of the ampoule using a completely dry pipette and transfer this solvent into the dropper bottle (R2).</p> <p>3. Fit the dropper to the bottle, carefully close it and shake the bottle containing the dehydrated active ingredient.</p> <p>4. Wait 10 minutes until it is completely dissolved.</p>
Notes	<p>These reagents are VERY sensitive to light: check the appearance of the reagents before transferring them into the dropper bottles and before each use. After transferring the contents of the ampoules into the dropper-bottles, put the bottles into a used Vitek card foil package (or other aluminum foil).</p>	
	<p>The Zym B reagent is normally yellow to amber in colour. Dispose of the reagent if any tint of pink is observed (sign of deterioration). Exposure to laboratory lighting for a short period of time (approx. 1 hour) will damage the reagent.</p>	<p>The James reagent may only be used if it is pale yellow. If a pink colour appears when the reagent is reconstituted with the solvent, wait until this pink colour has completely disappeared before using the reagent.</p>

REFERENCES:

- bioMerieux. (2010, 02). apiNH.

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50230API-NHPRO.doc	PRINT DATE: 19 April 2016

TITLE: API-NH	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50230	Status: Approved	
Distribution: Microbiology Test Manual	Page: 7 of 7	

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Dec2013	Initial Release	A.Darrach
2.0	31Mar16	Update of "Special Safety Precautions" to reflect risk assessment recommendations.	C. Russell