


STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

TITLE: Germ Tube	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50900	Status: Approved	
Distribution: Microbiology Test Manual	Page: 1 of 5	
Approved by: S. Asmussen, Manager of Diagnostic Services	Signed by: 	

INTRODUCTION:

This test is the most valuable rapid test for presumptive identification of *Candida albicans*. Yeast infections are among the most common fungal infections affecting humans and they are often seen as secondary invaders in immuno-compromised patients and those on chemotherapy.

PURPOSE:

The germ tube production test is used to differentiate *Candida albicans* from other *Candida species*.

PRINCIPLE:

Formation of germ tubes is associated with increased synthesis of protein and ribonucleic acid. The essential requirements for this synthesis are contained in tryptic soy broth and fetal bovine sera. The solution is lyophilized for stability.

SAMPLE INFORMATION:

Type	Well isolated presumptive yeast colonies
-------------	--

REAGENTS and/or MEDIA:

Type	Remel CAT#21068
Storage Requirements	Store lyophilized product in its original container at 2-8°C until used. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.
Stability	Store hydrated aliquots at -20°C for up to 4 months.

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50900GermTubePRO.doc	PRINT DATE: 19 April 2016

TITLE: Germ Tube	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50900	Status: Approved	
Distribution: Microbiology Test Manual	Page: 2 of 5	

SUPPLIES:

- Test tube
- Glass slides and cover slips
- Pasteur pipettes
- 35°C incubator

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

QUALITY CONTROL:

Quality Control is set up each day the test is performed using the following control organisms.

Positive: *Candida albicans* ATCC 10231
Negative: *Candida glabrata* ATCC 15126

- A TQC order is automatically generated to record the QC results

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC50900GermTubePRO.doc	PRINT DATE: 19 April 2016



TITLE: Germ Tube	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50900	Status: Approved	
Distribution: Microbiology Test Manual	Page: 3 of 5	

PROCEDURE INSTRUCTIONS:

Step	Action
Performing a Germ Tube Test	
1	In the plate log – Order ^GT
2	Rehydrate Germ Tube solution with sterile water according to the volume size indicated on the vial.
3	Aliquot 0.5 ml of solution into clean test tubes. Unused aliquots may be capped tightly and frozen at -20°C for up to 4 months.
4	Make a dilute suspension of a single yeast colony by touching the tip of a Pasteur pipette to the colony and emulsifying the cells in the solution. DO NOT inoculate the solution heavily; excessive inoculum causes a significant decrease in the percentage of cells forming germ tubes. Positive and negative controls should be tested simultaneously.
5	Incubate aerobically at 35-37°C for 2-4 hours.
6	Examine microscopically under high (40X) magnification for the presence of germ tubes.

INTERPRETATION OF RESULTS:

Positive test: Germ tubes appearing as short, tube-like structures with no constriction at juncture with yeast cell wall.





NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC50900GermTubePRO.doc

PRINT DATE: 19 April 2016

TITLE: Germ Tube	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50900	Status: Approved	
Distribution: Microbiology Test Manual	Page: 4 of 5	

Negative test: Yeast cells only or with attachment of pseudohyphae.



NOTES AND PRECAUTIONS:

1. Very heavy concentration of inoculum causes a significant decrease in the percentage of cells forming germ tubes.
2. *Candida stellatoidea*, which also produces germ tubes, is no longer a species and has been combined with *C.albicans*.
3. Some oral strains of *Candida tropicalis* have been shown to produce germ tubes.
4. *Candida dubliensis* is germ tube positive; however, careful microscopic observation will reveal that, unlike *C.albicans*, the short hyphal initials are constricted at the junction of the blastoconidium and the germ tube.
5. This test is only part of the overall scheme for identification. Further biochemical testing may be necessary for identification.

REFERENCES:

- Germ Tube Solution package insert, 2008
- Murray Patrick, Baron Ellen Jo, Jorgensen James, Landry Marie Louise, Pfaller Michael, Manual of Clinical Microbiology, 9th edition, 2007, p.335

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC50900GermTubePRO.doc

PRINT DATE: 19 April 2016

TITLE: Germ Tube	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC50900	Status: Approved	
Distribution: Microbiology Test Manual	Page: 5 of 5	

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	14/11/12	Initial Release	A.Darrach
2.0	31Mar16	Update of "Special Safety Precautions" to reflect risk assessment recommendations.	C. Russell

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC50900GermTubePRO.doc

PRINT DATE: 19 April 2016