


# STANTON TERRITORIAL HEALTH AUTHORITY

## Yellowknife, Northwest Territories

<b>TITLE: KOH StringTest</b>	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number: MIC51100</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 1 of 4</b>	
<b>Approved by:</b> S. Asmussen, Manager of Diagnostic Services	<b>Signed by:</b> 	

### **PURPOSE:**

A visible loopful of cells from a single, well-isolated colony is emulsified into a drop of 3% KOH. If the mixture becomes viscous and “strings” within 60 seconds of mixing (KOH-positive) then the colony is considered gram-negative. The reaction depends on the lysis of the gram-negative cell in the dilute alkali solution releasing cellular DNA to turn the suspension viscous. The formation of a string (DNA) in 3% KOH indicates that the isolate is a gram negative organism although it can be non-reactive while testing anaerobes.

### **REAGENTS and/or MEDIA:**

- 40% KOH

### **SUPPLIES:**

- 50mL conical tube
- Sterile blue loop
- Pipette
- Eppendorf

### **SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC51100KOHStringTestPRO.doc</b>	<b>PRINT DATE: 19 April 2016</b>

<b>TITLE: KOH StringTest</b>	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number: MIC51100</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 2 of 4</b>	

- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

**QUALITY CONTROL:**

Performed each time the test is performed:

- Positive control: *P.aeruginosa* ATCC2783
- Negative control: *S.aureus* ATCC25923

A QC order is generated in the TQC system: Resulting Worklist → MICS → 3STR

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Preparing the 3% KOH Solution</b>	
<b>1</b>	Using the equation: <b>C1V1=C2V2</b> Calculate the amount of 40% KOH solution needed for a 3 mL solution with a final concentration of 3% $(40)X=(3)(3)$ $X=9/40$ $X=0.225\text{mLs of KOH needed}$ Water needed: $3\text{mLs}-0.225\text{mLs}=2.8\text{mLs}$ Using an Eppendorf pipettor set to 225µL – pipette 40% KOH into a conical tube
<b>2</b>	Add 2.8 mL of sterile water and mix



<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC51100KOHStringTestPRO.doc</b>	<b>PRINT DATE: 19 April 2016</b>

<b>TITLE: KOH StringTest</b>	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number: MIC51100</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 3 of 4</b>	

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Performing the String Test</b>	
<b>1</b>	Label 3 slides as follows: <ol style="list-style-type: none"> <li><b>1. TEST</b></li> <li><b>2. POS Control</b></li> <li><b>3. NEG Control</b></li> </ol>
<b>2</b>	Add 1 drop of 3% KOH to each slide
<b>3</b>	Emulsify a loopful of organism in the KOH
<b>4</b>	Stir for a maximum of 60 seconds and slowly lift the loop observing for string formation

**EXPECTED RESULTS:**

<b>Positive</b>	Formation of a string occurs 
<b>Negative</b>	No string formation occurs 

**PRECAUTIONARY NOTES:**

- False negatives can occur when using too light of an inoculums and/or testing anaerobes
- False positives can occur when testing mucoid organisms and using too heaving of an inoculum

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME: MIC51100KOHStringTestPRO.doc</b>	<b>PRINT DATE: 19 April 2016</b>

<b>TITLE: KOH StringTest</b>	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number: MIC51100</b>	<b>Status: <span style="color: red;">Approved</span></b>	
<b>Distribution: Microbiology Test Manual</b>	<b>Page: 4 of 4</b>	

**REFERENCES:**

- Mount Sinai Microbiology Manual. (n.d.). Retrieved October 23, 2013, from <http://microbiology.mtsinai.on.ca/manual/tech/tech22.pdf>
- Scott Sutton, P. (n.d.). *The Gram Stain*. Retrieved October 23, 2013, from The Microbiology Network: <http://www.microbiol.org/resources/monographswhite-papers/the-gram-stain/>

**REVISION HISTORY:**

<b>REVISION</b>	<b>DATE</b>	<b>Description of Change</b>	<b>REQUESTED BY</b>
1.0	31Dec2013	Initial Release	A.Darrach
2.0	31Mar2016	Update of "Special Safety Precautions" to reflect risk assessment recommendations.	C. Russell