


STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

TITLE: Serology – Salmonella Serology	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC51915	Status: Approved	
Distribution: Microbiology Test Manual	Page: 1 of 5	
Approved by: S. Asmussen, Manager of Diagnostic Services	Signed by: 	

PURPOSE:

The genus *Salmonella* contains a wide variety of human pathogens. Identification requires culture isolation, biochemical characterizations and serological identifications. Serological identification involves mixing the suspected organism with antiserum containing specific *Salmonella* antibodies and observing for agglutination. Polyvalent “O /Vi” antisera are intended to aid in initial serogrouping. Factor Vi antisera is intended for further identification of *Salmonella typhi*

SAMPLE INFORMATION:

Type	24 hour old culture on Blood Agar; NLF, TSI: K/A, VITEK ID: <i>Salmonella sp</i>
Source	Any clinical sample

REAGENT INFORMATION:

Source	Pro-Lab Diagnostics REF: PL6000 and PL.6040
Stability	Stable until expiry date
Storage Requirements	2-8°C

SUPPLIES:

- Glass slides
- Sterile loop
- Magnifying lamp
- Normal saline

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC51915SalmonellaSerologyPRO.doc | **PRINT DATE:** 19 April 2016

TITLE: Serology – Salmonella Serology	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC51915	Status: Approved	
Distribution: Microbiology Test Manual	Page: 2 of 5	

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

QUALITY CONTROL:

A QC order is generated in the TQC system:

- Resulting Worklist→MICS→POLYO

A saline control is run with every test, see following Procedure

A positive control is run with every **Poly O** test


- *Salmonella enterica* ATCC14028

PROCEDURE INSTRUCTIONS:

Step	Action
Performing the Poly O Test	
1	Performed on isolates previously identified biochemically and by VITEK as <i>Salmonella</i> sp Label the frosted end of a glass slide with: TEST POLYO

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC51915SalmonellaSerologyPRO.doc	PRINT DATE: 19 April 2016

TITLE: Serology – Salmonella Serology	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC51915	Status: Approved	
Distribution: Microbiology Test Manual	Page: 3 of 5	


2	Using a wax pencil make two square boxes on the glass slide <div style="text-align: center;">  </div>
3	Add one drop of sterile saline to each square
4	Using a sterile loop, emulsify half of one suspect <i>Salmonella</i> colony from a 24hr culture into one drop of saline. Repeat for the other drop Mix to get a smooth suspension
5	Add one drop of the “Poly O and Vi” anti-serum to one box Add one drop of sterile saline to the other box to use as a control
6	Mix the antisera into the suspension with a sterile loop
7	Using the magnifying lens, gently rock the slide for one minute to observe for agglutination
IF	
POLYO and Vi: Agglutination SALINE CONTROL: No Agglutination	<ul style="list-style-type: none"> • Test is Positive. • Probable <i>Salmonella sp</i>. • Proceed to Step 8 for FactorVi testing
POLYO and Vi: No Agglutination SALINE CONTROL: No Agglutination	<ul style="list-style-type: none"> • Test is Negative. • <i>Salmonella sp</i> NOT isolated
POLYO and Vi: Agglutination SALINE CONTROL: Agglutination	<ul style="list-style-type: none"> • Organism is auto-agglutinating. • Reaction is not specific. Test is invalid

PROCEDURE INSTRUCTIONS:

Step	Action
Performing the Vi Test	
1	Label the frosted end of a glass slide with: TEST POLYVi

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC51915SalmonellaSerologyPRO.doc	PRINT DATE: 19 April 2016

TITLE: Serology – Salmonella Serology	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC51915	Status: Approved	
Distribution: Microbiology Test Manual	Page: 4 of 5	

2	<p>Using a wax pencil make two square boxes on the glass slide</p> <div style="text-align: center;">  </div>		
3	Add one drop of sterile saline to each square		
4	<p>Using a sterile loop, emulsify half of one suspect <i>Salmonella</i> colony from a 24hr culture into one drop of saline. Repeat for the other drop</p> <p>Mix to get a smooth suspension</p>		
5	<p>Add one drop of the Factor Vi anti-serum to one box</p> <p>Add one drop of sterile saline to the other box to use as a control</p>		
6	Mix the antisera into the suspension with a sterile loop		
7	Using the magnifying lens, gently rock the slide for one minute to observe for agglutination		
IF		THEN	
<p>Factor Vi: Agglutination SALINE CONTROL: No Agglutination</p>		<ul style="list-style-type: none"> Test is Positive. Probable <i>Salmonella typhi</i> Send to ProvLab Calgary for confirmation See Procedure MIC52815 Microbiology Isolate Referrals 	
<p>Factor Vi: No Agglutination SALINE CONTROL: No Agglutination</p>		<ul style="list-style-type: none"> Test is Negative for <i>Salmonella typhi</i> Send to ProvLab Calgary for speciation See Procedure MIC52815 Microbiology Isolate Referrals 	
<p>Factor Vi: Agglutination SALINE CONTROL: Agglutination</p>		<ul style="list-style-type: none"> Organism is auto-agglutinating. Reaction is not specific. Test is invalid. Send to ProvLab Calgary for speciation See Procedure MIC52815 Microbiology Isolate Referrals 	

TITLE: Serology – Salmonella Serology	Revision Date: 20-April-2018	Issue Date: 20-April-2016
Document Number: MIC51915	Status: Approved	
Distribution: Microbiology Test Manual	Page: 5 of 5	

LIMITATIONS OF PROCEDURE:

1. Use on cultures previously characterized biochemically as *Salmonella*, as other organisms have similar antigens and can give false positives
2. Sensitivity of the slide may be reduced if volumes of greater than 10µL are used

REFERENCES:

- Pro-Lab Diagnostics. (2008, 06). Salmonella Antisera.

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Dec2013	Initial Release	A.Darrach
2.0	31 Mar 2016	Update of "Special Safety Precautions" to reflect risk assessment recommendations.	C. Russell

NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

FILENAME: MIC51915SalmonellaSerologyPRO.doc | **PRINT DATE:** 19 April 2016