


# STANTON TERRITORIAL HEALTH AUTHORITY

## Yellowknife, Northwest Territories

<b>TITLE:</b> Water Testing – HPC Unit Dose SimPlate Method	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number:</b> MIC52615	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Microbiology Test Manual	<b>Page:</b> 1 of 5	
<b>Approved by:</b> S. Asmussen, Manager of Diagnostic Services	<b>Signed by:</b> 	

### PURPOSE:

The presence of heterotrophic bacteria in dialysate waters can lead to the development of a Gram negative toxin mediated pyrogenic reaction, bacteremia, and chronic inflammatory response syndrome. Hot tub waters are screened for the presence of *Pseudomonas aeruginosa*, the causative agent of a superficial skin infection known as “Hot tub folliculitis”. The HPC SimPlate method utilizes enzyme technology to target the most common enzymes of waterborne bacteria. The byproduct of the enzymatic reactions can be seen as fluorescence using a UV light. Additionally, hot tub waters have a MacConkey plate inoculated to aid in the detection of *Pseudomonas aeruginosa*.

### SAMPLE INFORMATION:

<b>Type</b>	Hot tub water, Dialysate water
<b>Volume</b>	10mL +/- 0.2mL
<b>Stability</b>	48hrs, refrigerated
<b>Storage Requirements</b>	2-8°C (refrigerated)
<b>Criteria for rejection and follow up action</b>	Waters: <10mL and/or >48hrs old

### REAGENTS and/or MEDIA:

- Sterile SimPlate with lid
- Sterile media tubes
- MacConkey agar (for Hot Tub waters only)

<b>NOTE:</b> This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
<b>FILENAME:</b> MIC52615HPCPRO.doc	<b>PRINT DATE:</b> 19 April 2016

<b>TITLE:</b> Water Testing – HPC Unit Dose SimPlate Method	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number:</b> MIC52615	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Microbiology Test Manual	<b>Page:</b> 2 of 5	

**SUPPLIES:**

- 10mL syringe
- UV light
- 35C incubator

**QUALITY CONTROL:**

Performed per shipment or Lot number – refer to Procedure **MIC60400WaterQC**

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Performing HPC Water Testing</b>	
<b>1</b>	Accession waters under the appropriate name – see LIS Water list. The LIS test code should already populate as waters are set up as Standing Orders. If this is a new water client – see Tech2  <b>LIS CODE: TSWAT, DIAST or DIAFS(Fort Smith only)</b>
<b>2</b>	Select Source(where required) – if Source is <b>HOT TUB</b> or <b>SPA WATER</b> – HPC will automatically order and a MacConkey (MAC) plate label will generate
<b>3</b>	If water is chlorinated – treat with sodium thiosulfate prior to analysis. Sample containers from Stanton contain sodium thiosulfate
<b>4</b>	Label requisition and samples
<b>5</b>	Label sterile media tube (green top), SimPlate and MAC plate (where required)
<b>6</b>	Mix water sample
<b>7</b>	Using 10mL syringe – aliquot off 10mL
<b>8</b>	Add to green top tube and shake
<b>9</b>	Allow powder in tube to dissolve
<b>10</b>	Remove SimPlate lid and pour in contents of green top tube
<b>11</b>	Replace lid and gently swirl to distribute the sample  <b>NOTE:</b> Air bubbles do not interfere with test

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME:</b> MIC52615HPCPRO.doc	<b>PRINT DATE:</b> 19 April 2016

<b>TITLE:</b> Water Testing – HPC Unit Dose SimPlate Method	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number:</b> MIC52615	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Microbiology Test Manual	<b>Page:</b> 3 of 5	

12	Tip the plate at a 90 degree angle so the excess water will drain into the absorbent at the bottom
13	Invert the plate onto the plastic lid and incubate at <b>35°C for 48hrs.</b> Results can be read from 45-72 hours from the start of incubation
14	While viewing under the <b>UV lamp</b> - count the number of wells showing fluorescence
15	Refer to the MPN table to determine the Most Probable Number of heterotrophic plate count bacteria in the original sample – see attached chart. THE LIS SYSTEM CALCULATES THIS FOR YOU – CHOOSE THE NUMBER OF POSITIVE WELLS IN THE KEYPAD
16	<b>HOT TUB/SPA WATERS</b> – Aliquot off an additional 1mL of water and flood a labeled MAC plate - incubate along with the SimPlate (35°C for 48hrs) <ul style="list-style-type: none"> <li>• Look for growth of <i>Pseudomonas aeruginosa</i> – Non-lactose fermenter, oxidase positive – send to VITEK for GNI</li> <li>• Report as: “Isolated: <i>Pseudomonas aeruginosa</i>”</li> <li>• If MAC plate is negative- Report as: “No <i>Pseudomonas aeruginosa</i> isolated”</li> </ul>

**RELATED DOCUMENTS:**

- MIC60400WaterQC Procedures

**REFERENCES:**

- Garcia, L. S. (2007). Bacterial Endotoxin Test (Limulus Amebocyte Lysate). In *Clinical Microbiology Procedures Handbook, volume 3* (p. 13.7.4).
- *SimPlate for HPC Water Testing Product Information*. (n.d.). Retrieved from [www.idexx.ca/view/xhtml/en\\_ca/water/simplate.jsf](http://www.idexx.ca/view/xhtml/en_ca/water/simplate.jsf)

<b>NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.</b>	
<b>FILENAME:</b> MIC52615HPCPRO.doc	<b>PRINT DATE:</b> 19 April 2016

<b>TITLE:</b> Water Testing – HPC Unit Dose SimPlate Method	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number:</b> MIC52615	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Microbiology Test Manual	<b>Page:</b> 4 of 5	

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	31Dec2013	Initial Release	A.Darrach
1.1	31March2016	Reviewed – No Changes	C. Russell

<b>TITLE:</b> Water Testing – HPC Unit Dose SimPlate Method	<b>Revision Date:</b> 20-April-2018	<b>Issue Date:</b> 20-April-2016
<b>Document Number:</b> MIC52615	<b>Status:</b> <b>Approved</b>	
<b>Distribution:</b> Microbiology Test Manual	<b>Page:</b> 5 of 5	

**APPENDIX:**

**Unit-Dose  
SimPlate For HPC  
Most Probable Number (MPN) Table**

# Positive Wells	MPN	95% confidence limits	
		lower	upper
0	<0.2	<0.03	<1.4
1	0.2	0.0	1.4
2	0.4	0.1	1.6
3	0.6	0.2	1.9
4	0.8	0.3	2.2
5	1.0	0.4	2.5
6	1.2	0.6	2.7
7	1.5	0.7	3.0
8	1.7	0.8	3.3
9	1.9	1.0	3.6
10	2.1	1.1	3.9
11	2.3	1.3	4.2
12	2.6	1.5	4.5
13	2.8	1.6	4.8
14	3.0	1.8	5.1
15	3.3	2.0	5.4
16	3.5	2.2	5.8
17	3.8	2.3	6.1
18	4.0	2.5	6.4
19	4.3	2.7	6.7
20	4.5	2.9	7.0
21	4.8	3.1	7.4
22	5.1	3.3	7.7
23	5.3	3.5	8.0
24	5.6	3.8	8.4
25	5.9	4.0	8.7
26	6.2	4.2	9.1
27	6.5	4.4	9.4
28	6.8	4.7	9.8
29	7.1	4.9	10.2
30	7.4	5.1	10.6
31	7.7	5.4	10.9
32	8.0	5.6	11.3
33	8.3	5.9	11.7
34	8.6	6.2	12.1
35	9.0	6.4	12.6
36	9.3	6.7	13.0
37	9.7	7.0	13.4
38	10.0	7.3	13.9
39	10.4	7.6	14.3
40	10.8	7.9	14.8
41	11.2	8.2	15.2
42	11.6	8.5	15.7

# Positive Wells	MPN	95% confidence limits	
		lower	upper
43	12.0	8.8	16.2
44	12.4	9.1	16.7
45	12.8	9.5	17.3
46	13.2	9.8	17.8
47	13.7	10.2	18.3
48	14.1	10.6	18.9
49	14.6	10.9	19.5
50	15.1	11.3	20.1
51	15.6	11.7	20.7
52	16.1	12.1	21.3
53	16.6	12.5	22.0
54	17.1	13.0	22.7
55	17.7	13.4	23.4
56	18.3	13.9	24.1
57	18.9	14.4	24.9
58	19.5	14.9	25.7
59	20.2	15.4	26.5
60	20.9	15.9	27.3
61	21.6	16.5	28.2
62	22.3	17.1	29.2
63	23.1	17.7	30.2
64	23.9	18.3	31.2
65	24.8	19.0	32.3
66	25.7	19.7	33.5
67	26.6	20.4	34.7
68	27.6	21.2	36.1
69	28.7	22.0	37.5
70	29.9	22.9	39.0
71	31.1	23.8	40.7
72	32.4	24.8	42.5
73	33.9	25.8	44.4
74	35.5	27.0	46.6
75	37.2	28.2	49.1
76	39.2	29.6	51.9
77	41.4	31.1	55.1
78	44.0	32.8	58.9
79	47.0	34.8	63.6
80	50.7	37.1	69.5
81	55.5	39.8	77.5
82	62.3	43.2	89.9
83	73.8	47.6	114.6
84	>73.8	>47.6	>114.6

MPN is per ml of the 10 ml sample added to the media tube (pour-off is accounted for).

**NOTE:** This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.