


STANTON TERRITORIAL HEALTH AUTHORITY

Yellowknife, Northwest Territories

| | | |
|--|---|-------------------------------------|
| TITLE: Bile Solubility | Revision Date: 20-April-2018 | Issue Date: 20-April-2016 |
| Document Number: MIC52900 | Status: Approved | |
| Distribution: Microbiology Test Manual | Page: 1 of 5 | |
| Approved by: S. Asmussen, Manager of Diagnostic Services | Signed by:  | |

PURPOSE:

Test is used for the presumptive identification and differentiation of *Streptococcus pneumoniae* from other alpha hemolytic streptococci.

INTRODUCTION:

Streptococcus pneumoniae possess an autolytic intracellular enzyme that lyses the organism's own cell wall during cell division. The bile solubility test is based on the observation that *S.pneumoniae* cells visibly lyse when 10% sodium desoxycholate is applied, while other alpha-hemolytic streptococci do not. Two methods can be utilized: tube test or the direct plate method.

SAMPLE INFORMATION:

| | |
|--|---|
| Type | Catalase Negative, alpha hemolytic Gram positive cocci |
| Source | Sterile sites |
| Stability | See expiration date |
| Storage Requirements | 20-25°C, RT |
| Criteria for rejection and follow up action | Visible signs of deterioration, color change, expired reagent |

REAGENTS:

- Bile Solubility Reagent: remel REF:R21206

| | |
|---|----------------------------------|
| NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use. | |
| FILENAME: MIC52900BileSolubilityPRO.doc | PRINT DATE: 19 April 2016 |

| | | |
|---|--|-------------------------------------|
| TITLE: Bile Solubility | Revision Date: 20-April-2018 | Issue Date: 20-April-2016 |
| Document Number: MIC52900 | Status: Approved | |
| Distribution: Microbiology Test Manual | Page: 2 of 5 | |

SUPPLIES:

- test tubes
- 37°C incubator
- Normal saline
- Pipette
- Sterile wooden stick

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

QUALITY CONTROL:

Test each new lot of reagent with positive and negative control:

- Positive control: *Streptococcus pneumoniae* ATCC49619
- Negative control: *Streptococcus salivarius* ATCC13419

| | |
|---|----------------------------------|
| NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use. | |
| FILENAME: MIC52900BileSolubilityPRO.doc | PRINT DATE: 19 April 2016 |

| | | |
|---|--|-------------------------------------|
| TITLE: Bile Solubility | Revision Date: 20-April-2018 | Issue Date: 20-April-2016 |
| Document Number: MIC52900 | Status: Approved | |
| Distribution: Microbiology Test Manual | Page: 3 of 5 | |

PROCEDURE INSTRUCTIONS: TUBE TEST

| Step | Action |
|--------------------------------------|---|
| Tube Test for Bile Solubility | |
| 1 | Prepare a 1.0mL saline suspension equivalent to a 1.0 McFarland of the test isolate from an 18-24hr culture into a glass test tube – label tube with the accession number and CONTROL. Example: 1234567CTRL |
| 2 | Aliquot 0.5mL of suspension from Step 1 into a second glass test tube – label with the accession number and TEST. Example: 1234567TEST |
| 3 | To the tube labeled CONTROL – add 0.5mL of sterile saline and gently mix |
| 4 | To the tube labeled TEST – add 0.5mL of 10% Bile Solubility reagent and gently mix |
| 5 | Incubate tubes at 37°C – examine hourly for clearing in the TEST suspension up to a maximum of 3hrs |
| 6 | The control suspension should remain turbid |


PROCEDURE INSTRUCTIONS: DIRECT PLATE METHOD



| Step | Action |
|--|---|
| Direct Plate Method for Bile Solubility | |
| 1 | Add 1 drop of 10% bile solubility reagent to a well-isolated, 18-24hr colony of the test isolate. To achieve well-isolated colonies, a sub-culture may have to be done. |
| 2 | Incubate plate aerobically - agar down and always level, in a 37°C incubator for 30mins. Plate MUST be kept level to prevent the possibility of non-pneumococcal colonies from washing away with the reagent (colonies carried away with reagent is not the same as auto-lysing). |
| 3 | Place the lid slightly ajar to allow for evaporation |
| 4 | Observe the colony for disintegration or solubility |

| | |
|---|----------------------------------|
| NOTE: This is a CONTROLLED document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use. | |
| FILENAME: MIC52900BileSolubilityPRO.doc | PRINT DATE: 19 April 2016 |


| | | |
|---|--|-------------------------------------|
| TITLE: Bile Solubility | Revision Date: 20-April-2018 | Issue Date: 20-April-2016 |
| Document Number: MIC52900 | Status: Approved | |
| Distribution: Microbiology Test Manual | Page: 4 of 5 | |

INTERPRETATION OF RESULTS:

| TUBE METHOD | |
|---|--|
| IF | THEN |
| TEST suspension: clear CONTROL suspension: turbid | Test is POSITIVE ; isolate is presumptive <i>S.pneumoniae</i> |
| TEST suspension: turbid CONTROL suspension: turbid | Test is NEGATIVE ; isolate is not <i>S.pneumoniae</i> |
| TEST suspension: turbid CONTROL suspension: clear | Test is invalid ; repeat testing – probable mislabeling of test tubes |
|  | |

| DIRECT PLATE METHOD | |
|---|--|
| IF | THEN |
| Colony has disintegrated | Test is POSITIVE ; isolate is presumptive <i>S.pneumoniae</i> |
| Colony remains intact | Test is NEGATIVE ; isolate is not <i>S.pneumoniae</i> |
| Before bile salts: | After bile salts: Colony has disintegrated – Test is POSITIVE |
|  |  |

| | | |
|---|--|-------------------------------------|
| TITLE: Bile Solubility | Revision Date: 20-April-2018 | Issue Date: 20-April-2016 |
| Document Number: MIC52900 | Status: Approved | |
| Distribution: Microbiology Test Manual | Page: 5 of 5 | |



REFERENCES:

- Remel. (2012, July). Bile Solubility Reagent (2% and 10% Sodium Desoxycholate).

REVISION HISTORY:

| REVISION | DATE | Description of Change | REQUESTED BY |
|----------|---------|--|--------------|
| 1.0 | 30Dec13 | Initial Release | Darrach (A) |
| 2.0 | 31Mar16 | Update of “Special Safety Precautions” to reflect risk assessment recommendations. | C. Russell |
| | | | |
| | | | |
| | | | |
| | | | |