

Document Name:

ALA (Porphyrin Production) Test

Approved By: JGD Bernier, A/Manager, Laboratory Services

Status: **APPROVED**

PURPOSE:

The ALA test is a rapid test used to determine the growth requirement for hemin (X factor) in the identification of *Haemophilus* species.

SAMPLE INFORMATION:

Type	Tiny Gram-negative rods or coccobacilli growing only on chocolate agar with the typical <i>Haemophilus</i> colonial morphology and which do not grow on BAP
Source	18-24h culture

REAGENTS and/or MEDIA:

Type	A.L.A. Reagent Disk
Source	Remel
Volume	1 disk
Stability	Stable until date of expiration indicated on the container
Storage Requirements	Store at 2-8°C Protect disks from moisture. Protect from light, as the substrate is highly light sensitive.
Criteria for rejection and follow up action	Do not use if: <ul style="list-style-type: none"> • The disk color has changed from white • The expiration date has passed • The desiccant has changed from blue to pink • There are other signs of deterioration

PLEASE NOTE:

- Any positive result from a sterile site for *Haemophilus influenza* must be sent immediately to the Provincial Lab Edmonton for typing as soon as identification is confirmed. Assure there is a purity plate made that can be used for this purpose and can be sent out the day the identification is confirmed. Refer to MIC10510.
- Any positive result from a sterile site for *Haemophilus influenza* must also be sent to NML for International Circumpolar Surveillance (ICS) program. Refer to MIC10520.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

Document Name: ALA (Porphyrin Production) Test	Document Number: MIC53300	
	Version No: 2.0	Page: 2 of 5
	Effective: 13 July, 2017	

SUPPLIES:

- Wooden applicator sticks
- Wood's lamp
- O₂ incubator

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potential infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used when there is a known or potential risk of exposure of splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes and other sharp objects should be strictly limited.

QUALITY CONTROL:

Quality control is set up each day the test is performed using the following control organisms:

- Positive Control: *Aggregatibacter aphrophilus* ATCC 7901
 Negative Control: *Haemophilus influenzae* ATCC 10211

A TQC order is automatically generated to record the QC results.

PROCEDURE INSTRUCTIONS:

Step	Action	
Performing an ALA test		
1	In the plate log - order ^ALA.	
2	Prior to inoculation, allow product to equilibrate to room temperature.	
3	Place ALA disk, with "A" side down, on the agar surface of the culture plate.	
4	Inoculate the disk with a heavy, visible inoculum removed from a pure, 18-24h culture of the test isolate.	
5	Incubate for up to 6 hours at 35°C in ambient air incubator.	
6	Examine disk at 1h under the ultraviolet light for reddish-orange fluorescence. If negative, re-incubate the test and examine periodically for up to 6h before reporting as negative.	
7	If:	Then:
	Reddish-orange fluorescence production	Positive
	No fluorescence production	Negative

INTERPRETATION OF RESULTS:

Result	Interpretation
Positive	The presence of red fluorescence indicating that the organism does not require X factor or hemin and that the ALA has been utilized NOT indicative of <i>Haemophilus influenzae</i>
Negative	The lack of red fluorescence indicating that the organism requires X factor and that the ALA has not been utilized INDICATIVE of <i>Haemophilus influenzae</i>

PLEASE NOTE:

- Any positive result from a sterile site for *Haemophilus influenzae* must be sent immediately to the Provincial Lab Edmonton for typing as soon as identification is confirmed. Assure there is a purity plate made that can be used for this purpose and can be sent out the day the identification is confirmed. Refer to MIC10510.
- Any positive result from a sterile site for *Haemophilus influenzae* must also be sent to NML for International Circumpolar Surveillance (ICS) program. Refer to MIC10520.

Document Name: ALA (Porphyrin Production) Test	Document Number: MIC53300	
	Version No: 2.0	Page: 4 of 5
	Effective: 13 July, 2017	

PROCEDURAL NOTES:

- Use for differentiating *Haemophilus species* only.
- Best results are obtained using a heavy inoculum.
- Examine for fluorescence in a darkened room.
- False negative reactions may occur if the inoculum is insufficient or if the culture is greater than 24hr old. Cultures being tested must not be older than 24hr.
- *Aggregatibacter aphrophilus* was previously known as either *Haemophilus paraphrophilus* or *Haemophilus aphrophilus*.

LIMITATIONS:

- Many organisms will give a positive reaction. If test is performed only on Gram-negative coccobacilli colonies that growth well on CHOC in 24h and not on BAP, results are for *Haemophilus* species.
- The ALA test will not separate *Haemophilus influenzae* from *Haemophilus haemolyticus*. The latter is rare and not pathogenic. It will sometimes grow on BAP without a “staph streak” if it is able to hemolyze the blood to supply it with V factor.
- Organisms that are strongly oxidase positive or catalase positive may give a false positive test. Such organisms make heme and its precursors from ALA in the process of synthesizing oxidase or catalase. Verify that the test organism resembles *Haemophilus* both by gram stain and colonial morphology before testing.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

REFERENCES:

- Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016
- Remel ALA Disk™ package insert, revised July 26, 2010

REVISION HISTORY:

REVISION	DATE	Description of Change	BY
1.0	28-Oct-2016	Initial Release – Test Implementation	L. Steven
2.0	13-Jul-2017	Updated to include note regarding the referral of H. influenza to Provincial Laboratory and NML	L. Steven

 Invalid signature

X 

Jennifer G. Daley Bernier, R.T. (CSMLS)
A/Manager, Laboratory Services
Signed by: Jennifer G. Daley Bernier