



Stanton Territorial Hospital
P.O. Box 10, 550 Byrne Road
YELLOWKNIFE NT X1A 2N1

Document Number: MIC31900

Version No: 1.0

Page: 1 of 14

Distribution:
Microbiology Culture Manual

Effective: 23 June, 2017

Date Reviewed: 23 June, 2017

Next Review: 23 June, 2019

Document Name:
Stool Culture

Approved By:
Jennifer G. Daley Bernier, A/ Manager, Laboratory Services

Status: APPROVED

PURPOSE:

To culture stool specimens for the detection of major enteric pathogens that commonly cause diarrheal illness: *Aeromonas spp.*, *Salmonella spp.*, *Shigella spp.*, *Campylobacter spp.*, *Escherichia coli 0157* and *Yersinia enterocolitica*.

Specific examination for other potential causes of gastroenteritis, such as *Plesiomonas shigelloides* and *Vibrio spp.* can be performed if specifically requested by the ordering physician and a history of travel or consumption of seafood is provided or in the investigation of an outbreak where standard procedures have failed to find the cause.

SAMPLE INFORMATION:

| | |
|--|---|
| Type | <ul style="list-style-type: none"> • Stool collected in enteric transport medium. • Stool in sterile container, if received within 2 hours of collection. • Rectal swab if feces cannot be obtained. |
| Source | Feces |
| Stability | Specimens should be received and processed with 24 hours of collection. |
| Storage Requirements | <ul style="list-style-type: none"> • Fresh stool, not in transport medium, needs to be processed within 2 hours after collection. • Fresh stool in transport medium may be refrigerated at 4°C and transported to the laboratory within 48 hours of collection. |
| Criteria for rejection and follow up action | <ol style="list-style-type: none"> 1. Unlabeled/mislabeled specimen. 2. Not in enteric transport media and more than 2 hours old. 3. If specimen in transport medium is delayed for more than 48 hours at 4°C or is delayed more than 24 hours at 25°C. 4. Fecal cultures received from adults and pediatric patients (3 years of age or older) hospitalized for more than 3 days, unless patient is known to be HIV positive or there is a cluster epidemic within the institution. 5. Do not reject stool samples from infants and toddlers until after the fourth day of hospitalization, since it may take longer to collect a stool sample from pediatric patients admitted with gastroenteritis. 6. Stool with barium. 7. Specimens submitted in Ova and Parasite collection containers. |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| | | |
|-------------------------------------|----------------------------------|----------------------|
| Document Name: Stool Culture | Document Number: MIC31900 | |
| | Version No: 1.0 | Page: 2 of 14 |
| | Effective: 23 June, 2017 | |

REAGENTS and/or MEDIA:

- Blood Agar (BAP), MacConkey Agar (MAC), Hektoen Agar (HEK), MacConkey Sorbitol Agar (MAC-S), Campylobacter Agar (CAMP), Yersinia Select Agar (CIN) and Selenite Broth
- Identification reagents: catalase, oxidase, Salmonella Vi and O antisera, E.coli 0157 serology, hippurate discs, TSI/Urea slants, etc.

SUPPLIES:

- Sterile pipettes and wooden sticks
- Glass test tubes
- Biosafety cabinet
- 35° O₂ incubator
- Vitek 2 Compact and supplies
- Campylobacter microaerophilic jar and pouch

SPECIAL SAFETY PRECAUTIONS:

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.


All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

QUALITY CONTROL:

- Refer to MIC60100 Non-Exempt Media Quality Control procedure
- Refer to Quality Control manual for reagent quality control procedures

| | |
|---|-----------------------------------|
| NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use. | |
| FILENAME: MIC31900StoolCulturePRO.docx | Print Date: 6/23/2017 11:59:00 AM |

PROCEDURE INSTRUCTIONS:

| Step | Action |
|---|---|
| Processing Feces for Stool Culture | |
| 1 | In the biosafety cabinet, mix specimen and use cotton-tipped swab to inoculate each plate. Lightly inoculate BA, MAC and MAC-S. Use a heavier inoculum for the HEK, CAMP and CIN as they are more selective. Break off swab into Selenite Broth and recap tube loosely. |
| 2 | <p>Streak for isolated growth using a disposable inoculation needle.</p>  <p>Streak out to cover the whole plate.</p> |
| 3 | <ul style="list-style-type: none"> • Incubate BA, MAC, MAC-S and HEK in 35° O₂ incubator. • Incubate CAMP, in campy jar with microaerophilic pouch, in 42° incubator for 72 hours. • Incubate CIN plate at room temperature for 48 hours. • Incubate selenite broth in O₂ incubator before leaving at 20:00. |
| 4 | <ul style="list-style-type: none"> • Subculture selenite broth first thing the next morning. Saturate a cotton-tipped swab from the top of the selenite broth and inoculate the Hektoen plate reserved for selenite subculture. • Streak for isolated growth using a disposable inoculation needle and incubate in 35° O₂ incubator for 18-24 hours. |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

INTERPRETATION OF RESULTS:

| Step | Action | |
|----------------------------------|---|---|
| Examination of Blood Agar | | |
| 1 | Predominate or pure growth <i>P.aeruginosa</i> , <i>S.aureus</i> or <i>Yeast species</i> : <ul style="list-style-type: none"> • Identify with GNI, WP, TC and RS as appropriate • Yeast species does not have to be identified | |
| 2 | Perform sweep oxidase to screen for <i>Aeromonas spp.</i> , which are oxidase positive and may be hemolytic. Note: If <i>Plesiomonas shigelloides</i> or <i>Vibrio species</i> was requested by ordering physician, these organisms are oxidase positive as well and will need to be investigated. | |
| 3 | If: | Then: |
| | Sweep oxidase negative | Discard plate, no <i>Aeromonas</i> , <i>Plesiomonas</i> or <i>Vibrio</i> isolated |
| | Sweep oxidase positive | Pick one representative colony of each morph type of oxidase positive colonies. Subculture each colony picked to a BA purity plate. |

| Oxidase Interpretation | |
|-------------------------------|--|
| Oxidase | Next Steps: |
| Negative | Discard purity plate. Go back to original BA plate and isolate oxidase positive colony again. |
| Positive | Vitek GNI to rule out <i>Aeromonas</i> , <i>Plesiomonas shigelloides</i> or <i>Vibrio species</i> <ul style="list-style-type: none"> • Perform GNS as per ASTM guidelines |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| Step | Action | |
|--------------------------------------|--|--|
| Examination of MacConkey Agar | | |
| 1 | Examine plate for Non-lactose fermenting colonies. | |
| | Record observations in LIS. | |
| | If no Non-lactose fermenting colonies present, record "All Lactose Fermenters" in LIS. | |
| 2 | If: | Then: |
| | Lactose fermenters | Discard plate, no enteric pathogens isolated. |
| | Non-lactose fermenters | Pick one representative colony of each morph type of colorless colonies. Subculture each colony picked to a TSI, urea slant and BA purity plate. |

| TSI/Urea Interpretation for MacConkey Growth | | | |
|--|----------------------|-------------|---|
| TSI | Urea | Other Tests | Next steps: |
| A/A Gas or H ₂ S+ | Positive | | Discard |
| A/A Gas and H ₂ S - | Positive | | Vitek GNI to rule out <i>Yersinia enterocolitica</i> Perform GNS as per ASTM guidelines |
| A/A Gas and H ₂ S - | Negative | Oxidase - | Discard |
| A/A Gas and H ₂ S - | Negative | Oxidase + | Vitek GNI to rule out <i>Vibrio species</i> Perform GNS as per ASTM guidelines |
| A/A Gas or H ₂ S+ | Negative | | Discard |
| K/K Gas or H ₂ S +/- | Positive or Negative | | Discard |
| K/A Gas or H ₂ S +/- | Positive | | Discard |
| K/A Gas and H ₂ S - | Negative | Oxidase + | Vitek GNI to rule out <i>Aeromonas</i> and <i>Plesiomonas shigelloides</i> Perform GNS as per ASTM guidelines |
| K/A Gas and H ₂ S - | Negative | Oxidase - | Vitek GN to rule out <i>Salmonella</i> and <i>Shigella</i> <ul style="list-style-type: none"> If identification is <i>Salmonella</i> species, perform <i>Salmonella</i> latex serology testing from BA purity plate Perform GNS as per ASTM guidelines |
| K/A Gas or H ₂ S + | Negative | | Vitek GN to rule out <i>Salmonella</i> <ul style="list-style-type: none"> If identification is <i>Salmonella</i> species, perform <i>Salmonella</i> latex serology testing from BA purity plate Perform GNS as per ASTM guidelines |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| Step | Action | |
|------------------------------------|---|---|
| Examination of Hektoen Agar | | |
| 1 | Examine plate for green (with or without black centers) or blue colonies. Record observations in LIS. If no green (with or without black centers) or blue colonies present, record "All yellow-orange or salmon pink" in LIS. | |
| 2 | If: | Then: |
| | Yellow-orange or salmon pink | Discard plate, no enteric pathogens isolated. |
| | Blue or green colonies with or without black centers | Pick one representative colony of each morph type of suspicious colonies. Subculture each colony picked to a TSI, urea slant and BA purity plate. |

| TSI/Urea Interpretation for Hektoen Growth | | | | |
|--|----------|-------------|---|--|
| TSI | Urea | Next steps: | | |
| Any reaction | Positive | Discard | | |
| K/K Gas or H ₂ S +/- | Negative | Discard | | |
| A/A Gas or H ₂ S +/- | Negative | Discard | | |
| K/A Gas and H ₂ S - | Negative | Vitek GNI | <u>Shigella:</u> <ul style="list-style-type: none"> Vitek 2 AST-N213 | <ul style="list-style-type: none"> Refer to ASTM for susceptibility reporting |
| K/A Gas or H ₂ S + | Negative | Vitek GNI | <u>Salmonella:</u> <ul style="list-style-type: none"> Perform Salmonella latex serology testing from BA. Vitek 2 AST-N213 | |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| Step | Action | |
|---|--|--|
| Examination of MacConkey Sorbitol Agar | | |
| 1 | Examine plate for Non-sorbitol fermenters. Record observations in LIS. If no non-sorbitol fermenters present, record "All Sorbitol Fermenters" in LIS. | |
| 2 | If: | Then: |
| | Sorbitol fermenters | Discard plate, no enteric pathogens isolated. |
| | Non-sorbitol fermenters | Pick one representative colony of each morph type of colorless colonies. Subculture colony to a BA purity plate From BA purity plate, perform indole and oxidase. |

| Indole/Oxidase Interpretation for Sorbitol MacConkey Growth | | |
|--|---|---|
| Indole | Oxidase | Next steps: |
| Negative | Any reaction | Discard |
| Positive | Positive | Discard |
| Positive | Negative | Perform E.coli 0157 latex agglutination test |
| | | <table border="0"> <tr> <td style="vertical-align: top;"> If Negative: <ul style="list-style-type: none"> Discard </td> <td style="vertical-align: top;"> If Positive: <ul style="list-style-type: none"> Perform GNI for identification to confirm E.coli Susceptibility testing is not routinely performed on stool isolates of <i>E.coli 0157:H7</i> </td> </tr> </table> |
| If Negative: <ul style="list-style-type: none"> Discard | If Positive: <ul style="list-style-type: none"> Perform GNI for identification to confirm E.coli Susceptibility testing is not routinely performed on stool isolates of <i>E.coli 0157:H7</i> | |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| Step | Action | |
|-------------------------------------|--|---|
| Examination of Yersinia Agar | | |
| 1 | Examine plate for Bull's Eye colonies. Record observations in LIS. If no bull's eye colonies present, record "No Bull's Eye Colonies at 48hrs" or if there is no growth present, record "No Growth at 48h" in LIS. | |
| 2 | IF | THEN |
| | No growth or no bulls 'eye colonies at 48hrs | Discard plate, <i>Yersinia</i> not isolated. |
| | Bull's eye colonies | Pick one representative colony of each morph type of bull's eye colonies. Subculture each colony picked to a TSI, urea and BA purity plate. |

| TSI/Urea Interpretation for CIN Growth | | |
|--|--------------|--|
| TSI | UREA | Next steps: |
| K/K Gas or H2S+/- | Any reaction | Discard |
| K/A Gas or H2S+/- | Any reaction | Discard |
| A/A Gas or H2S+/- | Negative | Discard |
| A/A Gas or H2S+ | Positive | Discard |
| A/A Gas and H2S- | Positive | <ul style="list-style-type: none"> Vitek GNI to rule out <i>Yersinia enterocolitica</i> Perform GNS as per ASTM guidelines |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| Step | Action | |
|--|--|---|
| Examination of Campylobacter Agar | | |
| 1 | Examine plate for grey, flat, irregular, spreading, sometimes mucoid colonies. Record observations in LIS. If no suspicious colonies present, record "No Growth at 72h" or "No Colonies Typical of Campy at 72h" in the LIS. | |
| 2 | If: | Then: |
| | No growth or no grey colonies | Discard plate, <i>Campylobacter</i> not isolated. |
| | Grey colonies | Pick one representative colony of each morph type of grey colonies colonies. Perform gram stain, catalase, oxidase and hippurate. |

| Interpretation for CAMP Growth | | | | |
|---|----------|----------|-----------|---------------------------------|
| Gram | Catalase | Oxidase | Hippurate | Next steps: |
| Small, curved Gram-negative bacilli | Positive | Positive | Positive | Report as <i>C.jejuni</i> |
| Small, curved Gram-negative bacilli | Positive | Positive | Negative | Send to Prov.Lab for further ID |
| Any other gram stain, catalase and oxidase reactions are not <i>Campylobacter</i> and can be discarded. | | | | |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

REPORTING RESULTS:

| Results | Action |
|---|--|
| No growth after 3 days | <ul style="list-style-type: none"> Report "Reduced commensal flora present. No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas species</i> isolated" Report "No <i>Plesiomonas shigelloides</i> or <i>Vibrio species</i> isolated" if requested. |
| Normal enteric flora isolated | <ul style="list-style-type: none"> Report "No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas species</i> isolated" Report "No <i>Plesiomonas shigelloides</i> or <i>Vibrio species</i> isolated" if requested. |
| Overgrowth of <i>S.aureus</i> , <i>P.aeruginosa</i> or Yeast | <ul style="list-style-type: none"> Report "No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas species</i> isolated" Report "Predominant or pure growth of <i>Staphylococcus aureus</i>" or "Predominant or pure growth of <i>Pseudomonas aeruginosa</i>" or "Predominant or pure growth of Yeast species" |
| <i>Aeromonas</i> isolated | <ul style="list-style-type: none"> Report "No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> isolated" Report "<i>Aeromonas species</i>" List quantitation as "Isolated" Report AST as per ASTM. |
| <i>Plesiomonas shigelloides</i> or <i>Vibrio species</i> isolated | <ul style="list-style-type: none"> If organisms are found incidentally, only report it if no other pathogens were isolated. Report "Report "No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas species</i>" Report "<i>Plesiomonas shigelloides</i> or <i>Vibrio species</i>" List quantitation as "Isolated" Report AST as per ASTM. |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| | |
|---------------------------|---|
| <p>Suspected Shigella</p> | <ul style="list-style-type: none"> • Report “No <i>Salmonella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas</i> species isolated” • Report “<i>Shigella group</i>” with isolate comment &SHIG to state: “Isolate has been referred to Alberta Provincial Laboratory for Public Health (Tel: 1 780 407 7121) for confirmation and typing” • List quantitation as “Possible” • Report AST as per ASTM. • Freeze specimen and enter into specimen isolate log. • Go to Order Entry; copy report to Chief Medical Officer of Health (HPU) and Infection Control Nurse (SOHS) if in-patient. • Add test ?REFE. Send to Prov. Lab for confirmation and serotyping. |
| <p>Salmonella</p> | <ul style="list-style-type: none"> • Report “No <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157, <i>Yersinia</i> or <i>Aeromonas</i> species isolated” • Report “<i>Salmonella spp</i>” with isolate comment &SALM to state: “Isolate has been referred to Alberta Provincial Laboratory for Public Health (Tel: 1 780 407 7121) for further identification and typing” • List quantitation as “Isolated” • Report AST as per ASTM. • Freeze specimen and enter into specimen isolate log. • Phone results to ordering location and Chief Medical Officer of Health (HPU1) and document in call log. Results must be phoned to HPU1 within 24 hours of isolation. • Go to Order Entry; copy report to Chief Medical Officer of Health (HPU1) and Infection Control Nurse (SOHS) if in-patient. • Add test ?REFE. Send to Prov. Lab for further identification and serotyping. |

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

| | |
|---------------|---|
| E.coli 0157 | <ul style="list-style-type: none"> • Report “No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>Yersinia</i> or <i>Aeromonas</i> species isolated” • Report “Escherichia coli 0157” with isolate comment &0157 to state: “Isolate has been referred to Alberta Provincial Laboratory for Public Health (Tel: 1 780 407 7121) for confirmation and typing” • List quantitation as “Isolated” • AST not performed on this organism. • Freeze specimen and enter into specimen isolate log. • Go to Order Entry; copy report to Chief Medical Officer of Health (HPU) and Infection Control Nurse (SOHS) if in-patient. • Add test ?REFE. Send to Prov. Lab for confirmation and serotyping. |
| Campylobacter | <ul style="list-style-type: none"> • Report “No <i>Salmonella</i>, <i>Shigella</i>, <i>E.coli 0157</i>, <i>Yersinia</i> or <i>Aeromonas</i> species isolated” • Report organism according to hippurate result. <ul style="list-style-type: none"> • Hippurate positive: Report “C.jejuni” • Hippurate negative: Report “Campylobacter spp”. with isolate comment &CAMP to state: “Isolate has been referred to Alberta Provincial Laboratory for Public Health (Tel: 1 780 407 7121) for further identification” • List quantitation as “Isolated” • Report AST as per ASTM. • Freeze specimen and enter into specimen isolate log. • Go to Order Entry; copy report to Chief Medical Officer of Health (HPU) and Infection Control Nurse (SOHS) if in-patient. • Add test ?REFE. Send to Prov. Lab for further identification if hippurate result is negative. |

| | | |
|-------------------------------------|----------------------------------|----------------|
| Document Name: Stool Culture | Document Number: MIC31900 | |
| | Version No: 1.0 | Page: 13 of 14 |
| | Effective: 23 June, 2017 | |

| | |
|----------|---|
| Yersinia | <ul style="list-style-type: none"> • Report “No <i>Salmonella</i>, <i>Shigella</i>, <i>Campylobacter</i>, <i>E.coli</i> 0157 or <i>Aeromonas</i> species isolated”. • Report “Yersinia spp” with isolate comment &YERS to state: “Isolate has been referred to Alberta Provincial Laboratory for Public Health (Tel: 1 780 407 7121) for confirmation and typing” • List quantitation as “Isolated” • Report AST as per ASTM. • Freeze specimen and enter into specimen isolate log. • Go to Order Entry; copy report to Chief Medical Officer of Health (HPU) and Infection Control Nurse (SOHS) if in-patient. • Add test ?REFE. Send to Prov. Lab for confirmation and serotyping. |
|----------|---|

NOTES AND PRECAUTIONS:

All cultured organisms sent to referral lab must be sent according to current transportation of dangerous goods protocols.

PROCEDURE NOTES:

- Direct examination of fecal cultures is not indicated
- Fecal cultures should not be performed for patients being treated with broad-spectrum antimicrobial agents because it is likely that this therapy is responsible for the diarrhea. These cultures may show overgrowth of other organisms including *Pseudomonas aeruginosa* and *Candida* species – the role of which in disease production is not clear

| | |
|---|------------------------------------|
| NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use. | |
| FILENAME: Stool Culture | Print Date: 23/06/2017 11:59:00 AM |

REFERENCES:

Clinical Microbiology Procedures Handbook, 4th edition, ASM Press, 2016

Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. 2015. Manual of Clinical Microbiology, 11th edition, ASM Press, Washington, D.C.

Northwest Territories Health and Social Services. NWT Reportable Diseases as per 2009 Public Health Act.

http://www.hss.gov.nt.ca/sites/default/files/nwt_communicable_disease_report_form.pdf

REVISION HISTORY:

| REVISION | DATE | Description of Change | REQUESTED BY |
|----------|-------------|-----------------------|--------------|
| 1.0 | 23 Jun 2017 | Initial Release | L. Steven |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

 Invalid signature

X 

Jennifer G. Daley Bernier, R.T. (CSMLS)
A/ Manager, Laboratory Services
Signed by: Jennifer G. Daley Bernier

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.