

Document Name:

DensiCHEK plus

Approved By:

Jennifer G. Daley Bernier, A/ Manager, Laboratory Services

Status: **APPROVED**

PURPOSE:

The DensiCHEK plus measures the organism optical density by two measurements

- Measurement of the optical density of the air before each actual reading, to compensate for variations in ambient parameters.
- Continuous measurement of the optical density of the solution within the tube.

Values are in McFarland units, proportional to bacterial concentrations of organisms isolated from clinical specimens.

This procedure will cover the following topics:



1. Replacing the batteries
2. Selecting plastic or glass test tube setting
3. Zeroing the instrument with saline blank
4. Preparing patient samples
5. Monthly maintenance
6. Cleaning



NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

1. Replacing the Batteries:




The instrument runs on 4 alkaline AAA batteries. The batteries should be replaced as soon as the low battery icon starts to display.

Step	Action
Replacing DensiCHEK plus batteries	
1	<p>Notice the low battery icon at the bottom left of the instrument window.</p> 
2	<p>Obtain 4 <u>alkaline</u> AAA Batteries. These can be ordered from:</p> <p style="text-align: center;">Stanton Territorial Authority Stores Stocked Items Catalogue Number: 9100906 Description: Batteries, AAA-<i>ALKALINE</i> Department Code: 1.71410</p> <p>WARNING: DO NOT USE HIGH ENERGY TYPE BATTERIES SUCH AS LITHIUM OR NiCad</p>
3	<p>Insert batteries respecting the correct polarities indicated inside the battery compartment.</p> 
4	<p>Turn on the instrument and ensure the instrument is set to the proper tube type needed. (refer to topic 2)</p>



2. Selecting a Tube Type (Glass or Plastic):

WARNING:

Vitek 2 demands the use of polystyrene plastic tubes whereas McFarland standards come in glass tubes. Users **MUST** ensure the instrument is set to the appropriate tube type prior to use. Failure to do this can lead to severe patient result errors.



Step	Action
Selecting Tube Type	
1	<p>Turn-on the instrument  and look for a black triangle on the screen to determine the currently selected tube type. The black triangle will be pointing towards the current setting (either glass or plastic)</p> <p>In the example below the instrument is on the plastic setting:</p> 
2	<p>To begin changing the setting press the menu key ONCE.</p>  <p>The screen will display the words SEL and the black triangle will be flashing below the current setting.</p>

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

3	<p>Press the read/enter key to toggle the triangle between settings.</p> 
4	<p>When you have the black triangle below the desired heading press the menu key once.</p> 
5	<p>Check the main screen to ensure the tube type setting is correct.</p>

3. Zeroing the Instrument with Saline

The instrument should be zeroed at the beginning of each set-up run.




Step	Action
Zeroing the DensiCHEK plus	
1	Turn the power on 
2	ENSURE THE TUBE TYPE SETTING IS CORRECT. If the setting is wrong change it as per topic 2 and start again.
3	Choose a plastic test tube that is free from scratches and add 3 mL sterile saline.
4	Insert test tube into instrument.
5	<p>Press the ZERO/SCROLL key and slowly rotate the test tube.</p>  <p>Ensure one full rotation is completed before the reading is displayed. The instrument will display a series of dashes followed by 0.00.</p>
6	Once the instrument has been zeroed it can be used to measure patient suspensions.

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

4. Preparing Patient Samples

Special Safety Precautions:

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

Step	Action
Preparing Patient Samples	
1	<p>Turn the power on . The instrument should be set to PLASTIC.</p>  <p>WARNING: FAILURE TO SET THE TUBE TYPE SETTING CORRECTLY WILL AFFECT INOCULUM DENSITY AND CAN CAUSE INCORRECT ID AND AST CARD RESULTS!!!! (refer to topic 2)</p>
2	<p>Ensure the instrument has been zeroed for the plastic setting. (refer to topic 3)</p>
3	<p>Select a clean plastic test tube free from scratches and ensure that no labeling is placed on the tube within the instruments' reading zone.</p> 
4	<p>Fill the test tube with 3.0mL of sterile saline and inoculate with colonies as per established procedures for making suspensions for Vitek 2.</p>


NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

Document Name: DensiCHEK plus	Document Number: MIC70100	
	Version No: 3.0	Page: 7 of 10
	Effective: 22 November, 2017	

5	<p>With the instrument ON, place tube in instrument and rotate slowly.</p> <p>Ensure one full rotation is completed before the reading is displayed.</p> <p>The instrument will display a series of dashes followed by a reading.</p>
6	Check that the McFarland value is within the acceptable range for card type.
7	<p>Adjust suspension if necessary.</p> <p><i>Note: If the instrument flashes 0.00 or 4.00, the suspension is outside the readable range of the instrument.</i></p>
8	Repeat steps for each new patient suspension.

5. Monthly Maintenance:

The DensiCHEK plus calibration must be verified monthly using 0.0, 0.5, 2.0 and 3.0 McFarland Standards.

Step	Action												
Monthly DensiCHEK plus Calibration													
1	Set the tube type to GLASS.												
2	Gently invert the 0.0 McFarland Standard several time then insert into instrument. DO NOT USE VORTEXER: AIR BUBBLES WILL AFFECT READING.												
3	<div style="text-align: center;">  </div> <p>Zero the instrument using the Zeroing key.</p> <p><i>Note: same steps as zeroing the saline blank but using the McFarland Standard instead.</i></p>												
4	<p>Read each standard by gently inverting to mix, cleaning outside of tube with kim wipe, inserting into instrument and rotating one full turn.</p> <p>Ensure the values obtained are within acceptable range.</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th style="width: 30%;">Standard</th> <th colspan="2" style="width: 70%;">Acceptable Range</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">0.5 McF</td> <td style="text-align: center;">0.44</td> <td style="text-align: center;">0.56</td> </tr> <tr> <td style="text-align: center;">2.0 McF</td> <td style="text-align: center;">1.85</td> <td style="text-align: center;">2.15</td> </tr> <tr> <td style="text-align: center;">3.0 McF</td> <td style="text-align: center;">2.79</td> <td style="text-align: center;">3.21</td> </tr> </tbody> </table>	Standard	Acceptable Range		0.5 McF	0.44	0.56	2.0 McF	1.85	2.15	3.0 McF	2.79	3.21
Standard	Acceptable Range												
0.5 McF	0.44	0.56											
2.0 McF	1.85	2.15											
3.0 McF	2.79	3.21											
5	Enter results onto MIC70110.1 Maintenance Record – Vitek 2. If results are out of range notify the Tech II.												

NOTE: This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

6. Cleaning:

The instrument should be cleaned any time a spill occurs.

Step	Action
DensiCHEK plus Cleaning	
1	Prepare a 10% bleach solution. **DO NOT USE ALCOHOL**
2	Wipe the surface of the DensiCHEK plus with the bleach solution.
3	Remove the adaptor.
4	Use a swab dipped in 10% bleach solution to clean the reading chamber surfaces.
5	Soak the adaptor in the 10% bleach solution.
6	Rinse the adaptor in plain water.
7	Fully dry the adaptor and re-insert into the instrument.
8	Perform calibration check.

References:

- DensiCHEK plus Instrument User Manual; Biomérieux EN 01/2010
- DensCHEK plus Standard Kit Package Insert Ref 21255-P1ML1, Biomérieux

REVISION HISTORY:

REVISION	DATE	Description of Change	REQUESTED BY
1.0	17-Sep-2012	Initial Release	M-L Dufresne
2.0	15-Feb-2017	Update format	L. Steven
3.0	22-Nov-2017	Update format and new instrumentation	L. Steven