

<b>NORTHWEST TERRITORIES</b> <b>Health and Social Services Authority</b>	<b>Laboratory</b> <b>Stanton Territorial Hospital</b> P.O. Box 10, 550 Byrne Road YELLOWKNIFE NT X1A 2N1	<b>Document Number:</b> MIC30800	
		<b>Version No:</b> 2.0	Page: 1 of 7
<b>Document Name:</b> CSF Culture		<b>Distribution:</b> <b>Microbiology Culture Manual</b>	
<b>Approved By:</b> Jennifer G. Daley Bernier, A/Manager, Laboratory Services		<b>Effective:</b> 11 January, 2017 <b>Date Reviewed:</b> 11 January, 2017 <b>Next Review:</b> 11 January, 2019	
		<b>Status:</b> <b>APPROVED</b>	

**PURPOSE:** To determine the presence or absence of bacterial pathogens in CSF specimens

**SAMPLE INFORMATION:**

<b>Special Precautions</b>	Refer to Policy B-0160: Specimens Containing Suspected Risk Group 3 Pathogens for Primary Specimen Handling Flow Chart
<b>Type</b>	CSF collected into clean, sterile, leak-proof centrifuge tubes to be transported to the laboratory immediately.
<b>Source</b>	<ul style="list-style-type: none"> <li>• Central nervous system shunt fluid</li> <li>• Fluid from Ommaya reservoirs</li> <li>• External ventricular drainage fluid</li> <li>• CSF from lumbar puncture</li> </ul>
<b>Volume</b>	Generally, 1-3 mL of CSF is required for the Microbiology Laboratory. If viral, fungal or mycobacterial testing is required, then at least 3-4 mL should be sent.
<b>Stability</b>	Transport to the laboratory immediately
<b>Storage Requirements</b>	If a delay in processing is anticipated, hold specimens at room temperature, do <b>NOT</b> refrigerate.
<b>Criteria for rejection and follow up action</b>	<ul style="list-style-type: none"> <li>• Insufficient volume for tests requested: contact the physician to prioritize requests.</li> <li>• Leaking specimens should be processed, but alert the physician of the possibility of contamination.</li> <li>• Improperly collected, labeled, transported or handled specimens should be processed. Waiver of responsibility form SCM40110 needs to be filled out by the responsible nurse.</li> </ul>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

Document Name: CSF Culture	Document Number: MIC30800	
	Version No: 2.0	Page: 2 of 7
	Effective: 11 January, 2017	

**REAGENTS and/or MEDIA:**

- Blood agar (BAP), Chocolate agar (CHOC) and MacConkey agar (MAC)
- Identification reagents: catalase, oxidase, rapid Staph, rapid Strep, etc.

**SUPPLIES:**

- Wooden sticks
- Disposable inoculation needles
- Microscope slides
- Biosafety cabinet
- 35° ambient air and 37° CO<sub>2</sub> incubators
- Vitek 2 and supplies

**SPECIAL SAFETY PRECAUTIONS:**

Containment Level 2 facilities, equipment, and operational practices for work involving infectious or potentially infectious materials or cultures.

- Lab gown must be worn when performing activities with potential pathogens.
- Gloves must be worn when direct skin contact with infected materials is unavoidable.
- Eye protection must be used where there is a known or potential risk of exposure to splashes.
- All procedures that may produce aerosols, or involve high concentrations or large volumes should be conducted in a biological safety cabinet (BSC).
- The use of needles, syringes, and other sharp objects should be strictly limited.

All patient specimens are assumed to be potentially infectious. Universal precautions must be followed. Since viable micro-organisms are used, all cultures must be handled with appropriate precautions. All equipment in contact with cultures should be decontaminated by appropriate methods.

**QUALITY CONTROL:**

Refer to Quality Control manual for reagent quality control procedures.

<b>NOTE:</b> This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.	
FILENAME: MIC30800CSFCulture	Print Date: 6/25/2018 1:23:00 PM

**COMMON BACTERIAL CAUSES OF ACUTE MENINGITIS:**


Age Group	Organism(s)
Neonate	<i>E.coli, S.agalactiae, Listeria monocytogenes</i>
Infants and Children	<i>S.agalactiae, H.influenzae, S.pneumoniae, N.meningitidis</i>
Adolescents and Young Adults	<i>S.pneumoniae, N.meningitidis</i>
Older Adults	<i>N.meningitidis, S.pneumoniae, H.influenzae, S.agalactiae, L.monocytogenes</i>
Ventriculoperitoneal Shunt Infections	<i>CNS, S.aureus, Streptococcus spp., Gram-negative bacilli, Candida albicans, Corynebacterium spp, Propionibacterium acnes</i>

**PROCEDURE INSTRUCTIONS:**

Step	Action
<b>Processing CSF Tube #4</b>	
1	Do not centrifuge, regardless of amount received.
2	Using aseptic technique, transfer entire specimen into a labelled, sterile red top tube. Para film for transport.
3	Forward specimen to the Provincial Laboratory for any requested viral tests as per Prov. Lab protocol. If transport time will be > 24h freeze at -70 for transport on dry ice.
4	If no viral testing is ordered, then freeze specimen in the -70 freezer.

Step	Action
<b>Processing CSF Tube #2</b>	
1	If >1 mL received <ul style="list-style-type: none"> <li>Centrifuge tube #2 at 3500 rpm for 10 minutes.</li> <li>Aseptically transfer the supernatant to a labeled sterile glass tube and set aside.</li> <li>Label 1 sterile ringed slide and plates with LIS plate labels.</li> </ul>
	If < 1 mL received <ul style="list-style-type: none"> <li>Do not centrifuge</li> <li>Label 1 sterile ringed slide and plates with LIS plate labels.</li> </ul>

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

<b>2</b>	<p>In the biosafety cabinet, using a sterile pipette:</p> <ul style="list-style-type: none"> <li>Aspirate fluid from the bottom of the collection tube.</li> <li>Place 1 - 2 drops each onto BAP, CHOC and MAC. Streak for isolated growth using a disposable inoculation needle. Streak out to cover the whole plate.</li> </ul>  <ul style="list-style-type: none"> <li>Prepare smear by placing 1 or 2 drops of CSF on a microscope slide. Allow the drop(s) to form one large drop. Do not spread the fluid.</li> <li>Shunt fluids should be also planted to THIO broth and held for 14 days.</li> </ul>
<b>3</b>	Place specimen sediment tube, supernatant tube and MAC in the O <sub>2</sub> incubator. Place BAP and CHOC plates in the CO <sub>2</sub> incubator in the designated tray.
<b>4</b>	Allow smear to dry and perform Gram Stain. Gram stain must be read before culture plates. Refer to MIC20115 – Gram Stain Procedure.
<b>5</b>	Interpret CSF stains immediately. During the regular Microbiology lab hours of 07:00 to 20:00, turnaround time for these gram stains is <1 hour. Outside the regular Microbiology lab hours, Microbiology Technologist may be called in if ordering physician determines the stain must be read immediately.
<b>6</b>	Immediately phone results of any positive stain results for microorganisms to ordering location and document the conversation within the LIS. Positive gram stain results need to be copied to Chief Medical Officer of Health (HPU1) and Infection Control (SOHS) if inpatient.
<b>7</b>	Send a pending report for all gram stains.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

**INTERPRETATION OF RESULTS:**

Step	Action
<b>Interpretation of CSF specimens</b>	
1	Examine plates after 24 hour incubation. Record observations in the LIS.
2	Re-incubate CO <sub>2</sub> plate(s) for an additional 48 hours. Discard O <sub>2</sub> plate.
3	If the Gram-stained smear is positive and there is no growth on the plates, re-incubate CO <sub>2</sub> plates for total of 7 days.
4	At 48 hours, record observations in the LIS for re-incubated plates.
5	At 72 hours, record observations in the LIS for re-incubated plates.
6	Identify and perform susceptibility testing as per DynaLIFE ASTM.
7	Freeze organism in glycerol and record in patient isolate log.
8	A copy of all positive reports must be sent to Chief Medical Officer of Health (HPU1).
9	A copy of all positive reports on inpatients must be sent to Infection Control (SOHS).
10	Any positive CSF for Group A Streptococcus or Streptococcus pneumoniae must be sent to Streptococcus Unit at NML Winnipeg for surveillance testing.
11	Any positive CSF for Haemophilus influenza or Neisseria meningitis must be sent immediately to the Provincial Lab Edmonton for typing as soon as identification is confirmed. Assure there is a purity plate made that can be used for this purpose and can be sent out the day the identification is confirmed. Provincial Lab will then forward the specimen to NML for surveillance testing.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

**REPORTING OF RESULTS:**

IF:	ACTION:
No growth after 1 day	<p><b>PRELIM:</b></p> <ul style="list-style-type: none"> <li>Report: <b>“No Growth After 1 Day. Further report to follow”</b></li> </ul>
No growth after 3 days	<p><b>FINAL:</b></p> <ul style="list-style-type: none"> <li>Report: <b>“No growth aerobically after 3 days”</b></li> </ul>
Growth of organism	<ul style="list-style-type: none"> <li>Report organism identification under the isolates tab.</li> <li>List quantitation.</li> <li>Report susceptibility results as per DynaLIFE ASTM.</li> <li>Refer to Reportable Diseases – Public Health Act as of September 2009 for reporting to HPU1.</li> <li>Refer to MIC35100 – Nosocomial Infection Notification Job Aid to determine if organism needs to be sent to Infection Control.</li> <li>Refer to L-0910-Laboratory: Critical Values for results that need to be phoned to ordering location.</li> </ul>

**PROCEDURE NOTES:**

1. A positive culture generally indicates infection with the organism.
2. Lack of pus cells in CSF does not rule out infection, especially in Listeriosis.
3. The most common cause of community acquired bacterial meningitis is Streptococcus pneumoniae.
4. Direct bacterial antigen testing is not recommended.
5. Since THIO is mainly a broth for anaerobes and does not support the growth of the most common pathogens in CSF, it is not recommended for routine CSF culture but should be used when Shunt Infection is suspected.

**NOTE:** This is a controlled document for internal use only. Any documents appearing in paper form are not controlled and should be checked against electronic version prior to use.

**REFERENCES:**

- Clinical Microbiology Procedures Handbook, 4<sup>th</sup> edition, ASM Press, 2016
- Jorgensen J.H., Pfaller M.A., Carroll K.C., Funke G., Landry M.L., Richter S.S., Warnock D.W. 2015. Manual of Clinical Microbiology, 11<sup>th</sup> edition, ASM Press, Washington, D.C.

**REVISION HISTORY:**

REVISION	DATE	Description of Change	REQUESTED BY
1.0	11 Jan 2017	Initial Release	L. Steven
2.0	24 Apr 2018	Change to reflect new Vitek 2 instrument	L. Steven